

MODEL 1700 LIQUID LEVEL INSTRUMENT

(HELIUM & NITROGEN VERSION)

**INSTALLATION, OPERATION, AND
MAINTENANCE INSTRUCTIONS**

American Magnetics, Inc.

P.O. Box 2509, 112 Flint Road, Oak Ridge, TN 37831-2509, Tel: 865-482-1056, Fax: 865-482-5472

Table of Contents

Cover Page.....	i
Table of Contents	i
List of Figures.....	v
List of Tables	ix
Foreword	xi
Purpose and Scope	xi
Contents of this Manual	xi
General Precautions	xii
Cryogen Safety Summary	xiii
Safety Legend	xiv
Equipment Warnings.....	xv
Other Manual Conventions	xv
Instrument Configuration	xv
Introduction	1
1.1 Model 1700 Instrument.....	1
1.1.1 Cryogenic (Capacitance-Based) Liquid Level Sensors	1
1.1.2 Superconducting Level Sensors	3
1.1.3 Digitally-Controlled	3
1.1.4 System Flexibility	3
1.1.5 Display	4
1.1.6 Intuitive Human-Interface Design	4
1.1.7 Measurement Flexibility	4
1.1.8 Real Time Clock	4
1.1.9 Valve Control Output.....	4
1.1.10 Analog Outputs	4
1.1.11 Signal Relays.....	4
1.1.12 Connectivity	5
1.2 Model 1700 Front Panel Layout	6
1.3 Model 1700 Rear Panel Layout	6
1.4 Model 1700 Specifications @ 25°C.....	7
Installation	11
2.1 Unpacking and Inspecting the Instrument	11

Table of Contents

2.2	Mounting the Model 1700 Instrument	11
2.3	Installing Capacitance-Based Liquid Level Sensors.....	11
2.4	Installing the Liquid Helium Level Sensor	12
2.5	Connecting the Sensor(s) to the Instrument	14
2.5.1	Connecting a Capacitance Sensor	14
2.5.2	Connecting a Liquid Helium Level Sensor.....	15
2.6	Setting Up an Autofill System	16
2.6.1	Autofill System Description.....	16
2.6.2	Autofill System Setup	18
2.7	Power Requirements	19
Operation		21
3.1	Energizing the Model 1700 Instrument.....	21
3.2	Screen Navigation	22
3.2.1	Home Screen Footer.....	22
3.3	Navigating the Instrument Menus	22
3.3.1	Menu Overview.....	23
3.3.2	Editing a Field.....	29
3.3.3	Menu Navigation.....	30
3.4	Capacitance (Liquid nitrogen) Level	30
3.4.1	Configure the instrument to display nitrogen level.....	30
3.5	Superconducting (Liquid Helium) Level	31
3.5.1	Configure the instrument to display helium level	31
3.5.2	Sample and Hold Operation	34
3.5.3	Other Liquid Helium Functions	34
3.5.4	Other Liquid Helium Functions	34
3.6	Alarms and Relays	35
3.6.1	Overview	35
3.6.1.1	Level-Based Alarms.....	35
3.6.1.2	Time-Based Fill Alarm (nitrogen only)	35
3.6.1.3	Multiple alarms	36
3.6.2	Configuring Alarm Setpoints	36
3.6.3	Acknowledging an Alarm	37
3.6.4	Muting an Alarm	38
3.6.5	Resetting the Autofill Timeout Alarm	38
3.7	Configure the Autofill Function.....	39
3.7.1	Autofill overview	39
3.7.2	Setting the autofill parameters	39
3.7.3	Enable the Autofill function.....	40
3.8	Select the appropriate units on the display.....	40

Table of Contents

3.9	Analog output signals.....	41
3.9.1	Connecting to the Aux Connector.....	41
3.9.2	Configuring the Analog Outputs.....	41
3.10	Ethernet Connectivity	41
3.10.1	IP Addressing Scheme	41
3.11	Abnormal Operation	42
3.11.1	Capacitance Sensor Contamination	42
3.11.2	Dirty Helium Sensor Operational Mode	43
3.11.3	Resetting the Instrument to Factory Defaults	43
3.12	Shutting the Instrument Down	44
Calibration.....		45
4.1	Setting the System Date and Time	45
4.2	Capacitance-based Level Calibration.....	46
4.2.1	Understanding the Sensor Active Length	46
4.2.2	Relationship between Calibration and Sensor Length	47
4.2.3	Variations in the Dielectric with Changing Density	48
4.2.4	Capacitance-based Sensor Calibration Methods.....	49
4.2.4.1	Selection of Capacitance Sensor Calibration Methods.....	50
4.2.5	Pre-Calibration Procedure.....	52
4.2.5.1	Enter Capacitance Sensor Information	52
4.2.5.2	Performing Loss of Sensor Calibration	54
4.2.6	Performing an Open Dewar Calibration	56
4.2.7	Approximate Calibration.....	57
4.2.8	Closed Dewar Calibration.....	61
4.2.8.1	Presetting the maximum and minimum calibration points	61
4.2.8.2	Completing the closed dewar calibration procedure.....	63
4.3	Superconductivity-Based (Liquid Helium) Level Calibration.....	67
4.3.1	Verify the Liquid Helium Sensor Type.....	67
4.3.2	Sensor Sample Interval	68
4.3.3	Sensor Active Length.....	69
4.3.4	Continuous Measure Time Limit	69
4.3.5	Sensor Name	70
Remote Interface Reference		71
5.1	SCPI Command Summary	71
5.2	SCPI Ethernet Communication	71
5.3	SCPI Serial (RS-232) Communication	72
5.3.1	Serial port connector and cabling.....	72
5.4	Command/return termination characters.....	72
5.5	Command Set Reference.....	73

Table of Contents

5.5.1	Commands for determining the instrument configuration	73
5.5.2	Commands for setting the units of measurement	78
5.5.3	Commands for configuring setpoints	79
5.5.4	Commands for setting the channel identifiers	82
5.5.5	Commands for making liquid level measurements	83
5.5.6	Commands for calibrating level sensors	85
5.5.7	System Commands	86
5.6	Error Codes	87
Service and Repair		89
6.1	Cleaning	89
6.2	User Replaceable Parts	89
6.3	Battery Replacement	90
6.3.1	Tools Required	90
6.3.2	Procedure	90
6.4	Fuse Replacement	91
6.4.1	Tools Required	91
6.4.2	Procedure	91
Appendix		93
A.1	Connector Wiring	93
A.1.1	Serial (RS-232) Connector	93
A.1.2	Ethernet Connector	94
A.1.3	Liquid Helium Connector J1 Wiring	95
A.1.4	Aux I/O Connector	96
A.2	Troubleshooting	97
A.3	Firmware Licenses	97
Glossary		105
Abbreviations and Acronyms Used in This Manual		105
Index		107

List of Figures

Figure 1-1	Model 1700 Instrument using internal oscillator/transmitter	2
Figure 1-2	Model 1700 Instrument using external oscillator/transmitter	3
Figure 1-3	Model 1700 (Dual Display) Via Web Browser	5
Figure 1-4	Model 1700 Front View; Dual Channel Instrument Shown ...	6
Figure 2-1	Model 1700 Instrument in an Autofill configuration	17
Figure 2-1	Typical Autofill Setup	17
Figure 2-2	Loss Of nitrogen Level Sensor	19
Figure 3-1	LN2 Home Screen	21
Figure 3-2	LHe Home Screen.....	21
Figure 3-3	Both LN2 and LHe Home Screen	22
Figure 3-4	Home Screen Footer	22
Figure 3-5	MENU button on HOME screen.....	22
Figure 3-6	MENU Selection Screen.....	23
Figure 3-7	HOME button on MENU screen.....	23
Figure 3-8	BACK button	23
Figure 3-9	Field Editing Footer	29
Figure 3-10	Model 1700 Menu Structure	30
Figure 3-11	INTERNAL Oscillator Indicated	31
Figure 3-12	Nitrogen Level Displayed on Home Screen	31
Figure 3-13	Type of LHe Level Sensor	32
Figure 3-14	Helium Level On Home Screen	32
Figure 3-15	Helium Level On Home Screen	32
Figure 3-16	Calibrate Helium Screen	32
Figure 3-17	Numeric Keypad.....	33
Figure 3-18	Calibrate Helium Screen	33
Figure 3-19	Setting the Sample Interval	33
Figure 3-20	Helium Level Measurement Mode	34
Figure 3-21	Helium Sensor Voltage	34
Figure 3-22	Dirty Sensor Mode.....	35
Figure 3-23	Fill Timeout Alarm Screen	36
Figure 3-24	Output Configuration Screen, Page 1	36
Figure 3-25	Relay Source Configuration	37
Figure 3-26	Relay Setpoints.....	37
Figure 3-27	Alarm Condition Above or Below Setpoints.....	37
Figure 3-28	Alarm Annunciator	38
Figure 3-29	LO LEVEL Alarm Condition Footer Displays	38
Figure 3-30	Muted Alarm Condition	38
Figure 3-31	Autofill Level Control Selection.....	39
Figure 3-32	Autofill Level Start Setting	39
Figure 3-33	Autofill Timeout Setting	39
Figure 3-34	LN2 Home Screen	40
Figure 3-35	Analog Outputs Source Selection	41

List of Figures

Figure 3-36 Network Selection	41
Figure 3-37 Editing Network Selection.....	41
Figure 3-38 Editing Network Settings	42
Figure 3-39 Saving Network Settings	42
Figure 3-40 System Menu, Page 2.....	43
Figure 3-41 Resetting Instrument to Factory Defaults	43
Figure 3-42 Invoking Instrument Shut Down	44
Figure 3-43 Confirming Instrument Shut Down.....	44
Figure 4-1 Menu Selection From Home Screen.....	45
Figure 4-2 Menu Selection Screen.....	45
Figure 4-3 System Menu, Page 1.....	45
Figure 4-4 System Menu, Page 2.....	46
Figure 4-5 Typical Capacitance-based Liquid Level Sensor	47
Figure 4-6 Dielectric vs. pressure for nitrogen under saturated conditions. 49	
Figure 4-7 Calibration method selection diagram.....	51
Figure 4-8 MENU Selection Button	52
Figure 4-9 SENSORS Selection Button	52
Figure 4-10 NITROGEN OSCILLATOR Selection.....	52
Figure 4-11 CALIBRATE NITROGEN Selection Button.....	53
Figure 4-12 SENSOR ACTIVE LENGTH field	53
Figure 4-13 Footer BACK Button Selection.....	53
Figure 4-14 SYSTEM Menu Selection	53
Figure 4-15 PERFORM MAX CAL. Selection Button	54
Figure 4-16 Home Selection Button	54
Figure 4-17 MENU Selection Button	54
Figure 4-18 SENSORS Selection Button	55
Figure 4-19 CALIBRATE NITROGEN Selection Button.....	55
Figure 4-20 NO SENSOR CAL button.....	55
Figure 4-21 MENU Selection Button	56
Figure 4-22 SENSORS Selection Button	56
Figure 4-23 CALIBRATE NITROGEN Selection Button.....	56
Figure 4-24 PERFORM MAX CAL. Selection Button	56
Figure 4-25 Updated MAX CAL Frequency.....	57
Figure 4-26 PERFORM MIN CAL. Selection Button.....	57
Figure 4-27 MENU Selection Button	58
Figure 4-28 SENSORS Selection Button	58
Figure 4-29 CALIBRATE NITROGEN Selection Button.....	58
Figure 4-30 PERFORM MAX CAL. Selection Button	58
Figure 4-31 Updated MAX CAL Frequency.....	59
Figure 4-32 PERFORM MIN CAL. Selection Button.....	59
Figure 4-33 APPROX. CAL. VALUE field	60

List of Figures

Figure 4-34 MAX CAL values edited.....	60
Figure 4-35 MENU Selection Button	62
Figure 4-36 SENSORS Selection Button	62
Figure 4-37 CALIBRATE NITROGEN Selection Button.....	62
Figure 4-38 PERFORM MIN CAL. Selection Button	62
Figure 4-39 PERFORM MAX CAL. Selection Button	63
Figure 4-40 PERFORM MAX CAL. Selection Button	63
Figure 4-41 PERFORM MAX CAL. Selection Button	64
Figure 4-42 PERFORM MAX CAL. Selection Button	65
Figure 4-43 Updated MAX CAL Frequency.....	66
Figure 4-44 Home Screen.....	67
Figure 4-45 MENU Selection Button	67
Figure 4-46 SENSORS Selection Button	67
Figure 4-47 Home Screen.....	68
Figure 4-48 MENU Selection Button	68
Figure 4-49 SENSORS Selection Button	68
Figure 4-50 Enter Liquid Helium Sensor Active Length	68
Figure 4-51 Enter Liquid Helium Sensor Active Length	69
Figure 4-52 Enter Liquid Helium Sensor Active Length	69
Figure 4-53 Footer BACK Button Selection.....	69
Figure 4-54 SYSTEM Menu Selection.....	70
Figure 4-55 HELIUM SENSOR NAME: Field.....	70
Figure 4-56 Home Selection Button	70
Figure A-1 Serial (RS-232) Pin Out	93
Figure A-2 Ethernet Connector Socket Pin out.....	94
Figure A-3 Liquid Helium (LHe) Level Sensor Connector	95
Figure A-4 Aux I/O Connector	96

List of Figures

List of Tables

Table 1-1	Model 1700 Rear Panel Description.....	6
Table 2-1	LHe Level Sensor Wire Identification.....	13
Table 2-2	Minimum recommended wire gauge for copper lead wire	16
Table 2-3	Standard Autofill Setup Description.....	17
Table 3-1	Model 1700 Instrument Home Screen Footer.....	22
Table 3-2	Model 1700 Dual Level Instrument Menu Structure.....	24
Table 3-3	Model 1700 Instrument Footer during editing a field.....	29
Table 3-4	Autofill Settings	40
Table A-1	Serial (RS-232) Pin Definitions	93
Table A-2	Ethernet RJ-45 Connector Pin Definitions	94
Table A-3	LHe Level Connector Pin Definitions	95
Table A-4	Aux I/O Pin Definitions.....	96

List of Tables

Foreword

Purpose and Scope

This manual contains the operation and maintenance instructions for the American Magnetics, Inc. Model 1700 Liquid Level Control Instrument and outlines applications for various system configurations. Since it is not possible to cover all equipment combinations for all magnet systems, only the most common configurations are discussed. The user is encouraged to contact an authorized AMI Technical Support Representative for information regarding specific configurations not explicitly covered in this manual. This manual refers to the instrument as configured for both liquid nitrogen and liquid helium operation.

Contents of this Manual

Introduction describes the functions, specifications, and characteristics of the Model 1700 Instrument. It provides illustrations of the front and rear panel layouts as well as documenting the performance specifications. Additional information is provided in the form of system block diagrams.

Installation describes how the Model 1700 Instrument is unpacked and installed in conjunction with ancillary equipment in typical systems. Block-level diagrams document the interconnects for various system configurations.

Operation describes how the Model 1700 Instrument is used to monitor and automatically control liquid nitrogen and liquid helium levels.

Calibration describes the various calibration techniques for liquid helium and liquid nitrogen level sensors.

Remote Interface Reference documents all remote commands and queries available through the Model 1700 Instrument RS-232 and Ethernet interfaces. A quick-reference summary of commands is provided as well as a detailed description of each.

The ***Appendix*** and ***Glossary*** sections support the information in the sections listed above. See the Appendix section when referenced from other sections. See the Glossary for any words or acronyms presented in the above sections, requiring a more complete understanding.

Foreword

General Precautions

General Precautions

Cryogen Safety

The two most common cryogenic liquids used in superconducting magnet systems are nitrogen and helium. Both of these cryogens are extremely cold at atmospheric pressure (-321°F and -452°F , respectively). The following paragraphs outline safe handling precautions for these liquids.

Personnel handling cryogenic liquids should be thoroughly instructed and trained as to the nature of the liquids. Training is essential to minimize accidental spilling. Due to the low temperature of these materials, a cryogen spilled on many objects or surfaces may damage the surface or cause the object to shatter, often in an explosive manner.

Inert gases released into a confined or inadequately ventilated space can displace sufficient oxygen to make the local atmosphere incapable of sustaining life. Liquefied gases are potentially extreme suffocation hazards since a small amount of liquid will vaporize and yield a very large volume of oxygen-displacing gas. Always ensure the location where the cryogen is used is well ventilated. Breathing air with insufficient oxygen content may cause unconsciousness without warning. If a space is suspect, purge the space completely with air and test before entry. If this is not possible, wear a forced-air respirator and enter only with a co-worker standing by wearing a forced-air respirator.

Cryogenic liquids, due to their extremely low temperatures, will also burn the skin in a similar manner as would hot liquids. Never permit cryogenic liquids to come into contact with the skin or allow liquid nitrogen to soak clothing. Serious burns may result from careless handling. Never touch uninsulated pipes or vessels containing cryogenic liquids. Flesh will stick to extremely cold materials. Even nonmetallic materials are dangerous to touch at low temperatures. The vapors expelled during the venting process are sufficiently cold to burn flesh or freeze optic tissues. Insulated gloves should be used to prevent frost-bite when operating valves on cryogenic tanks. Be cautious with valves on cryogenic systems; the temperature extremes they are typically subjected to cause seals to fail frequently.

In the event a person is burned by a cryogen or material cooled to cryogenic temperatures, the following first aid treatment should be given pending the arrival and treatment of a physician or other medical care worker:

1. If any cryogenic liquid contacts the skin or eyes, immediately flush the affected area gently with tepid water (102°F – 105°F , 38.9°C – 40.5°C) and then apply cold compresses.

2. Do not apply heat. Loosen any clothing that may restrict circulation. Apply a sterile protective dressing to the affected area.
3. If the skin is blistered or there is any chance that the eyes have been affected, get the patient immediately to a physician for treatment.

Containers of cryogenic liquids are self pressurizing (as the liquid boils off, vapor pressure increases). Hoses or lines used to transfer these liquids should never be sealed at both ends (i.e. by closing valves at both ends).

When pouring cryogenic liquids from one container to another, the receiving container should be cooled gradually to prevent damage by thermal shock. The liquid should be poured slowly to avoid spattering due to rapid boil off. The receiving vessel should be vented during the transfer.

Introduction of a substance at or near room temperature into a cryogenic liquid should be done with great caution. There may be a violent gas boil-off and a considerable amount of splashing as a result of this rapid boiling. There is also a chance that the material may crack or catastrophically fail due to forces caused by large differences in thermal contraction of different regions of the material. Personnel engaged in this type of activity should be instructed concerning this hazard and should always wear a full face shield and protective clothing. If severe spraying or splashing could occur, safety glasses or chemical goggles along with body length protective aprons will provide additional protection.

The properties of many materials at extremely low temperatures may be quite different from the properties that these same materials exhibit at room temperatures. Exercise extreme care when handling materials cooled to cryogenic temperatures until the properties of these materials under these conditions are known.

Metals to be used for use in cryogenic equipment application must possess sufficient physical properties at these low temperatures. Since ordinary carbon steels, and to somewhat a lesser extent, alloy steels, lose much of their ductility at low temperatures, they are considered unsatisfactory and sometimes unsafe for these applications. The austenitic Ni-Cr alloys exhibit good ductility at these low temperatures and the most widely used is 18-8 stainless steel. Copper, Monel[®], brass and aluminum are also considered satisfactory materials for cryogenic service.

Cryogen Safety Summary

Cryogenic systems are complex systems with the potential to seriously injure personnel or equipment if not operated according to procedures. The

Foreword

Safety Summary

use of safety mechanisms (pressure relief valves, rupture disks, etc.) in cryogenic systems is usually necessary.

Recommended Safety Equipment

The use of proper safety equipment is necessary. Such equipment may include, but not limited to, the following items:

- First Aid kit
- Fire extinguisher rated for class C fires
- Cryogenic gloves
- Face shield
- Signs to indicate that there are potentially dangerous cryogens in use in the area.

Safety Legend



Instruction manual symbol: the product is marked with this symbol when it is necessary for you to refer to the instruction manual in order to protect against damage to the product or personal injury.



Hazardous voltage symbol.



Alternating Current (Refer to IEC 417, No. 5032).



Off (Supply) (Refer to IEC 417, No. 5008).



On (Supply) (Refer to IEC 417, No. 5007).

Warning

The Warning sign denotes a hazard. It calls attention to a procedure or practice, which if not correctly adhered to, could result in personal injury. Do not proceed beyond a Warning sign until the indicated conditions are fully understood and met.

Caution

The Caution sign denotes a hazard. It calls attention to an operating procedure or practice, which if not adhered to, could cause damage or destruction of a part or all of the product. Do not proceed beyond a Caution sign until the indicated conditions are fully understood and met.

Equipment Warnings

Before energizing the equipment, the earth ground of the power receptacle must be verified to be at earth potential and able to carry the rated current of the power circuit. Using extension cords should be avoided. However, if an extension cord must be used, insure the ground conductor is intact and the cord is capable of carrying the rated current without excessive voltage drop.

In the event that the ground path becomes less than sufficient to carry the rated current of the power circuit, the equipment should be disconnected from power, labeled as unsafe, and removed from place of operation.

Do not operate this equipment in the presence of flammable gases. Doing so could result in a life-threatening explosion.

Do not modify this equipment in any way. If component replacement is required, return the equipment to AMI facilities as described in the troubleshooting section of this manual.

If used in a manner not specified in this manual, the protection provided by the design, manufacture and documentation of the Model 1700 Instrument may be impaired.

Other Manual Conventions

This manual refers to measuring liquid nitrogen (LN₂) when referring to capacitance-based level measurement since nitrogen is by far the most common cryogenic fluid used by this sort of instrument. This instrument can be used to measure most any cryogenic liquid.

Instrument Configuration

The Model 1700 Instrument is configured at time of purchase in several ways:

- As a capacitance-based (typically liquid nitrogen) level instrument/controller.
- As a helium level instrument for level sensors with active length up to 40 inches (102 cm) for either 4.2K or 2K LHe.
- As a helium level instrument for level sensors with active length up to 80 inches (203 cm) for either 4.2K or 2K LHe.
- As a combination nitrogen and helium instrument/controller. Note that there is only one control channel even though the instrument can simultaneously display both liquid levels.

Every configuration may be further customized by the following options:

- Table top, single rack mounting, dual rack mounting.
- Line cord: North American, European Schuko, Australia/NZ, China, UK, pigtailed ends.

The instrument part number, shown on the serialization label located on the underside of the instrument, as well as in a field in the instrument Menu, identifies the configuration according to the following key:

1700-A-B-C-D-E where

A indicates the mounting method:

Tbl = tabletop

SR19L = single rack mounted, 19" wide rack standard, instrument on left side

SR19R = single rack mounted, 19" wide rack standard, instrument on right side

SR10L = single rack mounted, 10" wide rack standard, instrument on left side

SR10R = single rack mounted, 10" wide rack standard, instrument on right side

DRL = dual rack mount, 19" wide rack standard, instrument on left side

DRR = dual rack mount, 19" wide rack standard, instrument on right side

B indicates the line cord shipped with the instrument:

N = North American

E = European, Schuko

A = Australian/New Zealand

C = Chinese

U = United Kingdom

P = India/pigtailed

C indicates the capacitance-based level configuration:

CAP = capacitance-based (typically nitrogen)

N = not configured

The sensor active length and length units are appended to the configuration code, ie ...CAP-10.4CM-...

D indicates the helium level configuration:

LVHe2K = helium for sensors with active length ≤ 40 inches

LVHe4K = helium for sensors with active length ≤ 40 inches

HVHe2K = helium for sensors with active length ≤ 80 inches

HVHe4K = helium for sensors with active length ≤ 80 inches

N = not configured

The sensor active length and length units are appended to the configuration code, ie ...LVHe2K-40IN-...

E is used to denote any instrument customization:

S = standard (no customization)

C = instrument modified.

1 Introduction

1.1 Model 1700 Instrument

The AMI Model 1700 Liquid Level Instrument is a sophisticated measurement and control instrument which provides monitoring liquid helium and/or capacitance-based level sensors as inputs and provides for automatic level control based on user set parameters.

At time of purchase, the Model 1700 will be configured as:

- A combination capacitance-based and helium level instrument/controller. Note that there is only one valve control channel (typically used for nitrogen auto fill systems) but both liquid levels can be displayed simultaneously.

The Model 1700 instrument will measure liquid helium level using a superconducting level sensor. The instrument can be configured for standard (1 - 40") or long sensors (up to 80"). The instrument can be used with either 4.2K or 2K liquid helium level sensors.

1.1.1 Cryogenic (Capacitance-Based) Liquid Level Sensors

The Model 1700 Instrument will be used with a capacitance-based sensor for all cryogenic liquids except for helium. Special insulated capacitance sensors are required for liquids with electrical conductivity, i.e. water.

The Model 1700 Instrument must use an oscillator device to measure cryogenic liquid levels. This instrument has an internal oscillator if the distance between the sensor and the instrument is less than or equal to 15 feet (457 cm) so no external oscillator is necessary. If the distance between the sensor and the instrument is greater than 15 feet, an external oscillator/transmitter unit is necessary. The instrument will automatically

Introduction

configure itself to use the internal oscillator unless the external oscillator/transmitter is connected to the BNC connector on the rear panel.




Figure 1-1. Model 1700 Instrument using internal oscillator/transmitter

Introduction

Model 1700 Liquid Level Instrument




Figure 1-2. Model 1700 Instrument using external oscillator/transmitter

1.1.2 Superconducting Level Sensors

The instrument can be used with a superconducting level sensor to measure liquid helium levels. The instrument will be configured at the factory to measure liquid helium levels. The instrument will be configured for either standard active length level sensors (1 - 40 inches) or extra long active length level sensors (1 - 80 inches). The instrument can be used with either 4.2K or 2K level sensor types.

1.1.3 Digitally-Controlled

The Model 1700 contains a microcomputer which controls analog data conversion, display/keypad functions, communications I/O, dry contact closures, generation of analog output signals and relay control of a mains power outlet for solenoid valve autofill applications.

1.1.4 System Flexibility

The Model 1700 instrument incorporates data converters to translate signals between the analog and digital domains. Precision instrumentation techniques and potentiometer-free designs are employed throughout the Model 1700 Instrument to ensure long term stability and accurate signal translation for a wide range of conditions.

Introduction

Model 1700 Liquid Level Instrument

1.1.5 Display

The Model 1700 Instrument has a 4.3" diagonal measure TFT (Thin Film Transistor) color liquid crystal display of 480 x 272 pixels. The display has a 4-wire resistive touch overlay for easy operator input.

1.1.6 Intuitive Human-Interface Design

The Model 1700 instrument is designed to simplify the touch-screen based user interface. All functions were analyzed and subsequently programmed so that the most commonly used functions are addressed with the least number of keystrokes. The menus are presented in a logical fashion so that the operation of the Model 1700 is intuitive to the user.

1.1.7 Measurement Flexibility

Depending on the instrument version purchased, the Model 1700 can be configured to monitor and display one capacitance-based liquid level sensor (LN₂ version) or one liquid helium level sensor (LHe version) or both simultaneously. Set points can be assigned to control two dry-contact relay outputs. The make or break function of each relay can be set independently. A solid-state relay allows mains power to be switched according to set points to operate a solenoid-operated valve or other load.

1.1.8 Real Time Clock

The Model 1700 Instrument incorporates a real time clock with automatic daylight savings time adjustment (if desired).

1.1.9 Valve Control Output

The Model 1700 Instrument has a switched 2 ampere at mains voltage output for energizing a solenoid operated flow valve or other loads. This outputs is controlled by a zero-crossing solid state relay.

The fill start is triggered by a level reading. The fill is stopped by reaching a user-set level, or exceeding a maximum fill time.

1.1.10 Analog Outputs

The Model 1700 Instrument has two analog outputs, a 0-10 V_{DC} voltage output and a 4-20 mA_{DC} current loop output. The 4-20 mA_{DC} loop output has 1500 V_{PK} circuit isolation. The outputs can be used simultaneously. The sources for either of the analog outputs can be either one of the two configured sensors. Both outputs can be driven from the same source, if desired.

1.1.11 Signal Relays

The Model 1700 Instrument has two signal relays that change state based on a measured input. The set points of these relays are user-selectable as is their function, alarm on level above or below set point. For the dual

Introduction

Model 1700 Liquid Level Instrument

channel instrument, both relays can be assigned to either channel. For the single channel instrument, both relays are assigned to the only channel.

1.1.12 Connectivity

The Model 1700 Instrument has a 10Base-T Ethernet connection as well as a 115,200 baud RS-232 port for connecting to other equipment. The instrument communicates with a SCPI-based command set. The command set is 100% backward compatible with the AMI Model 135, 136, 185, and 186 instruments when configured as a single channel instrument. When configured as a dual channel instrument, additional commands are incorporated into the command set.

The Model 1700 Instrument allows for remote operation with an external browser via TCP/IP connection. All commands that are available by touching the local screen are available via the web browser




Figure 1-3. Model 1700 (Dual Display) Via Web Browser

Note

Using a web browser to connect to the instrument allows different browser sessions to display different information. For instance, one browser window can display helium level only and a second browser window can display nitrogen level only. The instrument can be configured independently as well since its display is also a browser. Refer to “Configure the instrument to display nitrogen level” on page 30. Refer to “Configure the instrument to display helium level” on page 31.

1.2 Model 1700 Front Panel Layout



Figure 1-4. Model 1700 Front View; Dual Channel Instrument Shown

1.3 Model 1700 Rear Panel Layout




Table 1-1. Model 1700 Rear Panel Description

1	Computer Network Connector
2	Aux I/O Connector
3	RS-232 Serial Connector
4	Switched Valve Outlet Socket
5	Capacitive Sensor Input Connector
6	Mains Power Entry Connector
7	LHe Level Sensor Connector

Introduction

Model 1700 System Specifications

1.4 Model 1700 Specifications @ 25°C

System Architecture	
Display:	4.3" 24-bit color TFT display, 480x272 pixel with resistive touch screen
Sensor types:	Capacitance-based liquid level Superconductivity-based liquid level
Maximum length readout:	Capacitance-based liquid level - 999 in Superconductivity-based liquid level (LHe) - 80 in
Superconducting (LHe) sensor excitation:	Continuous reading or Sample and Hold mode
System operating firmware storage:	Flash memory
System clock:	Real time clock with automatic DST adjustment
Display measurement units:	liquid level in cm, in or percent
Level Measurement	
Resolution:	0.1%, 0.1 cm, 0.1 in
Accuracy:	±0.5% of active sensor length
Linearity:	±0.1% or 1 mm (whichever is greater)
Capacitance Sensor Excitation Voltage:	5 V _{DC}
Capacitance Transmitter Measurement Resolution:	0.7 pF
Superconducting Sensor Current:	4.2K LHe Temperature: 75 mA _{DC} nominal 2K LHe Temperature: 57 mA _{DC} nominal
Dirty Sensor Mode:	Approximately twice normal current for 1 second prior to normal measurement excitation
Superconducting Sensor Voltage:	4.2K Sensor: approximately 0.87 V _{DC} per inch of sensor active length @ 10K 2K Sensor: approximately 0.66 V _{DC} per inch of sensor active length @ 10K
Maximum Open Circuit Voltage:	48 or 96 V _{DC} , galvanically isolated
Operating Parameters	
Alarm Set points:	0% to 100%, adjustable; Alarm condition settable to above or below set point; Assignable to any channel as high or low alarm for dual instrument
Controller Output:	Line voltage @ 2 A _{AC} (maximum)
Sample and Hold Period:	1 second to 86,400 seconds (24 hrs)
Audible alarm:	3500 ± 500 Hz, 73 to 86 dB(A)
Analog Outputs	

Introduction

Model 1700 System Specifications

Output Types:	0-10 V _{DC} and simultaneous 4 - 20 mA _{DC} ; Each assignable to either level channel for dual instrument
4-20 mA Current Loop Power Supply Voltage:	12-32 V _{DC}
0-10 V _{DC} Recorder Output Output Load:	5k ohms or greater
Converter Resolution:	12 bits
Integral Non-linearity:	±1LSB
Differential Non-linearity ^a :	±1LSB

Relays	
№1 and/or №2 (W171DIP-7, or equivalent):	Contact Form: 1 Form A (SPST-NO) Maximum Switched Power: 0.29 W Maximum Switched Current: 3 A Switching Voltage: 60 V _{AC} / 100 V _{DC}

Auto Fill	
Controller output socket:	IEC 60320-13 socket on rear panel
Controller output power:	2 A _{AC} at line voltage
Autofill start/stop triggering:	Level-based; 0% to 100%, adjustable; Assignable to either nitrogen or helium channel
Fill Timeout Period:	1 minute to 99 hours, 59 minutes
Fill error alarm:	Fill time out

Communication Protocol	
Host computer network protocol:	10Base-T TCP/IP and RS-232 115,200 baud
IP Addressing:	DHCP or static, IPv4
Network connectivity and traffic indication:	Link and Activity LEDs on instrument rear panel
RS-232 connector specifications:	9-pin D-sub female connector to connect standard DTE 9-pin D-sub male connector using a standard straight cable
Communication command set:	SCPI-based.

Power Requirements	
Primary:	100-240 ±10% V _{AC} , 50-60 Hz, 2.2 A maximum (200 VA plus sum of controller output)
Backup battery for RTC:	CR2032

Physical	
Dimensions ^b :	Table top configuration: 3.8" H x 8.4" W x 11.4" D [97 mm H x 213 mm W x 290 mm D] Single rack mount configuration: 3.5" H x 19.0" W x 11.4" D [89 mm H x 483 mm W x 290 mm D]
Weight:	table-top configuration: 3.3 lbm [1.5 kG]; single rack-mount configuration: 4.0 lbm [1.8 kG]

Introduction

Model 1700 System Specifications

Environmental Limits	
Ambient Temperature:	Operating: 0°C to 40°C [32°F to 104°F]; Non-operating: -20 °C to 60 °C [-4 °F to 140 °F]
Relative Humidity:	0 to 95%; non-condensing
Maximum Instrument Background Field:	Gauss (TBD)

Standards	
Test Standards	Testing of Equipment for Measurement, Control, and Laboratory Use (IEC 61326-1:2012, EN 61326- 1) Electrostatic Discharge (ESD) (EN 61000-4-2) Radiated Immunity (EN 61000-4-3) Fast Transient Burst (EN 61000-4-4) Surges (EN 61000-4-5) Conducted Immunity (EN 61000-4-6) Power Frequency Magnetic Field (EN 61000-4-8) Voltage Dips and Interrupts (EN 61000-4-11) Harmonics (EN 61000-3-2) Flicker (EN 6100-3-3) Conducted Emissions (EN 55011/IEC/CISPR 11) Radiated Emissions (EN 55011/IEC/CISPR 11) Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use (IEC 61010-1)

- a. Guaranteed monotonic over operating temperature range
- b. H = height; W = width; D = depth

Introduction

Model 1700 System Specifications

2 Installation

Warning

Equipment warnings apply to all system installation configurations. Refer to “Equipment Warnings” on page xv, in the Forward to be familiar with the safety requirements for a system installation.

2.1 Unpacking and Inspecting the Instrument

Carefully remove the equipment, interconnecting cabling, and documentation CD (and/or printed material) from the shipping carton, and remove all packaging material.

Note

If there is any shipping damage, save all packing material and contact the shipping company representative to file a damage claim. Do not return to AMI unless prior authorization has been received.

2.2 Mounting the Model 1700 Instrument

If the Model 1700 Instrument is to be used as a table top model, place the equipment on a flat, secure surface.

If the Model 1700 Instrument is to be rack mounted, install it in a 19" wide instrument rack using the mounting hardware supplied by the rack cabinet manufacturer. The feet on the bottom of the instrument may be removed to facilitate rack mounting. Secure the front panel to the rack rail in each of the four corners.

2.3 Installing Capacitance-Based Liquid Level Sensors

Refer to the installation instructions provided with the level sensor(s). The following steps are general installation notes and should be used to supplement the installation instructions provided with the sensor.

Exercise care when installing the capacitance sensor since dents, crimps, bends or other physical distortions in the cylindrical capacitor will change electrical characteristics, possibly causing calibration errors and/or disruption of proper instrument operation. Before installing the sensor, review “Calibration” on page 45 to determine what, if any, calibration procedures may be necessary prior to operation.

Installation

Power Requirements

Note

The coaxial interconnecting cables and the transmitter should be mounted in such a manner as to avoid large temperature changes such as those encountered in the path of dewar vents.

1. Carefully remove the sensor from the shipping container and remove all packaging material.

Note

If there is any shipping damage, save all packing material and contact the shipping representative to file a damage claim. Do not return the instrument to AMI unless prior authorization has been received.

2. Install the sensor in the vessel using the specified fitting of the sensor.

Caution

- Ensure the sensor is mounted with the top vent hole located inside of the cryostat.

- Avoid installing in a location where icing may occur. Ice formations or moisture buildup on the BNC connector may cause the sensor to short out indicating a higher liquid level than actually exists.

- Moisture or contaminants in any of the BNC coaxial connectors can short out the sensor and cause measurement errors. A pack of non-conductive electrical connection lubricant (ECL or "Dielectric Tune-up Grease") has been included with the liquid level sensor packaging to reduce the possibility of this occurring. If desired, apply a small amount of ECL to any of the BNC connectors that may be exposed to moisture. Mate the ECL-coated connectors then remove any excess ECL from the outside of the connector. Added protection can be achieved by covering the ECL-coated connections with a short section of heat-shrink tubing.

- MSDS sheets for the ECL are available upon request.

- Exercise care when installing the sensor since dents, crimps, bends or other physical distortions in the thin wall capacitor will change electrical characteristics possibly causing calibration errors and/or disruption of proper instrument operation.

2.4 Installing the Liquid Helium Level Sensor

1. The sensor must be mounted with the electrical leads at the top.

Installation

Power Requirements

2. For minimum losses, mount the liquid helium sensor so that warm helium gas rising from the sensor can pass directly out of the dewar without contacting surfaces at 4.2K.
3. Do not mount the sensor in restricted areas (tubes, etc.) where the liquid level around the sensor might be depressed by pressure differences in the gas. Do not cover the holes in the sensor.
4. The sensor may be mounted by taping or clipping it to an appropriate support structure. Do not exert excess pressure on the sensor with the mounting device to avoid crushing the tube. Avoid constraining both ends of the sensor and allow for contraction of the sensor during cooldown.

Caution

Do not operate the sensor in a vacuum. Operating the sensor in a vacuum may cause thermal damage and/or destruction of the superconducting filament sensor. Do not inadvertently turn the instrument on with the sensor in an evacuated chamber. Operation in pumped liquid helium environments is acceptable to 1K as long as liquid helium is present.

5. Avoid installing in a location where icing (frozen water or gas) may occur since ice formations may cause erratic operation. Ice formation on the NbTi filament may stop the propagation of the normal (resistive) zone before it actually reaches the liquid/gas interface. This will give an indication of a higher helium level than actually exists.
6. Connect the sensor to the Model 1700 LHe Level Sensor connector on the instrument rear panel (refer to “Liquid Helium Connector J1 Wiring” on page 95). The liquid helium level sensor leads are color coded:

Table 2-1. LHe Level Sensor Wire Identification

Wire Function	Teflon Insulation Color	Formvar Insulation Color	Instrument Connector Pin
I+	Red	Red	1
V+	Blue	Green	8
V-	Yellow	Natural	6
I-	Black	Blue	7

Installation

Connecting the Transmitters to the Sensors

2.5 Connecting the Sensor(s) to the Instrument

2.5.1 Connecting a Capacitance Sensor

The sensor may be connected directly to the BNC connector on the instrument rear panel if the length of the coaxial cable is 15 feet or less. If the sensor is greater than 15 feet from the instrument, an external oscillator/transmitter unit must be used. Refer to figures “Model 1700 Instrument using internal oscillator/transmitter” on page 2 and “Model 1700 Instrument using external oscillator/transmitter” on page 3 as appropriate.

If the transmitter is connected to the sensor with a length of coaxial cable, the a 15 ft standard length cable, with part number of EH2362, is available from AMI. Contact the factory for details. Speak to an AMI Sales Engineer before using cables longer than 15 feet.

Note

Regarding the coaxial cable is used to connect the capacitance level sensor to the instrument or oscillator/transmitter, in order to maintain system performance and accuracy, the cable must be Trompeter TCC-75-2 and should not be longer than 15 feet [4.57 m]. If a longer length section of coaxial cable is necessary, please discuss with an Authorized AMI Technical Representative.

Note

If an oscillator/transmitter is used, the length of coaxial cable between the oscillator/transmitter and the instrument may be up to 500 feet in length.

Caution

Moisture or contaminants in any of the BNC coaxial connectors can short out the sensor and cause an erroneous readings or transmitter failure. A pack of non-conductive electrical connection lubricant (ECL or “Dielectric Tune-up Grease”) has been included with the liquid level sensor packaging to reduce the possibility of this occurring.

To connect the coaxial cable to the BNC connector on the capacitance level sensor:

1. Apply a small amount of ECL to any of the BNC connectors that may be exposed to moisture.
2. Mate the ECL-coated connectors then remove any excess ECL from the outside of the connector.

Installation

Connecting the Transmitters to the Sensors

3. Cover the ECL-coated connections with a short section of heat-shrink tubing, also included, for added moisture protection.

Note

MSDS sheets for the ECL are available upon request.

2.5.2 Connecting a Liquid Helium Level Sensor

The instrument is connected to the level sensor with a 4-conductor cable which has a 9-pin D-sub male connector on one end that mates with the connector used at the instrumentation feed through connector on the cryostat. This connector is typically a multi-pin circular type connector.

Prepare the sensor to be connected to the instrument by soldering the sensor leads to a male 9-pin D-Sub connector which will connect to the female 9-pin D-Sub connector on the transmitter. Refer to the *Appendix* of this manual and the AMI sensor manual for the proper pin out and wire color connections. Connect the sensor to the connector on the transmitter.

Warning



Although the sensor connector terminals are isolated from earth ground and therefore touching one terminal is not hazardous, the voltage between terminals is at a hazardous potential. The sensor connector is for use with an AMI LHe sensor and the wiring for the sensor is to have no live parts which are accessible. Conductors connected to its terminals must be insulated from user contact by basic insulation rated for 150 V_{AC} (Category I).

The lead wire for the sensor may be sized by the following equation:

$$R = 420 - 5.21L$$

where R is the maximum allowable resistance (in ohms) for each lead wire from the instrument to the sensor, and L is the active length of the

Installation

Setting Up an Autofill System

connected helium level sensor in inches. Tables for active sensor length vs. lead wire distance are provided below.

	$R=367$	$R=315$	$R=263$	$R=211$	$R=107$	$R=3.2$
Distance	$L=10''$	$L=20''$	$L=30''$	$L=40''$	$L=60''$	$L=80''$
10 ft.	36 AWG				36 AWG	34 AWG
20 ft.						30 AWG
30 ft.						28 AWG
40 ft.						27 AWG
50 ft.						24 AWG
100 ft.						22 AWG
200 ft.						16 AWG
500 ft.					32 AWG	

Table 2-2. Minimum recommended wire gauge for copper lead wire

Note

If the system is an Autofill system, proceed to section “Setting Up an Autofill System” on page 16. If the installed capacitance or liquid helium level sensor(s) are used for indication only (not autofill), proceed to “Energizing the Model 1700 Instrument” on page 21.

2.6 Setting Up an Autofill System

Autofill systems are useful in applications where liquid level in a dewar must be maintained without operator intervention.

To create an autofill system, a cryogenic liquid source must be attached to the target or level controlled dewar via a solenoid-controlled valve. The valve will be controlled by the Model 1700 Instrument.

2.6.1 Autofill System Description

For autofill, the system consists of a Model 1700 Instrument with a liquid level sensor, and a solenoid-operated flow valve. The instrument makes continuous level measurements and based on level, energizes the valve to

Installation

Setting Up an Autofill System

begin liquid transfer. The transfer is stopped when the measured level reaches a user-determined point.




Figure 2-1. Typical Autofill Setup

Table 2-3. Standard Autofill Setup Description

Number	Item
1	Model 1700 Liquid Level Instrument (Level Controller)
2	IEC60320 C13 socket labeled VALVE CONTROL
3	BNC connector labeled O/T
4	Solenoid-operated flow control valve line cord with IEC60320 C14 plug
5	Instrument IEC60329 C14 Power cord socket
6	Coaxial cable connecting the oscillator / transmitter and the instrument

Installation

Setting Up an Autofill System

Table 2-3. Standard Autofill Setup Description

Number	Item
7	Oscillator / Transmitter (optional; refer to Figure 2.5.1 on page 14)
8	Level Sensor for level controlled or target dewar
9	Coaxial cable connecting the Oscillator / Transmitter and the liquid level sensor (optional; refer to Figure 2.5.1 on page 14)
10	Fill port on target dewar
11	Transfer line attached to the target dewar and the fill solenoid valve on the source dewar
12	Solenoid-operated fill valve
13	Supply dewar relief valve
14	Supply dewar

2.6.2 Autofill System Setup

Caution

A relief valve must be used in autofill systems to ensure no cryogenic liquid can be trapped in a transfer line volume where expansion can cause damaging pressure. This can occur if the solenoid operated fill valve and the supply dewar isolation valve are closed, trapping a cryogenic liquid in a confined volume. All AMI transfer line systems include a relief valve to preclude this sort of event.

1. Mount the level sensor (**8**) in the target dewar.
2. Connect the transfer line (**11**) and fill solenoid valve (**12**) or supply manifold to the source dewar.
3. Connect the other end of the transfer line to the fill port (**10**) on the valve/manifold of the target dewar.
4. Connect the sensor to the instrument.
 - a. For distances of 6 feet and less, connect the coaxial cable (**6**) between the BNC connector on the liquid level sensor and the BNC connector on the back of the instrument labeled **O/T (3)**.
 - b. For distances greater than 6 feet, connect the coaxial cable (**9**) between the BNC connector on the liquid level sensor and the BNC connector on the oscillator / transmitter (**7**). Use a second length of coaxial cable (**6**) to connect between the oscillator /

Installation

Power Requirements

transmitter **(7)** and the BNC connector on the back of the instrument labeled **O/T (3)**.

5. Connect the solenoid valve **(12)** to the IEC60320 C13 valve socket¹ on the Model 1700 Instrument rear panel labeled **VALVE CONTROL (2)**.

Connect the IEC320 C14 inlet connector and instrument power cord **(5)** to an appropriate power source (see below).

Note

*Should the nitrogen level sensor become disconnected from the instrument, a **LOSS OF SENSOR** message will be displayed and the autofill valve will be shut (**M-CLOSED**). When the sensor connection has been restored, the instrument will display the level but the autofill state will have to be manually changed back to **AUTOFILL**.*




Figure 2-2. Loss Of nitrogen Level Sensor

2.7 Power Requirements



Warning

The Model 1700 Instrument operates on 50-60 Hz power and may be powered from 100-240 Vac. Insure that the input ground terminal is connected securely to an external earth ground.

Insure the detachable mains supply cord is of suitable rating, i.e. 10 A (min) at 125 Vac for North America.

Insure the power switch is in the OFF (**O**) position. Plug the Model 1700 Instrument line cord into the power entry module on the instrument rear panel and into the appropriate power receptacle.

1. The valve socket must be IEC60320 C13 type.

3 Operation

This section describes the operation of the Model 1700 Instrument.

3.1 Energizing the Model 1700 Instrument

1. Turn the power switch on the front panel of the instrument to the On (I) position. The display will briefly show the AMI logo and then boot information.

The boot process takes approximately 30 seconds. This time can be longer (approximately 2 minutes) if the instrument has been configured for a network connection and then is booted without the network present.

Note

If the instrument was purchased with a level sensor, the instrument will be configured and calibrated at the factory.

2. When the boot process is complete, the instrument will display the home (level) screen.
 - a. If the instrument is configured for only liquid nitrogen level, the display will look similar to the following:
 - b. If the instrument is configured for only liquid helium level, the display will look similar to the following:




Figure 3-1. LN₂ Home Screen




Figure 3-2. LHe Home Screen

Operation

Model 1700 Screen Navigation

- c. If the instrument is configured for both liquid nitrogen and liquid helium levels, the display will look similar to the following:
- d. If the instrument requires calibration¹, refer to the following chapter to calibrate the instrument with an AMI level sensor.




Figure 3-3. Both LN2 and LHe Home Screen

3.2 Screen Navigation

3.2.1 Home Screen Footer

Every screen has a footer. The level home screen (refer to Figure 3.2.1 as an example) displays both level information and buttons in the footer to navigate to other screens.

The home screen footer appears as follows:



Figure 3-4. Home Screen Footer

Table 3-1. Model 1700 Instrument Home Screen Footer

No.	Name	Function
①	MENU	Takes the user to the menu screen
②	ALARM	When illuminated, displays an alarm condition
③	AUTOFILL	Indicates the condition of the autofill function

3.3 Navigating the Instrument Menus

The menu system is invoked by pressing the Menu button in the lower left corner of the instrument level display screen.



Figure 3-5. MENU button on HOME screen

1. If the instrument was purchased with level sensor(s), the instrument will be shipped set up and calibrated.

Operation

Using the Model 1700 Menus

When invoked, the MENU screen will be displayed:




Figure 3-6. MENU Selection Screen

Pressing the Home icon in the lower left corner of the screen will display the Home screen.



Figure 3-7. HOME button on MENU screen

Once a MENU has been chosen, a back button will be presented to the right of the Home icon. This button will return the user to the previously displayed screen.



Figure 3-8. BACK button

3.3.1 Menu Overview

The Model 1700 Instrument displays menus on the graphic display to the left of the keypad. Press **MENU** on the menu keypad to display options on the graphic display. Menu options are listed in the following Table.

Operation

Using the Model 1700 Menus

Note

The following table shows all menu choices, some of which will not be shown based on the instrument configuration.

Table 3-2. Model 1700 Dual Level Instrument Menu Structure

Menu Label	Function	Field Type
SENSORS	NITROGEN OSCILLATOR:	Information: INTERNAL, EXTERNAL ^a
	N2 LEVEL ON HOME SCREEN?	Toggles between: YES, NO
	HELIUM SENSOR TYPE:	Toggles between: 2K, 4.2K
	He LEVEL ON HOME SCREEN?	Toggles between: YES, NO
	CALIBRATE NITROGEN	Transfer to another screen
	CALIBRATE HELIUM	Transfer to another screen

Operation

Using the Model 1700 Menus

Table 3-2. Model 1700 Dual Level Instrument Menu Structure

Menu Label	Function	Field Type
OUTPUTS	ALARM RELAY No.1 SOURCE:	Toggles between: DISABLED, NITROGEN, HELIUM
	(ALARM RELAY No. 1) SETPOINT:	Data entry <value> %
	(ALARM RELAY No. 1) ALARM WHEN LEVEL <state>	Toggles between: ≤, ≥ SETPOINT
	ALARM RELAY No.2 SOURCE:	Toggles between: DISABLED, NITROGEN, HELIUM
	(ALARM RELAY No. 2) SETPOINT:	Data entry <value> %
	(ALARM RELAY No. 2) ALARM WHEN LEVEL <state>	Toggles between: ≤, ≥ SETPOINT
	AUTOFILL SOURCE:	Toggles between: DISABLED, NITROGEN, HELIUM ^b
	AUTOFILL CONTROL: START:	Data entry <value> % ^b
	AUTOFILL CONTROL: STOP:	Data entry <value> % ^b
	FILL TIMEOUT (N2 ONLY):	Data entry <value> MINUTES ^b
	0-10 Vdc SOURCE:	Toggles between: DISABLED, NITROGEN, HELIUM ^b
	4-20 mA SOURCE:	Toggles between: DISABLED, NITROGEN, HELIUM ^b

Operation

Using the Model 1700 Menus

Table 3-2. Model 1700 Dual Level Instrument Menu Structure

Menu Label	Function	Field Type
NETWORK	ADDRESS: <value>	Data entry or Information
	NETMASK: <value>	Data entry or Information
	GATEWAY: <value>	Data entry or Information
	ADDRESSING:	Toggles between: DISABLED, STATIC, DYNAMIC
	MAC ADDRESS:	Information: <value>
SYSTEM (page 1)	NITROGEN SENSOR NAME:	Data entry: <value>
	HELIUM SENSOR NAME:	Data entry: <value>
	SYSTEM DATE & TIME:	Information: <values>^a
	(SYSTEM DATE & TIME) SET	Transfer to another screen
	(SYSTEM SETTINGS) PAGE 2	Transfer to another screen
SYSTEM: DATE & TIME (page 4)	SYSTEM DATE AND TIME: YEAR	Data entry
	SYSTEM DATE AND TIME: MONTH	Data entry
	SYSTEM DATE AND TIME: DAY	Data entry
	SYSTEM DATE AND TIME: HOUR	Data entry
	SYSTEM DATE AND TIME: MIN	Data entry

Operation

Using the Model 1700 Menus

Table 3-2. Model 1700 Dual Level Instrument Menu Structure

Menu Label	Function	Field Type
SYSTEM (page 2)	SERIAL NUMBER:	Information: <value> ^a
	HARDWARE VERSION:	Information: <value> ^a
	DATE OF MANUFACTURE:	Information: <value> ^a
	FIRMWARE VERSION:	Information: <value> ^a
	HELIUM POWER SUPPLY:	Information: NONE
	RESET INSTRUMENT TO FACTORY DEFAULTS	Transfer to another screen
SYSTEM (page 3)	RESET INSTRUMENT	Performs a function
CALIBRATE TOUCH SCREEN	Assists user in performing the instrument touch screen calibration	Transfer to another screen
SHUT DOWN SYSTEM	Shuts down the instrument in an orderly fashion which reduces boot time for the next power on.	Transfer to another screen

Operation

Using the Model 1700 Menus

Table 3-2. Model 1700 Dual Level Instrument Menu Structure

Menu Label	Function	Field Type
SENSORS: CALIBRATE NITROGEN	OSCILLATOR PERIOD:	Information: <value> ^a μ s
	SENSOR ACTIVE LENGTH:	Data entry: <value> cm
	PERFORM MAX CAL	Transfer to another screen and Information: <value> ^a μ s
	PERFORM MIN CAL	Transfer to another screen and Information: <value> ^a μ s
	NO SENSOR CAL	Transfer to another screen and Information: <value> ^a μ s
	APPROX CAL. VALUE:	Data entry: <value>
	(APPROX CAL. VALUE:) APPLY	Performs a function
SENSORS: CALIBRATE HELIUM^b	SENSOR SAMPLE INTERVAL:	Data entry: <value> min
	MEASURE SENSOR VOLTAGE	Performs a function: <value> V ^a
	SENSOR ACTIVE LENGTH:	Data entry: <value> cm
	CONTIN. MEASURE TIME LIMIT:	Data entry: <value> min
	DIRTY SENSOR MODE:	Toggles between: DISABLED, ENABLED)

a. Displays the state or value (display only).

b. OUTPUTS Menu page 2

3.3.2 Editing a Field

Once a field on a screen has been selected for editing, the footer changes to appear as follows:



Table 3-3. Model 1700 Instrument Footer during editing a field

Button No.	Name	Function	Reference Paragraph
①	SAVE	Saves the entries made on the screen. The footer changes to the footer shown in Table 3-1, above.	
②	CANCEL	Exits out of the screen, not saving entries. The footer changes to the footer shown in Table 3-1, above.	
③	ALARM Status	Refer to Table 3-1, above.	
④	AUTOFILL Status		

Operation

Using the Model 1700 Menus

3.3.3 Menu Navigation




Figure 3-10. Model 1700 Menu Structure

3.4 Capacitance (Liquid nitrogen) Level

3.4.1 Configure the instrument to display nitrogen level

Note

Operation

Using the Model 1700 Menus

If the instrument was purchased with a capacitance-based level sensor, Steps 1 through 4, below have already been performed.

1. From the main screen, choose the following: MENU > SENSORS.

2. The instrument has a built-in sensor oscillator which will be used if the level sensor is connected to the BNC connector on the rear panel of the instrument with a maximum of 6 foot of RG59/U coaxial cable (refer to Figure 1-1 on page 2). Ensure the “NITROGEN OSCILLATOR:” field displays “INTERNAL”.




Figure 3-11. INTERNAL Oscillator Indicated

3. If the sensor is greater than 15 feet from the instrument, an external oscillator/transmitter will have to be used and the “NITROGEN OSCILLATOR:” will display “EXTERNAL” (refer to Figure 1-2 on page 3).
4. Ensure that N2 LEVEL ON HOME SCREEN? is set to YES.
5. Press the Home icon at the bottom of the screen.




Figure 3-12. Nitrogen Level Displayed on Home Screen

3.5 Superconducting (Liquid Helium) Level

3.5.1 Configure the instrument to display helium level

Note

If the instrument was purchased with an AMI LHe level sensor, Steps 1 through 4, below have already been performed.

1. From the main screen, choose the following: MENU > SENSORS.

Operation

Using the Model 1700 Menus

2. Choose the type of LHe level sensor, 4.2K or 2K.




Figure 3-13. Type of LHe Level Sensor

3. Ensure that He LEVEL ON HOME SCREEN? is set to YES.




Figure 3-14. Helium Level On Home Screen

4. Press the **CALIBRATE HELIUM** button.




Figure 3-15. Helium Level On Home Screen

5. Press the **SENSOR ACTIVE LENGTH** field.




Figure 3-16. Calibrate Helium Screen

Operation

Using the Model 1700 Menus

6. In the numeric pop-up keypad, enter the sensor active length in centimeters. Press Enter when finished.




Figure 3-17. Numeric Keypad

7. Press in the **CONTIN. MEASURE TIME LIMIT** field.
8. Using the pop-up numeric keypad, enter the maximum amount of time that the sensor should remain energized in the **MEASURE CONTINUOUSLY** mode before automatically changing to the **SAMPLE AND HOLD** mode. This limit prevents inadvertent sensor energization for long periods of time which will cause excessive liquid helium boil off. When the sensor is energized (**MEASURE CONTINUOUSLY**) a timer is started and after the **CONTIN. MEASURE TIME LIMIT** is reached, the sensor is de-energized (switches back to **SAMPLE AND HOLD** mode).




Figure 3-18. Calibrate Helium Screen

9. Set the **SENSOR SAMPLE INTERVAL** to an appropriate value. This value is used by the **SAMPLE AND HOLD** timer to determine how often the reading is updated. The timer is started when the instrument is powered up or when this value is changed.




Figure 3-19. Setting the Sample Interval

10. Press the Home icon in the screen footer.

Operation

Using the Model 1700 Menus

3.5.2 Sample and Hold Operation

11. While viewing the home screen, toggle between **SAMPLE AND HOLD** and **MEASURE CONTINUALLY** modes by pressing the button.




Figure 3-20. Helium Level Measurement Mode

3.5.3 Other Liquid Helium Functions

3.5.4 Other Liquid Helium Functions

1. The instrument displays the LHe level sensor voltage. If the voltage is shown in light blue, it is the actual (real-time) voltage as the instrument is in **MEASURE CONTINUALLY** mode. If it is displayed in gray, it is the voltage measured the last time the sensor was energized (**SAMPLE AND HOLD** mode). To update the reading, press the **MEASURE SENSOR VOLTAGE** button and the level reading will be updated as will the displayed sensor voltage.




Figure 3-21. Helium Sensor Voltage

2. AMI expects the helium level sensor to be reasonably clean and free from oil, water, ice, etc. for proper operation. However, it is recognized that some experiments might result in some material being deposited on the sensor wire. Ice formation at some point on the sensor is a typical occurrence. Therefore, the Model 135/136 has the capability of increasing the current for a short period of time at the beginning of the measurement cycle (in the **SAMPLE** mode only) to try and drive the resistive zone of the sensor wire past the dirty region. This operation may or may not be successful depending on the degree of sensor contamination. This mode should be viewed as a stopgap measure only. If correct readings cannot be reestablished, the only choice is to warm the sensor or remove for cleaning or replacement.

Note

Operation in the dirty sensor mode increases liquid helium losses. Consequently, operation in this mode should not be used unless the sensor is known or anticipated to become dirty or the helium level measurement is in question due to unclean operation.

Operation

Using the Model 1700 Menus

To enable the **DIRTY SENSOR MODE**, Press the **DISABLED** button until **ENABLED** appears and press **SAVE**.

3. Press the home icon in the footer to return to the home screen.




Figure 3-22. Dirty Sensor Mode

3.6 Alarms and Relays

3.6.1 Overview

The Model 1700 Instrument has two types of alarms, level-based and time-based alarms.

3.6.1.1 Level-Based Alarms

The Model 1700 Instrument has two user-configurable level alarms. Each alarm can be triggered by either level measurement (for dual level configured instruments). Each alarm can be configured to be active above or below a user-defined setpoint.

When an alarm condition occurs, an audible alert will sound and visual indication on the front panel.

Alarm №1 and №2 have relays associated with them. These relays have Normally Open (NO) contacts. The alarm/relays can be assigned to either the helium or nitrogen channel and the alarm/relay can each be configured to have the alarm active when the reading is either \leq or \geq the setpoint. As an example of this setup flexibility, a level channel can be configured to have an alarm condition when the level is outside a normal operating band, or in a dual level configuration, each channel can have a low level alarm.

3.6.1.2 Time-Based Fill Alarm (nitrogen only)

The Model 1700 Instrument has an alarm to indicate that there is a problem with the LN_2 autofill function. If enabled, the instrument will start a timer when an autofill condition is initiated, and if the level has not

Operation

Using the Model 1700 Menus

reached the fill stop level within the user-set period of time, an Autofill Timeout alarm will occur.

This alarm will cause three things to occur:

1. the de-energizing of the fill valve socket on the instrument rear panel,
2. an audible alarm will sound,
3. a visual indication in the footer of the instrument front panel which shows repeatedly **ALARM > NITROGEN > TIMEOUT > MUTE?**¹ > **ALARM** > etc.
4. The blue AUTOFILL text in the right side of the screen footer will turn to red T-CLOSED indicating that the fill timeout has automatically closed the fill valve by deenergizing the power socket on the instrument rear panel.




Figure 3-23. Fill Timeout Alarm Screen

This Autofill Timeout function can be disabled by setting the interval to 0 minutes.

3.6.1.3 Multiple alarms

More than one alarm condition can occur at the same time. The footer will display the cause(s) of the alarm condition(s).

3.6.2 Configuring Alarm Setpoints

1. From the MENU choice, select OUTPUTS and the first page of the Output Configuration screen will be displayed.




Figure 3-24. Output Configuration Screen, Page 1

1. The **MUTE?** function is not applicable to the screen on a remote browser since audible alarms are not supported remotely

Operation

Using the Model 1700 Menus

2. If the instrument is configured to display both nitrogen and helium, toggle the Alarm Relay Source fields to choose the appropriate level source. If the instrument is configured for either nitrogen or helium, the source will automatically be assigned to the configured level.




Figure 3-25. Relay Source Configuration

3. Set the levels at which the alarm will be triggered in the Setpoint fields.




Figure 3-26. Relay Setpoints

4. Use the \leq or \geq button to toggle between the two states of alarm, either alarm when the indicated level is less than or equal to the setpoint or alarm when the indicated level is greater than or equal to the alarm setpoint.




Figure 3-27. Alarm Condition Above or Below Setpoints

3.6.3 Acknowledging an Alarm

Note

The alarms are not “latched” so if the alarm condition clears itself, the instrument will remove the alarm condition.

1. When an alarm is initiated, several things will occur:

Operation

Using the Model 1700 Menus

- a. The bar-graph level display that is causing the alarm condition as well as the sensor name will flash red.
- b. The **ALARM** button in the footer will indicate what the alarm condition is, either by showing LO LEVEL, HI LEVEL, or TIMEOUT.




Figure 3-28. Alarm Annunciator

- c. An audible alarm will be energized.
2. For example, a helium low level alarm will flash the following three displays in a repeating fashion:




Figure 3-29. LO LEVEL Alarm Condition Footer Displays

3.6.4 Muting an Alarm

The audible alarm can be muted by pressing the **ALARM** button in the footer. As long as the alarm condition occurs with muting enabled, the **ALARM** button in the footer will alternate between the alarm conditions (Figure 3-29) and **MUTED**.




Figure 3-30. Muted Alarm Condition

3.6.5 Resetting the Autofill Timeout Alarm

1. Press the red **T-CLOSED** text in the right side of the footer twice until **M-CLOSED** is displayed.
2. Press **SAVE**.
3. To restart the autofill process, Press the **M-CLOSED** annunciator once and the **AUTOFILL** annunciator will be displayed. Press **SAVE** and the autofill sequence will be enabled.

3.7 Configure the Autofill Function

3.7.1 Autofill overview

There are several variables that must be addressed to set up an autofill system. These include the level indication that will be used to control autofill, the Fill Start level (A), the Fill Stop Level (B), and the fill timeout interval.

3.7.2 Setting the autofill parameters

1. From the Menu screen, choose **OUTPUTS** and then choose **PAGE 2**.

2. Toggle the **AUTOFILL CONTROL** button until **NITROGEN** is displayed.




Figure 3-31. Autofill Level Control Selection

3. Enter the Fill **START** and **STOP** levels.




Figure 3-32. Autofill Level Start Setting

4. Enter the **FILL TIMEOUT (N2 ONLY)** interval. Refer to section 3.6.1.2 on page 35 for a description of the Autofill timeout function.
5. Press the home icon in the footer to return back to the level display.




Figure 3-33. Autofill Timeout Setting

Operation

Using the Model 1700 Menus

3.7.3 Enable the Autofill function

The autofill function must be enabled from the front panel of the instrument. After the autofill control loop has been configured (steps in section 3.7.2) the instrument will be left in the valve manually-closed state.

To enable the autofill function:

1. Press on the **M-CLOSED** icon in the display footer until it reads **AUTOFILL**.
2. Press **SAVE** to enable the autofill function.

The AUTOFILL button has three states:

Table 3-4. Autofill Settings

Function	Operation	Overrides
AUTOFILL	Maintains level between fill START and STOP setpoints.	Autofill will alarm and cease if fill valve stays open for \geq the FILL TIMEOUT setting (nitrogen AUTOFILL CONTROL only).
M-OPEN^a	Energizes the valve control socket on the rear panel.	None
M-CLOSED^b	De-energizes the valve control socket on the rear panel.	

- a. Manual Open
b. Manual Closed

3.8 Select the appropriate units on the display

Touch the units on the display to change the units. The available units are percent (**%**), inches (**in**), and centimeters (**cm**).




Figure 3-34. LN₂ Home Screen

Operation

Using the Model 1700 Menus

3.9 Analog output signals

3.9.1 Connecting to the Aux Connector

Refer to “Aux I/O Connector” on page 96 of the Appendix for a connector pin-out.

3.9.2 Configuring the Analog Outputs

1. From the **MENU** screen, choose **OUTPUTS**, then **Page 2**.
2. If necessary, choose the source for the **0-10 VDC** output and **4-20 mA** output.
3. Press the **SAVE** button to save the choice (or **CANCEL** to quit without making a change).
4. Press the home icon to go back to the home screen.




Figure 3-35. Analog Outputs Source Selection

3.10 Ethernet Connectivity

3.10.1 IP Addressing Scheme

1. From the **MENU** screen, choose **NETWORK**. The current settings will be displayed.
2. To change the settings, choose **EDIT**.




Figure 3-36. Network Selection




Figure 3-37. Editing Network Selection

Operation

Using the Model 1700 Menus

3. In the **ADDRESSING:** button, choose **STATIC**, **DYNAMIC**, or **DISABLED** as appropriate.
4. If **STATIC** is chosen for the addressing scheme, enter **IP ADDRESS**, **NETMASK**, and **GATEWAY** addresses that are appropriate for the connected network. Once an address field is touched, the pop-up keyboard will be presented for data entry. Press **SAVE** when done.
5. If the **DYNAMIC** addressing scheme is chosen, the **IP ADDRESS**, **NETMASK**, and **GATEWAY** addresses will automatically be assigned from a network DHCP server. Press **SAVE** when done.
6. If changes are made, select the **SAVE** button and then the instrument will reboot and reconfigure itself with the chosen settings.




Figure 3-38. Editing Network Settings

Note

Until the instrument is rebooted, the IP Addressing changes will not be applied.




Figure 3-39. Saving Network Settings

3.11 Abnormal Operation

3.11.1 Capacitance Sensor Contamination

To ensure proper instrument calibration and operation, care must be taken to ensure the sensor is kept free of contaminants and not subjected to any force which would physically distort the sensor. Water or other electrically conducting substances in the sensor will disturb the measured capacitance and the instrument's response. Physically distorting the sensor in any way will also cause abnormal instrument operation by introducing variations in the sensor capacitance not due to liquid level. The absolute calibration of the instrument can be inaccurate if care is not taken to ensure the sensor is in the proper environment.

Operation

Using the Model 1700 Menus

Cold sensors exposed to humidified air can show erroneous high level readings because the air contains moisture that can condense between the cold sensing tubes. A minute film of water can cause a shorted or partially shorted condition, which results in false level readings. As the sensor warms, the moisture may evaporate and the sensor will again read correctly. This is a physical phenomenon and does not indicate any problem with your AMI level equipment. Limit or eliminate exposure of cold sensors to humidified air to avoid this condition.

If a sensor should require cleaning and the sensor is for use with liquids other than liquid oxygen, flushing with pure alcohol is recommended. The sensor cannot be used again until all the alcohol has been evaporated. Under no circumstances should the sensor be disassembled.

For sensors to be used with liquid oxygen (LOX), although measures are taken to minimize oils and greases during manufacture, no special cleaning required for LOX service is provided by AMI. Certified LOX cleaning is the responsibility of the customer.

3.11.2 Dirty Helium Sensor Operational Mode

Refer to section 3.5.3 on page 34.

3.11.3 Resetting the Instrument to Factory Defaults

1. Press **MENU**, then **SYSTEM**, then **PAGE 2**.

2. Press **RESET INSTRUMENT TO FACTORY DEFAULTS**.




Figure 3-40. System Menu, Page 2

3. Press **RESET** and the instrument will be reset to factory defaults.




Figure 3-41. Resetting Instrument to Factory Defaults

Operation

Using the Model 1700 Menus

3.12 Shutting the Instrument Down

1. The instrument should be shut down by using the menu function **SHUT DOWN SYSTEM**.

Note

*The Model 1700 Instrument is a Linux-based computer system and in order to ensure the file system is properly unmounted, the **SHUT DOWN SYSTEM** function should be invoked. If it is not, i.e. the instrument is shut down by removing power via the front panel power switch, the next time the instrument boots up, it will have to scan it's memory system to ensure everything is in order.*

2. Choose **YES** to confirm shutdown.
3. When prompted, turn off the front panel power switch.




Figure 3-42. Invoking Instrument Shut Down



Figure 3-43. Confirming Instrument Shut Down

4 Calibration

Model 1700 instrument is calibrated at the factory for a specific length sensor(s) for use in a specific liquid(s). The calibration length(s) and calibration liquid(s) are listed on the calibration sticker on the bottom of the instrument. For capacitance sensors, if the factory calibration method utilized was approximate, the calibration length will be noted as an approximate value.

4.1 Setting the System Date and Time

1. From the home screen, choose **MENU**.




Figure 4-1. Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.




Figure 4-2. Menu Selection Screen

3. From the **SYSTEM CONFIGURATION, PAGE 1** screen, choose **SET**.




Figure 4-3. System Menu, Page 1

Calibration

Capacitance-Based Level Calibration

4. Edit the **YEAR**, **MONTH**, **DAY**, **HOURL**, and **MIN** fields as necessary. Touching in a field will launch the keyboard on the screen. Edit the information in the field as necessary and choose **Enter** to enter the data in the field and close the pop up keyboard.




Figure 4-4. System Menu, Page 2

Note

The clock is set to GMT at the factory and is battery backed. There is no provision in the instrument for automatic Daylight Savings Time correction.

5. Choose **SAVE** in the footer after all the fields have been edited as necessary.

LN2 VERSION

4.2 Capacitance-based Level Calibration

4.2.1 Understanding the Sensor Active Length

American Magnetics, Inc. fabricates the liquid level sensor with two vent holes; a lower vent hole in the side wall near the bottom which is typically the minimum liquid level calibration point and the upper vent hole in the sensor side typically near the top of the sensor. The liquid level location approximately 2.5 cm (1 in) below the upper vent hole is typically the 100% calibration point.

The Model 1700 Instrument requires the user to enter the calibrated, or active length, (physical distance between the Min and Max calibration locations on the sensor) in order for the absolute units function (inches, cm) to be displayed if desired.

Note

Without an active length entered, the instrument will not be able to read out in units other than percent.

Calibration

Capacitance-Based Level Calibration : Relationship Between Cal and Sensor Length




Figure 4-5. Typical Capacitance-based Liquid Level Sensor

The user must enter the sensor length in centimeters. Use the Active Length value noted on the level sensor documentation or measure the distance between the lower vent hole on the sensor and 1.0 inch (2.5 cm) below the upper vent hole on the sensor.

4.2.2 Relationship between Calibration and Sensor Length

The capacitance-based method of measuring the liquid level operates by measuring the frequency of an oscillator, which is contained in the

Calibration

Capacitance-Based Level Calibration : Variation in Dielectric with Density

oscillator/transmitter unit. As the liquid level varies, the value of the capacitance varies proportionally. Since the dielectric properties of liquids vary and the component tolerances for the sensor and oscillator introduce variations, a calibration is required to assure maximum accuracy for a specific sensor immersed in the target liquid. The calibration minimum and maximum settings correspond to the maximum and minimum oscillation frequencies, respectively, for a given sensor and target liquid configuration.

The length setting of the instrument is only provided as a means of scaling the 0% (minimum calibration) to 100% (maximum calibration) range of the measurement to meaningful units of length. During the calibration it is important to accurately measure the distance between the physical locations on the sensor corresponding to the maximum and minimum calibration points. The measured value for the length will be used in configuring the instrument for operation.

4.2.3 Variations in the Dielectric with Changing Density

For cryogenic liquids, the dielectric of the liquid will change with a change in density. The amount of change is dependent on the properties of the specific liquid. Figure 4-6 illustrates the variations in dielectric for nitrogen vs. pressure under *saturated* conditions.¹ Since the instrument uses a capacitance-based method for determining liquid level, such a change in the dielectric of the liquid will result in a shift in the level reading of the instrument. The calibration procedures described herein are most accurate when applied in situations where the operating conditions of the cryo-vessel are relatively constant, i.e. the operating pressure and temperature of the cryo-vessel are relatively constant.

To minimize the effects of shifts in the dielectric of the target liquid, perform a closed dewar calibration (see page 61) at the expected operating condition of the cryo-vessel. If this is not feasible, then calibrate the sensor at atmospheric pressure and use the approximate calibration method to compensate for the shift of the dielectric when the cryogenic liquid is under pressure. For this type of approximate calibration, the *reference liquid* will be the *target liquid* at atmospheric pressure — see page 57 for a detailed discussion of the approximate calibration method. If any questions exist in regard to calibration issues, contact AMI for assistance in determining the optimal calibration strategy.

Note

All references to “dielectric constant” herein refer to the unitless relative dielectric to ϵ_0 (ϵ_0 is the dielectric constant of a vacuum).

1. Data obtained from NIST Standard Reference Database 12.




Figure 4-6. Dielectric vs. pressure for nitrogen under saturated conditions.

4.2.4 Capacitance-based Sensor Calibration Methods

The most straightforward calibration method is the *Open Dewar Calibration* which requires the customer to have access to a filled dewar where the full active length of the sensor can be dipped. The *Closed Dewar Calibration* method can be performed in situations where it is not feasible for the customer to dip the sensor into an open dewar, such as situations where the target liquid is under pressure. The closed dewar calibration is more complex and may require initial preparations to insure success.

Occasionally customers ask AMI to calibrate an instrument and sensor for a liquid which is not available at AMI for calibration purposes and/or for a sensor which is too long to be calibrated at our facilities.

For the case of the target liquid being unavailable, AMI uses liquid nitrogen as the reference liquid and an *Approximate Calibration* is performed using mathematical manipulation of the ratio of the dielectric constants between liquid nitrogen and the desired liquid. This procedure is outlined in the *Approximate Calibration* section beginning on page 57. The technique is intended to provide the instrument with an approximate calibration so that it can be used immediately by the customer. However, the customer is still expected to perform a more accurate calibration where feasible, such as the open dewar or closed dewar calibration, with the target liquid.

Calibration

Open dewar calibration

For the case where a sensor is too long to be calibrated in AMI facilities, AMI will perform a partial length open dewar calibration in liquid nitrogen, and then calculate the MAX calibration point. A dielectric ratio may also be subsequently utilized to adjust for a target liquid other than liquid nitrogen. The customer is expected to perform a more accurate open dewar or closed dewar calibration if feasible.

As a quick guide for selection of the best calibration method available, a calibration selection diagram is presented below. If the instrument and sensor are purchased as a unit from AMI, then the factory calibration will be adequate in most cases. However, for the exceptions noted in the previous paragraphs (which are *approximate* calibrations), the customer should perform a more accurate open dewar or closed dewar calibration. A customer performed calibration is also required for sensors that are purchased as a separate item from the instrument, since the instrument and sensor were not both available for calibration at AMI facilities.

Each Model 1700 Instrument must be calibrated with a sensor before use. If the instrument was purchased with a sensor to be used in liquid nitrogen, the instrument has been calibrated at the factory.

4.2.4.1 Selection of Capacitance Sensor Calibration Methods

As a quick guide for selection of the best calibration method available, a calibration selection diagram is presented in 4-7. If the instrument and sensor are purchased as a unit from AMI, then the factory calibration, including sensor serial number and sensor physical parameter information entered at the factory, will be adequate in most cases. However, for the exceptions noted in the following paragraphs (which are *approximate* calibrations), the customer should perform a more accurate open dewar or closed dewar calibration if at all possible. A customer-performed calibration is also required for sensors that are purchased as a separate item from the instrument, since the instrument and sensor were not both available for calibration at AMI facilities.

The most straightforward calibration method is the *Open Dewar Calibration* which requires the customer to have access to a filled dewar where the full active length of the sensor can be dipped. The *Closed Dewar Calibration* method can be performed in situations where it is not feasible for the customer to dip the sensor into an open dewar, such as situations where the target liquid is under pressure. The closed dewar calibration is more complex and requires initial preparations to insure success.

Occasionally customers ask AMI to calibrate an instrument and Sensor Transmitter for a liquid which is not available at AMI for calibration purposes and/or for a sensor which is too long to be calibrated at our facilities. For the case of the target liquid being unavailable, AMI uses liquid nitrogen as the reference liquid and an *Approximate Calibration* is performed using mathematical manipulation of the ratio of the dielectric

Calibration

Capacitance-Based Level Calibration : Selection Flowchart




Figure 4-7. Calibration method selection diagram.

constants between liquid nitrogen and the desired liquid. The technique is intended to provide the instrument with an approximate calibration so that it can be used immediately by the customer. However, the customer is still expected to perform a more accurate calibration where feasible, such as the open dewar or closed dewar calibration, with the target liquid.

For the case where a sensor is too long to be calibrated in AMI facilities, AMI performs a partial length open dewar calibration in liquid nitrogen, and then calculates the maximum calibration point. A dielectric ratio (i.e. approximate calibration factor) may also be subsequently used to adjust

Calibration

Capacitance-Based Level Calibration : Pre-Calibration Procedure

for a target liquid other than liquid nitrogen. The customer is expected to perform a more accurate open dewar or closed dewar calibration if feasible.

4.2.5 Pre-Calibration Procedure

4.2.5.1 Enter Capacitance Sensor Information

1. Press the **MENU** button in the lower left corner of the display screen.




Figure 4-8. MENU Selection Button

2. Choose the **SENSORS** selection from the **MENU** screen.
3. If the capacitance sensor will be longer than 15 feet from the instrument, ensure an Oscillator/Transmitter is used between the instrument and the sensor.




Figure 4-9. SENSORS Selection Button

4. Verify that the oscillator source selection is correct, either **INTERNAL** or **EXTERNAL**. Note that this setting is auto-detected by the instrument,




Figure 4-10. NITROGEN OSCILLATOR Selection

Calibration

Capacitance-Based Level Calibration : Pre-Calibration Procedure

5. Press the **CALIBRATE NITROGEN** button




Figure 4-11. CALIBRATE NITROGEN Selection Button

6. Touch in the **SENSOR ACTIVE LENGTH** field and using the numerical keypad, enter the sensor active length in cm. Press **Enter** and then **SAVE** at the bottom of the screen when finished




Figure 4-12. SENSOR ACTIVE LENGTH field

7. Press the back button in the screen footer twice to revert back to the **MENU** screen.




Figure 4-13. Footer BACK Button Selection

8. Press the **SYSTEM** button.




Figure 4-14. SYSTEM Menu Selection

Calibration

Capacitance-Based Level Calibration : Pre-Calibration Procedure

9. Touch in the **NITROGEN SENSOR NAME:** field. The keyboard will be displayed. Edit the displayed name and press **ENTER** at the bottom of the screen.




Figure 4-15. PERFORM MAX CAL. Selection Button

10. Press the home icon button in the screen footer to return to the level display screen.
- 11.

4.2.5.2 Performing Loss of Sensor Calibration

The Model 1700 will detect when the liquid level sensor has been disconnected from the instrument and display a notification on the front panel. If the instrument has been configured for autofill, the fill valve will be shut, requiring operator intervention to resume autofill operation. This loss of sensor threshold must be calibrated as follows:




Figure 4-16. Home Selection Button

1. For a system where the internal oscillator is used, connect the coaxial cable to the BNC connector on the instrument rear panel but leave the cable disconnected from the sensor BNC connector. For a system where an external oscillator is used, connect the coaxial cable between the instrument and the oscillator/transmitter unit. Connect the second (6') coaxial cable to the input of the oscillator/transmitter unit but leave the cable disconnected from the sensor BNC connector.
2. Press the **MENU** button in the lower left corner of the display screen.




Figure 4-17. MENU Selection Button

Calibration

Capacitance-Based Level Calibration : Pre-Calibration Procedure

3. Choose the **SENSORS** selection from the **MENU** screen.




Figure 4-18. SENSORS Selection Button

5. From the **SENSORS** Menu, choose **CALIBRATE NITROGEN**.




Figure 4-19. CALIBRATE NITROGEN Selection Button

6. Select the **NO SENSOR CAL** button.
7. When the calibration procedure is completed, choose **SAVE** in the footer of the screen.
8. Press on the home icon to return to the home screen.
9. Connect the BNC cable to the liquid level sensor.




Figure 4-20. NO SENSOR CAL button

Calibration

Capacitance-Based Level Calibration : Open Dewar Cal Procedure

4.2.6 Performing an Open Dewar Calibration

1. Press the **MENU** button in the lower left corner of the display screen.




Figure 4-21. MENU Selection Button

2. Choose the **SENSORS** selection from the **MENU** screen.




Figure 4-22. SENSORS Selection

3. Press the **CALIBRATE NITROGEN** button




Figure 4-23. CALIBRATE NITROGEN Selection Button

4. Position the capacitance sensor in the target liquid at the 100% level. Hold the sensor at this location and press the **PERFORM MAX CAL** button.




Figure 4-24. PERFORM MAX CAL. Selection Button

Calibration

Capacitance-Based Level Calibration : Approximate Calibration Procedure

5. The instrument will display the following screen as it takes data for several seconds. Once the calibration measurement is completed, press the **SAVE** button.




Figure 4-25. Updated MAX CAL Frequency

6. Position the capacitance sensor in the target liquid at the 0% level. Hold the sensor at this level and press the **PERFORM MIN CAL** button.
7. After several seconds of displaying “**CALIBRATING...**”, the instrument will complete the calibration process. Press the **SAVE** button to save the new calibration set point.




Figure 4-26. PERFORM MIN CAL Selection Button

Note

*Note that the frequencies listed to the right of the **PERFORM MAX** and **MIN CAL** buttons are updated as well as the date and time stamps of the calibration points.*

Note

The Minimum and Maximum calibration can be performed in either order. Also, either the Minimum or Maximum calibration point can be updated without altering the other calibration point.

4.2.7 Approximate Calibration

This procedure is the least accurate form of calibration and should be used only when the aforementioned calibration procedures are not viable. The approximate calibration method can be used in cases where the sensor cannot be dipped into the target liquid, the full active length of the sensor cannot be dipped into an open dewar, or both. Approximate calibration may also be useful for situations where the sensor cannot be dipped into the target liquid under the expected operating pressure.

If the target liquid is not available for dipping, a substitute non-conducting reference liquid can be used. If the full length of the

Calibration

Capacitance-Based Level Calibration : Approximate Calibration Procedure

sensor cannot be dipped, then a partial length dip can be performed. If both situations are encountered, then a partial length dip can be performed in a substitute reference liquid.

1. Press the **MENU** button in the lower left corner of the display screen.




Figure 4-27. MENU Selection Button

2. Choose the **SENSORS** selection from the **MENU** screen.




Figure 4-28. SENSORS Selection

3. Press the **CALIBRATE NITROGEN** button




Figure 4-29. CALIBRATE NITROGEN Selection Button

4. Position the capacitance sensor in the target liquid at the 100% level. Hold the sensor at this location and press the **PERFORM MAX CAL** button.




Figure 4-30. PERFORM MAX CAL. Selection Button

Calibration

Capacitance-Based Level Calibration : Approximate Calibration Procedure

5. The instrument will display the following screen as it takes data for several seconds. Once the calibration measurement is completed, press the **SAVE** button.




Figure 4-31. Updated MAX CAL Frequency

6. Position the capacitance sensor in the target liquid at the 0% level. Hold the sensor at this level and press the **PERFORM MIN CAL** button.
7. After several seconds of displaying “**CALIBRATING...**”, the instrument will complete the calibration process. Press the **SAVE** button to save the new calibration set point.




Figure 4-32. PERFORM MIN CAL. Selection Button

Note

*Note that the frequencies listed to the right of the **PERFORM MAX** and **MIN CAL** buttons are updated as well as the date and time stamps of the calibration points.*

Note

The Minimum and Maximum calibration can be performed in either order. Also, either the Minimum or Maximum calibration point can be updated without altering the other calibration point.

8. Measure the distance between the bottom hole of the sensor and the location of the liquid level dipped for max calibration. This measured length is L_{dipped} .
9. The dielectric constant for the reference liquid, e_1 , and the target liquid, e_2 , must be known to complete the approximate calibration. These values must be placed in the equation:

$$\text{Approximate Calibration Factor} = \left[\frac{e_2 - 1}{e_1 - 1} \times 100 \right] \frac{L_{active}}{L_{dipped}}$$

Calibration

Capacitance-Based Level Calibration : Approximate Calibration Procedure

where L_{dipped} is the length of the sensor dipped in the reference liquid and L_{active} is the active sensor length.

Note

If the target liquid is available for dipping (i.e. the reference liquid and target liquid are the same), then the dielectric ratio, $(\epsilon_2-1)/(\epsilon_1-1)$, becomes 1. If the full active length of the sensor can be dipped, then the length ratio, L_{active}/L_{dipped} , becomes 1.

Note that $\epsilon_1 = 1.454$ for liquid nitrogen at -203°C at atmospheric pressure. Dielectric constants for several liquids are provided in the Appendix. The dielectric constant varies with temperature and pressure, therefore for best accuracy use the dielectric constant for the target liquid at the temperature and pressure maintained in the containing vessel.

10. Touch the **APPROX. CAL.** **VALUE:** field and using the numerical keypad that pops up, enter the Approx Cal Value to be applied.




Figure 4-33. APPROX. CAL. VALUE field

11. Press the **APPLY** button and note that the Approximate Calibration value will be used to scale the **MAX** Calibration frequency displayed adjacent to the **PERFORM MAX CAL** button and the value entered into the **APPROX. CAL. VALUE** field will vanish.




Figure 4-34. MAX CAL values edited

Note

The last approximate calibration factor is not retained in the instrument memory, therefore the effects of repeated approximate calibrations are cumulative.

Example: Purchased a 100" active length sensor for operation in liquid argon at atmospheric pressure, however only liquid nitrogen is available for calibration at a maximum depth of 30":

First, the sensor is dipped as far as possible into the liquid nitrogen and cooled. The minimum point is then set as outlined

Calibration

Capacitance-Based Level Calibration : Closed Dewar Calibration Procedure

in step 2. The maximum point is set as outlined in step 3 while the sensor is submerged 30" in liquid nitrogen. The dielectric constant for liquid nitrogen is 1.454 and for liquid argon is 1.53. Substituting all values into the approximate calibration factor equation yields:

$$\text{Approximate Calibration Factor} = \left[\frac{1.53 - 1}{1.454 - 1} \times 100 \right] \frac{100}{30} = 389.1$$

A value of 389.1 would be entered as the approximate calibration factor as outlined in steps 6 and 7. The sensor is now approximately calibrated for 100" active length operation in liquid argon.

12. The sensor can now be installed in the dewar containing the target liquid. The approximate calibration can be used until an open dewar or closed dewar calibration can be performed with the target liquid.

4.2.8 Closed Dewar Calibration

A calibration can be performed in a closed dewar system by monitoring the liquid level while transferring the target liquid to an initially empty (or near empty) dewar at a constant rate. In order to insure success with the closed dewar technique, it is necessary to prepare the instrument by presetting the calibration minimum and maximum calibration points outside the estimated level range. If the instrument is not prepared in this manner before the calibration procedure, it is possible to reach the maximum calibration point of the instrument before the target vessel is at the desired maximum level point. If minimum and maximum liquid level indication is available via some other means (e.g. flow calculation, visual determination, point sensors, etc.), then the presetting of the instrument is not necessary.

4.2.8.1 Presetting the maximum and minimum calibration points

The following procedure should be performed before installation of the sensor in the target cryo-vessel.

1. Connect the sensor coaxial cable to the BNC connector on the rear panel of the instrument (Refer to "Connecting the Sensor(s) to the Instrument" on page 14). *Do not connect the sensor.* Energize the instrument.

Calibration

Capacitance-Based Level Calibration : Closed Dewar Calibration Procedure

2. Press the **MENU** button in the lower left corner of the display screen.




Figure 4-35. MENU Selection Button

3. Choose the **SENSORS** selection from the **MENU** screen.




Figure 4-36. SENSORS Selection

4. Press the **CALIBRATE NITROGEN** button




Figure 4-37. CALIBRATE NITROGEN Selection Button

5. Press the **PERFORM MIN CAL** button.
6. After several seconds of displaying "**CALIBRATING...**", the instrument will complete the calibration process. Press the **SAVE** button to save the new calibration set point.




Figure 4-38. PERFORM MIN CAL. Selection Button

Note

*Note that the frequency listed to the right of the **PERFORM MIN CAL** button is updated as well as the date and time stamp of the calibration point.*

Calibration

Capacitance-Based Level Calibration : Closed Dewar Calibration Procedure

7. Connect the sensor to the oscillator coaxial cable that is connected to the instrument.
8. Perform the maximum level calibration by invoking the maximum calibration function by pressing **MENU > SENSORS > CALIBRATE NITROGEN > PERFORM MAX CAL.**




Figure 4-39. PERFORM MAX CAL. Selection Button

9. Save the calibration value.
10. Calculate the factor C_{adj} using the following equation:

$$C_{adj} = 120 \left[1 + \frac{2.1(L_{active})}{5.2(L_{total})} \right] \left[\frac{e - 1}{0.454} \right]$$

where L_{total} is the total sensor length in inches, L_{active} is the active sensor length in inches, and e is the dielectric constant of the target liquid.

11. Enter C_{adj} into the instrument by touching the **APPROX CAL VALUE:** field.
12. Using the pop up numeric keypad, enter the C_{adj} value and press the **APPLY** button.
13. With the sensor connected, again press the **PERFORM MIN CAL** button. The presetting procedure is complete. Proceed to the remainder of the closed dewar calibration procedure as presented below.




Figure 4-40. PERFORM MAX CAL. Selection Button

4.2.8.2 Completing the closed dewar calibration procedure

1. Install the sensor in the dewar and energize the instrument with the sensor connected to the instrument via the oscillator (if required) and extension cable(s) (see the system diagram on page 2).
2. Connect a strip chart recorder or graphical data logging application to the recorder output terminals on the rear panel of the instrument. If the recorder output is not available, the 4-20 mA current loop output may be used if installed, or an installed

Calibration

Capacitance-Based Level Calibration : Closed Dewar Calibration Procedure

communications option can be used to query the instrument for the liquid level at regular time intervals during the calibration procedure. If no remote monitoring or communication option is installed, the level display must be manually plotted vs. time during the procedure.

3. Refer to “Analog output signals” on page 41 to configure the recorder output or current loop output.
4. Commence filling the dewar. While the sensor is cooling down, there may be a slow drift in the displayed liquid level. However, when the liquid actually touches the bottom of the sensor, contact with the liquid surface may become apparent by virtue of more random and frequent fluctuations in the displayed liquid level. The liquid level trace will also start to show an increasing profile with positive slope.

Once the indications of the contact between the sensor and liquid become readily apparent, use the **PERFORM MIN CAL** procedure below to reset the new minimum calibration point. This point is the 0% level of the sensor when the **PERFORM MIN CAL** procedure is finished and saved becomes the 0% level.

Note

If the sensor is installed in the dewar with some small amount of liquid already in contact with the sensor, then the final minimum calibration point can be set before filling begins but after any thermally induced fluctuations in the observed output have diminished. However, note that the measured span of the liquid level is reduced by the initial level of liquid in contact with the sensor.

5. Perform the minimum level calibration by invoking the minimum calibration function by pressing **MENU > SENSORS > CALIBRATE NITROGEN > PERFORM MIN CAL**.
6. After several seconds of displaying “**CALIBRATING...**”, the instrument will complete the calibration process. Press the **SAVE** button to save the new calibration set point.




Figure 4-41. PERFORM MAX CAL.
Selection Button

Calibration

Capacitance-Based Level Calibration : Closed Dewar Calibration Procedure

Note

*Note that the frequency listed to the right of the **PERFORM MIN CAL** button is updated as well as the date and time stamp of the calibration point.*

- Continue the transfer while observing the liquid level trace on the strip chart recorder or computer display, whose slope is proportional to the transfer rate. The slope of the liquid level trace should decrease significantly when the liquid reaches the hole in the top of the sensor.

When the break in the slope of the level trace occurs (i.e. the slope of the level trace becomes 0 or horizontal), perform a **PERFORM MAX CAL** procedure below. The level on the sensor when the **PERFORM MAX CAL** procedure is finished and saved becomes the 100% level.

Note

If the instrument displayed a 100% reading before a break is observed in the slope of the level trace, then the maximum calibration point set prior to the current procedure has interfered. If this occurs, the customer has two options: 1) stop the procedure, repeatedly enter a value of 120 for C_{adj} (see steps 4 and 5 of the presetting procedure) until the current liquid level display falls below 100%, and then continue the procedure; or 2) continue the liquid transfer until the liquid level is determined to be 100% by means other than feedback from the instrument and then performing the maximum calibration procedure.

- Perform the maximum level calibration by invoking the maximum calibration function by pressing **MENU > SENSORS > CALIBRATE NITROGEN > PERFORM MAX CAL**.




Figure 4-42. PERFORM MAX CAL.
Selection Button

Calibration

Capacitance-Based Level Calibration : Closed Dewar Calibration Procedure

9. The instrument will display the following screen as it takes data for several seconds. Once the calibration measurement is completed, press the **SAVE** button.

Note

*Note that the frequency listed to the right of the **PERFORM MAX CAL** button is updated as well as the date and time stamp of the calibration point.*




Figure 4-43. Updated MAX CAL Frequency

10. To achieve a standard calibration of the sensor with the active region located from the lower hole to one inch below the upper hole, use the level data from the instrument to recalibrate the maximum point when the percent level corresponds to one inch below the upper hole. Use the following equation to determine the percent level at which to reset the maximum calibration point:

$$MAX_{percent} = 100 - 100 \left[\frac{1}{L_{active}} \right]$$

where L_{active} is the active length of the sensor in inches. This technique can be used assuming the sensor was built as a standard sensor. If the sensor was made in a custom configuration, refer to the sensor documentation and/or drawing or contact AMI.

Example: 20" active length sensor:

When the sensor is calibrated by the closed dewar procedure, the actual length of calibration will be 21" (distance between the bottom and top holes in the sensor). When the liquid is 1" below the upper hole, the display will show 95.2% [e.g. $100\% - (1"/21" \times 100\%)$]. When the liquid level reaches this point during usage, perform the **PERFORM MAX CAL** operation. The instrument and sensor are now calibrated with a standard active region of 20". The length setting of the sensor in the instrument should also be configured for 50.8 cm (20").

Proceed to the *Operation* section for directions for configuring the instrument.

Calibration

Superconductivity-Based (LHe) Calibration

4.3 Superconductivity-Based (Liquid Helium) Level Calibration

The instrument has been calibrated for Liquid helium Sensors at the AMI facility. No further calibration is needed.

However, to have the Model 1700 Instrument work correctly, at a minimum, the user must enter the correct active length for the sensor.

Note

If the instrument was purchased with a helium level sensor, the active length will be set prior to shipping.

4.3.1 Verify the Liquid Helium Sensor Type

1. Press the **MENU** button in the lower left corner of the display screen.




Figure 4-44. Home Screen

2. Choose the **SENSORS** button from the **MENU** screen.




Figure 4-45. MENU Selection Button

3. Verify that the **HELIUM SENSOR TYPE:** (either **4.2K** or **2K**) is correct for the sensor to be used with the instrument. If not, contact the factory for assistance.




Figure 4-46. SENSORS Selection

Calibration

Superconductivity-Based (LHe) Calibration : Sample Interval

4.3.2 Sensor Sample Interval

1. Press the **MENU** button in the lower left corner of the display screen.




Figure 4-47. Home Screen

2. Choose the **SENSORS** button from the **MENU** screen.




Figure 4-48. MENU Selection Button

3. Choose the **CALIBRATE HELIUM** selection.




Figure 4-49. SENSORS Selection

4. Press in the **SENSOR SAMPLE INTERVAL** field. A pop up numeric keypad will be launched.
5. Enter the desired sample interval time in minutes and press **Enter** and then **SAVE**.




Figure 4-50. Enter Liquid Helium Sensor Active Length

Calibration

Superconductivity-Based (LHe) Calibration : Sensor Active Length

4.3.3 Sensor Active Length

1. Press in the **SENSOR ACTIVE LENGTH** field. A pop up numeric keypad will be launched.
2. Enter the **SENSOR ACTIVE LENGTH** in centimeters and press **Enter** and then **SAVE**.




Figure 4-51. Enter Liquid Helium Sensor Active Length

4.3.4 Continuous Measure Time Limit

The Continuous measure time limit feature sets a maximum time that the instrument will keep the liquid helium level sensor energized. If the instrument is left in the **CONTINUOUS** mode for the **CONTIN. MEASURE TIME LIMIT** interval, the instrument will revert back to the Sample and Hold mode to prevent excessive liquid helium boil off.

1. Press in the **CONTIN. MEASURE TIME LIMIT** field. A pop up numeric keypad will be launched.
2. Enter the maximum time in minutes that the sensor should remain energized in the continuous mode.
3. Press **ENTER** on the keypad and then **SAVE** at the bottom of the screen.
4. Press the back button in the screen footer twice to revert back to the **MENU** screen.




Figure 4-52. Enter Liquid Helium Sensor Active Length




Figure 4-53. Footer BACK Button Selection

Calibration

Superconductivity-Based (LHe) Calibration : Sensor Name

4.3.5 Sensor Name

1. Press the **SYSTEM** button.




Figure 4-54. SYSTEM Menu Selection

2. Touch in the **HELIUM SENSOR NAME:** field. The keyboard will be displayed. Edit the displayed name and press **ENTER** and then **SAVE** at the bottom of the screen.




Figure 4-55. HELIUM SENSOR NAME: Field

3. Press the home icon button in the screen footer to return to the level display screen.




Figure 4-56. Home Selection Button

5 Remote Interface Reference

The Model 1700 Instrument provides both serial (RS-232) and Ethernet interfaces as standard features. The serial and Ethernet interfaces can be operated simultaneously. Separate output buffers are also provided for the serial and Ethernet return data. However, for optimal performance and simplicity of programming, AMI recommends limiting normal operation to one interface. An exception to this recommendation would be using the serial port as a debugging aid during programming of the Ethernet port, or vice-versa, which can prove to be a useful resource.

The Model 1700 also allows a browser connection via TCP/IP. The instrument's IP address can be ascertained by referring to the section titled "IP Addressing Scheme" on page 41. By using a browser to connect to the instrument, all functionality of the Model 1700 can be controlled via the browser.

5.1 SCPI Command Summary

The following manual conventions are used for SCPI (*Standard Commands for Programmable Instruments*) syntax for the remote interface commands:

- Braces { } enclose valid parameter choices.
- A vertical bar | separates multiple choices for each parameter.

For example, the command `CONFigure:TIME:DST {0|1}` indicates that the command `CONFigure:TIME:DST` has two parameter options: 0 or 1. The following section is a reference list of SCPI commands. Refer to the detailed description of each command for information regarding specific parameter choices and their meanings. Capitalized portions of the commands indicate acceptable abbreviations. Default settings are shown in bold.

5.2 SCPI Ethernet Communication

The Ethernet port via an RJ-45 connector on the rear of the instrument allows a computer to communicate with the instrument using the SCPI commands described in "Command Set Reference" on page 73. The host computer must run a Telnet program, such as PuTTY¹, and connect to port 7180.

1. <http://www.chiark.greenend.org.uk/~sgtatham/putty/>

Remote Interface Reference

5.3 SCPI Serial (RS-232) Communication

An RS-232 serial communication port is available as a 9-pin D-type connector on the rear panel of the instrument for serial communication function.

5.3.1 Serial port connector and cabling

An PC-compatible computer's serial port can be directly connected to the Model 1700 via a standard cable. Refer to the computer's documentation to determine which serial ports are available on a computer and the required connector type.

The Model 1700 uses three wires of the rear-panel DB25 connector: pin 2 (transmit), pin 3 (receive), and pin 7 (common). There is no software or hardware handshaking. The Model 1700 is classified as a DCE (Data Communication Equipment) device since it transmits data on pin 3 and receives data on pin 2. The instrument to which the Model 1700 is attached must do the opposite, i.e., transmit on pin 2 and receive on pin 3 (the requirements for a DTE, or Data Terminal Equipment device). If a serial-to-parallel converter is used, it must be capable of receiving data on pin 3 or the cable connected to the Model 1700 must interchange the wires between pins 2 and 3.

5.4 Command/return termination characters

All commands are transmitted and received as ASCII values and are case insensitive. The Model 1700 always transmits `<CR><LF>` (i.e. a *carriage return* followed by a *linefeed*) at the end of an serial transmission. The Model 1700 can accept `<CR>`, `<LF>`, `<CR><LF>`, or `<LF><CR>` as termination characters from an external computer.

The simplest method for communicating with the Model 1700 via RS-232 is by using the interactive mode of a commercially available terminal emulation program. The Model 1700 transmits and receives information at a baud rate of 115,200 and uses 8 data bits, no parity, and 1 stop bit. When the Model 1700 receives a terminated ASCII string, it always sends back a reply as soon as the string is processed. *When sending commands to the Model 1700, you must wait for the reply from the Model 1700 before sending another command even if the reply consists of only termination characters.* Otherwise, the shared input/output command buffer of the Model 1700 may become corrupted.

The host device can be a standard dot matrix printer connected via a serial-to-parallel converter, or connected directly with a printer capable of receiving serial data. Presumably, any serial-to-parallel converter which can be properly configured is acceptable. AMI has tested the Model 1700

Remote Interface Reference

Serial Command Set Reference

with a standard, low cost converter configured as a DTE device, 115,200 baud, 8 data bits, no parity, and 1 stop bit. In order to communicate with the host device, it is necessary to set the terminal program to the identical baud rate of the host device.

5.5 Command Set Reference

All commands sent to the Model 1700 are processed and the Model 1700 responds with a return value (if applicable) and termination. If the command is invalid, the Model 1700 will respond with an error code (see the *Error Codes* section). All return values including error codes are terminated with `<CR><LF>` (i.e. a *carriage return* followed by a *linefeed*). For those commands that do not return a value, the Model 1700 will return the `<CR><LF>` termination only.

The remote units settings are saved in non-volatile memory and are restored at power-up.

The Model 1700 instrument may be configured for reading liquid nitrogen, liquid helium, or both. Some commands will not be applicable if the instrument is not configured for certain level measurement.

5.5.1 Commands for determining the instrument configuration

Command: N2?

Function: Returns a 0 if the instrument is not configured to read liquid nitrogen level, a 1 if it is with the internal oscillator, and a 2 if it is with an external oscillator/transmitter.

Returns: 0, 1 or 2<CR><LF>

Default: N/A

Command: HE?

Function: Returns a 0 if the instrument is not configured to read liquid helium level, a 1 if instrument is configured to read 4.2K liquid helium level for sensors of active length ≤ 40 inches, 2 if instrument is configured to read 4.2K liquid helium level for sensors of active length ≤ 80 inches, a 3 if the instrument is configured to read 2K liquid helium level for sensors of active length ≤ 40 inches, 5 if instrument is configured to read 2K liquid helium level for sensors of active length ≤ 80 inches.

Returns: 0, 1, 2, 3, 4 or 5<CR><LF>

Default: N/A

Remote Interface Reference

Command Set Reference : Instrument Configuration Queries

Command: DISPLAY:N2?

Function: Returns a 0 if the instrument is not configured to display liquid nitrogen level on the home screen and a 1 if it is.

Returns: 0 or 1<CR><LF>

Default: N/A

Command: DISPLAY:HE?

Function: Returns a 0 if the instrument is not configured to display liquid helium level on the home screen and a 1 if it is.

Returns: 0 or 1<CR><LF>

Default: : N/A

Command: N2:UNIT?

Function: Returns the current liquid nitrogen level units in use.

Returns: C, I, or P<CR><LF>

Default: P

Command: HE:UNIT?

Function: Returns the current liquid helium level units in use.

Returns: C, I, or P<CR><LF>

Default: P

Command: RELAY1:CHannel?

Function: Returns a 0 if relay №1 is disabled, a 1 if the relay is assigned to the nitrogen channel, and a 2 if the relay is assigned to the helium channel.

Returns: 0, 1, or 2<CR><LF>

Default: 0

Command: RELAY2:CHannel?

Function: Returns a 0 if relay №2 is disabled, a 1 if the relay is assigned to the nitrogen channel, and a 2 if the relay is assigned to the helium channel.

Returns: 0, 1, or 2<CR><LF>

Default: 0

Command: FILL:CHannel?

Function: Returns a 0 if the auto fill relay is disabled, a 1 if the relay is assigned to the nitrogen channel, and a 2 if the relay is assigned to the helium channel.

Returns: 0, 1, or 2<CR><LF>

Default: 0

Remote Interface Reference

Command Set Reference : Instrument Configuration Queries

Command: RELAY1:Operation?

Function: Returns a 0 if relay №1 closes (alarms) when the level is above the setpoint and a 1 if the relay closes (alarms) when the relay is below the setpoint. By default, relay №1 is configured as the high level relay with alarm condition when level is greater than the setpoint.

Returns: 0 or 1<CR><LF>

Default: 0

Command: RELAY2:Operation?

Function: Returns a 0 if relay №2 closes (alarms) when the level is below the setpoint and a 1 if the relay closes (alarms) when the relay is above the setpoint. By default, relay №2 is configured as the low level relay with alarm condition when level is less than the setpoint.

Returns: 0 or 1<CR><LF>

Default: 0

Command: RELAY1:SETpoint?

Function: Returns the Relay №1 setpoint in the current units

Returns: <value><CR><LF>

Default: 90%

Command: RELAY2:SETpoint?

Function: Returns the Relay №2 setpoint in the current units

Returns: <value><CR><LF>

Default: 10%

Command: A

Function: Returns the A setpoint limit (auto fill stop level) in the current units

Returns: <value><CR><LF>

Default: 80%

Command: B

Function: Returns the A setpoint limit (auto fill start level) in the current units

Returns: <value><CR><LF>

Default: 20%

Remote Interface Reference

Command Set Reference : Instrument Configuration Queries

Command: N2:INTERVAL?

Function: Returns the fill timer setting in minutes if the instrument is configured for the nitrogen channel.

Returns: <value><CR><LF>

Default: 15

Command: HE:INTERVAL?

Function: Returns the sampling interval in minutes if the instrument is configured for the helium channel.

Returns: <value><CR><LF>

Default: 15

Command: DIRTY_SEN_MODE?

Function: Returns a 0 if the instrument's dirty sensor mode is disabled for the helium channel and a 1 if it is enabled.

Returns: <value><CR><LF>

Default: 0

Command: SOURCE:REC_OUT?

Function: Returns a 0 if the 0-10 V_{DC} Recorder Output is disabled, a 1 if it is configured for the nitrogen channel, and a 2 if it is configured for the helium channel.

Returns: 0, 1, or 2<CR><LF>

Default: 0

Command: SOURCE:CURREnt_LOOP?

Function: Returns a 0 if the 4-20 mA Current Loop Output is disabled, a 1 if it is configured for the nitrogen channel, and a 2 if it is configured for the helium channel.

Returns: 0, 1, or 2<CR><LF>

Default: 0

Command: NAME:SENSor:N2?

Function: Returns the name of the nitrogen level sensor.

Returns: <string><CR><LF>

Default: Nitrogen Level

Command: NAME:SENSor:HE?

Function: Returns the name of the helium level sensor.

Returns: <string><CR><LF>

Default: Helium Level

Remote Interface Reference

Command Set Reference : Instrument Configuration Queries

Command: SERIAL_NUMber?

Function: Returns the serial number of the instrument.

Returns: <string><CR><LF>

Default: N/A

Command: HardWare_VERSION?

Function: Returns the hardware version of the instrument.

Returns: <string><CR><LF>

Default: N/A

Command: DATE_MANUFACTure?

Function: Returns the date of manufacture of the instrument.

Returns: <string><CR><LF>

Default: N/A

Command: FirmWare_VERSION?

Function: Returns the firmware version of the instrument.

Returns: <string><CR><LF>

Default: N/A

Command: HElium_PowerSupply?

Function: Returns the helium power supply configuration of the instrument where 0 is not installed, 1 is standard, and 2 is XL version.

Returns: 0, 1, or 2<CR><LF>

Default: N/A

Remote Interface Reference

Command Set Reference : Setting Measurement Units

5.5.2 Commands for setting the units of measurement

Command: `CONFigure:N2:UNIT <value>`

Function: Sets the liquid nitrogen level units of measurement to percent (0, or PERCENT), centimeters (1, or CM), or inches (2, or INCH)

Returns: `<CR><LF>`

Default: N/A

Command: `CONFigure:HE:UNIT <value>`

Function: Sets the liquid helium level units of measurement to percent (0, or PERCENT), centimeters (1, or CM), or inches (2, or INCH)

Returns: `<CR><LF>`

Default: N/A

Remote Interface Reference

Command Set Reference : Configuring Setpoints

5.5.3 Commands for configuring setpoints

Command: `CONFigure:RELAY1:CHannel <value>`

Function: Assigns relay 1 to either no channel (disabled) (0), nitrogen (1), or helium (2).

Returns: `<CR><LF>`

Default: 0

Command: `CONFigure:RELAY2:CHannel <value>`

Function: Assigns relay 1 to either no channel (disabled) (0), nitrogen (1), or helium (2).

Returns: `<CR><LF>`

Default: 0

Command: `CONFigure:RELAY1:OPeration <value>`

Function: Configures relay №1 such that it closes (alarms) when the level is \leq the setpoint (0) or \geq the setpoint (1).

Returns: `<CR><LF>`

Default: 1

Command: `CONFigure:RELAY2:OPeration <value>`

Function: Configures relay №2 such that it closes (alarms) when the level is \leq the setpoint (0) or \geq the setpoint (1).

Returns: `<CR><LF>`

Default: 0

Command: `CONFigure:FILL:CHannel <value>`

Function: Assigns the auto fill control relay to either no channel (disabled) (0), nitrogen (1), or helium (2).

Returns: `<CR><LF>`

Default: 0

Note: If the relay is not assigned to a channel, the units are assumed to be percent.

Command: `CONFigure:RELAY1:SETpoint <value>`

Function: Configures the relay №1 trip setpoint in the current channel's units.

Returns: `<CR><LF>`

Default: 20

Note: If the relay is not assigned to a channel, the units are assumed to be percent.

Remote Interface Reference

Command Set Reference : Configuring Setpoints

Command: `CONFigure:RELAY2:SETpoint <value>`
Function: Configures the relay №2 trip setpoint in the current channel's units.
Returns: `<CR><LF>`
Default: 80

Command: `CONFigure:FILL:A <value>`
Function: Sets the A setpoint (control band upper limit).
Returns: `<CR><LF>`
Default: 60
Note: The A setpoint (autofill stop) must be greater than the B setpoint (fill start) and must also be between 0 and 100%.

Command: `CONFigure:FILL:B <value>`
Function: Sets the B setpoint (control band lower limit).
Returns: `<CR><LF>`
Default: 40
Note: The B setpoint (autofill start) must be less than the A setpoint (fill stop) and must also be between 0 and 100%.

Command: `CONFigure:INTerval:FILL <value>`
Function: Sets the fill timer in minutes.
Returns: `<CR><LF>`
Default: 15
Note: Setting the value of FILL to 0 disables the fill timer function.

Command: `CONFigure:INTerval:SAMPle <value>`
Function: Sets the liquid helium sampling interval in minutes.
Returns: `<CR><LF>`
Default: 15

Command: `CONFigure:HE:TIME_LIMIT <value>`
Function: Sets the continuous measurement time limit for liquid helium measurements in minutes.
Returns: `<CR><LF>`
Default: 15

Command: `CONFigure:SOURCE:REC_OUT <value>`
Function: Configures the 0-10 Vdc Recorder Output's source to disabled (0), assigned to the nitrogen channel (1), or the helium channel (2).
Returns: `<CR><LF>`
Default: 0

Remote Interface Reference

Command Set Reference : Configuring Setpoints

Command: `CONFigure:SOURCE:CURRENT_LOOP <value>`

Function: Configures the 4-20 mA Current Loop's output source to disabled (0), assigned to the nitrogen channel (1), or the helium channel (2).

Returns: `<CR><LF>`

Default: 0

Remote Interface Reference

Command Set Reference : Instrument Channel Identifiers

5.5.4 Commands for setting the channel identifiers

Command: `CONFigure:NAME:SENSOR:N2=<"string">`

Function: Sets the name of the nitrogen level sensor.

Returns: `<CR><LF>`

Default: Nitrogen Level

Command: `CONFigure:NAME:SENSOR:HE=<"string">`

Function: Sets the name of the helium level sensor.

Returns: `<CR><LF>`

Default: Helium Level

Remote Interface Reference

Command Set Reference : Measuring Level

5.5.5 Commands for making liquid level measurements

Command: MEASure:N2:LEVel?

Function: Returns the liquid nitrogen level in the current units.

Returns: <value><CR><LF>

Default: N/A

Command: MEASure:N2:PERIod?

Function: Returns the liquid nitrogen level measurement period in microseconds.

Returns: <value><CR><LF>

Default: N/A

Command: MEASure:HE:LEVel?

Function: Returns the liquid helium level in the current units.

Returns: <value><CR><LF>

Default: N/A

Note: If in sample and hold mode, the value returned will be the last sample taken, not a current reading.

Command: MEASure:HE:HOLD

Function: Changes liquid helium level measurement operation from Continuous to Sample and Hold.

Returns: HE:HOLDING<CR><LF>

Default: N/A

Command: MEASure:HE:CONTInuous

Function: Changes liquid helium level measurement operation from Sample and Hold to Continuous mode.

Returns: HE:SAMPLING<CR><LF>

Default: N/A

Command: MEASure:HE:SAMPlE

Function: Energizes the liquid helium level sensor, makes a reading and returns to Sample and Hold mode.

Returns: HE:SAMPLED<CR><LF>

Default: N/A

Command: MEASure:ADC0?

Function: Returns the liquid helium level sensor voltage in volts.

Returns: <value><CR><LF>

Default: N/A

Note: Will return a value of 0 if the helium level sensor is not energized when the command is issued.

Remote Interface Reference

Command Set Reference : Measuring Level

Command: MEASure:ADC1?

Function: Returns the liquid helium sensor power supply voltage in volts.

Returns: <value><CR><LF>

Default: N/A

Command: MEASure:ADC2?

Function: Returns the liquid helium sensor excitation current in milliamperes.

Returns: <value><CR><LF>

Default: N/A

Remote Interface Reference

Command Set Reference : Calibration Functions

5.5.6 Commands for calibrating level sensors

Command: MINCAL

Function: Performs a minimum calibration point calibration.

Returns: <CR><LF>

Default: N/A

Command: MAXCAL

Function: Performs a maximum calibration point calibration.

Returns: <CR><LF>

Default: N/A

Command: CONFigure:N2:LENGTH=<value>

Function: Configures the liquid nitrogen sensor active length in current units.

Returns: <CR><LF>

Default: N/A

Note: Returns -5 if the current units are percent.

Command: CONFigure:HE:LENGTH=<value>

Function: Configures the liquid helium sensor active length in current units.

Returns: <CR><LF>

Default: N/A

Note: Returns -5 if the current units are percent.

Command: NOSENSCAL

Function: Calibrates the loss of sensor point for the LN₂ measurement function.

Returns: <CR><LF>

Default: N/A

Remote Interface Reference

Command Set Reference : System Commands

5.5.7 System Commands

Command: `SYStem:BEEPer:IMMediate,<time>`

Function: The receipt of this command causes an audible tone to be generated by the instrument. Note that this command generates an event and therefore it has no associated `SYS:REBOOT` state or query form. The duration time parameter is specified in seconds.

Returns: `<CR><LF>`

Default: N/A

Command: `SYStem:BEEPer:STATe <Boolean>`

Function: Enables/disables the beeper. When `STATe 0` is selected, no instrument condition, except the `:SYStem:BEEPer:IMMediate` command, shall cause an audible beep to be emitted. At `SYS:REBOOT`, this value is reset to 1.

Returns: `<CR><LF>`

Default: 1

Command: `SYStem:DATE <year>,<month>,<day>`

Function: `<year>` as a four-digit number; `<month>` Range is 1 to 12 inclusive; `<day>` Number of days in the month.

Returns: `<CR><LF>`

Default: N/A

Command: `SYStem:DATE?`

Function: Returns the system date.

Returns: `<year>,<month>,<day><CR><LF>`

Default: N/A

Command: `SYStem:KLOCK<Boolean>`

Function: This command locks the local controls of an instrument. This includes any front panel, keyboard, or other local interfaces. This value cannot be reset to OFF, unless `SYStem:SECurity:STATe` is OFF. If `SYStem:SECurity:STATe` is OFF, the `KLOCK` value is set to OFF at `SYS:REBOOT`. If `SYStem:SECurity:STATe` is ON, `SYS:REBOOT` has no effect.

Returns: `<CR><LF>`

Default: OFF

Command: `SYStem:TIME <hour>,<minute>,<second>`

Function: This command is used to set the instrument's clock:

`<hour>` Range is 0 to 23 inclusive.

`<minute>` Range is 0 to 59 inclusive.

Remote Interface Reference

Command Set Reference : Calibration Functions

<second> Range is 0 to 60.

The query response message shall consist of three fields separated by commas: <hour>,<minute>,<second>

Returns: <CR><LF>

Default: N/A

Command: SYSTem:TIME?

Function: This command reads the instrument's clock.

Returns: <hour>,<minute>,<second><CR><LF>

Default: N/A

Command: SYSTem:REBOOT

Function: This command reboots the instrument.

Returns: <CR><LF>

Default: N/A

Command: SYSTem:RESTORE

Function: This command reboots the instrument and sets all parameters back to factory defaults.

Returns: <CR><LF>

Default: N/A

5.6 Error Codes

The Model 1700 returns specific error codes for invalid commands and/or arguments. If an error condition is returned, the command is not processed and the configuration of the instrument is not modified. The table below provides a list of error codes, their meaning, and any associated limits.

Error Code	Meaning	Valid Range
-1	LO (or relay №2 ^a) setpoint out of range	$0 \leq \text{LO (or relay №2}^a) \leq \text{LENGTH}$
-2	Fill B setpoint (fill start) out of range	$0 \leq B < A$
-3	Fill A setpoint (fill stop) out of range	$B < A \leq \text{LENGTH}$
-4	HI (or relay №1 ^a) setpoint out of range	$0 \leq \text{HI (or relay №1}^a) \leq \text{LENGTH}$
-5	Attempted to set or query for LENGTH in PERCENT units mode	

Remote Interface Reference

Command Set Reference : Calibration Functions

Error Code	Meaning	Valid Range
-6	Invalid argument, value out of maximum calibration range	$1 \text{ cm} \leq \text{value} \leq 650 \text{ cm}$
-7	INTERVAL setting out of range	$0 \leq \text{INTERVAL} \leq 600 \text{ min}$
-8	Unrecognized command	
-9	Invalid argument, value was negative or non-numeric	
-10	Approximate calibration <i>factor</i> out of range	$0.1 \leq \text{factor} \leq 999.9$
-11	Command exceeds SCPI input buffer limit	256 characters, including spaces, etc.

- a. Applies to dual instrument configuration

6 Service and Repair

6.1 Cleaning

To prevent electrical shock, disconnect the instrument from AC mains power and disconnect all connected wiring before cleaning. Clean the outside of the instrument using a soft, lint-free, cloth slightly dampened with water.

Do not use detergent or solvents.

Do not attempt internal cleaning.

6.2 User Replaceable Parts

Replacement parts for the instrument are listed in the table below.

AMI Part Number	Description
HG0128	Instrument foot
SA 1045	Single Rack Mount Kit
SA 1046	Dual Rack Mount Kit
EF1700	Fuse, 3 A, 250 Vac, 5x20 mm, fast acting, UL/CSA recognized.
HG0005	Battery, 3V lithium, 20mm x 3.2 mm coin cell; CR2032.

6.3 Battery Replacement

This section describes the procedure for replacing the battery on the instrument's main circuit board.

Warning

This procedure should only be performed by a technician who is familiar with electronic instrumentation and trained in electrical safety and ESD precautions. Always disconnect the power cord and any external wiring before removing the instrument cover.

Always disconnect all inputs, cords, and cables before disassembling the instrument.

6.3.1 Tools Required

- Torx Plus (T.M.) size 10 driver (Wera 028034 or equivalent)
- Torx Plus (T.M.) size 15 driver (Wera 028035 or equivalent)
- Small, flat-blade screw driver (for prying)

6.3.2 Procedure

1. Unplug the instrument from the AC power source.
2. Using the T-15 driver, remove the four 8-32 machine screws on the sides of the instrument cover. Set these screws aside as they will be re-used.
3. Using the T-10 driver, remove the four 6-32 machine screws on the rear of the instrument cover. Set these screws aside as they will be re-used.
4. Lift the instrument cover off of the instrument chassis and set aside.
5. Using the small, flat-blade screwdriver, carefully pry the battery from the holder BH1.
6. Install the new battery into the battery holder BH1.
7. Replace the top cover and secure using the eight machine screws which were removed previously.

6.4 Fuse Replacement

This section describes the procedure for replacing the two fuses on the instrument's main circuit board.

Warning

This procedure should only be performed by a technician who is familiar with electronic instrumentation and trained in electrical safety and ESD precautions. Always disconnect the power cord and any external wiring before removing the instrument cover.

Always disconnect all inputs, cords, and cables before disassembling the instrument.

6.4.1 Tools Required

- Torx Plus (T.M.) size 10 driver (Wera 028034 or equivalent)
- Torx Plus (T.M.) size 15 driver (Wera 028035 or equivalent)
- Small, flat-blade screw driver (for prying)

6.4.2 Procedure

1. Unplug the instrument from the AC power source.
2. Using the T-15 driver, remove the four 8-32 machine screws on the sides of the instrument cover. Set these screws aside as they will be re-used.
3. Using the T-10 driver, remove the four 6-32 machine screws on the rear of the instrument cover. Set these screws aside as they will be re-used.
4. Lift the instrument cover off of the instrument chassis and set aside.
5. Using the small, flat-blade screwdriver, carefully pry the fuse(s) from the fuse holders F1 and/or F2.
6. Install the new fuse(s) into the fuse holder(s) F1 and/or F2.
7. Replace the top cover and secure using the eight machine screws which were removed previously.

Appendix

A.1 Connector Wiring

The following sections document the connector pin outs and pin definitions.

A.1.1 Serial (RS-232) Connector




Figure A-1 Serial (RS-232) Pin Out

The RS-232 connector is a 9-pin D-sub female connector to connect standard DTE 9-pin D-sub male connector using a standard straight (not NULL) cable.

Table A-1 Serial (RS-232) Pin Definitions

Pin	Mnemonic	Function
1	N/C	
2	TXD	Transmit Data
3	RXD	Receive Data
4	N/C	
5	GND	Signal Ground
6	N/C	
7	CTS	Clear to Send
8	RTS	Request to Send
9	N/C	

A.1.2 Ethernet Connector




Figure A-2 Ethernet Connector Socket Pin out

Table A-2 Ethernet RJ-45 Connector Pin Definitions

Pin	Mnemonic	Function
1	TXD+	Transmit differential output +
2	TXD-	Transmit differential output -
3	RXD+	Transmit differential input +
4	not used	
5		
6	RXD	Transmit differential input -
7	not used	
8		

Appendix

Liquid Helium Connector J1 Wiring

A.1.3 Liquid Helium Connector J1 Wiring




Figure A-3 Liquid Helium (LHe) Level Sensor Connector

Table A-3 LHe Level Connector Pin Definitions

Pin	Function
1	Sensor I+ (Red)
2	Not used
3	Not used
4	Not used
5	Not used
6	Sensor V- (Yellow)
7	Sensor I- (Black)
8	Sensor V+ (Blue)
9	Not used

Appendix

Aux I/O Connector Pin Outs

A.1.4 Aux I/O Connector




Figure A-4 Aux I/O Connector

Table A-4 Aux I/O Pin Definitions

Pin	Function	Polarity
1	4-20 mA Current Loop	+
2		-
3	0-10 VDC Output	+
4		-
5	Relay № 2 Dry Contact	N/A
6		
7	Relay № 1 Dry Contact	
8		
9	External Reset	

A.2 Troubleshooting

The following paragraphs serve as an aid to assist the user in troubleshooting a potential problem with the Model 1700 Instrument. If the user is not comfortable in troubleshooting the system, contact an AMI Technical Support.

If the cause of the problem cannot be located, contact an AMI Technical Support Representative at +1 (865) 482-1056 for assistance. The AMI technical support group may also be reached by internet e-mail at:

support@americanmagnetics.com

A.3 Firmware Licenses

The Model 1700 firmware is based on a distribution of Debian Linux, with modifications to the Linux kernel by Technologic Systems and AMI, and additional user interface components by AMI. Some components of this firmware are licensed under agreements that require AMI to make source code available to interested parties. Other components require explicit acknowledgment of the authorship/ownership of the firmware and/or the terms under which it is licensed. In particular:

The Linux kernel version 2.6.34 is licensed under the GNU Public License, version 2. Source code for the version of the Linux kernel used in the Model 1700 is available from the Technologic Systems github repository at <https://github.com/embeddedarm/linux-2.6.34-ts471x.git>

AMI's modifications to that Linux kernel sources are available from: <http://firmware.americanmagnetics.com/1700/kernel-patches.tar>

The Model 1700 uses the "lighttpd" web server, available in source code form from <https://www.lighttpd.net>. It is made available under the following license:

Copyright (c) 2004, Jan Kneschke, incremental
All rights reserved.

Redistribution and use in source and binary forms, with or without modification, are permitted provided that the following conditions are met:

- Redistributions of source code must retain the above copyright notice, this list of conditions and the following disclaimer.
- Redistributions in binary form must reproduce the above copyright notice, this list of conditions and the following disclaimer in the documentation and/or other materials provided with the distribution.

- Neither the name of the 'incremental' nor the names of its contributors may be used to endorse or promote products derived from this software without specific prior written permission.

THIS SOFTWARE IS PROVIDED BY THE COPYRIGHT HOLDERS AND CONTRIBUTORS "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED. IN NO EVENT SHALL THE COPYRIGHT OWNER OR CONTRIBUTORS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.

The Model 1700 firmware includes an OpenSSH client and server and also uses OpenSSH libraries in other components. The individual components of OpenSSH are subject to different licenses, all of which permit free distribution of compiled versions of the software, though some require explicit notice. Those required notices are reproduced below:

NO WARRANTY

BECAUSE THE PROGRAM IS LICENSED FREE OF CHARGE, THERE IS NO WARRANTY FOR THE PROGRAM, TO THE EXTENT PERMITTED BY APPLICABLE LAW. EXCEPT WHEN OTHERWISE STATED IN WRITING THE COPYRIGHT HOLDERS AND/OR OTHER PARTIES PROVIDE THE PROGRAM "AS IS" WITHOUT WARRANTY OF ANY KIND, EITHER EXPRESSED OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. THE ENTIRE RISK AS TO THE QUALITY AND PERFORMANCE OF THE PROGRAM IS WITH YOU. SHOULD THE PROGRAM PROVE DEFECTIVE, YOU ASSUME THE COST OF ALL NECESSARY SERVICING, REPAIR OR CORRECTION.

IN NO EVENT UNLESS REQUIRED BY APPLICABLE LAW OR AGREED TO IN WRITING WILL ANY COPYRIGHT HOLDER, OR ANY OTHER PARTY WHO MAY MODIFY AND/OR REDISTRIBUTE THE PROGRAM AS PERMITTED ABOVE, BE LIABLE TO YOU FOR DAMAGES, INCLUDING ANY GENERAL, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES ARISING OUT OF THE USE OR INABILITY TO USE THE PROGRAM (INCLUDING BUT NOT LIMITED TO LOSS OF DATA OR DATA BEING RENDERED INACCURATE OR LOSSES SUSTAINED BY YOU OR THIRD PARTIES OR A FAILURE OF THE PROGRAM TO OPERATE WITH ANY OTHER PROGRAMS), EVEN IF SUCH HOLDER OR OTHER PARTY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.

The 32-bit CRC compensation attack detector in deattack.c was contributed by CORE SDI S.A. under a BSD-style license.

```
* Cryptographic attack detector for ssh - source code
*
* Copyright (c) 1998 CORE SDI S.A., Buenos Aires, Argentina.
*
* All rights reserved. Redistribution and use in source and binary
* forms, with or without modification, are permitted provided that
* this copyright notice is retained.
*
* THIS SOFTWARE IS PROVIDED ``AS IS'' AND ANY EXPRESS OR IMPLIED
* WARRANTIES ARE DISCLAIMED. IN NO EVENT SHALL CORE SDI S.A. BE
* LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY OR
* CONSEQUENTIAL DAMAGES RESULTING FROM THE USE OR MISUSE OF THIS
* SOFTWARE.
*
* Ariel Futoransky <futo@core-sdi.com>
* <http://www.core-sdi.com>
```

Appendix

Troubleshooting

ssh-keyscan was contributed by David Mazieres under a BSD-style license.

```
* Copyright 1995, 1996 by David Mazieres <dm@lcs.mit.edu>.
*
* Modification and redistribution in source and binary forms is
* permitted provided that due credit is given to the author and the
* OpenBSD project by leaving this copyright notice intact.
```

The Rijndael implementation by Vincent Rijmen, Antoon Bosselaers and Paulo Barreto is in the public domain and distributed with the following license:

```
* @version 3.0 (December 2000)
*
* Optimised ANSI C code for the Rijndael cipher (now AES)
*
* @author Vincent Rijmen <vincent.rijmen@esat.kuleuven.ac.be>
* @author Antoon Bosselaers <antoon.bosselaers@esat.kuleuven.ac.be>
* @author Paulo Barreto <paulo.barreto@terra.com.br>
*
* This code is hereby placed in the public domain.
*
* THIS SOFTWARE IS PROVIDED BY THE AUTHORS 'AS IS' AND ANY EXPRESS
* OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED
* WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE
* ARE DISCLAIMED. IN NO EVENT SHALL THE AUTHORS OR CONTRIBUTORS BE
* LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR
* CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF
* SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR
* BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY,
* WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE
* OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE,
* EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
```

One component of the ssh source code is under a 3-clause BSD license, held by the University of California, since we pulled these parts from original Berkeley code.

```
* Copyright (c) 1983, 1990, 1992, 1993, 1995
*   The Regents of the University of California. All rights reserved.
*
* Redistribution and use in source and binary forms, with or without
* modification, are permitted provided that the following conditions
* are met:
* 1. Redistributions of source code must retain the above copyright
*   notice, this list of conditions and the following disclaimer.
* 2. Redistributions in binary form must reproduce the above copyright
*   notice, this list of conditions and the following disclaimer in the
*   documentation and/or other materials provided with the distribution.
* 3. Neither the name of the University nor the names of its contributors
*   may be used to endorse or promote products derived from this software
*   without specific prior written permission.
*
* THIS SOFTWARE IS PROVIDED BY THE REGENTS AND CONTRIBUTORS 'AS IS' AND
* ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE
* IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE
* ARE DISCLAIMED. IN NO EVENT SHALL THE REGENTS OR CONTRIBUTORS BE LIABLE
* FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL
* DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS
* OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION)
* HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT
* LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY
* OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF
* SUCH DAMAGE.
```

Appendix

Troubleshooting

Remaining components of the software are provided under a standard 2-term BSD license with the following names as copyright holders:

Markus Friedl
Theo de Raadt
Niels Provos
Dug Song
Aaron Campbell
Damien Miller
Kevin Steves
Daniel Kouril
Wesley Griffin
Per Allansson
Nils Nordman
Simon Wilkinson

Portable OpenSSH additionally includes code from the following copyright holders, also under the 2-term BSD license:

Ben Lindstrom
Tim Rice
Andre Lucas
Chris Adams
Corinna Vinschen
Cray Inc.
Denis Parker
Gert Doering
Jakob Schlyter
Jason Downs
Juha Yrj?l
Michael Stone
Networks Associates Technology, Inc.
Solar Designer
Todd C. Miller
Wayne Schroeder
William Jones
Darren Tucker
Sun Microsystems
The SCO Group
Daniel Walsh
Red Hat, Inc
Simon Vallet / Genoscope

```
* Redistribution and use in source and binary forms, with or without
* modification, are permitted provided that the following conditions
* are met:
* 1. Redistributions of source code must retain the above copyright
*   notice, this list of conditions and the following disclaimer.
* 2. Redistributions in binary form must reproduce the above copyright
*   notice, this list of conditions and the following disclaimer in the
*   documentation and/or other materials provided with the distribution.
*
```

```
* THIS SOFTWARE IS PROVIDED BY THE AUTHOR ``AS IS'' AND ANY EXPRESS OR
* IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES
* OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE DISCLAIMED.
* IN NO EVENT SHALL THE AUTHOR BE LIABLE FOR ANY DIRECT, INDIRECT,
* INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT
* NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE,
* DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY
* THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT
* (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF
* THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE.
```


Appendix

Troubleshooting

Portable OpenSSH contains the following additional licenses:

```
a) md5crypt.c, md5crypt.h
* "THE BEER-WARE LICENSE" (Revision 42):
* <phk@login.dknet.dk> wrote this file. As long as you retain this
* notice you can do whatever you want with this stuff. If we meet
* some day, and you think this stuff is worth it, you can buy me a
* beer in return. Poul-Henning Kamp
  b) snprintf replacement
* Copyright Patrick Powell 1995
* This code is based on code written by Patrick Powell
* (papowell@astart.com) It may be used for any purpose as long as this
* notice remains intact on all source code distributions
  c) Compatibility code (openbsd-compat)
    Apart from the previously mentioned licenses, various pieces of code
    in the openbsd-compat/ subdirectory are licensed as follows:
    Some code is licensed under a 3-term BSD license, to the following
    copyright holders:
Todd C. Miller
Theo de Raadt
Damien Miller
Eric P. Allman
The Regents of the University of California
Constantin S. Svintsoff
* Redistribution and use in source and binary forms, with or without
* modification, are permitted provided that the following conditions
* are met:
* 1. Redistributions of source code must retain the above copyright
* notice, this list of conditions and the following disclaimer.
* 2. Redistributions in binary form must reproduce the above copyright
* notice, this list of conditions and the following disclaimer in the
* documentation and/or other materials provided with the distribution.
* 3. Neither the name of the University nor the names of its contributors
* may be used to endorse or promote products derived from this software
* without specific prior written permission.
*
* THIS SOFTWARE IS PROVIDED BY THE REGENTS AND CONTRIBUTORS ``AS IS'' AND
* ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE
* IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE
* ARE DISCLAIMED. IN NO EVENT SHALL THE REGENTS OR CONTRIBUTORS BE LIABLE
* FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL
* DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS
* OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION)
* HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT
* LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY
* OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF
* SUCH DAMAGE.

Some code is licensed under an ISC-style license, to the following copyright
holders:
Internet Software Consortium.
Todd C. Miller
Reyk Floeter
Chad Mynhier
* Permission to use, copy, modify, and distribute this software for any
* purpose with or without fee is hereby granted, provided that the above
* copyright notice and this permission notice appear in all copies.
*
* THE SOFTWARE IS PROVIDED "AS IS" AND TODD C. MILLER DISCLAIMS ALL
* WARRANTIES WITH REGARD TO THIS SOFTWARE INCLUDING ALL IMPLIED WARRANTIES
* OF MERCHANTABILITY AND FITNESS. IN NO EVENT SHALL TODD C. MILLER BE LIABLE
* FOR ANY SPECIAL, DIRECT, INDIRECT, OR CONSEQUENTIAL DAMAGES OR ANY DAMAGES
```

* WHATSOEVER RESULTING FROM LOSS OF USE, DATA OR PROFITS, WHETHER IN AN ACTION
* OF CONTRACT, NEGLIGENCE OR OTHER TORTIOUS ACTION, ARISING OUT OF OR IN
* CONNECTION WITH THE USE OR PERFORMANCE OF THIS SOFTWARE.

Some code is licensed under a MIT-style license to the following copyright holders:

Free Software Foundation, Inc.

```
* Permission is hereby granted, free of charge, to any person obtaining a *
* copy of this software and associated documentation files (the *
* "Software"), to deal in the Software without restriction, including *
* without limitation the rights to use, copy, modify, merge, publish, *
* distribute, distribute with modifications, sublicense, and/or sell *
* copies of the Software, and to permit persons to whom the Software is *
* furnished to do so, subject to the following conditions: *
* *
* The above copyright notice and this permission notice shall be included *
* in all copies or substantial portions of the Software. *
* *
* THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS *
* OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF *
* MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. *
* IN NO EVENT SHALL THE ABOVE COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, *
* DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR *
* OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR *
* THE USE OR OTHER DEALINGS IN THE SOFTWARE. *
* *
* Except as contained in this notice, the name(s) of the above copyright *
* holders shall not be used in advertising or otherwise to promote the *
* sale, use or other dealings in this Software without prior written *
* authorization. *
```

The Model 1700 firmware uses portions of the "tslib" touchscreen library, which are licensed under the GNU Public License, version 2.

The source code to tslib was obtained by AMI from github.com using the command:

```
git clone https://github.com/kergoth/tslib
```

The Model 1700's touchscreen browser was linked against the Qt libraries for The X Window System that were available from the Debian package repository. The source code for those libraries, as well as the compilers and other tools required to recompile those libraries and the browser, are available from the Debian repository using the normal Debian package manipulation commands, e.g. pkg-add or symantic. The source code for the browser is brief enough to be included here:

```
---begin file browser.cpp---
#include <QtGui>
#include <QtWebKit>
int main(int argc, char** argv) {
    QApplication app(argc, argv);
    QWebView view;
    view.setWindowFlags (Qt::FramelessWindowHint);
    view.showFullScreen ();
    view.setUrl(QUrl(argv[1]));
}
```

```
        return app.exec();
    }
---end file browser.cpp---
---begin file browser.pro---
QT += webkit
SOURCES = browser.cpp
---end file browser.pro---
```

The Model 1700 firmware uses the Jansson library for encoding and decoding messages in the JSON (JavaScript Object Notation) format. The Jansson library is subject to the following license:

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

The "home" icon was obtained from https://www.iconfinder.com/icons/40798/home_icon#size=128

The license terms stated on that page are: Free for commercial use (Do not redistribute)

The "back" icon was obtained from: https://www.iconfinder.com/icons/27865/play_right_arrow_icon#size=128

The license terms stated on that page are: Creative Commons Attribution-No Derivative Works 3.0 Unported

The X.org display server is licensed under the following conditions:

Permission is hereby granted, free of charge, to any person obtaining a copy of this software and associated documentation files (the "Software"), to deal in the Software without restriction, including without limitation the rights to use, copy, modify, merge, publish, distribute, sublicense, and/or sell copies of the Software, and to permit persons to whom the Software is furnished to do so, subject to the following conditions:

The above copyright notice and this permission notice shall be included in all copies or substantial portions of the Software.

THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM, OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE SOFTWARE.

Glossary

Abbreviations and Acronyms Used in This Manual

Term	Meaning
AC; ac	Alternating Current; strictly, electrical <i>current</i> that periodically reverses direction. Typically used also to describe an electrical power source in terms of the <i>voltage</i> . For example, 240 Vac.
ASCII	American Standard Code for Information Interchange; numerical representation of characters such as 'a' or '@' or an action (such as line-feed); 'plain' raw text with no formatting such as tabs, bold or underscoring
BNC connector	A miniature quick connect/disconnect RF connector used for coaxial cable, featuring two bayonet lugs on the female connector.
CR	Text Carriage-Return character
Cryogen	A substance for obtaining low temperatures. in the case of use with the Model 1700 instrument, a cryogen is a liquefied gas such as liquid nitrogen or liquid helium.
D-Sub	Term referring to the family of connectors containing an odd number of pins in two parallel rows with a 1-pin difference in pins-per-row (DB9, DB15, and DB25 are most common)
DB9	Type of electrical connector containing 9 pins arranged in two parallel rows of 4 pins and 5 pins each
DCE	Data Circuit-terminating Equipment - a device that sits between the Data Terminal Equipment (DTE) and a data transmission circuit.
DHCP	Dynamic Host Configuration Protocol; a computer networking protocol which dynamically distributes the IP address to networked devices
dt	Rate of change
DTE	Data Terminal Equipment - an end instrument that converts user information into signals or reconverts received signals. A DTE device communicates with the Data Circuit-terminating Equipment (DCE).
ECL	Electrical Connection Lubricant - also known as Dielectric Tune-up Grease, a protective lubricant that prevents corrosion.
E_o	Power supply output voltage
i, I	Electrical current flow
I_o	Power supply output current
IEC	International Electrotechnical Commission

Glossary

Term	Meaning
IEEE	Institute of Electrical and Electronics Engineers
I/O	Input/Output; The hardware and associated protocol that implement communication between information processing systems and/or devices. Inputs are the signals or data received by the system or device, and outputs are the signals or data sent from it.
IP	Internet Protocol; when used with "address", refers to a numerical Internet address
kG	kilogauss: a magnetic field unit of measurement
LED	Light-Emitting Diode; a semiconductor device that emits light when energized - used for visual status indication
LHe	Liquid Helium
Max	Maximum
Min	Minimum
MSDS	Material Safety Data Sheet - provides workers and emergency personnel with procedures for handling or working with a specific substance in a safe manner and includes information such as physical data, toxicity, health effects, first aid, reactivity, storage, disposal, protective equipment, and spill-handling procedures.
RG-59/U	A specific type of coaxial cable, often used for low-power video and RF signal connections, with a characteristic impedance of 75 ohms.
R_{lead}	Electrical circuit lead or wiring resistance
RS-232	RS-232 is a long-established standard and protocol for relatively low speed serial data communication between computers and related devices; originally established for teletypewriter communication.
SCPI	Standard Commands for Programmable Instruments
V	Volts
VA	Volt-amperes ($V \times I$); a unit of electrical reactive power
V_{lead}	Voltage ($I \times R$) developed across circuit lead or wiring resistance due to current flow
V_m	Magnet voltage
V_s	Power supply voltage

Index

A

- abbreviations and acronyms 105
- analog input
 - pin outs 95, 96
- approximate calibration 57
- autofill
 - standard 18
 - system setup 16

C

- calibrating
 - approximate 57
 - closed dewar 61
 - methods 50
- calibration
 - presetting MAX/MIN 61
- closed dewar calibration 61
- coaxial connectors, moisture protection 14
- command
 - conventions 71
- command summary 71–??
- Configuration, instrument xv
- connecting
 - liquid helium sensors 15
- cryogenic liquids xii
- cryogenic liquids, first aid xii

E

- ECL use 14
- Ethernet connection pin out 94

F

- first aid xii
- front panel layout 6

I

- installing
 - autofill system 16
 - capacitance sensors 11
 - earth ground xv
 - liquid helium sensors 13
 - mounting 11
 - power 19
 - unpacking 11
- instrument description 1
- Instrument part number definition xvi
- IP Address 105

L

- liquid helium sensor
 - connector J1 wiring 95
 - connector warning 15
 - lead wire sizing 15
 - sensor preparation 15
- liquid level system terminology 105

M

- menus 22
- method of measurement 47

O

- operation warnings 11

P

- power on/off 21
- power supply
 - operating characteristics 7

R

- rear panel layout 6
- remote interface reference - see command RS-232
 - connector pin definitions 94
 - pin out 93

S

- safety
 - cryogenic liquids xii
 - equipment xiv
 - legend xiv
 - liquid helium sensor connectors 15
 - warnings xv
- Safety, cryogenics xii
- sensor
 - active length 46
 - installation 11
- sensor connector J1 wiring 95, 96
- sensor physical parameters
 - explanation 46
- sensor transmitters
 - calibration methods 50
- sensor vent holes 46– 47
- serial communications
 - commands 73
 - interactive communication 72
 - remote calibration 86
 - terminators 72
- serial connector pin definitions 94
- serial pin out 93
- serial port connector/cables 72
- system features 1
- system specifications 7

V

- vent holes 46– 47

W

- Warnings, equipment xv

Index
