



MODEL 1700 LIQUID LEVEL INSTRUMENT

(LIQUID HELIUM VERSION)

**INSTALLATION, OPERATION, AND
MAINTENANCE INSTRUCTIONS**

American Magnetics, Inc.

P.O. Box 2509, 112 Flint Road, Oak Ridge, TN 37831-2509, Tel: 865-482-1056, Fax: 865-482-5472

Foreword	vii
Purpose and Scope	vii
Contents of this Manual	vii
General Precautions	viii
Cryogen Safety Summary	x
Safety Legend	x
Equipment Warnings	xi
Instrument Configuration	xi
 Introduction	 1
Model 1700 Instrument	1
Superconducting Level Sensors	1
Digitally-Controlled	1
System Flexibility	1
Display	2
Intuitive Human-Interface Design	2
Measurement Flexibility	2
Real Time Clock	2
Analog Outputs	2
Signal Relays	2
Connectivity	2
Model 1700 Front Panel Layout	4
Model 1700 Rear Panel Layout	5
Model 1700 Specifications @ 25°C	6

Installation.....	9
Unpacking and Inspecting the Instrument	9
Mounting the Model 1700 Instrument	9
Liquid Helium Level Sensor Installation	10
Connecting the Liquid Helium Level Sensor	11
Configuring an Autofill System	12
LHe Autofill System Setup	12
Power Requirements.	13
Connecting the Analog Outputs	13
0-10 VDC Recorder Output.	13
4-20 mA Current Loop Output.	13
Operation.....	15
Energizing the Model 1700 Instrument	15
Screen Navigation.....	15
Home Screen	15
Home Screen Footer.	16
Editing a Field	16
Navigating the Instrument Menus	17
Menu Structure	18
Screen Descriptions	19
Superconducting (Liquid Helium) Level	24
Configure the instrument to display helium level.	24
Sensor Burnout Protection	26
Sampled or Continuous Measurements	27
Other Liquid Helium Functions	27
Alarms and Relays	29
Alarm Status Screen	29
Level-Based Alarms	29
Multiple alarms	30
Configuring Alarm Setpoints	30
Acknowledging an Alarm.	31
Muting an Alarm	31
Configuring Relay Setpoints	31

Configure the Autofill Function	33
Setting the autofill parameters	33
Enable the Autofill function	34
Clearing the Autofill Timeout Alarm	34
Select the appropriate units on the display	35
Analog output signals	35
Configuring the Analog Outputs	35
Ethernet Connectivity	36
IP Addressing Scheme	36
Serial Connectivity	37
Configuring the RS-232 Settings	37
Abnormal Operation	38
Dirty Helium Sensor Operational Mode	38
Resetting the Instrument to Factory Defaults	39
Shutting the Instrument Down	39
Calibration	41
Setting the System Date and Time	41
Superconductivity-Based (Liquid Helium) Level Calibration	43
Verify the Liquid Helium Sensor Type	43
Sensor Sample Interval	44
Sensor Active Length	45
Continuous Measure Time Limit	45
Sensor Name	46
Remote Interface Reference	47
SCPI Command Summary	47
RS-232 Configuration	52
Serial port connector and cabling	52
Command/return termination characters	53
Ethernet Configuration	53
Ethernet Connector	53
Termination Characters	53

Port Assignment	54
Command Reference	55
System Related Commands	55
Display Configuration Commands and Queries	58
Relay Configuration Commands and Queries	58
Alarm Configuration Commands and Queries	60
Measurement Commands and Queries	62
Fill Control and Queries	62
HE Channel Sampling Commands and Queries	63
Assignment Commands and Queries	64
HE Channel Calibration Commands and Queries	65
Remote Units Commands and Queries	65
Error Codes	66
Legacy Command Reference	67
Commands for Controlling the Units of Measurement	67
Commands for Configuring Permanent Memory	68
Querying the Configuration	69
Returning a Level Measurement	69
Service and Repair	71
Cleaning	71
User Replaceable Parts	71
Battery Replacement	72
Low Battery Indication	72
Tools Required	73
Procedure	73
Fuse Replacement	74
Tools Required	74
Procedure	74
Firmware Upgrade Via Ethernet	75
Upgrade via SCP	76
Appendix	79
Serial (RS-232) Connector	79

Ethernet Connector80

Liquid Helium Connector J1 Wiring81

Aux I/O Connector82

Troubleshooting83

 Instrument Displays “LOSS OF SENSOR” Condition for LHe Level Measurement 83

 Instrument Displays “SENSOR PROTECTED” Condition for LHe Measurement 83

 System Test Screens 83

 System Logs 86

Glossary93

Abbreviations and Acronyms93

Index95

FOREWORD

PURPOSE AND SCOPE

This manual contains the operation and maintenance instructions for the American Magnetics, Inc. Model 1700 Liquid Level Control Instrument and outlines applications for various system configurations. Since it is not possible to cover all equipment combinations for all magnet systems, only the most common configurations are discussed. The user is encouraged to contact an authorized AMI Technical Support Representative for information regarding specific configurations not explicitly covered in this manual. This manual refers to the instrument as configured for liquid helium operation.

CONTENTS OF THIS MANUAL

Introduction describes the functions, specifications, and characteristics of the Model 1700 Instrument. It provides illustrations of the front and rear panel layouts as well as documenting the performance specifications. Additional information is provided in the form of system block diagrams.

Installation describes how the Model 1700 Instrument is unpacked and installed in conjunction with ancillary equipment in typical systems. Block-level diagrams document the interconnects for various system configurations.

Operation describes how the Model 1700 Instrument is used to monitor liquid helium levels.

Calibration describes the calibration technique for liquid helium level sensors.

Remote Interface Reference documents all remote commands and queries available through the Model 1700 Instrument RS-232 and Ethernet interfaces. A quick-reference summary of commands is provided as well as a detailed description of each.

The **Appendix** and **Glossary** sections support the information in the sections listed above. See the Appendix section when referenced from other sections. See the Glossary for any words or acronyms presented in the above sections, requiring a more complete understanding.

GENERAL PRECAUTIONS

Cryogen Safety

The two most common cryogenic liquids used in superconducting magnet systems are nitrogen and helium. Both of these cryogens are extremely cold at atmospheric pressure (-321°F and -452°F , respectively). The following paragraphs outline safe handling precautions for these liquids.

Personnel handling cryogenic liquids should be thoroughly instructed and trained as to the nature of the liquids. Training is essential to minimize accidental spilling. Due to the low temperature of these materials, a cryogen spilled on many objects or surfaces may damage the surface or cause the object to shatter, often in an explosive manner.

Inert gases released into a confined or inadequately ventilated space can displace sufficient oxygen to make the local atmosphere incapable of sustaining life. Liquefied gases are potentially extreme suffocation hazards since a small amount of liquid will vaporize and yield a very large volume of oxygen-displacing gas. Always ensure the location where the cryogen is used is well ventilated. Breathing air with insufficient oxygen content may cause unconsciousness without warning. If a space is suspect, purge the space completely with air and test before entry. If this is not possible, wear a forced-air respirator and enter only with a co-worker standing by wearing a forced-air respirator.

Cryogenic liquids, due to their extremely low temperatures, will also burn the skin in a similar manner as would hot liquids. Never permit cryogenic liquids to come into contact with the skin or allow liquid nitrogen to soak clothing. Serious burns may result from careless handling. Never touch uninsulated pipes or vessels containing cryogenic liquids. Flesh will stick to extremely cold materials. Even nonmetallic materials are dangerous to touch at low temperatures. The vapors expelled during the venting process are sufficiently cold to burn flesh or freeze optic tissues. Insulated gloves should be used to prevent frost-bite when operating valves on cryogenic tanks. Be cautious with valves on cryogenic

systems; the temperature extremes they are typically subjected to cause seals to fail frequently.

In the event a person is burned by a cryogen or material cooled to cryogenic temperatures, the following first aid treatment should be given pending the arrival and treatment of a physician or other medical care worker:

1. If any cryogenic liquid contacts the skin or eyes, immediately flush the affected area gently with tepid water (102°F – 105°F, 38.9°C – 40.5°C) and then apply cold compresses.
2. Do not apply heat. Loosen any clothing that may restrict circulation. Apply a sterile protective dressing to the affected area.
3. If the skin is blistered or there is any chance that the eyes have been affected, get the patient immediately to a physician for treatment.

Containers of cryogenic liquids are self pressurizing (as the liquid boils off, vapor pressure increases). Hoses or lines used to transfer these liquids should never be sealed at both ends (i.e. by closing valves at both ends).

When pouring cryogenic liquids from one container to another, the receiving container should be cooled gradually to prevent damage by thermal shock. The liquid should be poured slowly to avoid spattering due to rapid boil off. The receiving vessel should be vented during the transfer.

Introduction of a substance at or near room temperature into a cryogenic liquid should be done with great caution. There may be a violent gas boil-off and a considerable amount of splashing as a result of this rapid boiling. There is also a chance that the material may crack or catastrophically fail due to forces caused by large differences in thermal contraction of different regions of the material. Personnel engaged in this type of activity should be instructed concerning this hazard and should always wear a full face shield and protective clothing. If severe spraying or splashing could occur, safety glasses or chemical goggles along with body length protective aprons will provide additional protection.

The properties of many materials at extremely low temperatures may be quite different from the properties that these same materials exhibit at room temperatures. Exercise extreme care when handling materials cooled to cryogenic temperatures until the properties of these materials under these conditions are known.

Metals to be used for use in cryogenic equipment application must possess sufficient physical properties at these low temperatures. Since ordinary carbon steels, and to somewhat a lesser extent, alloy steels, lose much of their ductility at low temperatures, they are considered unsatisfactory and sometimes unsafe for these applications. The

austenitic Ni-Cr alloys exhibit good ductility at these low temperatures and the most widely used is 18-8 stainless steel. Copper, Monel[®], brass and aluminum are also considered satisfactory materials for cryogenic service.

CRYOGEN SAFETY SUMMARY

Cryogenic systems are complex systems with the potential to seriously injure personnel or equipment if not operated according to procedures. The use of safety mechanisms (pressure relief valves, rupture disks, etc.) in cryogenic systems is usually necessary.

Recommended Safety Equipment

The use of proper safety equipment is necessary. Such equipment may include, but not limited to, the following items:

- First Aid kit
- Fire extinguisher rated for class C fires
- Cryogenic gloves
- Face shield
- Signs to indicate that there are potentially dangerous cryogenics in use in the area.

SAFETY LEGEND



Instruction manual symbol: the product is marked with this symbol when it is necessary to refer to the instruction manual in order to protect against damage to the product or personal injury.



Hazardous voltage symbol.



Alternating Current (Refer to IEC 417, No. 5032).



Off (Supply) (Refer to IEC 417, No. 5008).



On (Supply) (Refer to IEC 417, No. 5007).

WARNING

The Warning sign denotes a hazard. It calls attention to a procedure or practice, which if not correctly adhered to, could result in personal injury. Do not proceed beyond a Warning sign until the indicated conditions are fully understood and met.

CAUTION

The Caution sign denotes a hazard. It calls attention to an operating procedure or practice, which if not adhered to, could cause damage or destruction of a part or all of the product. Do not proceed beyond a Caution sign until the indicated conditions are fully understood and met.

EQUIPMENT WARNINGS

Before energizing the equipment, the earth ground of the power receptacle must be verified to be at earth potential and able to carry the rated current of the power circuit. Using extension cords should be avoided. However, if an extension cord must be used, insure the ground conductor is intact and the cord is capable of carrying the rated current without excessive voltage drop.

In the event that the ground path becomes less than sufficient to carry the rated current of the power circuit, the equipment should be disconnected from power, labeled as unsafe, and removed from place of operation.

Do not operate this equipment in the presence of flammable gases. Doing so could result in a life-threatening explosion.

Do not modify this equipment in any way. If component replacement is required, return the equipment to AMI facilities as described in the troubleshooting section of this manual.

If used in a manner not specified in this manual, the protection provided by the design, manufacture and documentation of the Model 1700 Instrument may be impaired.

INSTRUMENT CONFIGURATION

The Model 1700 Instrument is configured at time of purchase as one in several ways:

-
- As a helium level instrument for level sensors with active length up to 80 inches (203 cm) for either 4.2K or 2K LHe.

Every configuration may be further customized by the following options:

- Table top, single rack mounting, dual rack mounting.

- Line cord: North American, European Schuko, Australia/NZ, China, UK, pig-tailed ends.

The instrument part number, shown on the serialization label located on the underside of the instrument, as well as in a field in the instrument Menu, identifies the configuration according to the following key:

1700-A-B-C-D-E where

A indicates the mounting method:

Tbl = tabletop

SR19L = single rack mounted, 19" wide rack standard, instrument on left side

SR19R = single rack mounted, 19" wide rack standard, instrument on right side

SR10L = single rack mounted, 10" wide rack standard, instrument on left side

SR10R = single rack mounted, 10" wide rack standard, instrument on right side

DRL = dual rack mount, 19" wide rack standard, instrument on left side

DRR = dual rack mount, 19" wide rack standard, instrument on right side

B indicates the line cord shipped with the instrument:

N = North American

E = European, Schuko

A = Australian/New Zealand

C = Chinese

U = United Kingdom

P = India/pigtailed

C indicates the capacitance-based level configuration:

N = not configured

D indicates the helium level configuration:

He2K = 2K helium for sensors

He4K = 4.2K helium for sensors

N = not configured

The sensor active length and length units are appended to the configuration code, ie ...He2K-40IN-...

E is used to denote any instrument customization:

S = standard (no customization)

C = instrument modified.

INTRODUCTION

MODEL 1700 INSTRUMENT

The AMI Model 1700 Liquid Level Instrument is a sophisticated measurement and control instrument which provides monitoring liquid helium level sensors as inputs and provides for automatic level control based on user set parameters.

At time of purchase, the Model 1700 will be configured as:

- A liquid helium level instrument/controller.

The Model 1700 instrument will measure liquid helium level using a superconducting level sensor. The instrument will be factory configured for either 4.2K or 2K liquid helium level sensors at time of purchase.

SUPERCONDUCTING LEVEL SENSORS

The instrument can be used with a superconducting level sensor to measure liquid helium levels. The instrument will be configured at the factory to measure liquid helium levels. The instrument is configured for level sensors of active lengths of 1 to 80 inches. The instrument can be used with either 4.2K or 2K level sensor types but must be configured and calibrated at time of purchase for the LHe level sensor type.

DIGITALLY- CONTROLLED

The Model 1700 contains a microcomputer which controls analog data conversion, display/keypad functions, communications I/O, dry contact closures, generation of analog output signals and relay control of a mains power outlet for solenoid valve autofill applications.

SYSTEM FLEXIBILITY

The Model 1700 instrument incorporates data converters to translate signals between the analog and digital domains. Precision instrumentation techniques and potentiometer-free designs are

employed throughout the Model 1700 Instrument to ensure long term stability and accurate signal translation for a wide range of conditions.

DISPLAY

The Model 1700 Instrument has a 4.3" diagonal measure TFT (Thin Film Transistor) color liquid crystal display of 480 x 272 pixels. The display has a 4-wire resistive touch overlay for easy operator input.

INTUITIVE HUMAN- INTERFACE DESIGN

The Model 1700 instrument is designed to simplify the touch-screen based user interface. All functions were analyzed and subsequently programmed so that the most commonly used functions are addressed with the least number of keystrokes. The menus are presented in a logical fashion so that the operation of the Model 1700 is intuitive to the user. Context-sensitive Help screens are also provided.

MEASUREMENT FLEXIBILITY

Set points can be assigned to control two dry-contact relay outputs. The make or break function of each relay can be set independently. A solid-state relay allows mains power to be switched according to set points to operate a solenoid-operated valve or other load.

REAL TIME CLOCK

The Model 1700 Instrument incorporates a real time clock, which can be manually set or synchronized via NTP, with support for worldwide timezones including automatic DST adjustment.

ANALOG OUTPUTS

The Model 1700 Instrument has two analog outputs, a 0-10 V_{DC} voltage output and a 4-20 mA_{DC} current loop output. The 4-20 mA_{DC} loop output has 1500 V_{PK} circuit isolation. The outputs can be used simultaneously.


SIGNAL RELAYS

The Model 1700 Instrument has two signal relays that change state based on a measured input. The set points of these relays are user-selectable as is their function, alarm on level above or below set point.

CONNECTIVITY

The Model 1700 Instrument has a 10Base-T Ethernet connection as well as an RS-232 port for connecting to other equipment. The instrument communicates with a SCPI-based command set. The command set is 100% backward compatible with the AMI Model 135, 136, 185, and 186 instruments when configured as a single channel instrument. When configured as a dual channel instrument, additional commands are incorporated into the command set.

The Model 1700 Instrument allows for remote operation with an external browser via TCP/IP connection. All commands that are available by touching the local screen are available via the web browser



Model 1700 Helium Instrument Via Web Browser


NOTE Using a web browser to connect to the instrument allows different browser sessions to display different information. For instance, one browser window can display helium level only and a second browser window can display nitrogen level only. The instrument can be configured independently as well since it's display is also a browser. Refer to "Configure the instrument to display helium level" on page 24.

MODEL 1700 FRONT PANEL LAYOUT



Model 1700 Front View; Dual Channel Instrument Shown

MODEL 1700 REAR PANEL LAYOUT



Model 1700 Rear Panel Description

1	Computer Network Connector
2	Aux I/O Connector
3	RS-232 Serial Connector
4	Switched Valve Outlet Socket (IEC 60320 C13)
5	Capacitive Sensor Input Connector
6	Mains Power Entry Connector (IEC 60320 C14)
7	LHe Level Sensor Connector

LHe OPTION

MODEL 1700 SPECIFICATIONS @ 25°C

System Architecture

Display:	4.3" 24-bit color TFT display, 480x272 pixel with resistive touch screen
Sensor Types:	Superconductivity-based liquid level
Maximum Length Readout:	Superconductivity-based liquid level (LHe) up to 80 in
Superconducting (LHe) Sensor Excitation:	Continuous reading or Sample and Hold mode
System Operating Firmware Storage:	microSD card
System Clock:	Real time clock with timezone support, automatic DST adjustment, and NTP synchronization
Display Measurement Units:	Liquid level in cm, in or percent

Level Measurement

Resolution:	0.1%, 0.1 cm, 0.1 in
Accuracy:	±0.5% of active sensor length
Linearity:	±0.1% or 1 mm (whichever is greater)
Superconducting LHe Sensor Current:	4.2K LHe Temperature: 75 mA _{DC} nominal 2K LHe Temperature: 57 mA _{DC} nominal
Dirty Sensor Mode:	Approximately twice normal current for 1 second prior to normal measurement excitation
Superconducting LHe Sensor Voltage:	4.2K Sensor: approximately 0.87 V _{DC} per inch of sensor active length @ 10K 2K Sensor: approximately 0.66 V _{DC} per inch of sensor active length @ 10K
Maximum Open Circuit Voltage:	96 V _{DC} , galvanically isolated

Operating Parameters

Alarm Set Points:	0% to 100%, adjustable; Alarm condition settable to above or below set point
Sample and Hold Period:	1 second to 86,400 seconds (24 hrs)
Audible alarm:	3500 ± 500 Hz, 73 to 86 dB(A)

Analog Outputs

Output Types:	0-10 V _{DC} and simultaneous 4 - 20 mA _{DC}
4-20 mA Current Loop Power Supply Voltage:	12-32 V _{DC}
0-10 V _{DC} Recorder Output Load:	50k ohms or greater
0-10 V _{DC} Recorder Output 0% Max Offset:	+20 mV
0-10 V _{DC} Recorder Output 100% Max Error:	± 80 mV
Voltage and Current Output Converter Resolution:	12 bits
Integral Non-linearity:	±1LSB
Differential Non-linearity ^a :	±1LSB

Relays

No1 and/or No2 (W171DIP-7, or equivalent):	Contact Form: 1 Form A (SPST-NO) Maximum Switched Current: 3 A Switching Voltage: 60 V _{AC} / 100 V _{DC}
Level Control (Solid State):	Rated Load Voltage: 2 A at 100 to 240 V _{AC} Load Voltage Range: 75 to 264 V _{AC} (Line Dependent) Load Current: 0.1 to 2 A Surge Current: 30 A (60 Hz, 1 cycle)

Communication Protocol

Host Computer Network Protocol:	Ethernet 10Base-T TCP/IP and RS-232
IP Addressing:	DHCP or static, IPv4
Network Connectivity and Traffic Indication:	Link and Activity LEDs on instrument rear panel
RS-232 Connector Specifications:	9-pin D-sub female connector to connect standard DTE 9-pin D-sub male connector using a standard straight cable
Communication Command Set:	SCPI-based. 100% backward compatible with the Model 134, 135, 136 instruments.

Power Requirements

Primary:	100-240 ±10% V _{AC} , 50-60 Hz, 2.2 A maximum (200 VA plus sum of controller output)
Backup Battery for RTC:	CR2032

Physical

Dimensions ^b :	Table top configuration: 3.8" H x 8.4" W x 11.4" D [97 mm H x 213 mm W x 290 mm D] Single rack mount configuration: 3.5" H x 19.0" W x 11.4" D [89 mm H x 483 mm W x 290 mm D]
Weight:	Table-top configuration: 3.3 lbm [1.5 kG]; Single rack-mount configuration: 4.0 lbm [1.8 kG]

Environmental Limits

Ambient Temperature:	Operating: 0°C to 40°C [32°F to 104°F]; Non-operating: -20 °C to 60 °C [-4 °F to 140 °F]
Relative Humidity:	0 to 95%; non-condensing
Maximum Instrument Background Field:	Gauss (TBD)

Standards

Test Standards	<p>Testing of Equipment for Measurement, Control, and Laboratory Use (IEC 61326-1:2012, EN 61326-1)</p> <p>Electrostatic Discharge (ESD) (EN 61000-4-2)</p> <p>Radiated Immunity (EN 61000-4-3)</p> <p>Fast Transient Burst (EN 61000-4-4)</p> <p>Surges (EN 61000-4-5)</p> <p>Conducted Immunity (EN 61000-4-6)</p> <p>Power Frequency Magnetic Field (EN 61000-4-8)</p> <p>Voltage Dips and Interrupts (EN 61000-4-11)</p> <p>Harmonics (EN 61000-3-2)</p> <p>Flicker (EN 6100-3-3)</p> <p>Conducted Emissions (EN 55011/IEC/CISPR 11)</p> <p>Radiated Emissions (EN 55011/IEC/CISPR 11)</p> <p>Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use (IEC 61010-1)</p>
----------------	---

- a. Guaranteed monotonic over operating temperature range
- b. H = height; W = width; D = depth

INSTALLATION

WARNING

Equipment warnings apply to all system installation configurations. Refer to “Equipment Warnings” on page xi, in the Foreword to be familiar with the safety requirements for a system installation.

UNPACKING AND INSPECTING THE INSTRUMENT

Carefully remove the equipment, interconnecting cabling, and documentation from the shipping carton, and remove all packaging material.

NOTE

If there is any shipping damage, save all packing material and contact the shipping company representative to file a damage claim. Do not return to AMI unless prior authorization has been received.

MOUNTING THE MODEL 1700 INSTRUMENT

If the Model 1700 Instrument is to be used as a table top model, place the equipment on a flat, secure surface.

If the Model 1700 Instrument is to be rack mounted, install it in a 19" wide instrument rack using the mounting hardware supplied by the rack cabinet manufacturer. The feet on the bottom of the instrument may be removed to facilitate rack mounting by using a prying device to remove the glossy black cap from the center of the foot and then prying the rubber foot itself from the chassis. Secure the front panel to the rack rail in each of the four corners.

LIQUID HELIUM LEVEL SENSOR INSTALLATION

1. The sensor must be mounted with the electrical leads at the top.
2. For minimum losses, mount the liquid helium sensor so that warm helium gas rising from the sensor can pass directly out of the dewar without contacting surfaces at 4.2K.
3. Do not mount the sensor in restricted areas (tubes, etc.) where the liquid level around the sensor might be depressed by pressure differences in the gas. Do not cover the holes in the sensor.
4. The sensor may be mounted by taping or clipping it to an appropriate support structure. Do not exert excess pressure on the sensor with the mounting device to avoid crushing the tube. Avoid constraining both ends of the sensor and allow for contraction of the sensor during cooldown.

CAUTION

Do not operate the sensor in a vacuum. Operating the sensor in a vacuum may cause thermal damage and/or destruction of the superconducting filament sensor. Do not inadvertently turn the instrument on with the sensor in an evacuated chamber. Operation in pumped liquid helium environments is acceptable to 1K as long as liquid helium is present.

5. Avoid installing in a location where icing (frozen water or gas) may occur since ice formations may cause erratic operation. Ice formation on the NbTi filament may stop the propagation of the normal (resistive) zone before it actually reaches the liquid/gas interface. This will give an indication of a higher helium level than actually exists.
6. Connect the sensor to the Model 1700 LHe Level Sensor connector on the instrument rear panel (refer to "Liquid Helium Connector J1 Wiring" on page 81). The liquid helium level sensor leads are color coded:

LHe Level Sensor Wire Identification

Wire Function	Teflon Insulation Color	Formvar Insulation Color	Instrument Connector Pin
I+	Red	Red	1
V+	Blue	Green	8
V-	Yellow	Natural	6
I-	Black	Blue	7

CONNECTING THE LIQUID HELIUM LEVEL SENSOR

The instrument is connected to the liquid helium level sensor with a 4-conductor cable which has a 9-pin D-sub male connector on one end that mates with the connector used at the instrumentation feed through connector on the cryostat. This connector is typically a multi-pin circular type connector.

Prepare the sensor to be connected to the instrument by soldering the sensor leads to a male 9-pin D-Sub connector which will connect to the female 9-pin D-Sub connector on the transmitter. Refer to the table on page 10 and the *Appendix* of this manual as well as the AMI sensor manual for the proper pin out and wire color connections. Connect the sensor to the **LHe** connector on the rear panel.

WARNING



Although the sensor connector terminals are isolated from earth ground and therefore touching one terminal is not hazardous, the voltage between terminals is at a hazardous potential. The sensor connector is for use with an AMI LHe sensor and the wiring for the sensor is to have no live parts which are accessible. Conductors connected to its terminals must be insulated from user contact by basic insulation rated for 150 V_{AC} (Category I).

The lead wire for the sensor may be sized by the following equation:

$$R = 420 - 5.21(L_A)$$

where R is the maximum allowable resistance (in ohms) for each lead wire from the instrument to the sensor, and L_A is the *active length* of the connected helium level sensor in inches. Values for active sensor length vs. lead wire distance are provided in the table below.

Minimum recommended wire gauge for copper lead wire

Distance	$R=367$	$R=315$	$R=263$	$R=211$	$R=107$	$R=3.2$
	$L_A=10''$	$L_A=20''$	$L_A=30''$	$L_A=40''$	$L_A=60''$	$L_A=80''$
10 ft.	36 AWG				36 AWG	34 AWG
20 ft.						30 AWG
30 ft.						28 AWG
40 ft.						27 AWG
50 ft.						24 AWG
100 ft.						22 AWG
200 ft.						16 AWG
500 ft.					32 AWG	

CONFIGURING AN AUTOFILL SYSTEM

Autofill systems are useful in applications where liquid level in a dewar must be maintained without operator intervention.

To create an autofill system, a cryogenic liquid source must be attached to the target or level controlled dewar via a customized valving system for LHe application. The valve should be controlled by the Model 1700 Instrument.

LHe AUTOFILL SYSTEM SETUP

An autofill system for LHe application is a more complex issue than simply selecting a solenoid-controlled valve that can be used for other cryogenics such as LN2. The valving control must be highly-insulated to avoid vaporizing the transferred LHe. Contact an AMI Technical Support Representative to inquire about available LHe autofill solutions as this is a continually changing marketplace.

POWER REQUIREMENTS

WARNING

The Model 1700 Instrument operates on 50-60 Hz power and may be powered from 100-240 VAC. Insure that the input ground terminal is connected securely to an external earth ground.

Ensure the detachable mains supply cord is of suitable rating, i.e. 10 A (min) at 125 VAC for North America.

Ensure the power switch is in the OFF (●) position. Plug the Model 1700 Instrument line cord into the power entry module on the instrument rear panel and into the appropriate power receptacle.

CONNECTING THE ANALOG OUTPUTS

The Model 1700 instrument provides 0-10 V_{DC} and 4-20 mA_{DC} outputs as standard. These outputs may be optionally connected to external equipment per the instructions below.

NOTE

Both analog outputs can be connected to the same or different level measurement source as desired, and both will operate simultaneously.

0-10 V_{DC} RECORDER OUTPUT

The Model 1700 instrument provides a 0-10 V_{DC} output on the AUX I/O connector (refer to “Aux I/O Pin Definitions” on page 82) corresponding to 0-100% of liquid level. This output can be assigned to either level measurement channel (refer to “Analog output signals” on page 35). Connect the voltage output on the AUX I/O connector to a suitable receiver being careful not to exceed the impedance restrictions (refer to “Model 1700 Specifications @ 25°C” on page 6).


4-20 mA CURRENT LOOP OUTPUT

The Model 1700 instrument provides a 4-20 mA_{DC} output on the AUX I/O connector (refer to “Aux I/O Pin Definitions” on page 82) corresponding to 0-100% of liquid level. This output can be assigned to either level measurement channel (refer to “Analog output signals” on page 35).

The instrument requires an external power supply for current loop operation. Connect a suitable power supply as shown below.

CAUTION

It is extremely important to observe all polarities and to not exceed +32VDC for the loop power supply in order to prevent damage to the 4-20mA driver circuit.



OPERATION

This section describes the operation of the Model 1700 Instrument.

ENERGIZING THE MODEL 1700 INSTRUMENT

1. Turn the power switch on the front panel of the instrument to the On (I) position. The display will briefly show a boot image and then display the home screen showing level(s).
2. The boot process takes approximately 30 seconds. Boot time can be longer (approximately 2 minutes) if the instrument has been configured for a network connection and then is booted without the network present.
3. When the boot process is complete, the instrument will display the home (level) screen.


SCREEN NAVIGATION

HOME SCREEN

NOTE If the instrument was purchased with a level sensor, the instrument will be pre-configured and calibrated at the factory.

The display will look similar to the image shown at right:

If the instrument requires calibration¹, refer to the following chapter to calibrate the instrument with an AMI level sensor.



LHe Home Screen

HOME SCREEN FOOTER

Every screen has a footer. The level home screen displays both status information and buttons in the footer to navigate to other screens.

The home screen footer appears as shown at right:



Home Screen Footer

Model 1700 Instrument Home Screen Footer

No.	Name	Function
①	MENU	Takes the user to the main menu screen.
②	ALARM	When illuminated, indicates an active alarm condition. Touch to enter the ALARM STATUS screen.
③	AUTOFILL	When illuminated, indicates the AUTOFILL function is enabled. Touch to enter the AUTOFILL configuration screen.
④	HELP	Takes the user to a Help screen where available.

EDITING A FIELD

Once a field or multi-choice button on a screen has been edited, the footer changes to appear as follows:



Field Editing Footer

Model 1700 Instrument Footer during editing a field

No.	Name	Function
①	SAVE	Saves the entries made on the screen. The footer then reverts to the default footer as shown in the previous section.
②	CANCEL	Cancels a change while <i>not saving entries</i> . The footer then reverts to the default footer as shown in the previous section.

The functions of the **ALARM**, **AUTOFILL**, and **Help** footer items are as described in the previous section.

1. If the instrument was purchased with level sensor(s), the instrument will be shipped set up and calibrated.


NAVIGATING THE INSTRUMENT MENUS

The menu system is invoked by pressing the **MENU** button in the lower left corner of the instrument level display screen.



MENU button on HOME screen

When invoked, the **MENU** screen will be displayed:



MENU Selection Screen

Pressing **HOME** in the lower left corner of the screen will display the level display screen.




HOME button on MENU screen

Once a **MENU** has been chosen, a back button will be presented to the right **HOME**. This button will return the user to the previously displayed screen.



BACK button

MENU STRUCTURE






Model 1700 Menu Structure



SCREEN DESCRIPTIONS

The Model 1700 Instrument displays various screens on the graphic display when a MENU item is chosen. The items for each screen are described in the following table.




Model 1700 Helium Level Instrument Screen Descriptions

Screen ^a	Item Label	Field Type or Function
SENSOR CONFIGURATION(S)  MENU > SENSORS	HELIUM SENSOR TYPE:	Information: 2K, 4.2K ^a
	HELIUM LEVEL ON HOME SCREEN?	Toggles between: YES, NO
	SENSOR NAME(S)	Move to another screen
	CAL He	Move to another screen
	CAL He	Move to another screen
SENSOR NAME CONFIGURATION  MENU > SENSORS > SENSOR NAME(S)	HELIUM SENSOR NAME:	Data entry: <value>
LIQUID HELIUM SENSOR CALIBRATION  MENU > SENSORS > CAL He	SENSOR SAMPLE INTERVAL:	Data entry: <value> min
	MEASURE SENSOR V/A	Performs a function: <values> V mA ^a
	SENSOR ACTIVE LENGTH:	Data entry: <value> cm or in (tap units to change)
	CONTIN. MEASURE TIME LIMIT:	Data entry: <value> min
	DIRTY SENSOR MODE:	Toggles between: DISABLED, ENABLED
	SENSOR BURNOUT PROTECTION LIMIT:	Data entry: <value> %





Model 1700 Helium Level Instrument Screen Descriptions

Screen ^a	Item Label	Field Type or Function
OUTPUT CONFIGURATION: ALARMS  MENU > OUTPUTS > ALARMS	ALARM 1 SOURCE:	Toggles between: DISABLED, HELIUM
	ALARM 1 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	ALARM 1 WHEN LEVEL <state>	Toggles between: ≤, ≥
	ALARM 2 SOURCE:	Toggles between: DISABLED, HELIUM
	ALARM 2 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	ALARM 2 WHEN LEVEL <state>	Toggles between: ≤, ≥
	RELAYS	Move to another screen
	AUTOFILL	Move to another screen
	ANALOG OUTPUTS	Move to another screen
OUTPUT CONFIGURATION: RELAYS  MENU > OUTPUTS > RELAYS	RELAY 1 SOURCE:	Toggles between: DISABLED, HELIUM
	RELAY 1 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	RELAY 1 CLOSES WHEN LEVEL <state>	Toggles between: ≤, ≥
	RELAY 2 SOURCE:	Toggles between: DISABLED, HELIUM
	RELAY 2 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	RELAY 2 CLOSES WHEN LEVEL <state>	Toggles between: ≤, ≥
	ALARMS	Move to another screen
	AUTOFILL	Move to another screen
	ANALOG OUTPUTS	Move to another screen




Model 1700 Helium Level Instrument Screen Descriptions

Screen ^a	Item Label	Field Type or Function
OUTPUT CONFIGURATION: AUTOFILL  MENU > OUTPUTS > AUTOFILL	AutoFILL CONTROL:	Toggles between: DISABLED, HELIUM
	START:	Data entry <value> %, in, or cm (tap units to change)
	STOP:	Data entry <value> %, in, or cm (tap units to change)
	FILL TIMEOUT:	Data entry <value> MINUTES
	RELAYS	Move to another screen
	ALARMS	Move to another screen
	ANALOG OUTPUTS	Move to another screen
OUTPUT CONFIG.: ANALOG OUTPUTS  MENU > OUTPUTS > ANALOG OUTPUTS	OUTPUT 0-10 VDC SOURCE:	Toggles between: DISABLED, HELIUM
	OUTPUT 4-20 mA SOURCE:	Toggles between: DISABLED, HELIUM
	RELAYS	Move to another screen
	ALARMS	Move to another screen
	AUTOFILL	Move to another screen
NETWORK CONFIGURATION  MENU > NETWORK	Addressing Scheme:	Toggles between: DYNAMIC, DISABLED, STATIC
	Address:	Information: <value> ^b (DYNAMIC) Data entry: <value> (STATIC)
	Netmask:	Data entry: <value> ^b
	Gateway:	Data entry: <value> ^b
	Hostname:	Data entry: <value>
	MAC Address:	Information: <value> ^b

Model 1700 Helium Level Instrument Screen Descriptions

Screen ^a	Item Label	Field Type or Function
SYSTEM CONFIGURATION  MENU > SYSTEM	Instrument Serial Number:	Information: <value> ^a
	HW: FPGA: FW:	Information: <values> ^a
	Date of Manufacture:	Information: <value> ^a
	Helium Power Supply:	Information: <value> ^a
	SET TIME	Move to another screen
	UPDATE	Move to another screen
	RS-232 SETUP	Move to another screen
	SYSTEM TEST	Move to another screen
SYSTEM DATE/TIME  MENU > SYSTEM > SET TIME	Current System Date and Time:	Information ^a
	System Timezone:	Information ^a
	SET MANUALLY	Move to another screen
	SET FROM NTP SERVER	Action with status screen
FIRMWARE UPDATE INSTALL  MENU > SYSTEM > UPDATE	Select Firmware Version:	Selection list: <values> and then moves to INSTALL verification on selection
RS-232 SERIAL PORT CONFIGURATION  MENU > SYSTEM > RS-232 SETUP	Current Configuration:	Information ^b
	Configuration Choice:	Choose alternate function by touching
	Echo Setting:	Choose ECHO OFF or ECHO ON
	Line Ending:	Choose NO LF or ADD LF
	Baud rate:	Choose desired baud rate

Model 1700 Helium Level Instrument Screen Descriptions

Screen ^a	Item Label	Field Type or Function
SYSTEM TEST  MENU > SYSTEM > SYSTST	TURN ON ALL LCD PIXELS	Move to multi-color test screen
	TURN ON RED LCD PIXELS	Move to red test screen
	TURN ON GREEN LCD PIXELS	Move to green test screen
	TURN ON BLUE LCD PIXELS	Move to blue test screen
	HARDWARE CONTROL	Move to hardware test screen
	TOUCH SCREEN ALIBRATE	Move to touch calibration screen
	TOUCH SCREEN TEST	Move to touch test screen
	CURSOR ON/CURSOR OFF	Toggles visibility of cursor
HISTORY  MENU > HISTORY	LHe	Enables/disables helium trace
	<	Move graph back in time
	IN	Zoom in (less time per division)
	OUT	Zoom out (more time per division)
	>	Move graph forward in time
	>>	Shift graph to most recent levels
SHUTDOWN  MENU > SHUTDOWN	SHUTDOWN	Shuts down the instrument in an orderly fashion which reduces boot time for the next power on.

a. Displays the state or value (display only).

b. Value displayed is chosen by the DHCP server for the network in DYNAMIC mode.


SUPERCONDUCTING (LIQUID HELIUM) LEVEL

CONFIGURE THE INSTRUMENT TO DISPLAY HELIUM LEVEL


NOTE If the instrument was purchased with an AMI LHe level sensor, Steps 1 through 6, below have already been performed.

1. From the main screen, choose the following: **MENU > SENSORS**.


2. The instrument will display the instrument configuration vis-à-vis the LHe level sensor, 4.2K or 2K. Note that this configuration is set at time of manufacture and is not alterable by the user.




3. Ensure that **HELIUM LEVEL ON HOME SCREEN?** is set to **YES**.



4. Press the **CAL He** button.




5. Press the **SENSOR ACTIVE LENGTH** field.



Calibrate Helium Screen


NOTE The entered sensor active length *must* match the actual active length of the connected LHe sensor for accurate results.

6. In the numeric pop-up keypad, enter the sensor active length in centimeters. Press **ENTER** when finished.



Numeric Keypad


7. Touch in the **CONTIN. MEASURE TIME LIMIT** field.
8. Using the pop-up numeric keypad, enter the maximum amount of time that the sensor should remain energized in the **UPDATING** mode before automatically changing to the **HOLDING** mode. This limit prevents inadvertent sensor energization for long periods of time which will cause excessive liquid helium boil off. When the sensor is energized (**UPDATING**) a timer is started and after the **CONTIN. MEASURE TIME LIMIT** is reached, the sensor is de-energized (switches back to **HOLDING** mode).



Calibrate Helium Screen

NOTE Setting the **CONTIN. MEASURE TIME LIMIT** to 0 (zero) disables the timeout function and the sensor will remain energized indefinitely until the operator proactively exits to **HOLDING** mode.


9. Set the **SENSOR SAMPLE INTERVAL** to an appropriate value. This value is used by the **HOLDING** timer to determine how often the sensor is sampled and the display is updated. The timer is started when the instrument is powered up, when this value is changed, or when entering the **HOLDING** (sampled) mode.



Setting the Sample Interval

NOTE Setting the **SAMPLE INTERVAL** to 0 (zero) disables the sampling function and the sensor is energized at all times.

10. The **SENSOR BURNOUT PROTECTION LIMIT** is typically set at approximately 5%. This limit specifies the amount of excess sensed resistance that is allowed before the Model 1700 engages the burnout protection feature (see the next section) primarily intended to protect the sensor filament in a *vacuum* environment. If you wish to allow a higher excess resistance limit in a warm, gas-filled cryostat, increase the value by 2% increments until the burnout protection is disengaged.



Setting the Burnout Protection Limit

11. Press **HOME** in the footer to return to the home screen.


SENSOR BURNOUT PROTECTION

Operating the sensor in a *vacuum* without contact with liquid helium can cause thermal damage and/or destruction of the superconducting filament. The Model 1700 is equipped with circuitry that automatically prevents this type of damage from occurring (the sensor active length setting must be correctly configured). If excessive sensor resistance is detected, the Model 1700 display will indicate burnout protection is in effect and the sensor current will be turned off for a minimum of six seconds after which time the instrument will attempt to resume normal operation. When the sensor resistance has returned to within the expected range, normal operation will be restored. No relay or alarm states are affected.



Helium Level Sensor Protection Mode

If an attached sensor is installed in a *warm, gas-filled* cryostat (not in a vacuum), the sensor burnout protection may activate as soon as the Model 1700 is energized and a liquid helium level measurement is attempted. This is normal until the cryostat and sensor are cooled in the helium environment. However, you may adjust the **SENSOR BURNOUT PROTECTION LIMIT** setting in the **MENU > SENSORS > CAL He** screen (see page 26) to allow the sensing function to operate in the warm, gas-filled cryostat if desired.



Helium Level Sensor Protection Mode
in a Dual Display

SAMPLED OR CONTINUOUS MEASUREMENTS


While viewing the home screen, toggle between **HOLDING** (sampled) and **UPDATING** (continuous measurement) modes by pressing the button.



Helium Level Measurement Mode

OTHER LIQUID HELIUM FUNCTIONS

In the **MENU > SENSORS > CAL He** screen, the instrument displays the LHe level sensor voltage. If the voltage is shown in light blue, it is the actual (real-time) voltage as the instrument is in **UPDATING** mode. If it is displayed in gray, it is the voltage measured the last time the sensor was energized (**HOLDING** mode). To update the reading, press the **MEASURE SENSOR V/A** button and the level reading will be updated as will the measured sensor voltage and excitation current.



Helium Sensor Voltage


NOTE Nominal LHe sensor excitation current is 75 mA for 4K sensors, and 57 mA for 2K sensors.

AMI expects the helium level sensor to be reasonably clean and free from oil, water, ice, etc. for proper operation. However, it is recognized that some experiments might result in some material being deposited on the sensor wire. Ice formation at some point on the sensor is a typical occurrence. Therefore, the Model 1700 has the capability of increasing the current for a short period of time at the beginning of the measurement cycle (in the **HOLDING** mode only) to try and drive the resistive zone of the sensor wire past the dirty region. This is termed **DIRTY SENSOR MODE**. This operation may or may not be successful depending on the degree of sensor contamination. This mode should be viewed as a stopgap measure only. If correct readings cannot be reestablished, the only choice is to warm the sensor or remove for cleaning or replacement.

NOTE Operation in the dirty sensor mode increases liquid helium losses. Consequently, operation in this mode should not be used unless the sensor is known or anticipated to become dirty or the helium level measurement is in question due to unclean operation.

To enable the **DIRTY SENSOR MODE**, Press the **DISABLED** button until **ENABLED** appears and press **SAVE**.

Press the home icon in the footer to return to the home screen.




Dirty Sensor Mode

ALARMS AND RELAYS

The Model 1700 Instrument has two types of alarms, level-based and time-based alarms.

ALARM STATUS SCREEN

Prior revisions of the Model 1700 have attempted to indicate the alarm conditions in the main screen footer region, including cycling through multiple alarm texts in the footer. This has been replaced by an Alarm Status screen that summarizes all alarm states and provides a button for muting.



Alarm Status Screen

The Alarm Status screen can be accessed at any time by tapping on the **ALARM** text in the footer. All alarms are summarized in the Alarm Status screen. If the Alarm Status screen shows the “ALARMING” button (as illustrated above), then touching the button will mute the alarm(s). This mute function will remain in effect until a new alarm is activated, or the status of any active alarm changes. Shortcuts may also be provided on the Alarm Status screen to enter the Level Alarms and Autofill configuration screens.

LEVEL-BASED ALARMS

The Model 1700 Instrument has two user-configurable level alarms. Each alarm can be triggered by either level measurement (for dual level configured instruments). Each alarm can be configured to be active above or below a user-defined setpoint.

When an alarm condition occurs, an audible alert will sound and visual indication will appear on the front panel.

Two relays can also be configured with independent setpoints to actuate independently of the alarm states. These relays have Normally Open (NO) contacts. The relays can be assigned to either the helium or nitrogen channel and each relay can be configured to close when the reading is either \leq or \geq the setpoint.

As an example of this setup flexibility, a level channel can be configured to have an *alarm condition* when the level is outside a normal operating band. The *relay actions* can be set to act identically, or configured for other conditions, per the user-defined application.

MULTIPLE ALARMS


More than one alarm condition can occur at the same time. If that is the case, the Alarm Status screen summarizes all the alarm states. Access the Alarm Status screen by touching the **ALARM** text in the footer.

NOTE Prior versions of the instrument flashed a series of alarms in the screen footer. This has been superseded by the Alarm Status screen.

CONFIGURING ALARM SETPOINTS


1. From the **MENU** screen, select **OUTPUTS** and the **OUTPUT CONFIGURATION: ALARMS** screen will be displayed.

(If already in another **OUTPUT CONFIGURATION** screen, choose the **ALARMS** button.)



Output ALARMS Configuration Screen


2. Ensure the Alarm Source fields are set to **HELIUM** (or **DISABLED**).



Relay Source Configuration


3. Set the levels at which the alarm will be triggered in the Setpoint fields highlighted in red in the image at right.

NOTE Tapping the units symbol (% , in, or cm) allows selection of the available unit selections for the selected source liquid.



ALARMS Setpoints

4. Use the \leq or \geq button to toggle between the two states of alarm, either alarm when the indicated level is less than or equal to the setpoint or alarm when the indicated level is greater than or equal to the alarm setpoint.



Alarm Condition Above or Below Setpoints

ACKNOWLEDGING AN ALARM

NOTE The alarms are not “latched” so if the alarm condition clears itself, the instrument will remove the alarm condition.

1. When an alarm is initiated, several things will occur:


- a. The **ALARM** text in the footer will illuminate in bold, red text.
- b. An audible alarm will be energized.
- c. The Alarm Status screen, accessible by touching the **ALARM** text in the footer, will display a summary of all alarms.



Alarm Annunciator

MUTING AN ALARM

To mute the audible alarm, first touch the **ALARM** text in the footer. The Alarm Status screen will appear with the “ALARMING” status shown. Tap the “ALARMING” button to mute all active alarms. A change in any alarm state will also re-enable the audible alarm for any remaining active alarms. There is no permanent mute.




Muted Alarm Condition

CONFIGURING RELAY SETPOINTS


The relay outputs are available via the AUX I/O connector on the rear panel of the instrument per the pinout provided on page 82.

1. From the **MENU** screen, select **OUTPUTS** and then choose **RELAYS**.



OUTPUT CONFIGURATION: RELAYS
Screen


2. Toggle the Relay Source fields to choose **HELIUM** (or **DISABLED**).



RELAYS Source Configuration


3. Set the levels at which the relay will be closed in the Setpoint fields.

NOTE Tapping the units symbol (% , in, or cm) allows selection of the available unit selections for the selected source liquid.



RELAYS Setpoints

4. Use the \leq or \geq button to toggle between the two states of actuation. Either the relay will close when the indicated level is less than or equal to the setpoint, or when the indicated level is greater than or equal to the setpoint.




Relay Closed Above or Below Setpoints

CONFIGURE THE AUTOFILL FUNCTION

There are several variables that must be addressed to set up an autofill system. These include the level indication that will be used to control autofill, the Fill Start level (A), and the Fill Stop Level (B).

SETTING THE AUTOFILL PARAMETERS


- From the **MENU** screen, choose **OUTPUTS** and then choose **AUTOFILL**. Toggle the **Auto-FILL CONTROL** button until **HELIUM** is displayed.



- Enter the Fill **START** and **STOP** levels.


NOTE If the **START** value is set to 0.0% the autofill function will be effectively disabled.

NOTE Tapping the units symbol (% , in, or cm) allows selection of the available unit selections for the selected source liquid.



- Enter the **FILL TIMEOUT** interval. Refer to page 34 for a description of the Autofill timeout function.

NOTE If the **FILL TIMEOUT** value is set to 0 minutes, the timeout function is disabled.



- Press **HOME** in the footer to return back to the level display.

ENABLE THE AUTOFILL FUNCTION

The autofill function must first be enabled. After the autofill control loop has been configured (the parameters in the prior section) the instrument will enter the **AUTO-OFF** state.

To enable the autonomous autofill function:

1. Touch the **AUTO-OFF** text in the display footer until it reads **AUTO-ON**.
2. Touch **SAVE** to enable the autonomous autofill function.

The AUTOFILL function has four manually-selectable states by touching the fill state indicator either in the **HOME** screen or in the footer when not in the **HOME** screen:

AUTOFILL States

Function	Operation	Overrides
AUTO-ON	Automatically maintains level between fill START and STOP set-point (i.e. autonomous mode).	Autofill will alarm and cease if fill valve stays open for \geq the FILL TIMEOUT setting .
M-OPEN	Manual Open: Energizes the valve control socket on the rear panel.	Any active ALARM with \geq setpoint defined for the fill liquid will close valve.
M-CLOSED	Manual Close: De-energizes the valve control socket on the rear panel.	None
AUTO-OFF	Disables the AUTOFILL function.	None


CLEARING THE AUTOFILL TIMEOUT ALARM

1. Touch the **TIMEOUT** text in the **HOME** screen or in the footer until **AUTO-ON** is displayed.
2. Touch **SAVE** in the footer. This will restart the autonomous autofill function.
3. To clear the **TIMEOUT** *without restarting* the autofill function, choose the **M-CLOSED** state and then choose **SAVE** in the footer.

SELECT THE APPROPRIATE UNITS ON THE DISPLAY

Touch the units on the display to change the units. The available units are percent (%), inches (in), and centimeters (cm).

NOTE Tapping the units symbol (% , in, or cm) in a configuration screen also allows selection of the available unit selections. The units selection is used for all indication and configuration values.




LHe Home Screen

ANALOG OUTPUT SIGNALS

Refer to page 82 of the Appendix for the AUX connector pin-out.

CONFIGURING THE ANALOG OUTPUTS

1. From the **MENU** screen, choose **OUTPUTS**, then **ANALOG OUTPUTS**.
2. If necessary, choose the source for the **0-10 VDC** output and **4-20 mA** output.
3. Press the **SAVE** button to save the choice (or **CANCEL** to quit without making a change).
4. Press **HOME** to go back to the home screen.



ETHERNET CONNECTIVITY

IP ADDRESSING SCHEME

1. From the **MENU** screen, choose **NETWORK**. The current settings will be displayed.

NETWORK CONFIGURATION

Addressing Scheme: **DYNAMIC** Address: 192.168.1.101
 Netmask: 255.255.255.0
 Gateway: 192.168.1.1

Hostname: Model1700 Network is up
 MAC Address: 00:D0:69:50:F5:8F

SAVE CANCEL

Network Configuration Screen

2. Touch the addressing scheme button, to choose **STATIC**, **DYNAMIC**, or **DISABLED** as appropriate.

NETWORK CONFIGURATION

Addressing Scheme: **DYNAMIC** Address: 192.168.1.101
 Netmask: 255.255.255.0
 Gateway: 192.168.1.1

Hostname: Model1700 Network is up
 MAC Address: 00:D0:69:50:F5:8F

SAVE CANCEL

Editing the Networking Mode

3. If **STATIC** is chosen for the addressing scheme, enter **ADDRESS**, **NETMASK**, and **GATEWAY** values that are appropriate for the connected network. Once an address field is touched, the pop-up keyboard will be presented for data entry. Press **SAVE** when done.

NETWORK CONFIGURATION

Addressing Scheme: **STATIC** Address: 192.168.1.101
 Netmask: 255.255.255.0
 Gateway: 192.168.1.1

Hostname: Model1700 Network is up
 MAC Address: 00:D0:69:50:F5:8F

HOME BACK ALARM Help

Editing Static Network Parameters


4. Set the **HOSTNAME** field as desired.
5. If the **DYNAMIC** addressing scheme is chosen, the **IP ADDRESS**, **NETMASK**, and **GATEWAY** addresses will automatically be assigned from a network DHCP server. Press **SAVE** when done.
6. If changes are made, select the **SAVE** button and then the instrument will reconfigure itself with the chosen network settings.

SERIAL CONNECTIVITY

The serial (RS-232) connectivity can be customized per the following instructions.


CONFIGURING THE RS-232 SETTINGS

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.



MENU Selection Screen

3. From the **SYSTEM CONFIGURATION** screen choose **RS232 SETUP**.




SYSTEM CONFIGURATION Screen

4. The RS-232 Configuration screen shows the current function for the serial connection.



RS-232 SERIAL PORT CONFIGURATION
Screen

5. Choose the alternate function for the serial connection if desired.



RS-232 SERIAL PORT CONFIGURATION
Screen

6. Choose to customize the port behavior for character echo, line ending, and baud rate.



RS-232 SERIAL PORT CONFIGURATION
Screen

7. If changes are made, select the **SAVE** button in the footer.


ABNORMAL OPERATION

DIRTY HELIUM SENSOR OPERATIONAL MODE

Refer to “Other Liquid Helium Functions” on page 27.

RESETTING THE INSTRUMENT TO FACTORY DEFAULTS

1. Press **MENU**, then **SYSTEM**, then **FACTORY RESET**.



System Configuration Screen

2. Press **RESTORE**.
3. The instrument will be reset to factory defaults.




Resetting Instrument to Factory Defaults

SHUTTING THE INSTRUMENT DOWN

NOTE The Model 1700 Instrument is a Linux-based computer system and in order to ensure the file system is properly unmounted, the **SHUTDOWN** function should be invoked. If it is not, i.e. the instrument is shut down by removing power via the front panel power switch, the next time the instrument boots up, it will have to scan the memory system to ensure everything is in order.

1. The instrument should be shut down by using the menu function **SHUTDOWN**.



Invoking Instrument Shut Down

2. Choose **SHUTDOWN** again to confirm.
3. When prompted, turn off the front panel power switch.




Confirming Instrument Shutdown

CALIBRATION

Model 1700 instrument is calibrated at the factory for a specific length sensor(s) for use in a specific liquid(s). The calibration length(s) and calibration liquid(s) are listed on the calibration sticker on the bottom of the instrument.


SETTING THE SYSTEM DATE AND TIME

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.




MENU Selection Screen

- From the **SYSTEM CONFIGURATION** screen, choose **SET TIME**.




SYSTEM CONFIGURATION Screen

- In the **SYSTEM DATE/TIME** screen, choose to either set the date and time *manually*, or set it automatically from an NTP server if the instrument has Internet access.



SYSTEM DATE/TIME Screen

- For the **SET MANUALLY** selection, edit the **TIMEZONE**, **YEAR**, **MONTH**, **DAY**, **HOUR**, and **MIN** fields as necessary. Touching in a field will launch the keyboard on the screen. Edit the information in the field as necessary and choose **ENTER** to enter the data in the field and close the pop up keyboard.



SYSTEM DATE/TIME Manual Entry

NOTE The clock is set to GMT at the factory and is battery backed. The **TIMEZONE** field must be set to a value from the TZ database.¹ The list of **TIMEZONE** values can be found at:

https://en.wikipedia.org/wiki/List_of_tz_database_time_zones

The **TIMEZONE** entry must exactly match the entire field as listed in the above URL in the *TZ database name* column and it is case sensitive. An instrument reboot may be required for the **TIMEZONE** to take effect. If the entered **TIMEZONE** supports Daylight Savings Time (DST), then the local time will be automatically adjusted to DST.

- Choose **SAVE** in the footer after all the fields have been edited as necessary.

1. The tz database is a collaborative compilation of information about the world's time zones, primarily intended for use with computer programs and operating systems. The topic of worldwide timezones is a complex issue so please refer to the on-line resource noted above.


SUPERCONDUCTIVITY-BASED (LIQUID HELIUM) LEVEL CALIBRATION

The instrument has been internally calibrated for use with AMI liquid helium level sensors at the AMI facility. However, it is still necessary to verify, and if necessary enter, the *correct active length* for the connected sensor.

NOTE If the instrument was purchased *with* a helium level sensor, the active length will be set at AMI prior to shipping.


VERIFY THE LIQUID HELIUM SENSOR TYPE

1. Press the **MENU** button in the lower left corner of the display screen.




Home Screen for Helium Level Display

2. Choose the **SENSORS** button from the **MENU** screen.



SENSORS Selection


3. Verify that the **HELIUM SENSOR TYPE:** (either **4.2K** or **2K**) is correct for the sensor to be used with the instrument. If not, contact the factory for assistance.



SENSORS Selection Button


SENSOR SAMPLE INTERVAL

1. Press the **MENU** button in the lower left corner of the display screen.




Home Screen for Helium Level Display

2. Choose the **SENSORS** button from the **MENU** screen.




SENSORS Selection

3. Choose the **CAL He** selection.



SENSORS Selection Button

4. Press in the **SENSOR SAMPLE INTERVAL** field. A pop up numeric keypad will be launched.
5. Enter the desired sample interval time in minutes and press **ENTER** and then **SAVE**.




Enter the SENSOR SAMPLE INTERVAL

SENSOR ACTIVE LENGTH

1. Press in the **SENSOR ACTIVE LENGTH** field. A pop up numeric keypad will be launched.
2. Enter the **SENSOR ACTIVE LENGTH** in the selected units and press **ENTER** and then **SAVE**.

NOTE Tapping the units symbol (in or cm) allows selection of the available unit selections for that channel.




Enter **SENSOR ACTIVE LENGTH**


CONTINUOUS MEASURE TIME LIMIT

The Continuous measure time limit feature sets a maximum time that the instrument will keep the liquid helium level sensor energized. If the instrument is left in the **CONTINUOUS** mode for the **CONTIN. MEASURE TIME LIMIT** interval, the instrument will revert back to the sample and **HOLDING** mode to prevent excessive liquid helium boil off.

1. Press in the **CONTIN. MEASURE TIME LIMIT** field. A pop up numeric keypad will be launched.
2. Enter the maximum time in minutes that the sensor should remain energized in the continuous mode. Entering a value of zero disables the feature.
3. Press **ENTER** on the keypad and then **SAVE** at the bottom of the screen.
4. Press the **HOME** button in the screen footer to return to the home screen.



Enter **CONTIN. MEASURE TIME LIMIT**



Footer **HOME** Button Selection


SENSOR NAME

1. Press the **MENU > SENSORS > SENSOR NAME(S)** buttons.




SENSOR NAME(S) Selection

2. Touch in the **HELIUM SENSOR NAME:** field. The keyboard will be displayed. Edit the displayed name and press **ENTER** and then **SAVE** at the bottom of the screen.



HELIUM SENSOR NAME: Field

3. Press **HOME** in the screen footer to return to the level display screen.



HOME Selection Button

REMOTE INTERFACE REFERENCE

The Model 1700 Instrument provides both serial (RS-232) and Ethernet interfaces as standard features. The serial and Ethernet interfaces can be operated simultaneously. Separate output buffers are provided for the serial and Ethernet return data. The serial and Ethernet interfaces accept both legacy and SCPI syntax for commands and queries.

The Model 1700 provides a web browser connection (port 80) via TCP/IP and presents a mirror of the local display which has hot spots that can be clicked or tapped, and which also accepts keyboard entry for textual and numerical input. The instrument's IP address can be ascertained by referring to the section titled "IP Addressing Scheme" on page 36. By using a web browser to connect to the instrument, all functionality of the Model 1700 can be controlled via the remote connection.

Communication with the Model 1700 via Ethernet or RS-232 (serial) interface can also be accomplished via the interactive mode of a commercially available terminal emulation program, or by socket (port 7180) or COM port programming. When the Model 1700 receives a terminated ASCII string, it always sends back a reply as soon as the string is processed. *When sending commands to the Model 1700, you should wait for the reply from the Model 1700 before sending another command even if the reply consists of only termination characters.*

SCPI COMMAND SUMMARY

The following manual conventions are used for SCPI (*Standard Commands for Programmable Instruments*) syntax for the remote interface commands:

- Braces { } enclose valid parameter choices.
- A vertical bar | separates multiple choices for each parameter.

- Triangle brackets < > indicate that you must supply a value.
- Parentheses () within < > indicate alternative units are available.
- Capitalized portions of the commands indicate acceptable abbreviations.

For example, the command `CONFigure:RELAY1:CHannel {0|1|2}` indicates that the command `CONFigure:RELAY1:CHannel` has three parameter options: 0, 1, or 2, with 0 (disabled) being the default value.

The following section is a reference list of SCPI commands. Refer to the detailed description of each command for information regarding specific parameter choices and their meanings. Default settings are shown in bold.

NOTE The commands listed below are the standards-compliant form of the remote interface SCPI language. A set of legacy commands which do not conform to the SCPI standards are also available for backward compatibility with *existing installations* of the Model 13x series (see page 67).

System-Related Commands

(see page 55 for more information)

*IDN?
*RST
*TST?

N2?
HE?

SERial_NUMber?
DATE_MANUfacture?
HardWare_VERsion?
FirmWare_VERsion?
Scpi_VERsion?

SYStem:BEEPer:IMMediate <seconds>
SYStem:BEEPer:STATe {0|1}
SYStem:BEEPer:STATe?

SYStem:KLOCK {0|1}
SYStem:KLOCK?

SYStem:DATE <yyyy>,<mm>,<dd>
SYStem:DATE?

SYStem:TIME <hh>,<mm>,<ss>
SYStem:TIME?

SYStem:REBOOT

IPV4addr?
MACADDR?
GATEWAY?
HOSTname?

SCREENCAP
SYStem:LOcal
SYStem:REMOte

Display Configuration Commands and Queries

(see page 58 for more information)

DISPlay:HE?
CONFigure:NAME:SENSor:HE "<string>"
NAME:SENSor:HE?

Relay Configuration Commands and Queries

(see page 58 for more information)

CONFigure:RELAy1:CHannel {0|1|2}
RELAy1:CHannel?

CONFigure:RELAy1:SETpoint <level>
RELAy1:SETpoint?

CONFigure:RELAy1:OPeration {0|1}
RELAy1:OPeration?

CONFigure:RELAy2:CHannel {0|1|2}
RELAy2:CHannel?

CONFigure:RELAy2:SETpoint <level>
RELAy2:SETpoint?

CONFigure:RELAy2:OPeration {0|1}
RELAy2:OPeration?

{RELAy1|RELAy2}:STATUs?

Alarm Configuration Commands and Queries

(see page 60 for more information)

CONFigure:ALArm1:CHannel {0|1|2}
ALArm1:CHannel?

CONFigure:ALArm1:SETpoint <level>
ALArm1:SETpoint?

CONFigure:ALArm1:OPeration {0|1}
ALArm1:OPeration?

CONFigure:ALArm2:CHannel {0|1|2}
ALArm2:CHannel?

CONFigure:ALArm2:SETpoint <level>
ALArm2:SETpoint?

CONFigure:ALArm2:OPeration {0|1}
ALArm2:OPeration?

{ALArm1|ALArm2}:STATUs?

ALARm:MUTE {0|1} *or* {NO|YES}
ALARm:MUTE?

Measurement Commands and Queries

(see page 62 for more information)

MEASure:HE:LEVel?
MEASure:HE:HOLD
MEASure:HE:CONTInuous
MEASure:HE:SAMPle
MEASure:HE:VOLTage?

MEASure:ADC0?
MEASure:ADC1?
MEASure:ADC2?

Fill Control and Queries

(see page 62 for more information)

CONFigure:FILL:CHannel {0|1|2}
FILL:CHannel?

CONFigure:FILL:A <level>
FILL:A?

CONFigure:FILL:B <level>
FILL:B?

CONFigure:INTerval:FILL <minutes>
INTerval:FILL?

HE Channel Sampling Commands and Queries

(see page 65 for more information)

CONFigure:INTerval:SAMPlE <minutes>
INTerval:SAMPlE?

CONFigure:HE:TIME_limit <minutes>
HE:TIME_limit?

LOG_LHE {0|1}
LOG_LHE?

HElIum_PowerSupply?
DIRTy_sen_mode?

Assignment Commands and Queries

(see page 64 for more information)

CONFigure:SOURCE:REC_out {0|1|2}
SOURCE:REC_out?

CONFigure:SOURCE:CURrent_LOOP {0|1|2}
SOURCE:CURrent_LOOP?

HE Channel Calibration Commands and Queries

(see page 65 for more information)

```
CONFigure:HE:LENgth <value>  
HE:LENgth?
```

Remote Units Commands and Queries

(see page 65 for more information)

```
CONFigure:HE:UNIT {0|1|2} or {PERCENT|INCH|CM}  
HE:UNIT?
```

RS-232 CONFIGURATION

An RS-232 serial communication port is available as a 9-pin D-type connector on the rear panel of the instrument for serial communication function.

The Model 1700 uses the following *fixed* parameters related to the RS-232 interface:

- *Parity*: No Parity
- *Data Bits*: 8 Data Bits
- *Number of Start Bits*: 1 bit
- *Number of Stop Bits*: 1 bit
- *Flow Control*: None

The baud rate default is 115200, but is adjustable in the SYSTEM CONFIGURATION > RS232 SETUP menu (see page 37).

SERIAL PORT CONNECTOR AND CABLING

A PC-compatible computer's serial port can be directly connected to the Model 1700 via a standard USB-to-serial cable. Refer to the computer's documentation to determine which ports are available on a computer and the required connector type.

The Model 1700 uses only three wires of the rear-panel DB9 connector: pin 2 (transmit), pin 3 (receive), and pin 7 (common). There are no software or hardware handshaking options. The Model 1700 is classified as a DCE (Data Communication Equipment) device since it transmits data on pin 3 and receives data on pin 2. The instrument to which the Model 1700 is attached must do the opposite, i.e., transmit on pin 2 and

receive on pin 3 (the requirements for a DTE, or Data Terminal Equipment device).

COMMAND/RETURN TERMINATION CHARACTERS

All commands are transmitted and received as ASCII values and are case insensitive. The Model 1700 can be configured to return `<CR>` (i.e. a *carriage return*) or `<CR><LF>` (i.e. a *carriage return* followed by a *linefeed*) at the end of a serial transmission (see page 37). The Model 1700 can accept `<CR>`, `<LF>`, `<CR><LF>`, or `<LF><CR>` as termination characters from an external computer.

ETHERNET CONFIGURATION

The Model 1700 provides a 10/100Base-T Ethernet interface as a standard feature. It complies with the IEEE 802.3u 100Base-TX and 802.3 10Base-T standards.

The Model 1700 allows its IP address, subnet mask and gateway IP address to be assigned either statically or dynamically. To make these values static and assign them manually, set IP Address Assignment to Static (see page 36) and then set the values using the Edit parameter list. To enable the values to be dynamically assigned by a network DHCP server, set IP Address Assignment to DHCP.

The system name (also known as *host name* or *computer name*), can be set using the Network screen (see page 36) or via remote communications (either Ethernet or RS-232).

All network parameters (even those assigned by a DHCP server) can be viewed using the Network submenu (see page 36).

ETHERNET CONNECTOR

The Model 1700 uses a standard RJ-45 jack for Ethernet communications. The Ethernet jack pinout is fully documented on page 80 in the *Appendix*.

TERMINATION CHARACTERS

All commands and queries are transmitted and received as ASCII values and are case insensitive. The Model 1700 always transmits `<CR><LF>` (a *carriage return* followed by a *linefeed*) at the end of an Ethernet transmission. The Model 430 can accept `<CR>`, `<LF>`, `<CR><LF>`, or `<LF><CR>` as termination characters from an external computer.

PORT ASSIGNMENT

The Model 1700 accepts remote connections to port 7180. Multiple connections to port 7180 are allowed.

COMMAND REFERENCE

All commands sent to the Model 1700 are processed and the Model 1700 responds with a return value (if applicable) and termination. If the command is invalid, the Model 1700 will respond with an error code (see the *Error Codes* section). All return values including error codes are terminated with `<CR><LF>` (i.e. a *carriage return* followed by a *linefeed*). For those commands that do not return a value, the Model 1700 will return the `<CR><LF>` termination only.

The remote units settings are saved in non-volatile memory and are restored at power-up.

The Model 1700 instrument may be configured for reading liquid nitrogen, liquid helium, or both. Some commands will not be applicable if the instrument is not configured for certain level measurement.

SYSTEM RELATED COMMANDS

- *IDN?

Returns the identification string of the Model 1700. The identification string contains the manufacturer name, model number, serial number, and firmware revision code. Example output:

```
AMERICAN MAGNETICS INC.,MODEL 1700,1700-16-002,2.4.32015
```

- *RST

Performs a Factory Restore if a restore file is available. All prior settings are lost!

- *TST?

Returns a value incremented by “1” for each query to the requesting interface if unit is functioning. Return value does not indicate any operational status other than a functioning interface.

- N2?

Returns “0” if the instrument is not configured to read liquid nitrogen level.

- HE?

Returns “0” if the instrument is not configured to read liquid helium level, “1” if the instrument is configured to read 4.2K liquid helium level for sensors of active length ≤ 40 inches, “2” if the instrument is configured to read 4.2K liquid helium level for sensors of active length ≤ 80 inches, “3”

if the instrument is configured to read 2K liquid helium level for sensors of active length ≤ 40 inches, or "4" if the instrument is configured to read 2K liquid helium level for sensors of active length ≤ 80 inches.

- SERIAL_NUMber?

Returns the serial number of the instrument as a string, for example: 1700-16-002.

- DATE_MANUfacture?

Returns the date of manufacture as a string in the form *yyyy-mm-dd*, for example: 2016-11-30.

- HardWare_VERsion?

Returns the hardware version of the instrument.

- FirmWare_VERsion?

Returns the firmware version of the instrument, for example: 2.4.32015.

- Scpi_VERsion?

Returns the SCPI command version of the instrument.

- SYStem:BEEPer:IMMediate <seconds>

The receipt of this command with a valid argument causes an audible tone to be generated by the instrument. The duration time parameter is specified in seconds.

- SYStem:BEEPer:STATe {0|1}

Enables/disables the beeper. When STATe "0" (OFF) is selected, no instrument condition, except the SYStem:BEEPer:IMMediate command, shall cause an audible beep to be emitted. When the instrument restarts, this value is reset to a default of "1" (ON).

- SYStem:BEEPer:STATe?

Returns present beeper state. A value of "0" indicates the beeper is disabled. A value of "1" indicates it is enabled.

- SYStem:KLOCK {0|1}

This command locks the local controls of an instrument if set to a value of "1" (ON). This includes all front panel, keyboard, or other local

interfaces. The state of the lock defaults to “0” (OFF) when the instrument is restarted.

- `SYStem:KLOCK?`

Returns the state of the instrument local controls lock. A return value of “0” indicates the lock is OFF. A return value of “1” indicates the local interfaces, such as the touch display, of the instrument are locked.

- `SYStem:DATE <yyyy>,<mm>,<dd>`

Sets the date per *<yyyy>* as a four-digit number; *<mm>* month with range of 1 to 12 inclusive; *<dd>* day of the month. Do not include the `<>` characters in the command.

- `SYStem:DATE?`

Returns the system date if the form *<yyyy>,<mm>,<dd>*.

- `SYStem:TIME <hh>,<mm>,<ss>`

This command is used to set the instrument's clock:

<hh> Range of hours is 0 to 23 inclusive.

<mm> Range of minutes is 0 to 59 inclusive.

<ss> Range of seconds is 0 to 60.

Do not include the `<>` characters in the command.

- `SYStem:TIME?`

Returns the system time if the form *<hh>,<mm>,<ss>*.

- `SYStem:REBOOT`

Reboots the instrument.

- `SYStem:RESTORE`

Reboots the instrument and sets all parameters back to factory defaults. All prior settings are lost!

- `IPV4ADDR?`

Returns the presently assigned TCP/IP version 4 address in the form:
`xxx.xxx.xxx.xxx`

- `MACADDR?`

Returns the MAC address of the network interface.

- GATEWAY?

Returns the Gateway address of the network interface.

- HOSTname?

Returns the hostname of the network interface.

- SCREENCAP

Snaps a TGA-formatted (TARGA) image of the current local display screen of the instrument and places it in the “Log files” folder of the instrument’s network share.

- SYSTem:LOCal

Enables the front panel touchscreen. All front panels controls are enabled by default after a power-up or REBOOT command.

- SYSTem:REMOte

Disables front panel touchscreen.

DISPLAY CONFIGURATION COMMANDS AND QUERIES

- DISPLAY:N2?

Returns a “0” if the instrument is not configured to display liquid nitrogen level on the home screen and a “1” if it is.

- DISPLAY:HE?

Returns a “0” if the instrument is not configured to display liquid helium level on the home screen and a “1” if it is.

- CONFigure:NAME:SENSor:HE "<string>"

Sets the displayed name of the helium level sensor.

- NAME:SENSor:HE?

Returns the displayed name of the helium level sensor.

RELAY CONFIGURATION COMMANDS AND QUERIES

- CONFigure:RELAy1:CHannel {0|1|2}

Assigns relay №1 as disabled (0), or to nitrogen (1) or helium (2).

- RELay1:CHannel?

Returns a "0" if relay №1 is disabled, a "1" if the relay is assigned to the nitrogen channel, and a "2" if the relay is assigned to the helium channel. An error return code is generated if attempting to assign the relay to a non-existent measurement channel in the present instrument configuration.

- CONFigure:RELay1:SETpoint <level>

Configures the relay №1 trip setpoint in the currently assigned channel's units.

- RELay1:SETpoint?

Returns the relay №1 setpoint in the current units.

- CONFigure:RELay1:OPeration {0|1}

Configures relay №1 such that it closes when the level is \leq the setpoint (0) or \geq the setpoint (1).

- RELay1:OPeration?

Returns a "0" if relay №1 closes when the level is \leq the setpoint and a "1" if the relay closes when the level is \geq the setpoint.

- CONFigure:RELay2:CHannel {0|1|2}

Assigns relay №2 to disabled (0), or to nitrogen (1) or helium (2). An error return code is generated if attempting to assign the relay to a non-existent measurement channel in the present instrument configuration.

- RELay2:CHannel?

Returns a "0" if relay №2 is disabled, a "1" if the relay is assigned to the nitrogen channel, and a "2" if the relay is assigned to the helium channel.

- CONFigure:RELay2:SETpoint <level>

Configures the relay №2 trip setpoint in the currently assigned channel's units.

- RELay2:SETpoint?

Returns the Relay №2 setpoint in the current units.

- `CONFigure:RELAy2:OPeration {0|1}`

Configures relay №2 such that it closes when the level is \leq the setpoint (0) or \geq the setpoint (1)

- `RELAy2:OPeration?`

Returns a “0” if relay №2 closes when the level is \leq the setpoint and a “1” if the relay closes when the level is \geq the setpoint.

- `{RELAy1|RELAy2}:STATus?`

Returns the present status of the specified relay. A return value of “0” indicates the relay is not closed. A value of “1” indicates the relay is closed.

ALARM CONFIGURATION COMMANDS AND QUERIES

- `CONFigure:ALArml:CHannel {0|1|2}`

Assigns alarm №1 as disabled (0), or to nitrogen (1) or helium (2).

- `ALArml:CHannel?`

Returns a “0” if alarm №1 is disabled, a “1” if the alarm is assigned to the nitrogen channel, and a “2” if the alarm is assigned to the helium channel. An error return code is generated if attempting to assign the alarm to a non-existent measurement channel in the present instrument configuration.

- `CONFigure:ALArml:SETpoint <level>`

Configures the alarm №1 trip setpoint in the currently assigned channel's units.

- `ALArml:SETpoint?`

Returns the alarm №1 setpoint in the current units.

- `CONFigure:ALArml:OPeration {0|1}`

Configures alarm №1 such that it alarms when the level is \leq the setpoint (0) or \geq the setpoint (1).

- `ALArml:OPeration?`

Returns a “0” if alarm №1 alarms when the level is \leq the setpoint and a “1” if it alarms when the level is \geq the setpoint. Alarm №1 is considered as the “legacy” HI level alarm when the alarm condition is set to level \geq a setpoint.

- `CONFigure:ALArm2:CHannel {0|1|2}`

Assigns alarm №2 to disabled (0), or to nitrogen (1) or helium (2). An error return code is generated if attempting to assign the alarm to a non-existent measurement channel in the present instrument configuration.

- `ALArm2:CHannel?`

Returns a "0" if alarm №2 is disabled, a "1" if the alarm is assigned to the nitrogen channel, and a "2" if the alarm is assigned to the helium channel.

- `CONFigure:ALArm2:SETPoint <level>`

Configures the alarm №2 trip setpoint in the currently assigned channel's units.

- `ALArm2:SETPoint?`

Returns the alarm №2 setpoint in the current units.

- `CONFigure:ALArm2:OPeration {0|1}`

Configures alarm №2 such that it alarms when the level is \leq the setpoint (0) or \geq the setpoint (1)

- `ALArm2:OPeration?`

Returns a "0" if alarm №2 alarms when the level is \leq the setpoint and a "1" if it alarms when the level is \geq the setpoint. Alarm №2 is considered as the "legacy" LO level alarm when the alarm condition is set to level \leq a setpoint.

- `{ALArm1|ALArm2}:STATus?`

Returns the present status of the specified alarm. A return value of "0" indicates the alarm is inactive. A value of "1" indicates the alarm is active.

- `ALARm:MUTE {0|1} or {NO|YES}`

Mutes the audible alarm for any active alarm(s) when set to "1" or "YES". Unmutes the audible alarm if set to "0" or "NO".

- `ALARm:MUTE?`

Returns "0" if audible alarm is muted. Returns "1" if not muted.

MEASUREMENT COMMANDS AND QUERIES

- `MEASure:HE:LEVel?`

Returns the liquid helium level in the current units.

- `MEASure:HE:HOLD`

Changes liquid helium level measurement operation from Continuous to Sample and Hold.

- `MEASure:HE:CONTInuous`

Changes liquid helium level measurement operation from Sample and Hold to Continuous mode.

- `MEASure:HE:SAMPlE`

Energizes the liquid helium level sensor, makes a reading and returns to Sample and Hold mode.

- `MEASure:HE:VOLTage?`

Returns the last measured liquid helium level sensor voltage in volts.

- `MEASure:ADC0?`

Returns the last measured liquid helium level sensor voltage in volts.

- `MEASure:ADC1?`

Returns the liquid helium sensor power supply voltage in volts.

- `MEASure:ADC2?`

Returns the liquid helium sensor excitation current in milliamperes.

FILL CONTROL AND QUERIES

- `CONFigure:FILL:CHannel {0|1|2}`

Assigns the auto fill control relay to either no channel (disabled) (0), nitrogen (1), or helium (2). An error return code is generated if attempting to assign the relay to a non-existent measurement channel in the present instrument configuration.

- `FILL:CHannel?`

Returns a "0" if the auto fill relay is disabled, a "1" if the relay is assigned to the nitrogen channel, and a "2" if the relay is assigned to the helium channel.

- `CONFigure:FILL:A <level>`

Sets the A setpoint (control band upper limit) in the assigned channel's current units. The A setpoint (autofill stop) must be greater than the B setpoint (fill start) and must also be between 0% and 100%.

- `FILL:A?`

Returns the A setpoint (control band upper limit) in the assigned channel's current units.

- `CONFigure:FILL:B <level>`

Sets the B setpoint (control band lower limit) in the assigned channel's current units. The B setpoint (autofill start) must be less than the A setpoint (fill stop) and must also be between 0% and 100%.

- `FILL:B?`

Returns the B setpoint (control band lower limit) in the assigned channel's current units.

- `CONFigure:INTERval:FILL <minutes>`

Sets the fill timer in minutes. If the level does not reach or exceed the A setpoint within the fill time, the fill is terminated. Setting the value of FILL to "0" disables the fill timer function.

- `INTERval:FILL?`

Returns the fill timer setting in minutes.

HE CHANNEL SAMPLING COMMANDS AND QUERIES

- `CONFigure:INTERval:SAMPlE <minutes>`

Sets the liquid helium sampling interval in minutes.

- `INTERval:SAMPlE?`

Returns the sampling interval in minutes if the instrument is configured with a helium channel.

- `CONFigure:HE:TIME_limit <minutes>`

Sets the continuous measurement time limit for liquid helium measurements in minutes, after which time the instrument will return to sampling mode.

- `HE:TIME_limit?`

Returns the continuous measurement time limit in minutes.

- `LOG_LHE {0|1}`

Enables (1) or disables (0) liquid helium level logging.

- `LOG_LHE?`

Returns “0” if the liquid helium level logging is disabled, or “1” if enabled.

- `HElium_PowerSupply?`

Returns the helium power supply configuration of the instrument where “0” is not installed, “1” is standard, and “2” is XL version.

- `DIRTy_sen_mode?`

Returns a “0” if the instrument’s dirty sensor mode is disabled for the helium channel and a “1” if it is enabled.

ASSIGNMENT COMMANDS AND QUERIES

- `CONFigure:SOURCE:REC_out {0|1|2}`

Configures the 0-10 V_{DC} Recorder Output source to disabled (0), assigned to the nitrogen channel (1), or the helium channel (2). An error return code is generated if attempting to assign the output to a non-existent measurement channel in the present instrument configuration.

- `SOURCE:REC_out?`

Returns a “0” if the 0-10 V_{DC} Recorder Output is disabled, a “1” if it is configured for the nitrogen channel, and a “2” if it is configured for the helium channel.

- `CONFigure:SOURCE:CURrent_LOOP {0|1|2}`

Configures the 4-20 mA Current Loop output source to disabled (0), assigned to the nitrogen channel (1), or the helium channel (2). An error return code is generated if attempting to assign the output to a non-existent measurement channel in the present instrument configuration.

- `SOURCE:CURrent_LOOP?`

Returns a “0” if the 4-20 mA Current Loop output is disabled, a “1” if it is configured for the nitrogen channel, and a “2” if it is configured for the helium channel.

HE CHANNEL

CALIBRATION

COMMANDS AND

QUERIES

- `CONFigure:HE:LENgth <value>`

Configures the liquid helium sensor active length in current units. Returns an error if the current units are percent.

- `HE:LENgth?`

Returns the configured active length of the liquid helium sensor in current units. Returns an error code if the current units are percent.

REMOTE UNITS

COMMANDS AND

QUERIES

- `CONFigure:HE:UNIT {0|1|2} or {PERCENT|INCH|CM}`

Sets the liquid helium level units of measurement to percent (0 or `PERCENT`), inches (1 or `INCH`), or centimeters (2 or `CM`). The default is `PERCENT` units.

- `HE:UNIT?`

Returns the current liquid helium level units in use as either `C`, `I`, or `%`.

ERROR CODES

The Model 1700 returns specific error codes for invalid commands and/or arguments. If an error condition is returned, the command is not processed and the configuration of the instrument is not modified. The table below provides a list of error codes, their meaning, and any associated limits.

Error Code	Meaning	Valid Range
-1	LO (or relay №2 ^a) setpoint out of range	$0 \leq \text{LO (or relay №2}^{\text{a}}) \leq \text{LENGTH}$
-2	Fill B setpoint (fill start) out of range	$0 \leq B < A$
-3	Fill A setpoint (fill stop) out of range	$B < A \leq \text{LENGTH}$
-4	HI (or relay №1 ^a) setpoint out of range	$0 \leq \text{HI (or relay №1}^{\text{a}}) \leq \text{LENGTH}$
-5	Attempted to set or query for LENGTH in PERCENT units mode	
-6	Invalid argument, value out of maximum calibration range	$0 < \text{value} \leq 213.36 \text{ cm (HE)}$
-7	INTERVAL setting out of range	$0 \leq \text{INTERVAL} \leq 999 \text{ min}^{\text{b}}$ (or) $0 \leq \text{INTERVAL} \leq 99999 \text{ min}^{\text{c}}$
-8	Unrecognized command	
-9	Invalid argument, value was negative or non-numeric	
-11	Command exceeds SCPI input buffer limit	256 characters, including spaces, etc.
-12	Command invalid for selected channel or interface in present configuration	

a. Applies to dual instrument configuration

b. For Version HW = 4712

c. For Version HW = 4900

LEGACY COMMAND REFERENCE

The legacy command set is included for compatibility with existing Model 13x installations. New installations should consider using the modern command set that conforms more closely with the SCPI specification. The modern and legacy command sets can be mixed as desired.

All legacy command or queries sent to the Model 1700 are processed and the Model 1700 responds with a return value (if applicable) and termination. If the command is invalid, the Model 1700 will respond with an error code (see the *Error Codes* section). All return values including error codes are terminated with <CR><LF> (i.e. a *carriage return* followed by a *linefeed*). For those commands that do not return a value, the Model 1700 will return the <CR><LF> termination only.

COMMANDS FOR CONTROLLING THE UNITS OF MEASUREMENT

Command:	CM	Function:	Sets the units of measurement to centimeters	Returns:	<CR><LF>
Command:	INCH	Function:	Sets the units of measurement to inches	Returns:	<CR><LF>
Command:	PERCENT	Function:	Sets the measurement to % of sensor length	Returns:	<CR><LF>
Command:	UNIT	Function:	Returns the current units in use	Returns:	C, I, or % <CR><LF>

The CM command sets the units of measurement to centimeters and the INCH command selects inches. The PERCENT command sets the units of measurement to the percentage of sensor length that is immersed in liquid helium. **The units of measurement selected through the remote interface are controlled independently from the units shown on the front panel display.**

The remote units setting is *automatically* saved in permanent memory and is restored at power-up. The UNIT command returns a one character value (and termination) indicating the current units — C for centimeters, I for inches, or % for percentage.

COMMANDS FOR CONFIGURING PERMANENT MEMORY

Command:	HI=<value>	Function:	Configures the HI setpoint limit	Returns:	<CR><LF>
Command:	LO=<value>	Function:	Configures the LO setpoint limit	Returns:	<CR><LF>
Command:	A=<value>	Function:	Configures the A setpoint (control band upper limit)	Returns:	<CR><LF>
Command:	B=<value>	Function:	Configures the B setpoint (control band lower limit)	Returns:	<CR><LF>
Command:	INTERVAL=<value>	Function:	Configures the sampling interval in minutes	Returns:	<CR><LF>
Command:	LENGTH=<value>	Function:	Configures the active sensor length	Returns:	<CR><LF>
Command:	SAVE	Function:	None (for backward compatibility only).	Returns:	<CR><LF>

The HI and LO command configure the high and low setpoint limit values, respectively. For example, HI=90.0 would configure the high setpoint limit to 90.0 in whichever units of measurement last selected through the serial interface. The A and B commands configure the upper limit and lower limit of the control band, respectively. The HI, LO, A, and B commands are compatible with the percent units selection.

NOTE The HI setting is associated with hardware relay №1. The LO setting is associated with the hardware relay №2. This association is fixed and cannot be changed.

The LENGTH command configures the active sensor length setting in the current units. LENGTH=35.0 would configure the active sensor length to 35.0 units of centimeters or inches.

NOTE The LENGTH=<value> command will only function if CM or INCH are currently selected as the units of measurement. The LENGTH command does not configure the Model 1700 if the units of measurement are PERCENT.

The INTERVAL command sets the sampling interval used when the instrument is set for sampling mode operation (as opposed to continuous update). INTERVAL=0.1 sets the sampling interval to 0.1 minutes.

The HI, LO, A, B, INTERVAL, LENGTH, and current remote units settings are automatically saved to permanent memory. Saved settings are then recalled each time the power is turned off and then reapplied to the instrument on subsequent power on.

QUERYING THE CONFIGURATION

Command:	HI	Function:	Returns the HI setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	LO	Function:	Returns the LO setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	A	Function:	Returns the A setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	B	Function:	Returns the B setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	INTERVAL	Function:	Returns the sampling interval setting in minutes	Returns:	<value> <CR><LF>
Command:	LENGTH	Function:	Returns the active sensor length in the current units	Returns:	<value> <CR><LF>

The HI, LO, A, B, INTERVAL, and LENGTH queries return the current configuration of the instrument. Each return value is terminated with <CR><LF>.

RETURNING A LEVEL MEASUREMENT

Command:	LEVEL	Function:	Returns the liquid level in the current units	Returns:	<value> <CR><LF>
----------	-------	-----------	---	----------	---------------------

The LEVEL query returns the liquid helium level in the current units selected through the communication interface.

SERVICE AND REPAIR

CLEANING

To prevent electrical shock, disconnect the instrument from AC mains power and disconnect all connected wiring before cleaning. Clean the outside of the instrument using a soft, lint-free, cloth slightly dampened with water.

Do not use detergent or solvents.

Do not attempt internal cleaning.

USER REPLACEABLE PARTS

Replacement parts for the instrument are listed in the table below.

AMI Part Number	Description
HG0128	Instrument foot
SA 1045	Single Rack Mount Kit
SA 1046	Dual Rack Mount Kit
EF1700	Fuse, 3 A, 250 Vac, 5x20 mm, fast acting, UL/CSA recognized.
HG0005	Battery, 3V lithium, 20mm x 3.2 mm coin cell; CR2032.

BATTERY REPLACEMENT

This section describes the procedure for replacing the battery on the instrument's main circuit board.


WARNING This procedure should only be performed by a technician who is familiar with electronic instrumentation and trained in electrical safety and ESD precautions. Always disconnect the power cord and any external wiring before removing the instrument cover.

Always disconnect all inputs, cords, and cables before disassembling the instrument.

LOW BATTERY INDICATION


Battery status is displayed in the **SYSTEM DATE & TIME** dialog. To access the dialog:

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.




MENU Selection Screen

3. From the **SYSTEM CONFIGURATION** screen choose **SET TIME**.



SYSTEM CONFIGURATION Screen

4. The battery status is displayed in the **SYSTEM DATE & TIME** screen. If the battery is weak and should be replaced soon, the message will appear as shown at right.



Weak Battery Indication

5. If the battery is expired, the message will appear as shown at right.

NOTE The NOTIFY button state will cause a warning message to appear in the the HOME screen footer. NO NOTIFY supresses the HOME screen warning.



Expired Battery Indication

The following sections detail how to replace the battery.

TOOLS REQUIRED

- Torx Plus (T.M.) size 10 driver (Wera 028034 or equivalent)
- Torx Plus (T.M.) size 15 driver (Wera 028035 or equivalent)
- Small, flat-blade screw driver (for prying)

PROCEDURE

1. Unplug the instrument from the AC power source.
2. Using the T-15 driver, remove the four 8-32 machine screws on the sides of the instrument cover. Set these screws aside as they will be re-used.
3. Using the T-10 driver, remove the four 6-32 machine screws on the rear of the instrument cover. Set these screws aside as they will be re-used.
4. Lift the instrument cover off of the instrument chassis and set aside.

5. Using the small, flat-blade screwdriver, carefully pry the battery from the holder BH1.
6. Install the new battery into the battery holder BH1.
7. Replace the top cover and secure using the eight machine screws which were removed previously.

FUSE REPLACEMENT

This section describes the procedure for replacing the two fuses on the instrument's main circuit board.

WARNING This procedure should only be performed by a technician who is familiar with electronic instrumentation and trained in electrical safety and ESD precautions. Always disconnect the power cord and any external wiring before removing the instrument cover.

Always disconnect all inputs, cords, and cables before disassembling the instrument.

TOOLS REQUIRED

- Torx Plus (T.M.) size 10 driver (Wera 028034 or equivalent)
- Torx Plus (T.M.) size 15 driver (Wera 028035 or equivalent)
- Small, flat-blade screw driver (for prying)


PROCEDURE

1. Unplug the instrument from the AC power source.
2. Using the T-15 driver, remove the four 8-32 machine screws on the sides of the instrument cover. Set these screws aside as they will be re-used.
3. Using the T-10 driver, remove the four 6-32 machine screws on the rear of the instrument cover. Set these screws aside as they will be re-used.
4. Lift the instrument cover off of the instrument chassis and set aside.
5. Using the small, flat-blade screwdriver, carefully pry the fuse(s) from the fuse holders F1 and/or F2.
6. Install the new fuse(s) into the fuse holder(s) F1 and/or F2.
7. Replace the top cover and secure using the eight machine screws which were removed previously.

FIRMWARE UPGRADE VIA ETHERNET

The following procedure can be utilized by users to upgrade the Model 1700 by downloading firmware upgrade archives (.tar.gz) from the AMI support website.

1. Connect the instrument to a computer network via Ethernet cable and energize the instrument.
2. When the instrument has completed booting, note the IP Address: **MENU > NETWORK > Address**. There should be a message in green on the screen indicating "Network is up".




Network Configuration Screen


3. Determine the currently installed version of firmware. From the home page, press **MENU > SYSTEM > FW Version**.
4. Download the available new version of firmware and place on a PC on the same network as the instrument.

NOTE The following steps assume the computer is Windows-based.

5. Start Windows Explorer and enter the instrument's IP address noted from Step 2 above in the address field. Be sure to use the format such as \\192.168.1.101. Ensure the Explorer screen looks similar to the following.



- Open a second Windows Explorer and locate the previously downloaded firmware file:




- Drag the firmware file to the instrument's Application Update folder in the first Explorer screen.
- Install the new firmware by pressing **MENU > SYSTEM > UPDATE** on the instrument and choose the name of the firmware file copied to the instrument. It is likely the top choice of the files listed.
- Press **INSTALL** on the instrument when prompted.
- The instrument will install the new firmware and reboot when completed.

UPGRADE VIA SCP

If the instrument's file shares are not visible on the local network, it is possible to copy the firmware upgrade directly to the instrument using the `root` login and password using SCP (Secure Copy) over Ethernet.

- Connect the instrument to the local computer network via Ethernet cable and energize the instrument.
- When the instrument has completed booting, note the IP Address: **MENU > NETWORK > Address**. There should be a message in green on the screen indicating "Network is up".



Network Configuration Screen

- Determine the currently installed version of firmware. From the home page, press **MENU > SYSTEM > FW Version**.
- Download the available new version of firmware and place on a Windows, Linux, or Mac computer on the same network as the instrument.

5. Open a command line application and set the current directory to the folder where the newly download firmware from the prior step resides.
6. Enter the following command where the example IP address of 192.168.1.101 is replaced by the IP of the instrument noted in step 2 above, and the example firmware file of Ver-2.22719.15.tar.gz is replaced by the actual file downloaded from AMI:

```
scp ./Ver-2.220719.15.tar.gz root@192.168.1.101:/update
```

Note that you will be prompted to enter the root account password to complete the file transfer. If you do not know the password, contact an Authorized AMI Technical Support Representative for assistance.

7. Install the new firmware by pressing **MENU > SYSTEM > UPDATE** on the instrument and choose the name of the firmware file copied to the instrument via SCP. It is likely the top choice of the files listed.
8. Press **INSTALL** on the instrument when prompted.

The instrument will install the new firmware and reboot when completed.

APPENDIX

SERIAL (RS-232) CONNECTOR




Serial (RS-232) Pin Out

The RS-232 connector is a 9-pin D-sub female connector to connect standard DTE 9-pin D-sub male connector using a standard straight (not NULL) cable.

Serial (RS-232) Pin Definitions

Pin	Mnemonic	Function
1	N/C	
2	TXD	Transmit Data
3	RXD	Receive Data
4	N/C	
5	GND	Signal Ground
6	N/C	
7	N/C	
8	N/C	
9	N/C	

ETHERNET CONNECTOR




Ethernet Connector Socket Pin out

Ethernet RJ-45 Connector Pin Definitions

Pin	Mnemonic	Function
1	TXD+	Transmit differential output +
2	TXD-	Transmit differential output -
3	RXD+	Transmit differential input +
4	not used	
5		
6	RXD	Transmit differential input -
7	not used	
8		

LIQUID HELIUM CONNECTOR J1 WIRING



Liquid Helium (LHe) Level Sensor Connector

The LHe connector is a 9-pin D-sub female connector.

LHe Level Connector Pin Definitions

Pin	Function	AMI Standard Wiring Color
1	Sensor I+	Red
2	+48 VDC P/S Test Point	
3	+5 VDC P/S Test Point	
4	-12 VDC P/S Test Point	
5	P/S Test Points Common	
6	Sensor V-	Yellow
7	Sensor I-	Black
8	Sensor V+	Blue
9	+12 VDC P/S Test Point	

AUX I/O CONNECTOR



Aux I/O Connector

The AUX I/O connector is a 9-pin D-sub female connector.

Aux I/O Pin Definitions

Pin	Function	Polarity
1	4-20 mA Current Loop	+
2		—
3	0-10 VDC Output	+
4		—
5	Relay № 1 Dry Contact	N/A
6		
7	Relay № 2 Dry Contact	N/A
8		
9	N/A	

NOTE If the Model 1700 is being used to replace a Model 13x instrument, a M-M 9-pin gender changer will be necessary to change the female AUX I/O connector on the Model 1700 to a male type.

TROUBLESHOOTING

The following paragraphs serve as an aid to assist the user in troubleshooting a potential problem with the Model 1700 Instrument. If the user is not comfortable in troubleshooting the system, contact an AMI Technical Support.

If the cause of the problem cannot be located, contact an AMI Technical Support Representative at +1 (865) 482-1056 for assistance. The AMI technical support group may also be reached by internet e-mail at:

support@americanmagnetics.com

INSTRUMENT DISPLAYS “LOSS OF SENSOR” CONDITION FOR LHe LEVEL MEASUREMENT

This message indicates that the instrument cannot get sensor excitation current to flow in the sensor circuit¹. Check that the sensor is wired correctly and that the sensor is connected to the proper connector on the back of the instrument. When the wiring is corrected, place the instrument in MEASURE CONTINUALLY mode and the Loss of Sensor condition should be cleared at which time the instrument can be placed back in SAMPLE & HOLD mode if desired.

INSTRUMENT DISPLAYS “SENSOR PROTECTED” CONDITION FOR LHe MEASUREMENT

This message indicates that the LHe level sensor circuit is of a higher impedance than expected for the selected active sensor length.² The protection means that the normal sensor excitation has been stopped and the sensor is intermittently pulsing the current to see if the resistance is back in the expected range. This condition is usually seen when the instrument is set for an active length greater than the level sensor connected to the instrument. See “Sensor Active Length” on page 45. This condition will also occur if a LHe level sensor is energized in a vacuum since the level sensor’s sensing element will heat up to an abnormally high temperature and with its positive temperature coefficient, present a larger than expected resistance to the instrument.

SYSTEM TEST SCREENS


The instrument provides various system test functions that can be performed from the front panel touchscreen. The system test functions are not available to the web browser interface nor the SCPI remote

1. The threshold is approximately 1 mA.

2. The threshold is 5% greater resistance than the SENSOR ACTIVE LENGTH resistance.


interface over serial or Ethernet connections. To access the test functions:

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.




MENU Selection Screen

3. From the **SYSTEM CONFIGURATION** screen choose **SYSTEM TEST**.



SYSTEM CONFIGURATION Screen


4. In the **SYSTEM TEST** screen, several selections are available to test the touchscreen colors. Pressing a pixel test button will cause the screen to show a scrolling test pattern. When touched again, the display will then return to normal.



SYSTEM TEST Pixel Tests


5. To calibrate or test the touch function of the screen, use the **TOUCH SCREEN CALIBRATE** or **TOUCH SCREEN TEST** selections in the **SYSTEM TEST** screen.

The calibrate function will display a screen with touch targets. Touch all the targets as they are displayed to complete the calibration. The instrument will then automatically store the new touch calibration.




SYSTEM TEST Touch Cal/Test

6. The **HARDWARE CONTROL** selection in the **SYSTEM TEST** screen allows the user to observe in real time and/or manually test various functions of the instrument hardware. Press the **HARDWARE CONTROL** button to move to the detailed test screen.



Hardware Test Selection

7. The **HARDWARE CONTROL** screen includes several controls to manually exercise the hardware. However, to do so the operator must press the Control button to toggle it from **AUTO** to **MANUAL** mode.
8. Once in **MANUAL** mode, the user can touch any of the buttons to toggle the state of the hardware. The **Analog Outputs** button is, however, not a toggle but will move the operator to the **ANALOG OUTPUTS TEST** screen.




Hardware Control Screen

NOTE In the **AUTO** mode, the Test screen displays the real time measurement or state of each labeled function.

NOTE When changing to **MANUAL** test mode, please be aware that the automatic functions of the instrument will be overridden. Functions such as autofill will be interrupted, or a fill initiated, by the manual actions. Be sure the system is in a *safe condition* before exercising manual control.


9. The **ANALOG OUTPUTS TEST** screen requires the operator press the Control button to toggle from **AUTO** to the **MANUAL** mode to exercise control over the output values. The output values range from 0 to 4095 where 0 is the minimum output value and 4095 is the maximum output value. Tap a value to enter a manual override.



Analog Outputs Tests Screen

NOTE The hardware control state will return to **AUTO** after a few seconds if the operator does not explicitly return to automatic control. However, the state of the hardware or autofill functions are *not guaranteed to be the same* as when the **MANUAL** control state was entered.

10. Finally, the **CURSOR ON/OFF** button in the **SYSTEM TEST** screen allows the operator to turn the cursor off and on. The cursor display can be useful for verifying the touchscreen response if it appears to function incorrectly.



CURSOR Display Selection


11. Use the **HOME** button to return to the main level display.

SYSTEM LOGS

The Model 1700 provides a logging feature that saves detailed information regarding the liquid level and system events. The information can be very useful in troubleshooting suspected system errors.

The Model 17XX family of liquid level instruments keeps log files in a local Ethernet network accessible folder via Samba file shares. The log

folder is labeled “Log files” and can found on a network by opening a file explorer at the instrument IP address or hostname.



In the Log files folder, there will be a variety of logging files:

He.log: Contains a log of liquid helium measurement events, with each entry timestamped and a record of the hardware bits at that time.

N2-1.log: Contains a log of liquid nitrogen measurement events, with each entry timestamped and a record of the hardware bits at that time.

operations.csv: Contains a log of all operations/events, with each entry timestamped and appropriate data for the event in a comma separated file format.

LOG FILE FORMATS

The text based logging format is in a comma separated value format. It consist of a unix timestamp, followed by the current level, and then the status bits encoded in a hexadecimal number.

The status bits contains information for both the channels, even if the channels are disabled.

The status bits are a summation of the hardware flags, each flag has a unique bit position, and so the status byte displayed is a hexadecimal representation of those bits.

Refer to the Status Bits table on page 91 for the bit values assigned for each hardware flag.

A log entry is made when either the level changes by 0.1% or if the status bits change from last entry. Levels and bits are checked every second.

A Unix timestamp is the time in seconds since Jan 1, 1970 in Coordinated Universal Time.

DESCRIPTION OF THE HE.LOG FORMAT

Here is an example of a measurement, with the helium level changing from 26.3 to 45.1:

```

      timestamp, %level, status bits
#1  1626706702,  26.3,028
#2  1626706703,  45.1,038

```

Entry #1: timestamped at Mon Jul 19 2021 10:58:22 (Eastern Daylight Time) with a previously read level of 26.3. The status bits (028) indicated that the external/internal N2 sensor oscillator selection was internal, and that the helium sensor was active and dirty sensor mode was off. The level updated when the helium measurement cycle completed, as recorded in the next entry.

NOTE Note that even if an instrument is not configured for nitrogen measuring, there will be inactive status bits for that nitrogen channel in the status word.

Entry #2: timestamped at Mon Jul 19 2021 10:58:22 with a last read level of 45.1. The status bits (028) indicated that the external/internal N2 sensor oscillator selection was internal, and that the helium sensor was active and dirty sensor mode was off. The level had changed between entry #2 and #3, so when the unit stopped measuring (turning off the helium sensor) it updated the level, which changed to 45.1

Refer to the Status Bits table on page 91 for the bit values assigned for each hardware flag.

DESCRIPTION OF THE N2-1.LOG FORMAT

Here is an example of measurements of changing level, with a loss of sensor event for the liquid nitrogen measurements:

```

      timestamp, %level, status bits
#1  1626708988,  88.2,000020
#2  1626708989,  77.8,000020
#3  1626708990,  54.1,000020
#4  1626708991,  41.3,000020
#5  1626708992,  34.4,000020
#6  1626708994,  16.1,000020
#7  1626708996,   0.0,000221

```


#8 1626708999, 16.1,000020

Entry #1: timestamped at Mon Jul 19 2021 11:36:28 (Eastern Daylight Time) with a read level of 88.2%. The status bits (000020) indicated that the external/internal N2 sensor oscillator selection was internal.

Entries #2-#6: shows the level changing over a 5 second span, settling at 16.1%

Entry #7: timestamped at Mon Jul 19 2021 11:36:36 showed a loss of sensor event. The level will be 0.0%, and the status bits show that the sensor was disconnected and the alarm bit was set.

Entry #8: timestamped at Mon Jul 19 2021 11:36:36 showed back to normal operation, with a reading of 16.1%.

Refer to the Log Status Bits table on page 91 for information on the status bits.

DESCRIPTION OF THE OPERATIONS LOGGING FILE

The `operations.csv` file is a text file that records a log of all operations/events: remote, local, and SCPI.

This file contains operations log entries in a comma separated value format, as follows:

```
<unix timestamp>, <class>, <data field 1>, <data field 2>, <data field 3>
```

A data field can be empty, that is, no text.

An overview of the classes:

Class	Description of Event
AF	Master System Autofill
AL	Alarm(s)
BT	Battery Status
CE	System commands
N1	Liquid Nitrogen channel sensor
LE	GUI page load event
LF	GUI link (AML function link)
LL	GUI link (AML page link)
LS	GUI link (system variable)
RY	Hardware Relay Status
SC	SCPI command that modified the instrument
SF	State file operation
PU	Power Up
TG	GUI toggle

VM Modification of a instrument variable

More detailed information on each class is available from your Authorized AMI Technical Representative.

Here is a simple example of operations logging file entries for a shutdown initiation:

```

    timestamp, class,data1,data2,data3
#1 1626712173,LL,MENU,menu.aml,
#2 1626712191,LL,SHUTDOWN,../shutdown.aml,
#3 1626712200,LL,SHUTDOWN ,../shutdown.aml,
#4 1626712200,SF,shmem_dump,checksum,6A7E48E4
#5 1626712200,CE,shutdown,,

```

Line by line description:

```
1626712173,LL,MENU,menu.aml,
```

The operator pressed the MENU button on the home page.

Menu.aml indicated the screen page loaded was menu.aml, that is, the MENU page.

The class LL indicates that this was a link label button event, that is, the label on the button was MENU, and that the linked page loaded was menu.aml.

```
1626712191,LL,SHUTDOWN,../shutdown.aml,
```

The operator pressed the SHUTDOWN menu button which then loaded the ../shutdown.aml page.

```
1626712200,LL,SHUTDOWN ,../shutdown.aml,
```

The operator touched the SHUTDOWN command button which then executed the shutdown procedure.

```
1626712200,SF,shmem_dump,checksum,6A7E48E4
```

The instrument saved the variables to the state file. Class SF indicates this is a state file operation.

```
1626712200,CE,shutdown,,
```

The command executed (CE) was the shutdown command. The data fields #2 and #3 in this line are empty as shown by the two commas after the shutdown command.

LOG STATUS BITS TABLE**Status Bit Representation (in hexadecimal)**

STATUS_Alarm	0x00001
STATUS_relay1	0x00002
STATUS_relay2	0x00004
STATUS_Ext_1_osc	0x00020
STATUS_Valve_1_open	0x00040
STATUS_N2_1_LossOfSensor	0x00200
STATUS_N2_1_SensorShorted	0x00400
STATUS_PowerUp	0x00800
STATUS_He_1_Burnout_Timeout	0x04000
STATUS_He_1_LossOfSensor	0x08000

GLOSSARY

ABBREVIATIONS AND ACRONYMS

Term	Meaning
AC; ac	Alternating Current; strictly, electrical <i>current</i> that periodically reverses direction. Typically used also to describe an electrical power source in terms of the <i>voltage</i> . For example, 240 VAC.
ASCII	American Standard Code for Information Interchange; numerical representation of characters such as 'a' or '@' or an action (such as line-feed); 'plain' raw text with no formatting such as tabs, bold or underscoring
BNC connector	A miniature quick connect/disconnect RF connector used for coaxial cable, featuring two bayonet lugs on the female connector.
CR or <CR>	Text Carriage-Return character
Cryogen	A substance for obtaining low temperatures. In the case of use with the Model 1700 instrument, a cryogen is a liquefied gas such as liquid nitrogen or liquid helium.
D-Sub	Term referring to the family of connectors containing an odd number of pins in two parallel rows with a 1-pin difference in pins-per-row (DB9, DB15, and DB25 are most common)
DB9	Type of electrical connector containing 9 pins arranged in two parallel rows of 4 pins and 5 pins each
DCE	Data Circuit-terminating Equipment - a device that sits between the Data Terminal Equipment (DTE) and a data transmission circuit.
DHCP	Dynamic Host Configuration Protocol; a computer networking protocol which dynamically distributes the IP address to networked devices
dt	Rate of change
DTE	Data Terminal Equipment - an end instrument that converts user information into signals or reconverts received signals. A DTE device communicates with the Data Circuit-terminating Equipment (DCE).

Term	Meaning
ECL	Electrical Connection Lubricant - also known as Dielectric Tune-up Grease, a protective lubricant that prevents corrosion.
E_o	Power supply output voltage
i, I	Electrical current flow
I_o	Power supply output current
IEC	International Electrotechnical Commission
IEEE	Institute of Electrical and Electronics Engineers
I/O	Input/Output; The hardware and associated protocol that implement communication between information processing systems and/or devices. Inputs are the signals or data received by the system or device, and outputs are the signals or data sent from it.
IP	Internet Protocol; when used with "address", refers to a numerical Internet address
kG	kilogauss: a magnetic field unit of measurement
LED	Light-Emitting Diode; a semiconductor device that emits light when energized - used for visual status indication
LF or <LF>	Text Linefeed character
LHe	Liquid Helium
LN2	Liquid Nitrogen
Max	Maximum
Min	Minimum
MSDS	Material Safety Data Sheet - provides workers and emergency personnel with procedures for handling or working with a specific substance in a safe manner and includes information such as physical data, toxicity, health effects, first aid, reactivity, storage, disposal, protective equipment, and spill-handling procedures.
R_{lead}	Electrical circuit lead or wiring resistance
RS-232	RS-232 is a long-established standard and protocol for relatively low speed serial data communication between computers and related devices; originally established for teletypewriter communication.
SCPI	Standard Commands for Programmable Instruments
V	Volts
VA	Volt-amperes ($V \times I$); a unit of electrical reactive power
V_{lead}	Voltage ($I \times R$) developed across circuit lead or wiring resistance due to current flow
V_m	Magnet voltage
V_s	Power supply voltage

A

abbreviations and acronyms 93

abnormal operation 38

Alarm Status screen 29

alarms 29

- alarm status screen 29
- configure source and level 30
- indication and acknowledgement 31
- multiple 30
- muting 29, 31
- remote configuration 60
- states 31

AMI contacts 83

analog outputs

- connection 13
- connector 82
- current loop diagram 14
- current loop output 13
- recorder output 13
- remote configuration 64
- simultaneous usage 13
- source selection 35
- specifications 7

autofill

- AUTO-OFF 34
- AUTO-ON 34
- home screen indicator 16
- manual overrides 34
- M-CLOSED 34
- M-OPEN 34
- parameters 33
- remote configuration 62
- source selection 33
- start and stop levels 33
- states 34
- system setup 12
- timeout 33
- timeout reset 34

AUTO-OFF 34

AUTO-ON 34

B

battery backup 7

battery replacement 72

burnout protection 26

C

cables

- serial 52

calibration

- LHe active length 45

- liquid helium sensors 43
- remote LHe calibration 65
- touch screen 85

choosing units 35

command conventions 47

command summary 47–52

commands - see remote commands

communications

- commands - see remote commands
- error codes 66
- Ethernet 7
- Ethernet connector 7
- protocols 7
- remote, browser-based 2
- response format and termination 47
- serial 7
- serial connector 79

configuration

- alarm source and level 30
- alarm state 31
- alarms and relays 29
- analog outputs 35
- autofill 33
- capacitance xii
- custom xiii
- dirty sensor mode 28
- DST 42
- Ethernet 36
- helium xii
- instrument options xi
- level-based alarms 29
- line cord xii
- mounting xii
- muting alarms 31
- relay actuation 32
- relays source and level 31
- reset to defaults 39
- set date/time 41–42
- shutdown 39
- timezone 42

connections

- liquid helium sensors 11

connector

- analog outputs 82
- Ethernet 53, 80
- LHe sensor 81
- serial port 52

connectors

- RJ-45 - see Ethernet

contact AMI 83

cryogenic liquids viii

- first aid, first aid ix
- safety viii

current loop
 connection diagram 14
current loop specifications 7

D

dimensions 7
dirty sensor mode 28

display

- ALARM indicator/button 16
- autofill control 34
- AUTOFILL indicator/button 16
- BACK button 17
- field editing footer 16
- HELP button 16
- home screen 15
- home screen footer 16
- LHe level 24
- LHe sampling mode 27, 44
- MENU button 16, 17
- multiple alarm conditions 30
- numeric keypad 25
- remote configuration 58
- SAVE/CANCEL buttons 16
- size and resolution 2, 6
- units 35

DST support 42

E

environmental limits 8

equipment

- recommended for safety x

Ethernet

- configuration 36
- connector 7, 53
- connector pin out 80
- IP addressing 7
- ip parameters 53
- remote commands 55
- remote port 7180 54
- static ip parameters 36
- termination characters 53

Ethernet connector pin out 80

F

fill timeout 33
 disable 33

fill timeout for LN2
 reset 34

footer

- editing 16
- home screen 16

front panel layout 4
 power switch 15
fuse replacement 74

G

glossary of terms 93
 DHCP 93

H

home screen 15

- AUTOFILL button 34
- AUTO-OFF button 34
- AUTO-ON button 34
- footer 16
- M-CLOSED button 34
- menu button 17
- M-OPEN button 34
- show LHe level 24
- units selection 35

I

installation

- autofill system 12
- liquid helium sensors 10
- mounting 9
- power 13
- unpacking 9

installing

- earth ground xi

instrument

- battery replacement 72
- cleaning 71
- fuse replacement 74
- part number definition xii
- replacement parts 71
- troubleshooting 83

instrument description 1

ip address 53

L

level units 35

liquid helium sensor

- avoiding ice formation 10, 28
- burnout protection 26, 83
- calibration 43
- configuration 24
- connecting 11
- connector J1 wiring 81
- connector warning 11
- continuous mode time limit 25, 45
- lead wire sizing 11
- loss of sensor 83

- mounting 10
- sampling interval 44, 63
- sensor excitation values 27
- sensor name 46
- sensor preparation 11
- sensor sample interval 26
- setting active length 25, 43
- type 4.2K or 2K 1, 24, 43
- vacuum caution 10

liquid level system terminology 93

logging 86

- LHe level 63

loss of sensor

- LHe troubleshooting 83

M

M-CLOSED 34

measurement

- continuous mode time limit 25, 45
- LHe burnout protection 26
- remote query 62
- sampling interval for LHe 26
- sampling mode for LHe 27, 44

menu tree

- overview 18

menus 17

- descriptions 19
- overall structure 18

M-OPEN 34

mounting

- helium sensor 10
- removing bottom feet 9

muting alarms 31

N

network

- address 36
- configuration 36
- connector 53
- DHCP 53
- DHCP definition 93
- gateway 36
- hostname 36
- netmask 36
- port 7180 53
- static or DHCP mode 36

O

operation warnings x, 9, 11, 13, 72, 74

P

part number definition xii

power on/off 15

power requirements 7, 13

power supply

- operating parameters 6

R

rear panel layout 5

recorder output specifications 7

relays 29

- actuation states 32
- configure source and level 31
- ratings 7
- remote configuration 58
- type 7

remote commands

- alarms configuration 60
- analog output assignment 64
- display configuration 58
- error codes 66
- fill control 62
- legacy commands 67
- LHe calibration 65
- LHe sampling 63
- measurement 62
- relays configuration 58
- system related 55
- system-related commands 48
- units 65

remote control

- browser-based 2

remote error codes 66

remote interface reference

- command summary 50, 51, 52

remote interface reference - see remote commands

reset to factory defaults 39

RJ-45 connector - see Ethernet

RS-232 configuration

- parameters 52

RS-232 connector 79

RS-232 setup 37

S

safety

- cryogenic liquids viii
- equipment x

- legend x
 - liquid helium sensor connectors 11
 - warnings xi
 - sampling interval 26
 - screen descriptions
 - HISTORY 23
 - NETWORK 21
 - OUTPUTS
 - ALARMS 20
 - ANALOG OUTPUTS 21
 - AUTOFILL 21
 - RELAYS 20
 - SENSORS 19
 - CAL He 19
 - SENSOR NAME 19
 - SHUTDOWN 23
 - SYSTEM 22
 - RS232 SETUP 22
 - SET TIME 22
 - SYSTST 23
 - UPDATE 22
 - sensor
 - LHe active length 45
 - LHe calibration 43
 - LHe sensor name 46
 - LHe sensor test points 81
 - protected LHe troubleshooting 83
 - serial communications
 - baud rate 7, 38
 - commands 55
 - connector 7
 - connector pin out 79
 - factory reset 57
 - interactive communication 47
 - set date 57
 - set time 57
 - system reboot 57
 - terminators 53
 - serial port
 - baud rate 38
 - echo 38
 - function 38
 - setup 37
 - termination 38, 53
 - serial port connector/cables 52
 - service
 - AMI contact 83
 - service - see instrument
 - shutdown 39
 - specifications
 - analog output 7
 - display 6
 - environmental limits 8
 - level measurements 6
 - physical 7
 - power 7
 - relays 7
 - standards conformance 8
 - system
 - date and time 41–42
 - logging 86
 - remote configuration 55
 - shutdown 39
 - system features 1
 - system specifications 6
 - system test
 - analog outputs 86
 - hardware controls 85
 - screen pixels 84
 - touch screen cal 85
- T**
- test points for LHe sensor 81
 - TIMEZONE selection 42
 - troubleshooting common issues 83
 - system test screens 83
- U**
- units 35
 - remote configuration 65
 - user replaceable parts 71
- V**
- vacuum caution
 - helium sensor 10
- W**
- warnings
 - equipment xi
 - weight 7