



# **MODEL 1700 LIQUID LEVEL INSTRUMENT**

**(LIQUID NITROGEN VERSION)**

**INSTALLATION, OPERATION, AND  
MAINTENANCE INSTRUCTIONS**

***American Magnetics, Inc.***

---

P.O. Box 2509, 112 Flint Road, Oak Ridge, TN 37831-2509, Tel: 865-482-1056, Fax: 865-482-5472



<b>Foreword .....</b>	<b>vii</b>
<b>Purpose and Scope .....</b>	<b>vii</b>
<b>Contents of this Manual .....</b>	<b>vii</b>
<b>General Precautions .....</b>	<b>viii</b>
<b>Cryogen Safety Summary .....</b>	<b>x</b>
<b>Safety Legend .....</b>	<b>x</b>
<b>Equipment Warnings .....</b>	<b>xi</b>
<b>Other Manual Conventions .....</b>	<b>xi</b>
<b>Instrument Configuration .....</b>	<b>xii</b>
 <b>Introduction .....</b>	 <b>1</b>
<b>Model 1700 Instrument .....</b>	<b>1</b>
Cryogenic (Capacitance-Based) Liquid Level Sensors .....	1
Digitally-Controlled .....	3
System Flexibility .....	3
Display .....	3
Intuitive Human-Interface Design .....	3
Measurement Flexibility .....	4
Real Time Clock .....	4
Valve Control Output .....	4
Analog Outputs .....	4
Signal Relays .....	4
Connectivity .....	4
<b>Model 1700 Front Panel Layout .....</b>	<b>6</b>
<b>Model 1700 Rear Panel Layout .....</b>	<b>7</b>
<b>Model 1700 Specifications @ 25°C .....</b>	<b>8</b>
 <b>Installation .....</b>	 <b>11</b>
<b>Unpacking and Inspecting the Instrument .....</b>	<b>11</b>
<b>Mounting the Model 1700 Instrument .....</b>	<b>11</b>
<b>Capacitance-Based Liquid Level Sensor Installation .....</b>	<b>12</b>

Connecting the Capacitance Sensor . . . . .	13
<b>Configuring an Autofill System . . . . .</b>	<b>15</b>
LN2 Autofill System Description . . . . .	15
LN2 Autofill System Setup . . . . .	16
<b>Power Requirements. . . . .</b>	<b>18</b>
<b>Connecting the Analog Outputs . . . . .</b>	<b>18</b>
0-10 VDC Recorder Output. . . . .	18
4-20 mA Current Loop Output. . . . .	18
<b>Operation.....</b>	<b>21</b>
Energizing the Model 1700 Instrument . . . . .	21
<b>Screen Navigation. . . . .</b>	<b>21</b>
Home Screen . . . . .	21
Home Screen Footer. . . . .	22
Editing a Field . . . . .	22
<b>Navigating the Instrument Menus . . . . .</b>	<b>23</b>
Menu Structure . . . . .	24
Screen Descriptions . . . . .	25
<b>Capacitance (Liquid Nitrogen) Level. . . . .</b>	<b>30</b>
Configure the instrument to display nitrogen level . . . . .	30
<b>Alarms and Relays . . . . .</b>	<b>31</b>
Alarm Status Screen . . . . .	31
Level-Based Alarms . . . . .	31
Time-Based Fill Alarm. . . . .	32
Multiple alarms . . . . .	32
Configuring Alarm Setpoints . . . . .	33
Acknowledging an Alarm. . . . .	33
Muting an Alarm . . . . .	34
Configuring Relay Setpoints . . . . .	34
<b>Configure the Autofill Function . . . . .</b>	<b>36</b>
Setting the autofill parameters . . . . .	36
Enable the Autofill function . . . . .	37
Clearing the Autofill Timeout Alarm. . . . .	37
<b>Select the appropriate units on the display . . . . .</b>	<b>38</b>
<b>Damping Configuration . . . . .</b>	<b>38</b>

<b>Analog output signals</b> .....	<b>39</b>
Configuring the Analog Outputs .....	39
<b>Ethernet Connectivity</b> .....	<b>40</b>
IP Addressing Scheme .....	40
<b>Serial Connectivity</b> .....	<b>41</b>
Configuring the RS-232 Settings .....	41
<b>Abnormal Operation</b> .....	<b>42</b>
Capacitance Sensor Contamination .....	42
Resetting the Instrument to Factory Defaults .....	43
<b>Shutting the Instrument Down</b> .....	<b>44</b>
<b>Calibration</b> .....	<b>45</b>
<b>Setting the System Date and Time</b> .....	<b>45</b>
<b>Capacitance-based Level Calibration</b> .....	<b>47</b>
Understanding the Sensor Active Length .....	47
Relationship between Calibration and Sensor Length .....	48
Variations in the Dielectric with Changing Density .....	48
Capacitance-based Sensor Calibration Methods .....	49
Pre-Calibration Procedure .....	51
<b>Open Dewar Calibration</b> .....	<b>55</b>
<b>Closed Dewar Calibration</b> .....	<b>57</b>
Presetting the maximum and minimum calibration points .....	57
Completing the closed dewar calibration procedure .....	59
<b>Approximate Calibration</b> .....	<b>62</b>
<b>Remote Interface Reference</b> .....	<b>65</b>
<b>SCPI Command Summary</b> .....	<b>65</b>
<b>RS-232 Configuration</b> .....	<b>70</b>
Serial port connector and cabling .....	70
Command/return termination characters .....	70
<b>Ethernet Configuration</b> .....	<b>71</b>
Ethernet Connector .....	71
Termination Characters .....	71

Port Assignment .....	71
<b>Command Reference .....</b>	<b>72</b>
System Related Commands .....	72
Display Configuration Commands and Queries .....	75
Relay Configuration Commands and Queries .....	75
Alarm Configuration Commands and Queries .....	77
Measurement Commands and Queries .....	79
Fill Control and Queries .....	79
Assignment Commands and Queries .....	80
N2 Channel Calibration Commands and Queries .....	80
Remote Units Commands and Queries .....	81
<b>Error Codes .....</b>	<b>82</b>
<b>Legacy Command Reference .....</b>	<b>83</b>
Commands for Controlling the Units of Measurement .....	83
Commands for Configuring Permanent Memory .....	84
Querying the Configuration .....	85
Returning a Level Measurement .....	85
<b>Service and Repair .....</b>	<b>87</b>
<b>Cleaning .....</b>	<b>87</b>
<b>User Replaceable Parts .....</b>	<b>87</b>
<b>Battery Replacement. ....</b>	<b>88</b>
Low Battery Indication. ....	88
Tools Required .....	89
Procedure .....	89
<b>Fuse Replacement. ....</b>	<b>90</b>
Tools Required .....	90
Procedure .....	90
<b>Firmware Upgrade Via Ethernet. ....</b>	<b>91</b>
Upgrade via SCP .....	92
<b>Appendix .....</b>	<b>95</b>
<b>Serial (RS-232) Connector .....</b>	<b>95</b>
<b>Ethernet Connector .....</b>	<b>96</b>

.....97

**Aux I/O Connector .....98**

**Dielectric Constants for Cryogenic Liquids .....99**

**Troubleshooting.....100**

    Instrument Displays “LOSS OF SENSOR” Condition for LN2 Level Measurement ..... 100

    Instrument Displays “SENSOR SHORTED” Condition for LN2 Level Measurement ..... 100

    System Test Screens ..... 100

    System Logs..... 103

**Glossary .....109**

**Abbreviations and Acronyms .....109**

**Index .....113**





# FOREWORD

---

---

## PURPOSE AND SCOPE

---

This manual contains the operation and maintenance instructions for the American Magnetics, Inc. Model 1700 Liquid Level Control Instrument and outlines applications for various system configurations. Since it is not possible to cover all equipment combinations for all magnet systems, only the most common configurations are discussed. The user is encouraged to contact an authorized AMI Technical Support Representative for information regarding specific configurations not explicitly covered in this manual. This manual refers to the instrument as configured for liquid nitrogen operation.

---

## CONTENTS OF THIS MANUAL

---

**Introduction** describes the functions, specifications, and characteristics of the Model 1700 Instrument. It provides illustrations of the front and rear panel layouts as well as documenting the performance specifications. Additional information is provided in the form of system block diagrams.

**Installation** describes how the Model 1700 Instrument is unpacked and installed in conjunction with ancillary equipment in typical systems. Block-level diagrams document the interconnects for various system configurations.

**Operation** describes how the Model 1700 Instrument is used to monitor and automatically control liquid nitrogen levels.

**Calibration** describes the various calibration techniques for liquid nitrogen level sensors.

**Remote Interface Reference** documents all remote commands and queries available through the Model 1700 Instrument RS-232 and Ethernet interfaces. A quick-reference summary of commands is provided as well as a detailed description of each.

The **Appendix** and **Glossary** sections support the information in the sections listed above. See the Appendix section when referenced from other sections. See the Glossary for any words or acronyms presented in the above sections, requiring a more complete understanding.

---

## GENERAL PRECAUTIONS

---

### Cryogen Safety

The two most common cryogenic liquids used in superconducting magnet systems are nitrogen and helium. Both of these cryogens are extremely cold at atmospheric pressure ( $-321^{\circ}\text{F}$  and  $-452^{\circ}\text{F}$ , respectively). The following paragraphs outline safe handling precautions for these liquids.

Personnel handling cryogenic liquids should be thoroughly instructed and trained as to the nature of the liquids. Training is essential to minimize accidental spilling. Due to the low temperature of these materials, a cryogen spilled on many objects or surfaces may damage the surface or cause the object to shatter, often in an explosive manner.

Inert gases released into a confined or inadequately ventilated space can displace sufficient oxygen to make the local atmosphere incapable of sustaining life. Liquefied gases are potentially extreme suffocation hazards since a small amount of liquid will vaporize and yield a very large volume of oxygen-displacing gas. Always ensure the location where the cryogen is used is well ventilated. Breathing air with insufficient oxygen content may cause unconsciousness without warning. If a space is suspect, purge the space completely with air and test before entry. If this is not possible, wear a forced-air respirator and enter only with a co-worker standing by wearing a forced-air respirator.

Cryogenic liquids, due to their extremely low temperatures, will also burn the skin in a similar manner as would hot liquids. Never permit cryogenic liquids to come into contact with the skin or allow liquid nitrogen to soak clothing. Serious burns may result from careless handling. Never touch uninsulated pipes or vessels containing cryogenic liquids. Flesh will stick to extremely cold materials. Even nonmetallic materials are dangerous to touch at low temperatures. The vapors expelled during the venting process are sufficiently cold to burn flesh or freeze optic tissues. Insulated gloves should be used to prevent frost-bite when operating valves on cryogenic tanks. Be cautious with valves on cryogenic

systems; the temperature extremes they are typically subjected to cause seals to fail frequently.

In the event a person is burned by a cryogen or material cooled to cryogenic temperatures, the following first aid treatment should be given pending the arrival and treatment of a physician or other medical care worker:

1. If any cryogenic liquid contacts the skin or eyes, immediately flush the affected area gently with tepid water (102°F – 105°F, 38.9°C – 40.5°C) and then apply cold compresses.
2. Do not apply heat. Loosen any clothing that may restrict circulation. Apply a sterile protective dressing to the affected area.
3. If the skin is blistered or there is any chance that the eyes have been affected, get the patient immediately to a physician for treatment.

Containers of cryogenic liquids are self pressurizing (as the liquid boils off, vapor pressure increases). Hoses or lines used to transfer these liquids should never be sealed at both ends (i.e. by closing valves at both ends).

When pouring cryogenic liquids from one container to another, the receiving container should be cooled gradually to prevent damage by thermal shock. The liquid should be poured slowly to avoid spattering due to rapid boil off. The receiving vessel should be vented during the transfer.

Introduction of a substance at or near room temperature into a cryogenic liquid should be done with great caution. There may be a violent gas boil-off and a considerable amount of splashing as a result of this rapid boiling. There is also a chance that the material may crack or catastrophically fail due to forces caused by large differences in thermal contraction of different regions of the material. Personnel engaged in this type of activity should be instructed concerning this hazard and should always wear a full face shield and protective clothing. If severe spraying or splashing could occur, safety glasses or chemical goggles along with body length protective aprons will provide additional protection.

The properties of many materials at extremely low temperatures may be quite different from the properties that these same materials exhibit at room temperatures. Exercise extreme care when handling materials cooled to cryogenic temperatures until the properties of these materials under these conditions are known.

Metals to be used for use in cryogenic equipment application must possess sufficient physical properties at these low temperatures. Since ordinary carbon steels, and to somewhat a lesser extent, alloy steels, lose much of their ductility at low temperatures, they are considered unsatisfactory and sometimes unsafe for these applications. The

austenitic Ni-Cr alloys exhibit good ductility at these low temperatures and the most widely used is 18-8 stainless steel. Copper, Monel<sup>®</sup>, brass and aluminum are also considered satisfactory materials for cryogenic service.

---

## CRYOGEN SAFETY SUMMARY

---

Cryogenic systems are complex systems with the potential to seriously injure personnel or equipment if not operated according to procedures. The use of safety mechanisms (pressure relief valves, rupture disks, etc.) in cryogenic systems is usually necessary.

### Recommended Safety Equipment

The use of proper safety equipment is necessary. Such equipment may include, but not limited to, the following items:

- First Aid kit
- Fire extinguisher rated for class C fires
- Cryogenic gloves
- Face shield
- Signs to indicate that there are potentially dangerous cryogenics in use in the area.

---

## SAFETY LEGEND

---



Instruction manual symbol: the product is marked with this symbol when it is necessary to refer to the instruction manual in order to protect against damage to the product or personal injury.



Hazardous voltage symbol.



Alternating Current (Refer to IEC 417, No. 5032).



Off (Supply) (Refer to IEC 417, No. 5008).



On (Supply) (Refer to IEC 417, No. 5007).

### **WARNING**

The Warning sign denotes a hazard. It calls attention to a procedure or practice, which if not correctly adhered to, could result in personal injury. Do not proceed beyond a Warning sign until the indicated conditions are fully understood and met.

**CAUTION**

The Caution sign denotes a hazard. It calls attention to an operating procedure or practice, which if not adhered to, could cause damage or destruction of a part or all of the product. Do not proceed beyond a Caution sign until the indicated conditions are fully understood and met.

---

## EQUIPMENT WARNINGS

---

Before energizing the equipment, the earth ground of the power receptacle must be verified to be at earth potential and able to carry the rated current of the power circuit. Using extension cords should be avoided. However, if an extension cord must be used, insure the ground conductor is intact and the cord is capable of carrying the rated current without excessive voltage drop.

In the event that the ground path becomes less than sufficient to carry the rated current of the power circuit, the equipment should be disconnected from power, labeled as unsafe, and removed from place of operation.

Do not operate this equipment in the presence of flammable gases. Doing so could result in a life-threatening explosion.

Do not modify this equipment in any way. If component replacement is required, return the equipment to AMI facilities as described in the troubleshooting section of this manual.

If used in a manner not specified in this manual, the protection provided by the design, manufacture and documentation of the Model 1700 Instrument may be impaired.

---

## OTHER MANUAL CONVENTIONS

---

This manual refers to measuring liquid nitrogen (LN<sub>2</sub>) when referring to capacitance-based level measurement since nitrogen is by far the most common cryogenic liquid measured by a capacitance-based sensor. The capacitance sensor technology can be used to measure most any cryogenic liquid, with the exception of liquid helium (LHe) which requires a superconductor-based sensor.

---

## INSTRUMENT CONFIGURATION

---

The Model 1700 Instrument is configured at time of purchase as a capacitance-based (typically liquid nitrogen) level instrument/controller.

Every configuration may be further customized by the following options:

- Table top, single rack mounting, dual rack mounting.
- Line cord: North American, European Schuko, Australia/NZ, China, UK, pig-tailed ends.

The instrument part number, shown on the serialization label located on the underside of the instrument, as well as in a field in the instrument Menu, identifies the configuration according to the following key:

**1700-A-B-C-D-E** where

**A** indicates the mounting method:

Tbl = tabletop

SR19L = single rack mounted, 19" wide rack standard, instrument on left side

SR19R = single rack mounted, 19" wide rack standard, instrument on right side

SR10L = single rack mounted, 10" wide rack standard, instrument on left side

SR10R = single rack mounted, 10" wide rack standard, instrument on right side

DRL = dual rack mount, 19" wide rack standard, instrument on left side

DRR = dual rack mount, 19" wide rack standard, instrument on right side

**B** indicates the line cord shipped with the instrument:

N = North American

E = European, Schuko

A = Australian/New Zealand

C = Chinese

**U** = United Kingdom

**P** = India/pigtailed

**C** indicates the capacitance-based level configuration:

**CAP** = capacitance-based (typically nitrogen)

**N** = not configured

The sensor active length and length units are appended to the configuration code, ie ...CAP-10.4CM-...

**D** indicates the helium level configuration:

**N** = not configured

**E** is used to denote any instrument customization:

**S** = standard (no customization)

**C** = instrument modified.





---

# INTRODUCTION

---

---

## MODEL 1700 INSTRUMENT

---

The AMI Model 1700 Liquid Level Instrument is a sophisticated measurement and control instrument which provides monitoring capacitance-based level sensors as inputs and provides for automatic level control based on user set parameters.

At time of purchase, the Model 1700 will be configured as:

- A capacitance-based (typically liquid nitrogen)<sup>1</sup> level instrument/controller.

---


### CRYOGENIC (CAPACITANCE-BASED) LIQUID LEVEL SENSORS

The Model 1700 Instrument will be used with a capacitance-based sensor for all cryogenic liquids except for helium. Specially insulated capacitance sensors are required for liquids with electrical conductivity, e.g. water.


The Model 1700 Instrument must use an oscillator device to measure cryogenic liquid levels. This instrument has an internal oscillator if the distance between the sensor and the instrument is less than or equal to 15 feet (457 cm) so no external oscillator is necessary. If the distance between the sensor and the instrument is greater than 15 feet, an external oscillator/transmitter unit is necessary. The instrument will automatically configure itself at boot time to use the internal oscillator unless the external oscillator/transmitter is connected to the BNC connector on the rear panel.

---

1. The instrument can be used to measure/control any cryogenic liquid but the most common is nitrogen.



**Model 1700 Instrument using internal oscillator/transmitter**



**Model 1700 Instrument using external oscillator/transmitter**

## DIGITALLY- CONTROLLED

The Model 1700 contains a microcomputer which controls analog data conversion, display/keypad functions, communications I/O, dry contact closures, generation of analog output signals and relay control of a mains power outlet for solenoid valve autofill applications.

## SYSTEM FLEXIBILITY

The Model 1700 instrument incorporates data converters to translate signals between the analog and digital domains. Precision instrumentation techniques and potentiometer-free designs are employed throughout the Model 1700 Instrument to ensure long term stability and accurate signal translation for a wide range of conditions.

## DISPLAY

The Model 1700 Instrument has a 4.3" diagonal measure TFT (Thin Film Transistor) color liquid crystal display of 480 x 272 pixels. The display has a 4-wire resistive touch overlay for easy operator input.

## INTUITIVE HUMAN- INTERFACE DESIGN

The Model 1700 instrument is designed to simplify the touch-screen based user interface. All functions were analyzed and subsequently programmed so that the most commonly used functions are addressed

with the least number of keystrokes. The menus are presented in a logical fashion so that the operation of the Model 1700 is intuitive to the user. Context-sensitive Help screens are also provided.

---

## MEASUREMENT FLEXIBILITY

Set points can be assigned to control two dry-contact relay outputs. The make or break function of each relay can be set independently. A solid-state relay allows mains power to be switched according to set points to operate a solenoid-operated valve or other load.

---

## REAL TIME CLOCK

The Model 1700 Instrument incorporates a real time clock, which can be manually set or synchronized via NTP, with support for worldwide timezones including automatic DST adjustment.

---

## VALVE CONTROL OUTPUT

The Model 1700 Instrument has a switched 2 ampere at mains voltage output for energizing a solenoid operated flow valve or other loads. This output is controlled by a zero-crossing solid state relay.

The fill start is triggered by a level reading. The fill is stopped by reaching a user-set level, or exceeding a maximum fill time.

---

## ANALOG OUTPUTS

The Model 1700 Instrument has two analog outputs, a 0-10 V<sub>DC</sub> voltage output and a 4-20 mA<sub>DC</sub> current loop output. The 4-20 mA<sub>DC</sub> loop output has 1500 V<sub>PK</sub> circuit isolation. The outputs can be used simultaneously.

---

## SIGNAL RELAYS


The Model 1700 Instrument has two signal relays that change state based on a measured input. The set points of these relays are user-selectable as is their function, alarm on level above or below set point.

---

## CONNECTIVITY

The Model 1700 Instrument has a 10Base-T Ethernet connection as well as an RS-232 port for connecting to other equipment. The instrument communicates with a SCPI-based command set. The command set is 100% backward compatible with the AMI Model 135, 136, 185, and 186 instruments when configured as a single channel instrument. When configured as a dual channel instrument, additional commands are incorporated into the command set.

The Model 1700 Instrument allows for remote operation with an external browser via TCP/IP connection. All commands that are available by touching the local screen are available via the web browser



Model 1700 Nitrogen Instrument Via Web Browser

**NOTE** Using a web browser to connect to the instrument allows different browser sessions to display different information. For instance, one browser window can display helium level only and a second browser window can display nitrogen level only. The instrument can be configured independently as well since it's display is also a browser. Refer to "Configure the instrument to display nitrogen level" on page 30.

---


## MODEL 1700 FRONT PANEL LAYOUT

---



Model 1700 Front View; Dual Channel Instrument Shown

MODEL 1700 REAR PANEL LAYOUT



Model 1700 Rear Panel Description

LN2 VERSION

1	Computer Network Connector
2	Aux I/O Connector
3	RS-232 Serial Connector
4	Switched Valve Outlet Socket (IEC 60320 C13)
5	Capacitive Sensor Input Connector
6	Mains Power Entry Connector (IEC 60320 C14)
7	LHe Level Sensor Connector

## MODEL 1700 SPECIFICATIONS @ 25°C

### System Architecture

Display:	4.3" 24-bit color TFT display, 480x272 pixel with resistive touch screen
Sensor Types:	Capacitance-based liquid level
Maximum Length Readout:	Capacitance-based liquid level up to 999 in
System Operating Firmware Storage:	microSD card
System Clock:	Real time clock with timezone support, automatic DST adjustment, and NTP synchronization
Display Measurement Units:	Liquid level in cm, in or percent

### Level Measurement

Resolution:	0.1%, 0.1 cm, 0.1 in
Accuracy:	±0.5% of active sensor length
Linearity:	±0.1% or 1 mm (whichever is greater)
Capacitance Sensor Excitation Voltage:	5 V <sub>DC</sub>
Capacitance Transmitter Measurement Resolution:	0.7 pF
Extension Cable Limits with External Oscillator:	6-500 ft

### Operating Parameters

Alarm Set Points:	0% to 100%, adjustable; Alarm condition settable to above or below set point
Controller Output:	Line voltage @ 2 A <sub>AC</sub> (maximum)
Sample and Hold Period:	1 second to 86,400 seconds (24 hrs)
Audible alarm:	3500 ± 500 Hz, 73 to 86 dB(A)

### Analog Outputs

Output Types:	0-10 V <sub>DC</sub> and simultaneous 4 - 20 mA <sub>DC</sub>
4-20 mA Current Loop Power Supply Voltage:	12-32 V <sub>DC</sub>
0-10 V <sub>DC</sub> Recorder Output Load:	50k ohms or greater
0-10 V <sub>DC</sub> Recorder Output 0% Max Offset:	+20 mV
0-10 V <sub>DC</sub> Recorder Output 100% Max Error:	± 80 mV
Voltage and Current Output Converter Resolution:	12 bits
Integral Non-linearity:	±1LSB
Differential Non-linearity <sup>a</sup> :	±1LSB



**Relays**

№1 and/or №2 (W171DIP-7, or equivalent):	Contact Form: 1 Form A (SPST-NO) Maximum Switched Current: 3 A Switching Voltage: 60 V <sub>AC</sub> / 100 V <sub>DC</sub>
Level Control (Solid State):	Rated Load Voltage: 2 A at 100 to 240 V <sub>AC</sub> Load Voltage Range: 75 to 264 V <sub>AC</sub> (Line Dependent) Load Current: 0.1 to 2 A Surge Current: 30 A (60 Hz, 1 cycle)

**Auto Fill**

Controller Output Socket:	IEC 60320-13 socket on rear panel
Controller Output Power:	2 A <sub>AC</sub> at line voltage
Autofill Start/Stop Triggering:	Level-based; 0% to 100%, adjustable; Assignable to either nitrogen or helium channel
Fill Timeout Period:	1 minute to 99 hours, 59 minutes
Fill Error Alarm:	Fill time out

**Communication Protocol**

Host Computer Network Protocol:	Ethernet 10Base-T TCP/IP and RS-232
IP Addressing:	DHCP or static, IPv4
Network Connectivity and Traffic Indication:	Link and Activity LEDs on instrument rear panel
RS-232 Connector Specifications:	9-pin D-sub female connector to connect standard DTE 9-pin D-sub male connector using a standard straight cable
Communication Command Set:	SCPI-based. 100% backward compatible with the Model 184, 185, 186 instruments.

**Power Requirements**

Primary:	100-240 ±10% V <sub>AC</sub> , 50-60 Hz, 2.2 A maximum (200 VA plus sum of controller output)
Backup Battery for RTC:	CR2032

**Physical**

Dimensions <sup>b</sup> :	Table top configuration: 3.8" H x 8.4" W x 11.4" D [97 mm H x 213 mm W x 290 mm D] Single rack mount configuration: 3.5" H x 19.0" W x 11.4" D [89 mm H x 483 mm W x 290 mm D]
Weight:	Table-top configuration: 3.3 lbm [1.5 kG]; Single rack-mount configuration: 4.0 lbm [1.8 kG]

**Environmental Limits**

Ambient Temperature:	Operating: 0°C to 40°C [32°F to 104°F]; Non-operating: -20 °C to 60 °C [-4 °F to 140 °F]
Relative Humidity:	0 to 95%; non-condensing
Maximum Instrument Background Field:	Gauss (TBD)

**Standards**

Test Standards	<p>Testing of Equipment for Measurement, Control, and Laboratory Use (IEC 61326-1:2012, EN 61326-1)</p> <p>Electrostatic Discharge (ESD) (EN 61000-4-2)</p> <p>Radiated Immunity (EN 61000-4-3)</p> <p>Fast Transient Burst (EN 61000-4-4)</p> <p>Surges (EN 61000-4-5)</p> <p>Conducted Immunity (EN 61000-4-6)</p> <p>Power Frequency Magnetic Field (EN 61000-4-8)</p> <p>Voltage Dips and Interrupts (EN 61000-4-11)</p> <p>Harmonics (EN 61000-3-2)</p> <p>Flicker (EN 6100-3-3)</p> <p>Conducted Emissions (EN 55011/IEC/CISPR 11)</p> <p>Radiated Emissions (EN 55011/IEC/CISPR 11)</p> <p>Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use (IEC 61010-1)</p>
----------------	---

- a. Guaranteed monotonic over operating temperature range
- b. H = height; W = width; D = depth

---

# INSTALLATION

---

**WARNING**

Equipment warnings apply to all system installation configurations. Refer to “Equipment Warnings” on page xi, in the Foreword to be familiar with the safety requirements for a system installation.

---

## UNPACKING AND INSPECTING THE INSTRUMENT

---

Carefully remove the equipment, interconnecting cabling, and documentation from the shipping carton, and remove all packaging material.

**NOTE**

If there is any shipping damage, save all packing material and contact the shipping company representative to file a damage claim. Do not return to AMI unless prior authorization has been received.

---

## MOUNTING THE MODEL 1700 INSTRUMENT

---

If the Model 1700 Instrument is to be used as a table top model, place the equipment on a flat, secure surface.

If the Model 1700 Instrument is to be rack mounted, install it in a 19" wide instrument rack using the mounting hardware supplied by the rack cabinet manufacturer. The feet on the bottom of the instrument may be removed to facilitate rack mounting by using a prying device to remove the glossy black cap from the center of the foot and then prying the rubber foot itself from the chassis. Secure the front panel to the rack rail in each of the four corners.

---

## CAPACITANCE-BASED LIQUID LEVEL SENSOR INSTALLATION

---

Refer to the installation instructions provided with the level sensor(s). The following steps are general installation notes and should be used to supplement the installation instructions provided with the sensor.

Exercise care when installing the capacitance sensor since dents, crimps, bends or other physical distortions in the cylindrical capacitor will change electrical characteristics, possibly causing calibration errors and/or disruption of proper instrument operation. Before installing the sensor, review “Calibration” on page 45 to determine what, if any, calibration procedures may be necessary prior to operation.

**NOTE** The coaxial interconnecting cables and the transmitter should be mounted in such a manner as to avoid large temperature changes such as those encountered in the path of dewar vents.

1. Carefully remove the sensor from the shipping container and remove all packaging material.

**NOTE** If there is any shipping damage, save all packing material and contact the shipping representative to file a damage claim. Do not return the instrument to AMI unless prior authorization has been received.

2. Install the sensor in the vessel using the specified fitting of the sensor.

**CAUTION** Ensure the sensor is mounted with the top vent hole located inside of the cryostat.

**CAUTION** Avoid installing in a location where icing may occur. Ice formations or moisture buildup on the BNC connector may cause the sensor to short out indicating a higher liquid level than actually exists.

**CAUTION** Moisture or contaminants in any of the BNC coaxial connectors can short out the sensor and cause measurement errors. A pack of non-conductive electrical connection lubricant (ECL or “Dielectric Tune-up Grease”) has been included with the liquid level sensor packaging to reduce the possibility of this occurring. If desired, apply a small amount of ECL to any of the BNC connectors that may be exposed to moisture. Mate the ECL-coated connectors

then remove any excess ECL from the outside of the connector. Added protection can be achieved by covering the ECL-coated connections with a short section of heat-shrink tubing.

MSDS sheets for the ECL are available upon request.

**CAUTION**

Exercise care when installing the sensor since dents, crimps, bends or other physical distortions in the thin wall capacitor will change electrical characteristics possibly causing calibration errors and/or disruption of proper instrument operation.

---

**CONNECTING THE  
CAPACITANCE SENSOR**

The capacitance sensor may be connected directly to the BNC connector on the instrument rear panel if the length of the coaxial cable is 15 feet or less. If the sensor is greater than 15 feet from the instrument, an external oscillator/transmitter unit must be used. Refer to figures "Model 1700 Instrument using internal oscillator/transmitter" on page 2 and "Model 1700 Instrument using external oscillator/transmitter" on page 3 as appropriate.

If the transmitter is connected to the sensor with a length of coaxial cable, the a 15 ft standard length cable, with part number of EH2362, is available from AMI. Contact the factory for details. Speak to an AMI Sales Engineer before using cables longer than 15 feet.

**NOTE**

In order to maintain system performance and accuracy, the coaxial cable used to connect the capacitance level sensor to the instrument or oscillator/transmitter must be Trompeter TCC-75-2 or equivalent and should not be longer than 15 feet [4.57 m]. If a longer length section of coaxial cable is necessary, please discuss with an Authorized AMI Technical Representative.

**NOTE**

If an oscillator/transmitter is used, the length of coaxial cable between the oscillator/transmitter and the instrument may be up to 500 feet in length.

**CAUTION**

Moisture or contaminants in any of the BNC coaxial connectors can short out the sensor and cause an erroneous readings or transmitter failure. A pack of non-conductive electrical connection lubricant (ECL or "Dielectric Tune-up Grease") has been included with the liquid level sensor packaging to reduce the possibility of this occurring.

To connect the coaxial cable to the BNC connector on the capacitance level sensor:

1. Apply a small amount of ECL to any of the BNC connectors that may be exposed to moisture.
2. Mate the ECL-coated connectors then remove any excess ECL from the outside of the connector.
3. Cover the ECL-coated connections with a short section of heat-shrink tubing, also included, for added moisture protection.

**NOTE** MSDS sheets for the ECL are available upon request.


## CONFIGURING AN AUTOFILL SYSTEM

Autofill systems are useful in applications where liquid level in a dewar must be maintained without operator intervention.

To create an autofill system, a cryogenic liquid source must be attached to the target or level controlled dewar via a solenoid-controlled valve. The valve should be controlled by the Model 1700 Instrument.

### LN2 AUTOFILL SYSTEM DESCRIPTION

For LN2 autofill, the system consists of a Model 1700 Instrument with a liquid level sensor, and a solenoid-operated flow valve. The instrument makes continuous level measurements and based on level, energizes the valve to begin liquid transfer. The transfer is stopped when the measured level reaches a user-determined point.



Typical LN2 Autofill Setup

## Standard Autofill Setup Description

Number	Item
<b>1</b>	Model 1700 Liquid Level Instrument (Level Controller)
<b>2</b>	IEC60320 C13 socket labeled <b>VALVE CONTROL</b>
<b>3</b>	BNC connector labeled <b>O/T</b>
<b>4</b>	Solenoid-operated flow control valve line cord with IEC60320 C14 plug
<b>5</b>	Instrument IEC60329 C14 Power cord inlet
<b>6</b>	Coaxial cable connecting the oscillator / transmitter and the instrument
<b>7</b>	<b>Oscillator / Transmitter</b> (optional; refer to the figure on page 13)
<b>8</b>	Level Sensor for level controlled or target dewar
<b>9</b>	Coaxial cable connecting the Oscillator / Transmitter and the liquid level sensor (optional; refer to the figure on page 13)
<b>10</b>	Fill port on target dewar
<b>11</b>	Transfer line attached to the target dewar and the fill solenoid valve on the source dewar
<b>12</b>	Solenoid-operated fill valve
<b>13</b>	Supply dewar relief valve
<b>14</b>	Supply dewar

## LN2 AUTOFILL SYSTEM SETUP

**CAUTION**

A relief valve must be used in autofill systems to ensure no cryogenic liquid can be trapped in a transfer line volume where expansion can cause damaging pressure. This can occur if the solenoid operated fill valve and the supply dewar isolation valve are closed, trapping a cryogenic liquid in a confined volume. All AMI LN2 transfer line systems include a relief valve to preclude this sort of event.

**CAUTION**

The switched valve outlet socket switches the line voltage to which the instrument mains power is connected. Ensure any connected valve is rated to support the local VAC line voltage.


1. Mount the level sensor (**8**) in the target dewar.
2. Connect the transfer line (**11**) and fill solenoid valve (**12**) or supply manifold to the source dewar.
3. Connect the other end of the transfer line to the fill port (**10**) on the valve/ manifold of the target dewar.



4. Connect the capacitance sensor to the instrument.
  - a. For distances of 15 feet and less, connect the coaxial cable (6) between the BNC connector on the liquid level sensor and the BNC connector on the back of the instrument labeled **O/T (3)**.
  - b. For distances greater than 15 feet, connect the coaxial cable (9) between the BNC connector on the liquid level sensor and the BNC connector on the oscillator / transmitter (7). Use a second length of coaxial cable (6) to connect between the oscillator / transmitter (7) and the BNC connector on the back of the instrument labeled **O/T (3)**.
5. Connect the solenoid valve (12) to the IEC60320 C13 valve socket<sup>1</sup> on the Model 1700 Instrument rear panel labeled **VALVE CONTROL (2)**.

Connect the IEC320 C14 inlet connector and instrument power cord (5) to an appropriate power source (see below).

**NOTE** Should the nitrogen level sensor become disconnected from the instrument or an interconnecting cable fails, a **LOSS OF SENSOR** message will be displayed (as shown at right) and the autofill valve will be shut (**M-CLOSED**). When the sensor connection has been restored, the instrument will display the level but the autofill state will have to be manually changed back to **AUTOFILL**.



1. The valve socket must be IEC60320 C13 type.

---

## POWER REQUIREMENTS

---

**WARNING**

The Model 1700 Instrument operates on 50-60 Hz power and may be powered from 100-240 VAC. Insure that the input ground terminal is connected securely to an external earth ground.

Ensure the detachable mains supply cord is of suitable rating, i.e. 10 A (min) at 125 VAC for North America.

Ensure the power switch is in the OFF (●) position. Plug the Model 1700 Instrument line cord into the power entry module on the instrument rear panel and into the appropriate power receptacle.

---

## CONNECTING THE ANALOG OUTPUTS

---

The Model 1700 instrument provides 0-10 V<sub>DC</sub> and 4-20 mA<sub>DC</sub> outputs as standard. These outputs may be optionally connected to external equipment per the instructions below.

**NOTE**

Both analog outputs can be connected to the same or different level measurement source as desired, and both will operate simultaneously.

---

### 0-10 V<sub>DC</sub> RECORDER OUTPUT

The Model 1700 instrument provides a 0-10 V<sub>DC</sub> output on the AUX I/O connector (refer to “Aux I/O Pin Definitions” on page 98) corresponding to 0-100% of liquid level. This output can be assigned to either level measurement channel (refer to “Analog output signals” on page 39). Connect the voltage output on the AUX I/O connector to a suitable receiver being careful not to exceed the impedance restrictions (refer to “Model 1700 Specifications @ 25°C” on page 8).

---


### 4-20 mA CURRENT LOOP OUTPUT

The Model 1700 instrument provides a 4-20 mA<sub>DC</sub> output on the AUX I/O connector (refer to “Aux I/O Pin Definitions” on page 98) corresponding to 0-100% of liquid level. This output can be assigned to either level measurement channel (refer to “Analog output signals” on page 39).

The instrument requires an external power supply for current loop operation. Connect a suitable power supply as shown below.

**CAUTION**

It is extremely important to observe all polarities and to not exceed +32VDC for the loop power supply in order to prevent damage to the 4-20mA driver circuit.





# OPERATION

This section describes the operation of the Model 1700 Instrument.

## ENERGIZING THE MODEL 1700 INSTRUMENT

1. Turn the power switch on the front panel of the instrument to the On (I) position. The display will briefly show a boot image and then display the home screen showing level(s).
2. The boot process takes approximately 30 seconds. Boot time can be longer (approximately 2 minutes) if the instrument has been configured for a network connection and then is booted without the network present.
3. When the boot process is complete, the instrument will display the home (level) screen.

## SCREEN NAVIGATION

### HOME SCREEN

**NOTE** If the instrument was purchased with a level sensor, the instrument will be pre-configured and calibrated at the factory.

The display will look similar to the image shown at right:

If the instrument requires calibration<sup>1</sup>, refer to the following chapter to calibrate the instrument with an AMI level sensor.



LN<sub>2</sub> Home Screen

## HOME SCREEN FOOTER

Every screen has a footer. The level home screen displays both status information and buttons in the footer to navigate to other screens.

The home screen footer appears as shown at right:



Home Screen Footer

### Model 1700 Instrument Home Screen Footer

No.	Name	Function
①	MENU	Takes the user to the main menu screen.
②	ALARM	When illuminated, indicates an active alarm condition. Touch to enter the ALARM STATUS screen.
③	AUTOFILL	When illuminated, indicates the AUTOFILL function is enabled. Touch to enter the AUTOFILL configuration screen.
④	HELP	Takes the user to a Help screen where available.

## EDITING A FIELD

Once a field or multi-choice button on a screen has been edited, the footer changes to appear as follows:



Field Editing Footer

### Model 1700 Instrument Footer during editing a field

No.	Name	Function
①	SAVE	Saves the entries made on the screen. The footer then reverts to the default footer as shown in the previous section.
②	CANCEL	Cancels a change while <i>not saving entries</i> . The footer then reverts to the default footer as shown in the previous section.

The functions of the **ALARM**, **AUTOFILL**, and **Help** footer items are as described in the previous section.

1. If the instrument was purchased with level sensor(s), the instrument will be shipped set up and calibrated.

---

## NAVIGATING THE INSTRUMENT MENUS


---

The menu system is invoked by pressing the **MENU** button in the lower left corner of the instrument level display screen.



MENU button on HOME screen

When invoked, the **MENU** screen will be displayed:



MENU Selection Screen

Pressing **HOME** in the lower left corner of the screen will display the level display screen.



HOME button on MENU screen

Once a **MENU** has been chosen, a back button will be presented to the right **HOME**. This button will return the user to the previously displayed screen.



BACK button

```

graph TD
    MENU[MENU] --> SENSORS[SENSORS]
    MENU --> OUTPUTS[OUTPUTS]
    MENU --> NETWORK[NETWORK]
    MENU --> SYSTEM[SYSTEM]
    MENU --> HISTORY[HISTORY PLOT]
    MENU --> SHUTDOWN[SHUTDOWN]

    SENSORS --> S1[NITROGEN OSCILLATOR USED:  
<INTERNAL | EXTERNAL>]
    SENSORS --> S2[NITROGEN LEVEL ON HOME SCREEN?  
(YES | NO)]
    SENSORS --> S3[NITROGEN LEVEL DAMPING]
    S3 --> S3_1[Undamped O/T Period: <period> μs]
    S3 --> S3_2[Undamped LN2 Level: <level> %]
    S3 --> S3_3[Damping Time Period: <period> seconds]
    S3 --> S3_4[Damped O/T Period: <period> μs]
    S3 --> S3_5[Damping LN2 Level: <level> %]
    SENSORS --> S4[CAL N2]
    S4 --> S4_1[OSCILLATOR PERIOD:  
<value> μs]
    S4 --> S4_2[SENSOR ACTIVE LENGTH  
<min.n> cm]
    S4 --> S4_3[PERFORM MAX CAL]
    S4 --> S4_4[PERFORM MIN CAL]
    S4 --> S4_5[LOSS OF SENSOR CAL]
    S4 --> S4_6[APPROX. (MAX) CAL MULTIPLIER  
<value>]
    S4 --> S4_7[APPROX. (MAX) CAL MULTIPLIER  
(NOTIFY | NO NOTIFY)]
    SENSORS --> S5[SENSOR NAME]
    S5 --> S5_1[NITROGEN SENSOR NAME:  
<value>]

    OUTPUTS --> O1[ALARMS]
    O1 --> O1_1[ALARM No.1 SOURCE: (DISABLED | NITROGEN)]
    O1_1 --> O1_1_1[ALARM No.1 WHEN LEVEL (≥ | ≤) <level>%]
    O1 --> O1_2[ALARM No.2 SOURCE: (DISABLED | NITROGEN)]
    O1_2 --> O1_2_1[ALARM No.2 WHEN LEVEL (≥ | ≤) <level>%]
    OUTPUTS --> O2[RELAYS]
    O2 --> O2_1[RELAY No.1 SOURCE: (DISABLED | NITROGEN)]
    O2_1 --> O2_1_1[RELAY No.1 CLOSSES WHEN LEVEL (≥ | ≤) <level>%]
    O2 --> O2_2[RELAY No.2 SOURCE: (DISABLED | NITROGEN)]
    O2_2 --> O2_2_1[RELAY No.2 CLOSSES WHEN LEVEL (≥ | ≤) <level>%]
    OUTPUTS --> O3[AUTOFILL]
    O3 --> O3_1[AutoFill CONTROL: (DISABLED | NITROGEN)]
    O3_1 --> O3_1_1[START: <level> %]
    O3_1 --> O3_1_2[STOP: <level> %]
    O3_1 --> O3_1_3[FILL TIMEOUT: <m> MINUTES]
    OUTPUTS --> O4[ANALOG OUTPUTS]
    O4 --> O4_1[OUTPUT 0-10 SOURCE: (DISABLED | NITROGEN)]
    O4 --> O4_2[OUTPUT 4-20 mA SOURCE: (DISABLED | NITROGEN)]

    NETWORK --> N1[ADDRESSING SCHEME (STATIC | DYNAMIC | DISABLED)]
    N1 --> N1_1[ADDRESS <value>]
    N1 --> N1_2[NETMASK <value>]
    N1 --> N1_3[GATEWAY <value>]
    N1 --> N1_4[HOSTNAME <value>]
    N1 --> N1_5[MAC ADDRESS <value>]

    SYSTEM --> S1[Instrument Serial Number <value>]
    S1 --> S1_1[HW: <value>]
    S1 --> S1_2[FPGA: <value>]
    S1 --> S1_3[FW: <value>]
    S1 --> S1_4[Date of Manufacture: <value>]
    S1 --> S1_5[SET TIME]
    S1_5 --> S1_5_1[CURRENT SYSTEM DATE AND TIME: <value>]
    S1_5 --> S1_5_2[SYSTEM TIMEZONE: <value>]
    S1_5 --> S1_5_3[BATTERY GOOD <NOTIFY | NO NOTIFY>]
    S1_5 --> S1_5_4[SET MANUALLY]
    S1_5_4 --> S1_5_4_1[TIMEZONE: <value>]
    S1_5_4 --> S1_5_4_2[YEAR: <value>]
    S1_5_4 --> S1_5_4_3[MONTH: <value>]
    S1_5_4 --> S1_5_4_4[DAY: <value>]
    S1_5_4 --> S1_5_4_5[HOURL: <value>]
    S1_5_4 --> S1_5_4_6[MIN: <value>]
    S1_5_4 --> S1_5_4_7[SET FROM NTP SERVER]
    S1 --> S1_6[UPDATE]
    S1_6 --> S1_6_1[<Choose Firmware Version>]
    S1 --> S1_7[RS-232 SETUP]
    S1_7 --> S1_7_1[RS-232 PORT CONFIGURATION (LINUX LOGIN | SCPI)]
    S1_7 --> S1_7_2[ECHO (ON | OFF)]
    S1_7 --> S1_7_3[LINEFEED (ADD | NO)]
    S1_7 --> S1_7_4[Baud rate: <value>]
    S1 --> S1_8[SYSTEM TEST]
    S1_8 --> S1_8_1[TURN ON ALL LCD PIXELS]
    S1_8 --> S1_8_2[TURN ON RED LCD PIXELS]
    S1_8 --> S1_8_3[TURN ON GREEN LCD PIXELS]
    S1_8 --> S1_8_4[TURN ON BLUE LCD PIXELS]
    S1_8 --> S1_8_5[HARDWARE CONTROL]
    S1_8_5 --> S1_8_5_1[SYSTEM HARDWARE CONTROL: (AUTO | MANUAL)]
    S1_8_5_1 --> S1_8_5_1_1[VALVE (OPEN | CLOSED)]
    S1_8_5_1 --> S1_8_5_1_2[Osc (INTERN | EXTERN) <period> μs]
    S1_8_5_1 --> S1_8_5_1_3[RELAY 1 (OPEN | CLOSED)]
    S1_8_5_1 --> S1_8_5_1_4[RELAY 2 (OPEN | CLOSED)]
    S1_8_5_1 --> S1_8_5_1_5[SOUND (ON | OFF)]
    S1_8_5_1 --> S1_8_5_1_6[Analog Outputs]
    S1_8_5_1_6 --> S1_8_5_1_6_1[Analog Outputs Hardware Control: (AUTO | MANUAL)]
    S1_8_5_1_6_1 --> S1_8_5_1_6_1_1[Analog Outputs Hardware Control: AUTO]
    S1_8_5_1_6_1_1 --> S1_8_5_1_6_1_1_1[0-10 VDC Voltage DAC #1 value: <value>]
    S1_8_5_1_6_1_1 --> S1_8_5_1_6_1_1_2[4-20 mADC Current DAC #1 value: <value>]
    S1_8_5_1_6_1 --> S1_8_5_1_6_1_2[Analog Outputs Hardware Control: MANUAL]
    S1_8_5_1_6_1_2 --> S1_8_5_1_6_1_2_1[0-10 VDC Voltage DAC #1 value: (nnnn)]
    S1_8_5_1_6_1_2 --> S1_8_5_1_6_1_2_2[4-20 mADC Current DAC #1 value: (nnnn)]
    S1_8 --> S1_8_6[TOUCH SCREEN CALIBRATE]
    S1_8 --> S1_8_7[TOUCH SCREEN TEST]
    S1_8_7 --> S1_8_7_1[DRAG]
    S1_8_7 --> S1_8_7_2[DRAW]
    S1_8_7 --> S1_8_7_3[QUIT]
    S1_8 --> S1_8_8[CURSOR (ON | OFF)]
    S1 --> S1_9[FACTORY RESTORE]
    S1_9 --> S1_9_1[RESTORE]

    HISTORY --> H1[LN2 ON/OFF]
    H1 --> H1_1[<]
    H1 --> H1_2[IN]
    H1 --> H1_3[OUT]
    H1 --> H1_4[>]
    H1 --> H1_5[>>]




    SHUTDOWN --> SH1[SHUTDOWN]
  
```






## SCREEN DESCRIPTIONS

The Model 1700 Instrument displays various screens on the graphic display when a MENU item is chosen. The items for each screen are described in the following table.




### Model 1700 Level Instrument Screen Descriptions

Screen <sup>a</sup>	Item Label	Field Type or Function
<b>SENSOR CONFIGURATION(S)</b>  MENU > SENSORS	NITROGEN OSCILLATOR USED:	Information: INTERNAL, EXTERNAL <sup>b</sup>
	NITROGEN LEVEL ON HOME SCREEN?	Toggles between: YES, NO
	NITROGEN LEVEL DAMPING	Move to another screen
	CAL N2	Move to another screen
	SENSOR NAME(S)	Move to another screen
<b>LIQUID NITROGEN LEVEL DAMPING</b>  MENU > SENSORS > NITROGEN LEVEL DAMPING	Undamped O/T Period:	Information: <value> <sup>b</sup> μs (present value)
	Undamped LN2 Level:	Information: <value> <sup>b</sup> % (present value)
	Damping Time Period:	Data entry: <value> seconds
	Damped O/T Period:	Information: <value> <sup>b</sup> μs (present value)
	Damped LN2 Level:	Information: <value> <sup>b</sup> % (present value)
<b>LIQUID NITROGEN SENSOR CALIBRATION</b>  MENU > SENSORS > CAL N2	OSCILLATOR PERIOD:	Information: <value> <sup>b</sup> μs (present value)
	SENSOR ACTIVE LENGTH:	Data entry: <value> cm or in (tap units to change)
	PERFORM MAX CAL	Move to another screen and Information: <value> <sup>b</sup> in μs followed by the date/time of last entry
	PERFORM MIN CAL	Move to another screen and Information: <value> <sup>b</sup> in μs followed by the date/time of last entry
	NO SENSOR CAL	Move to another screen and Information: <value> <sup>b</sup> in μs followed by the date/time of last entry
	APPROX MAX CAL:	Data entry: <value> followed by a button that toggles between: NO NOTIFY, NOTIFY <sup>c</sup>





## Model 1700 Level Instrument Screen Descriptions

Screen <sup>a</sup>	Item Label	Field Type or Function
<b>SENSOR NAME CONFIGURATION</b>  <p>MENU &gt; SENSORS &gt; SENSOR NAME(S)</p>	NITROGEN SENSOR NAME:	Data entry: <value>
<b>OUTPUT CONFIGURATION: ALARMS</b>  <p>MENU &gt; OUTPUTS &gt; ALARMS</p>	ALARM 1 SOURCE:	Toggles between: DISABLED, NITROGEN
	ALARM 1 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	ALARM 1 WHEN LEVEL <state>	Toggles between: ≤, ≥
	ALARM 2 SOURCE:	Toggles between: DISABLED, NITROGEN
	ALARM 2 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	ALARM 2 WHEN LEVEL <state>	Toggles between: ≤, ≥
	RELAYS	Move to another screen
	AUTOFILL	Move to another screen
	ANALOG OUTPUTS	Move to another screen
<b>OUTPUT CONFIGURATION: RELAYS</b>  <p>MENU &gt; OUTPUTS &gt; RELAYS</p>	RELAY 1 SOURCE:	Toggles between: DISABLED, NITROGEN
	RELAY 1 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	RELAY 1 CLOSES WHEN LEVEL <state>	Toggles between: ≤, ≥
	RELAY 2 SOURCE:	Toggles between: DISABLED, NITROGEN
	RELAY 2 SETPOINT:	Data entry <value> %, in, or cm (tap units to change)
	RELAY 2 CLOSES WHEN LEVEL <state>	Toggles between: ≤, ≥
	ALARMS	Move to another screen
	AUTOFILL	Move to another screen
	ANALOG OUTPUTS	Move to another screen




## Model 1700 Level Instrument Screen Descriptions

Screen <sup>a</sup>	Item Label	Field Type or Function
<b>OUTPUT CONFIGURATION: AUTOFILL</b>  MENU > OUTPUTS > AUTOFILL	AutoFILL CONTROL:	Toggles between: DISABLED, NITROGEN
	START:	Data entry <value> %, in, or cm (tap units to change)
	STOP:	Data entry <value> %, in, or cm (tap units to change)
	FILL TIMEOUT:	Data entry <value> MINUTES
	RELAYS	Move to another screen
	ALARMS	Move to another screen
	ANALOG OUTPUTS	Move to another screen
<b>OUTPUT CONFIG.: ANALOG OUTPUTS</b>  MENU > OUTPUTS > ANALOG OUTPUTS	OUTPUT 0-10 VDC SOURCE:	Toggles between: DISABLED, NITROGEN
	OUTPUT 4-20 mA SOURCE:	Toggles between: DISABLED, NITROGEN
	RELAYS	Move to another screen
	ALARMS	Move to another screen
	AUTOFILL	Move to another screen
<b>NETWORK CONFIGURATION</b>  MENU > NETWORK	Addressing Scheme:	Toggles between: DYNAMIC, DISABLED, STATIC
	Address:	Information: <value> <sup>d</sup> (DYNAMIC) Data entry: <value> (STATIC)
	Netmask:	Data entry: <value> <sup>b</sup>
	Gateway:	Data entry: <value> <sup>b</sup>
	Hostname:	Data entry: <value>
	MAC Address:	Information: <value> <sup>b</sup>

## Model 1700 Level Instrument Screen Descriptions

Screen <sup>a</sup>	Item Label	Field Type or Function
<b>SYSTEM CONFIGURATION</b>  MENU > SYSTEM	Instrument Serial Number:	Information: <value> <sup>b</sup>
	HW: FPGA: FW:	Information: <values> <sup>b</sup>
	Date of Manufacture:	Information: <value> <sup>b</sup>
	SET TIME	Move to another screen
	UPDATE	Move to another screen
	RS-232 SETUP	Move to another screen
	SYSTEM TEST	Move to another screen
<b>SYSTEM DATE/TIME</b>  MENU > SYSTEM > SET TIME	Current System Date and Time:	Information <sup>b</sup>
	System Timezone:	Information <sup>b</sup>
	SET MANUALLY	Move to another screen
	SET FROM NTP SERVER	Action with status screen
<b>FIRMWARE UPDATE INSTALL</b>  MENU > SYSTEM > UPDATE	Select Firmware Version:	Selection list: <values> and then moves to INSTALL verification on selection
<b>RS-232 SERIAL PORT CONFIGURATION</b>  MENU > SYSTEM > RS-232 SETUP	Current Configuration:	Information <sup>b</sup>
	Configuration Choice:	Choose alternate function by touching
	Echo Setting:	Choose ECHO OFF or ECHO ON
	Line Ending:	Choose NO LF or ADD LF
	Baud rate:	Choose desired baud rate

## Model 1700 Level Instrument Screen Descriptions

Screen <sup>a</sup>	Item Label	Field Type or Function
<b>SYSTEM TEST</b>  MENU > SYSTEM > SYSTST	TURN ON ALL LCD PIXELS	Move to multi-color test screen
	TURN ON RED LCD PIXELS	Move to red test screen
	TURN ON GREEN LCD PIXELS	Move to green test screen
	TURN ON BLUE LCD PIXELS	Move to blue test screen
	HARDWARE CONTROL	Move to hardware test screen
	TOUCH SCREEN ALIBRATE	Move to touch calibration screen
	TOUCH SCREEN TEST	Move to touch test screen
	CURSOR ON/CURSOR OFF	Toggles visibility of cursor
<b>HISTORY</b>  MENU > HISTORY	LN2	Enables/disables nitrogen trace
	<	Move graph back in time
	IN	Zoom in (less time per division)
	OUT	Zoom out (more time per division)
	>	Move graph forward in time
	>>	Shift graph to most recent levels
<b>SHUTDOWN</b>  MENU > SHUTDOWN	SHUTDOWN	Shuts down the instrument in an orderly fashion which reduces boot time for the next power on.

- The path below each illustration indicates the necessary actions in the menu structure required to display the indicated screen.
- Displays the state or value (display only).
- If the NOTIFY state is selected and the Approx. Cal factor is not equal to 1.00, a brief message will be displayed indicating an approximate calibration is in effect when the instrument boots.
- Value displayed is chosen by the DHCP server for the network in DYNAMIC mode.


## CAPACITANCE (LIQUID NITROGEN) LEVEL

### CONFIGURE THE INSTRUMENT TO DISPLAY NITROGEN LEVEL

**NOTE** If the instrument was purchased with a capacitance-based level sensor, Steps 1 through 5, below have already been performed.


1. From the main screen, choose the following: **MENU > SENSORS**.

2. The instrument has a built-in sensor oscillator which will be used if the level sensor is connected to the BNC connector on the rear panel of the instrument with a maximum of 15 feet of RG59/U coaxial cable (refer to the figure on page 2). Ensure the **NITROGEN OSCILLATOR USED** field displays **INTERNAL**.



INTERNAL Oscillator Indicated

- a. If the sensor is greater than 15 feet from the instrument, an external oscillator/transmitter is required and the **NITROGEN OSCILLATOR USED** will display **EXTERNAL** (refer to the figure on page 3).
3. Ensure that **NITROGEN LEVEL ON HOME SCREEN?** is set to YES.
  4. Press the **HOME** icon at the bottom of the screen.



Nitrogen Level Displayed on Home Screen

---

## ALARMS AND RELAYS


---

The Model 1700 Instrument has two types of alarms, level-based and time-based alarms.

---

### ALARM STATUS SCREEN

Prior revisions of the Model 1700 have attempted to indicate the alarm conditions in the main screen footer region, including cycling through multiple alarm texts in the footer. This has been replaced by an Alarm Status screen that summarizes all alarm states and provides a button for muting.



Alarm Status Screen

The Alarm Status screen can be accessed at any time by tapping on the **ALARM** text in the footer. All alarms are summarized in the Alarm Status screen. If the Alarm Status screen shows the “ALARMING” button (as illustrated above), then touching the button will mute the alarm(s). This mute function will remain in effect until a new alarm is activated, or the status of any active alarm changes. Shortcuts may also be provided on the Alarm Status screen to enter the Level Alarms and Autofill configuration screens.

---

### LEVEL-BASED ALARMS

The Model 1700 Instrument has two user-configurable level alarms. Each alarm can be triggered by either level measurement (for dual level configured instruments). Each alarm can be configured to be active above or below a user-defined setpoint.

When an alarm condition occurs, an audible alert will sound and visual indication will appear on the front panel.

Two relays can also be configured with independent setpoints to actuate independently of the alarm states. These relays have Normally Open (NO) contacts. The relays can be assigned to either the helium or nitrogen channel and each relay can be configured to close when the reading is either  $\leq$  or  $\geq$  the setpoint.

As an example of this setup flexibility, a level channel can be configured to have an *alarm condition* when the level is outside a normal operating band. The *relay actions* can be set to act identically, or configured for other conditions, per the user-defined application.


## TIME-BASED FILL ALARM

The Model 1700 Instrument has an alarm to indicate that there is a problem with the autofill function. If enabled, the instrument will start a timer when an autofill condition is initiated, and if the level has not reached the fill stop level within the user-set period of time, an Autofill Timeout alarm will occur.

This alarm will cause three things to occur:

1. The de-energizing of the fill valve socket on the instrument rear panel,
2. An audible alarm will sound,
3. A **TIMEOUT** visual indication on the Home screen (as shown at right) or in the footer of the instrument front panel.

The Autofill Timeout function can be disabled by setting the Fill Timeout interval to 0 minutes.



Fill Timeout Alarm Indicators

## MULTIPLE ALARMS

More than one alarm condition can occur at the same time. If that is the case, the Alarm Status screen summarizes all the alarm states. Access the Alarm Status screen by touching the **ALARM** text in the footer.


**NOTE** Prior versions of the instrument flashed a series of alarms in the screen footer. This has been superseded by the Alarm Status screen.




## CONFIGURING ALARM SETPOINTS

1. From the **MENU** screen, select **OUTPUTS** and the **OUTPUT CONFIGURATION: ALARMS** screen will be displayed.

(If already in another **OUTPUT CONFIGURATION** screen, choose the **ALARMS** button.)




2. Ensure the Alarm Source fields are set to **NITROGEN** (or **DISABLED**).




3. Set the levels at which the alarm will be triggered in the Setpoint fields highlighted in red in the image at right.

**NOTE** Tapping the units symbol (% , in, or cm) allows selection of the available unit selections for the selected source liquid.



4. Use the  $\leq$  or  $\geq$  button to toggle between the two states of alarm, either alarm when the indicated level is less than or equal to the setpoint or alarm when the indicated level is greater than or equal to the alarm setpoint.



## ACKNOWLEDGING AN ALARM

**NOTE** The alarms are not “latched” so if the alarm condition clears itself, the instrument will remove the alarm condition.

1. When an alarm is initiated, several things will occur:


- a. The **ALARM** text in the footer will illuminate in bold, red text.
- b. An audible alarm will be energized.
- c. The Alarm Status screen, accessible by touching the **ALARM** text in the footer, will display a summary of all alarms.



Alarm Annunciator

## MUTING AN ALARM

To mute the audible alarm, first touch the **ALARM** text in the footer. The Alarm Status screen will appear with the “ALARMING” status shown. Tap the “ALARMING” button to mute all active alarms. A change in any alarm state will also re-enable the audible alarm for any remaining active alarms. There is no permanent mute.




Muted Alarm Condition


## CONFIGURING RELAY SETPOINTS

The relay outputs are available via the AUX I/O connector on the rear panel of the instrument per the pinout provided on page 98.

1. From the **MENU** screen, select **OUTPUTS** and then choose **RELAYS**.

OUTPUT CONFIGURATION: RELAYS  
Screen


- Toggle the Relay Source fields to choose **NITROGEN** (or **DISABLED**).



RELAYS Source Configuration


- Set the levels at which the relay will be closed in the Setpoint fields.

**NOTE** Tapping the units symbol (% , in, or cm) allows selection of the available unit selections for the selected source liquid.



RELAYS Setpoints

- Use the  $\leq$  or  $\geq$  button to toggle between the two states of actuation. Either the relay will close when the indicated level is less than or equal to the setpoint, or when the indicated level is greater than or equal to the setpoint.




Relay Closed Above or Below Setpoints

## CONFIGURE THE AUTOFILL FUNCTION

There are several variables that must be addressed to set up an autofill system. These include the level indication that will be used to control autofill, the Fill Start level (A), the Fill Stop Level (B), and the Fill Timeout interval.

### SETTING THE AUTOFILL PARAMETERS


1. From the **MENU** screen, choose **OUTPUTS** and then choose **AUTOFILL**.
2. Toggle the **AutoFILL CONTROL** button until **NITROGEN** is displayed.



3. Enter the Fill **START** and **STOP** levels.


**NOTE** If the **START** value is set to 0.0% the autofill function will be effectively disabled.

**NOTE** Tapping the units symbol (% , in, or cm) allows selection of the available unit selections for the selected source liquid.



4. Enter the **FILL TIMEOUT** interval. Refer to page 37 for a description of the Autofill time-out function.

**NOTE** If the **FILL TIMEOUT** value is set to 0 minutes, the timeout function is disabled.



5. Press **HOME** in the footer to return back to the level display.

## ENABLE THE AUTOFILL FUNCTION

The autofill function must first be enabled. After the autofill control loop has been configured (the parameters in the prior section) the instrument will enter the **AUTO-OFF** state.

To enable the autonomous autofill function:

1. Touch the **AUTO-OFF** text in the display footer until it reads **AUTO-ON**.
2. Touch **SAVE** to enable the autonomous autofill function.

The AUTOFILL function has four manually-selectable states by touching the fill state indicator either in the **HOME** screen or in the footer when not in the **HOME** screen:

**AUTOFILL States**

Function	Operation	Overrides
<b>AUTO-ON</b>	Automatically maintains level between fill START and STOP set-point (i.e. autonomous mode).	Autofill will alarm and cease if fill valve stays open for $\geq$ the FILL TIMEOUT setting .
<b>M-OPEN</b>	Manual Open: Energizes the valve control socket on the rear panel.	Any active ALARM with $\geq$ setpoint defined for the fill liquid will close valve.
<b>M-CLOSED</b>	Manual Close: De-energizes the valve control socket on the rear panel.	None
<b>AUTO-OFF</b>	Disables the AUTOFILL function.	None

## CLEARING THE AUTOFILL TIMEOUT ALARM

**NOTE** A Loss of Sensor error condition (see page 17) on the nitrogen channel will also halt an active autofill in addition to the autofill timeout feature.

1. Touch the **TIMEOUT** text in the **HOME** screen or in the footer until **AUTO-ON** is displayed.
2. Touch **SAVE** in the footer. This will restart the autonomous autofill function.
3. To clear the **TIMEOUT** *without restarting* the autofill function, choose the **M-CLOSED** state and then choose **SAVE** in the footer.

## SELECT THE APPROPRIATE UNITS ON THE DISPLAY

Touch the units on the display to change the units. The available units are percent (%), inches (in), and centimeters (cm).

**NOTE** Tapping the units symbol (% , in, or cm) in a configuration screen also allows selection of the available unit selections. The units selection is used for all indication and configuration values.



LN<sub>2</sub> Home Screen

## DAMPING CONFIGURATION

The level measurement function offers a damping feature which can be used to dampen fluctuations in the nitrogen level reading that arise from sloshing liquids or other surface disruptions.

To set the damping for the nitrogen measurement:


1. From the main screen, choose the following: **MENU > SENSORS**. Ensure the nitrogen level measurement is enabled and displayed on the **HOME** screen.
2. Touch the **NITROGEN LEVEL DAMPING** button to enter the damping configuration screen.



Choose Nitrogen Level Damping

3. In the **LEVEL DAMPING** screen, enter the desired damping time in seconds. The damping is a first-order low-pass filter function where the damping period is equal to *five times the filter time constant* (99% of final reading).

A damping period of 0.0 seconds means the damping function is disabled.



Damping Configuration


**NOTE** If the damping period is set to a non-zero value, the damping function will then be applied to the level measurement used for *all* indication and control functions for the liquid nitrogen level. This includes analog output and fill cycle functions.

## ANALOG OUTPUT SIGNALS

Refer to page 98 of the Appendix for the AUX connector pin-out.

### CONFIGURING THE ANALOG OUTPUTS

1. From the **MENU** screen, choose **OUTPUTS**, then **ANALOG OUTPUTS**.
2. If necessary, choose the source for the **0-10 VDC** output and **4-20 mA** output.
3. Press the **SAVE** button to save the choice (or **CANCEL** to quit without making a change).
4. Press **HOME** to go back to the home screen.




ANALOG OUTPUTS Source Selection

## ETHERNET CONNECTIVITY


### IP ADDRESSING SCHEME

1. From the **MENU** screen, choose **NETWORK**. The current settings will be displayed.




Network Configuration Screen

2. Touch the addressing scheme button, to choose **STATIC**, **DYNAMIC**, or **DISABLED** as appropriate.



Editing the Networking Mode

3. If **STATIC** is chosen for the addressing scheme, enter **ADDRESS**, **NETMASK**, and **GATEWAY** values that are appropriate for the connected network. Once an address field is touched, the pop-up keyboard will be presented for data entry. Press **SAVE** when done.



Editing Static Network Parameters

4. Set the **HOSTNAME** field as desired.
5. If the **DYNAMIC** addressing scheme is chosen, the **IP ADDRESS**, **NETMASK**, and **GATEWAY** addresses will automatically be assigned from a network DHCP server. Press **SAVE** when done.
6. If changes are made, select the **SAVE** button and then the instrument will reconfigure itself with the chosen network settings.




## SERIAL CONNECTIVITY

The serial (RS-232) connectivity can be customized per the following instructions.


### CONFIGURING THE RS-232 SETTINGS

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.



MENU Selection Screen

3. From the **SYSTEM CONFIGURATION** screen choose **RS232 SETUP**.



SYSTEM CONFIGURATION Screen

4. The RS-232 Configuration screen shows the current function for the serial connection.



5. Choose the alternate function for the serial connection if desired.



6. Choose to customize the port behavior for character echo, line ending, and baud rate.



7. If changes are made, select the **SAVE** button in the footer.

## ABNORMAL OPERATION

### CAPACITANCE SENSOR CONTAMINATION

To ensure proper instrument calibration and operation, care must be taken to ensure the sensor is kept free of contaminants and not subjected to any force which would physically distort the sensor. Water or other electrically conducting substances in the sensor will disturb the measured capacitance and the instrument's response. Physically distorting the sensor in any way will also cause abnormal instrument operation by introducing variations in the sensor capacitance not due to liquid level. The absolute calibration of the instrument can be inaccurate if care is not taken to ensure the sensor is in the proper environment.

Cold sensors exposed to humidified air can show erroneous high level readings because the air contains moisture that can condense between the cold sensing tubes. A minute film of water can cause a shorted or partially shorted condition, which results in false level readings. As the sensor warms, the moisture may evaporate and the sensor will again read correctly. This is a physical phenomenon and does not indicate any problem with your AMI level equipment. Limit or eliminate exposure of cold sensors to humidified air to avoid this condition.

If a sensor should require cleaning and the sensor is for use with liquids other than liquid oxygen, flushing with pure alcohol is recommended. The sensor cannot be used again until all the alcohol has been evaporated. Under no circumstances should the sensor be disassembled.

For sensors to be used with liquid oxygen (LOX), although measures are taken to minimize oils and greases during manufacture, no special cleaning required for LOX service is provided by AMI. Certified LOX cleaning is the responsibility of the customer.


## RESETTING THE INSTRUMENT TO FACTORY DEFAULTS

1. Press **MENU**, then **SYSTEM**, then **FACTORY RESET**.



System Configuration Screen

2. Press **RESTORE**.
3. The instrument will be reset to factory defaults.




Resetting Instrument to Factory Defaults

## SHUTTING THE INSTRUMENT DOWN

**NOTE** The Model 1700 Instrument is a Linux-based computer system and in order to ensure the file system is properly unmounted, the **SHUTDOWN** function should be invoked. If it is not, i.e. the instrument is shut down by removing power via the front panel power switch, the next time the instrument boots up, it will have to scan the memory system to ensure everything is in order.

1. The instrument should be shut down by using the menu function **SHUTDOWN**.



Invoking Instrument Shut Down

2. Choose **SHUTDOWN** again to confirm.
3. When prompted, turn off the front panel power switch.



Confirming Instrument Shutdown

# CALIBRATION

Model 1700 instrument is calibrated at the factory for a specific length sensor(s) for use in a specific liquid(s). The calibration length(s) and calibration liquid(s) are listed on the calibration sticker on the bottom of the instrument. For capacitance sensors, if the factory calibration method utilized was approximate, the calibration length will be noted as an approximate value.


## SETTING THE SYSTEM DATE AND TIME

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.




MENU Selection Screen

- From the **SYSTEM CONFIGURATION** screen, choose **SET TIME**.




SYSTEM CONFIGURATION Screen

- In the **SYSTEM DATE/TIME** screen, choose to either set the date and time *manually*, or set it automatically from an NTP server if the instrument has Internet access.



SYSTEM DATE/TIME Screen

- For the **SET MANUALLY** selection, edit the **TIMEZONE**, **YEAR**, **MONTH**, **DAY**, **HOURL**, and **MIN** fields as necessary. Touching in a field will launch the keyboard on the screen. Edit the information in the field as necessary and choose **ENTER** to enter the data in the field and close the pop up keyboard.



SYSTEM DATE/TIME Manual Entry

**NOTE** The clock is set to GMT at the factory and is battery backed. The **TIMEZONE** field must be set to a value from the TZ database.<sup>1</sup> The list of **TIMEZONE** values can be found at:

[https://en.wikipedia.org/wiki/List\\_of\\_tz\\_database\\_time\\_zones](https://en.wikipedia.org/wiki/List_of_tz_database_time_zones)

The **TIMEZONE** entry must exactly match the entire field as listed in the above URL in the *TZ database name* column and it is case sensitive. An instrument reboot may be required for the **TIMEZONE** to take effect. If the entered **TIMEZONE** supports Daylight Savings Time (DST), then the local time will be automatically adjusted to DST.

- Choose **SAVE** in the footer after all the fields have been edited as necessary.


1. The tz database is a collaborative compilation of information about the world's time zones, primarily intended for use with computer programs and operating systems. The topic of worldwide timezones is a complex issue so please refer to the on-line resource noted above.

## CAPACITANCE-BASED LEVEL CALIBRATION

### UNDERSTANDING THE SENSOR ACTIVE LENGTH

American Magnetics, Inc. fabricates the liquid level sensor with two vent holes; a lower vent hole in the side wall near the bottom which is typically the minimum liquid level calibration point and the upper vent hole in the sensor side typically near the top of the sensor. The liquid level location approximately 2.5 cm (1 in) below the upper vent hole is typically the 100% calibration point.

**NOTE** Without entry of the active length, the instrument will not be able to read out in units other than percent.



The Model 1700 Instrument requires the user to enter the calibrated, or active length, (physical distance between the MIN and MAX calibration locations on the sensor) in order for the absolute units function (inches, cm) to be displayed if desired.

The user must enter the sensor length in centimeters. Use the *Active Length* value noted on the level sensor documentation or measure the distance between the lower vent hole on the sensor and 1.0 inch (2.5 cm) below the upper vent hole on the sensor as illustrated at left.

Typical Capacitance-based Liquid Level Sensor

---

## RELATIONSHIP BETWEEN CALIBRATION AND SENSOR LENGTH

The capacitance-based method of measuring the liquid level operates by measuring the period of a signal from an oscillator, which can be internal or contained in an external oscillator/transmitter unit. As the liquid level varies, the value of the sensor capacitance varies proportionally. Since the dielectric properties of liquids vary and the component tolerances for the sensor and oscillator introduce variations, a calibration is required to assure maximum accuracy for a specific sensor immersed in the target liquid. The calibration MIN and MAX settings correspond to the minimum and maximum oscillation periods, respectively, for a given sensor and target liquid configuration.

The length setting of the instrument is only provided as a means of scaling the 0% (minimum calibration) to 100% (maximum calibration) range of the measurement to meaningful units of length. During the calibration it is important to accurately measure the distance between the physical locations on the sensor corresponding to the MAX and MIN calibration points. The measured value for the length will be used in configuring the instrument for operation.

---

## VARIATIONS IN THE DIELECTRIC WITH CHANGING DENSITY

For cryogenic liquids, the dielectric of the liquid will change with a change in density. The amount of change is dependent on the properties of the specific liquid. The figure below illustrates the variations in dielectric for nitrogen vs. pressure under *saturated* conditions.<sup>1</sup> Since the instrument uses a capacitance-based method for determining liquid level, such a change in the dielectric of the liquid will result in a shift in the level reading of the instrument. The calibration procedures described herein are most accurate when applied in situations where the operating conditions of the cryo-vessel are relatively constant, i.e. the operating pressure and temperature of the cryo-vessel are relatively constant.


To minimize the effects of shifts in the dielectric of the target liquid, perform a closed dewar calibration (see page 57) at the expected operating condition of the cryo-vessel. If this is not feasible, then calibrate the sensor at atmospheric pressure and use the approximate calibration method to compensate for the shift of the dielectric when the cryogenic liquid is under pressure. For this type of approximate calibration, the *reference liquid* will be the *target liquid* at atmospheric pressure — see page 62 for a detailed discussion of the approximate calibration method. If any questions exist in regard to calibration issues, contact AMI for assistance in determining the optimal calibration strategy.

**NOTE** All references to “dielectric constant” herein refer to the unitless relative dielectric to  $\epsilon_0$  ( $\epsilon_0$  is the dielectric constant of a vacuum).

---

1. Data obtained from NIST Standard Reference Database 12.





**Dielectric vs. pressure for nitrogen under saturated conditions.**

## CAPACITANCE-BASED SENSOR CALIBRATION METHODS

The most straightforward calibration method is the *Open Dewar Calibration* which requires the customer to have access to a filled dewar where the full active length of the sensor can be dipped. The *Closed Dewar Calibration* method can be performed in situations where it is not feasible for the customer to dip the sensor into an open dewar, such as situations where the target liquid is under pressure. The closed dewar calibration is more complex and may require initial preparations to insure success and avoid saturation of the signal during the calibration.


Occasionally customers ask AMI to calibrate an instrument and sensor for a liquid which is not available at AMI for calibration purposes and/or for a sensor which is too long to be calibrated at our facilities.

*For the case of the target liquid being unavailable*, AMI uses liquid nitrogen as the reference liquid and an *Approximate Calibration* is performed using mathematical manipulation of the ratio of the dielectric constants between liquid nitrogen and the desired liquid. This procedure is outlined in the *Approximate Calibration* section beginning on page 62. The technique is intended to provide the instrument with an approximate calibration so that it can be used immediately by the customer. However, the customer is still expected to perform a more accurate calibration where feasible, such as the open dewar or closed dewar calibration, with the target liquid.

For the case where a sensor is too long to be calibrated in AMI facilities, AMI will perform a partial length open dewar calibration in liquid nitrogen, and then calculate the MAX calibration point. A dielectric ratio may also be subsequently utilized to adjust for a target liquid other than liquid nitrogen. The customer is expected to perform a more accurate open dewar or closed dewar calibration if feasible.

### SELECTION OF CAPACITANCE SENSOR CALIBRATION METHODS

As a quick guide for selection of the best calibration method available, a calibration selection diagram is presented below.



Calibration method selection diagram.


If the instrument and sensor are purchased as a unit from AMI, then the factory calibration, including sensor serial number and sensor physical parameter information entered at the factory, will be adequate in most cases. However, for the exceptions noted in the previous paragraphs (which are *approximate* calibrations), the customer should perform a more accurate open dewar or closed dewar calibration if at all possible. A customer-performed calibration is also *required* for sensors that are purchased as a separate item from the instrument, since the instrument and sensor were not both available for calibration at AMI facilities.

Each Model 1700 Instrument must be calibrated with a sensor before use. If the instrument was purchased with a sensor to be used in liquid nitrogen, the instrument has been calibrated at the factory.


## PRE-CALIBRATION PROCEDURE

### ENTER CAPACITANCE SENSOR INFORMATION


1. Press the **MENU** button in the lower left corner of the display screen.




2. Choose the **SENSORS** selection from the **MENU** screen.
3. If the capacitance sensor will be longer than 15 feet from the instrument, ensure an Oscillator/Transmitter is used between the instrument and the sensor.



4. Verify that the oscillator source selection is correct, either **INTERNAL** or **EXTERNAL**. Note that this setting is auto-detected by the instrument at power-up or reboot,




5. Press the **CAL N2** button



CALIBRATE N2 Selection

6. Touch in the **SENSOR ACTIVE LENGTH** field and using the numerical keypad, enter the sensor active length in the selected units. Press **ENTER** and then **SAVE** at the bottom of the screen when finished



SENSOR ACTIVE LENGTH field


**NOTE** Tapping the units symbol (in or cm) allows selection of the available unit selections for that channel.

7. Press the **SENSOR NAME(S)** button.




SENSOR NAME(S) Selection

8. Touch in the **NITROGEN SENSOR NAME:** field. The keyboard will be displayed. Edit the displayed name and press **ENTER** at the bottom of the screen.



PERFORM MAX CAL. Selection Button

9. Press **HOME** in the screen footer to return to the level display screen.



Home Selection Button

### PERFORMING LOSS OF SENSOR CALIBRATION

The Model 1700 will detect when the liquid level sensor has been disconnected from the instrument and display a notification on the front panel. If the instrument has been configured for autofill, the fill valve will be shut, requiring operator intervention to resume autofill operation. This loss of sensor threshold must be calibrated as follows:


For a system where the internal oscillator is used, connect the coaxial cable to the BNC connector on the instrument rear panel but leave the cable disconnected from the sensor BNC connector. For a system where an external oscillator is used, connect the coaxial cable between the instrument and the oscillator/transmitter unit. Connect the second (6 ft) coaxial cable to the input of the oscillator/transmitter unit but leave the cable disconnected from the sensor BNC connector.

1. Press the **MENU** button in the lower left corner of the display screen.




MENU Selection Button

2. Choose the **SENSORS** selection from the **MENU** screen.




SENSORS Selection

3. From the **SENSOR CONFIGURATION(S)** screen, choose **CAL N2**.



CAL N2 Selection

4. Select the **NO SENSOR CAL** button.
5. When the calibration procedure is completed, choose **SAVE** in the footer of the screen.
6. Press **HOME** to return to the home screen.
7. Connect the BNC cable to the liquid level sensor.




NO SENSOR CAL button

## OPEN DEWAR CALIBRATION


The *Open Dewar Calibration* requires the customer to have access to a filled dewar where the full active length of the sensor can be dipped.

1. Press the **MENU** button in the lower left corner of the display screen.




MENU Selection Button

2. Choose the **SENSORS** selection from the **MENU** screen.




SENSORS Selection

3. Press the **CAL N2** button




CAL N2 Selection

4. Position the capacitance sensor in the target liquid at the 100% level. Hold the sensor at this location and press the **PERFORM MAX CAL** button.




PERFORM MAX CAL Selection

5. The instrument will display the following screen as it takes data. Once the calibration measurement stabilizes, press the **SAVE** button. The operator *must* press the **SAVE** button for the new cal point to be saved.



Updating MAX CAL Period

6. Position the capacitance sensor in the target liquid at the 0% level. Hold the sensor at this level and press the **PERFORM MIN CAL** button.



PERFORM MIN CAL Selection

7. While displaying “**CALIBRATING...**”, the instrument will display the **MIN LEVEL PERIOD**. Press the **SAVE** button to save the new minimum level calibration point.

**NOTE** Note that the periods listed to the right of the **PERFORM MAX CAL** and **PERFORM MIN CAL** buttons are updated as well as the date and time stamps of the calibration points.

**NOTE** The minimum and maximum calibration can be performed in either order. Also, either the minimum or maximum calibration point can be updated without altering the other calibration point.




## CLOSED DEWAR CALIBRATION

A calibration can be performed in a closed dewar system by monitoring the liquid level while transferring the target liquid to an initially empty (or near empty) dewar at a constant rate. In order to insure success with the closed dewar technique, it is necessary to prepare the instrument by presetting the calibration minimum and maximum calibration points outside the estimated level range. If the instrument is not prepared in this manner before the calibration procedure, it is possible to reach the maximum calibration point of the instrument before the target vessel is at the desired maximum level point. If minimum and maximum liquid level indication is available via some other means (e.g. flow calculation, visual determination, point sensors, etc.), then the presetting of the instrument is not necessary.


### PRESETTING THE MAXIMUM AND MINIMUM CALIBRATION POINTS

The following procedure should be performed before installation of the sensor in the target cryo-vessel.

1. Connect the sensor coaxial cable to the BNC connector on the rear panel of the instrument (see page 13). *Do not connect the sensor.* Energize the instrument.
2. Press the **MENU** button in the lower left corner of the display screen.
3. Choose the **SENSORS** selection from the **MENU** screen.




MENU Selection Button




SENSORS Selection

4. Press the **CAL N2** button



CAL N2 Selection


5. Press the **PERFORM MIN CAL** button.
6. While displaying “**CALIBRATING...**”, the instrument will display the **MIN LEVEL PERIOD**. Press the **SAVE** button to save the new minimum level calibration point.



PERFORM MIN CAL Selection

**NOTE** Note that the period listed to the right of the **PERFORM MIN CAL** button is updated as well as the date and time stamp of the calibration point.

7. Connect the sensor to the oscillator coaxial cable that is connected to the instrument.
8. Perform the maximum level calibration by pressing **PERFORM MAX CAL**.
9. While displaying “**CALIBRATING...**”, the instrument will display the **MAX LEVEL PERIOD**. Press the **SAVE** button to save the new maximum level calibration point.




PERFORM MAX CAL Selection

10. Calculate the factor  $C_{adj}$  using the following equation:

$$C_{adj} = 1.20 \left[ 1 + \frac{2.1(L_{active})}{5.2(L_{total})} \right] \left[ \frac{\epsilon - 1}{0.454} \right]$$

where  $L_{total}$  is the total sensor length in inches,  $L_{active}$  is the active sensor length in inches, and  $\epsilon$  is the dielectric constant of the *target liquid*.

11. Enter  $C_{adj}$  into the instrument by touching the **APPROX MAX CAL** field.
12. Using the pop up numeric keypad, enter the  $C_{adj}$  value and press the **APPLY** button.
13. With the sensor connected, again press the **PERFORM MIN CAL** button and **SAVE** the result. The presetting procedure is complete. Proceed to the remainder of the closed dewar calibration procedure as presented below.



APPROX MAX CAL Entry


## COMPLETING THE CLOSED DEWAR CALIBRATION PROCEDURE

1. Install the sensor in the dewar and energize the instrument with the sensor connected to the instrument via the oscillator (if required) and extension cable(s) (see the system diagram on page 2).
2. Connect a strip chart recorder or graphical data logging application to the recorder output terminals on the rear panel of the instrument. If the recorder output is not available, the 4-20 mA current loop output may be used if installed, or an installed communications option can be used to query the instrument for the liquid level at regular time intervals during the calibration procedure. If no remote monitoring or communication option is installed, the level display must be manually plotted vs. time during the procedure.
3. Refer to "Analog output signals" on page 39. to configure the recorder output or current loop output.
4. Commence filling the dewar. While the sensor is cooling down, there may be a slow drift in the displayed liquid level. However, when the liquid actually touches the bottom of the sensor, contact with the liquid surface may become apparent by virtue of more random and frequent fluctuations in the displayed liquid level. The liquid level trace will also start to show an increasing profile with positive slope.

Once the indications of the contact between the sensor and liquid become readily apparent, use the **PERFORM MIN CAL** procedure below to save the new minimum calibration point. This point is the 0% level of the sensor when the **PERFORM MIN CAL** procedure is finished and saved becomes the 0% level.

**NOTE** If the sensor is installed in the dewar with some small amount of liquid already in contact with the sensor, then the final minimum calibration point can be set before filling begins but after any thermally induced fluctuations in the observed output have diminished. However, note that the measured span of the liquid level is reduced by the initial level of liquid in contact with the sensor.

5. Press the **PERFORM MIN CAL** button.
6. While displaying “**CALIBRATING...**”, the instrument will display the **MIN LEVEL PERIOD**. Press the **SAVE** button to save the new minimum level calibration point.



PERFORM MIN CAL Selection


**NOTE** Note that the period listed to the right of the **PERFORM MIN CAL** button is updated as well as the date and time stamp of the calibration point.

7. Continue the transfer while observing the liquid level trace on the strip chart recorder or computer display, whose slope is proportional to the transfer rate. The slope of the liquid level trace should decrease significantly when the liquid reaches the hole in the top of the sensor.

When the break in the slope of the level trace occurs (i.e. the slope of the level trace becomes 0 or horizontal), perform a **PERFORM MAX CAL** procedure below. The level on the sensor when the **PERFORM MAX CAL** procedure is finished and saved becomes the 100% level.


**NOTE** If the instrument displayed a 100% reading before a break is observed in the slope of the level trace, then the maximum calibration point set prior to the current procedure has interfered. If this occurs, the customer has two options: 1) stop the procedure, repeatedly multiply  $C_{adj}$  by 1.20 (see steps 4 and 5 of the presetting procedure) and enter as the new APPROX MAX CAL value until the current liquid level display falls below 100%, and then continue the procedure; or 2) continue the liquid transfer until the liquid level is determined to be 100% by means other than feedback from the instrument and then performing the maximum calibration procedure.

8. Perform the maximum level calibration by pressing **PERFORM MAX CAL**.



PERFORM MAX CAL Selection

9. While displaying “**CALIBRATING...**”, the instrument will display the **MAX LEVEL PERIOD**. Press the **SAVE** button to save the new maximum level calibration point.



Updating MAX CAL Period

**NOTE** Note that the period listed to the right of the **PERFORM MAX CAL** button is updated as well as the date and time stamp of the calibration point.

10. To achieve a standard calibration of the sensor with the active region located from the lower hole to one inch below the upper hole, use the level data from the instrument to recalibrate the maximum point when the percent level corresponds to one inch below the upper hole. Use the following equation to determine the percent level at which to reset the maximum calibration point:

$$MAX_{percent} = 100 - 100 \left[ \frac{1}{L_{active}} \right]$$

where  $L_{active}$  is the active length of the sensor in inches. This technique can be used assuming the sensor was built as a standard sensor. If the sensor was made in a custom configuration, refer to the sensor documentation and/or drawing or contact AMI.

**Example:** 20" active length sensor:

When the sensor is calibrated by the closed dewar procedure, the actual length of calibration will be 21" (distance between the bottom and top holes in the sensor). When the liquid is 1" below the upper hole, the display will show 95.2% [e.g.  $100\% - (1"/21" \times 100\%)$ ]. When the liquid level reaches this point during usage, perform the **PERFORM MAX CAL** operation and **SAVE** the result. The instrument and sensor are now calibrated with a standard active region of 20". The length setting of the sensor in the instrument should also be configured for 50.8 cm (20").

## APPROXIMATE CALIBRATION

This procedure is the least accurate form of calibration and should be used only when the aforementioned calibration procedures are not viable. The approximate calibration method can be used in cases where the sensor cannot be dipped into the target liquid, the full active length of the sensor cannot be dipped into an open dewar, or both. Approximate calibration may also be useful for situations where the sensor cannot be dipped into the target liquid under the expected operating pressure.


If the target liquid is not available for dipping, a substitute non-conducting reference liquid can be used. If the full length of the sensor cannot be dipped, then a partial length dip can be performed. If both situations are encountered, then a partial length dip can be performed in a substitute reference liquid.

1. Press the **MENU** button in the lower left corner of the display screen.




MENU Selection Button

2. Choose the **SENSORS** selection from the **MENU** screen.




SENSORS Selection

3. Press the **CAL N2** button




CAL N2 Selection

4. While displaying “**CALIBRATING...**”, the instrument will display the **MAX LEVEL PERIOD**. Press the **SAVE** button to save the new maximum level calibration point.



Updating MAX CAL Period

5. Position the capacitance sensor in the target liquid at the 0% level. Hold the sensor at this level and press the **PERFORM MIN CAL** button.
6. While displaying “**CALIBRATING...**”, the instrument will display the **MIN LEVEL PERIOD**. Press the **SAVE** button to save the new minimum level calibration point.



PERFORM MIN CAL Selection

**NOTE** Note that the periods listed to the right of the **PERFORM MAX CAL** and **PERFORM MIN CAL** buttons are updated as well as the date and time stamps of the calibration points.

**NOTE** The minimum and maximum calibrations can be performed in either order. Also, either the minimum or maximum calibration point can be updated without altering the other calibration point.

7. Measure the distance between the bottom hole of the sensor and the location of the liquid level dipped for max calibration. This measured length is  $L_{dipped}$ .
8. The dielectric constant for the reference liquid,  $\epsilon_1$ , and the target liquid,  $\epsilon_2$ , must be known to complete the approximate calibration. These values must be placed in the equation:

$$\text{APPROX MAX CAL} = \left[ \frac{\epsilon_2 - 1}{\epsilon_1 - 1} \right] \frac{L_{active}}{L_{dipped}}$$

where  $L_{dipped}$  is the length of the sensor dipped in the reference liquid and  $L_{active}$  is the active sensor length.

**NOTE** If the target liquid is available for dipping (i.e. the reference liquid and target liquid are the same), then the dielectric ratio,  $(\epsilon_2 - 1) / (\epsilon_1 - 1)$ , becomes 1. If the full active length of the sensor can be dipped, then the length ratio,  $L_{active} / L_{dipped}$ , becomes 1.

Note that  $\epsilon_1 = 1.454$  for liquid nitrogen at  $-203^\circ\text{C}$  at atmospheric pressure. Dielectric constants for several liquids are provided in the Appendix. The dielectric constant varies with temperature and pressure, therefore for best accuracy use the dielectric constant for the target liquid at the temperature and pressure maintained in the containing vessel.

9. Touch the **APPROX MAX CAL** field and using the numerical keypad that pops up, enter the value to be applied and press **SAVE** in the footer.

LIQUID NITROGEN SENSOR CALIBRATION			
OSCILLATOR PERIOD:	96.207	μs	
SENSOR ACTIVE LENGTH:	80.0	cm	
PERFORM MAX CAL	100.991	μs	2020-11-05 09:56
PERFORM MIN CAL	95.695	μs	2020-11-05 09:56
NO SENSOR CAL	73.881	μs	2020-11-05 09:57
APPROX MAX CAL	1.000		NO NOTIFY
HOME BACK ALARM			Help

**APPROX MAX CAL Entry**

**NOTE** The calibration values are retained in the instrument memory, therefore it is possible to repeatedly adjust the **APPROX MAX CAL** value without losing the original MIN/MAX calibration points.

**Example:** Purchased a 100" active length sensor for operation in liquid argon at atmospheric pressure, however only liquid nitrogen is available for calibration at a maximum depth of 30":

First, the sensor is dipped as far as possible into the liquid nitrogen and cooled. The minimum point is then set as outlined in step 2. The maximum point is set as outlined in step 3 while the sensor is submerged 30" in liquid nitrogen. The dielectric constant for liquid nitrogen is 1.454 and for liquid argon is 1.53. Substituting all values into the approximate calibration factor equation yields:

$$\text{APPROX MAX CAL} = \left[ \frac{1.53 - 1}{1.454 - 1} \right] \frac{100}{30} = 3.891$$

A value of 3.891 would be entered as the **APPROX MAX CAL** factor as outlined in step 9 above. The sensor is now approximately calibrated for 100" active length operation in liquid argon.

10. The sensor can now be installed in the dewar containing the target liquid. The approximate calibration can be used until an open dewar or closed dewar calibration can be performed with the target liquid.

Proceed to the *Operation* section for directions for configuring the instrument.



# REMOTE INTERFACE REFERENCE

---

The Model 1700 Instrument provides both serial (RS-232) and Ethernet interfaces as standard features. The serial and Ethernet interfaces can be operated simultaneously. Separate output buffers are provided for the serial and Ethernet return data. The serial and Ethernet interfaces accept both legacy and SCPI syntax for commands and queries.

The Model 1700 provides a web browser connection (port 80) via TCP/IP and presents a mirror of the local display which has hot spots that can be clicked or tapped, and which also accepts keyboard entry for textual and numerical input. The instrument's IP address can be ascertained by referring to the section titled "IP Addressing Scheme" on page 40. By using a web browser to connect to the instrument, all functionality of the Model 1700 can be controlled via the remote connection.

Communication with the Model 1700 via Ethernet or RS-232 (serial) interface can also be accomplished via the interactive mode of a commercially available terminal emulation program, or by socket (port 7180) or COM port programming. When the Model 1700 receives a terminated ASCII string, it always sends back a reply as soon as the string is processed. *When sending commands to the Model 1700, you should wait for the reply from the Model 1700 before sending another command even if the reply consists of only termination characters.*

---

## SCPI COMMAND SUMMARY

---

The following manual conventions are used for SCPI (*Standard Commands for Programmable Instruments*) syntax for the remote interface commands:

- Braces { } enclose valid parameter choices.
- A vertical bar | separates multiple choices for each parameter.

- Triangle brackets < > indicate that you must supply a value.
- Parentheses ( ) within < > indicate alternative units are available.
- Capitalized portions of the commands indicate acceptable abbreviations.

For example, the command `CONFigure:RELAY1:CHannel {0|1|2}` indicates that the command `CONFigure:RELAY1:CHannel` has three parameter options: 0, 1, or 2, with 0 (disabled) being the default value.

The following section is a reference list of SCPI commands. Refer to the detailed description of each command for information regarding specific parameter choices and their meanings. Default settings are shown in bold.

**NOTE** The commands listed below are the standards-compliant form of the remote interface SCPI language. A set of legacy commands which do not conform to the SCPI standards are also available for backward compatibility with *existing installations* of the Model 18x series (see page 83).

### System-Related Commands

(see page 72 for more information)

\*IDN?  
\*RST  
\*TST?

N2?  
HE?

SERial\_NUMber?  
DATE\_MANUfacture?  
HardWare\_VERsion?  
FirmWare\_VERsion?  
Scpi\_VERsion?

SYStem:BEEPer:IMMediate <seconds>  
SYStem:BEEPer:STATe {0|1}  
SYStem:BEEPer:STATe?

SYStem:KLOCK {0|1}  
SYStem:KLOCK?

SYStem:DATE <yyyy>,<mm>,<dd>  
SYStem:DATE?

SYStem:TIME <hh>,<mm>,<ss>  
SYStem:TIME?

SYStem:REBOOT

IPV4addr?  
MACADDR?  
GATEWAY?  
HOSTname?

SCREENCAP  
SYStem:LOcal  
SYStem:REMOte

#### Display Configuration Commands and Queries

*(see page 75 for more information)*

DISPlay:N2?  
CONFigure:NAME:SENSor:N2 "<string>"  
NAME:SENSor:N2?

#### Relay Configuration Commands and Queries

*(see page 75 for more information)*

CONFigure:RELAy1:CHannel {0|1|2}  
RELAy1:CHannel?

CONFigure:RELAy1:SETpoint <level>  
RELAy1:SETpoint?

CONFigure:RELAy1:OPeration {0|1}  
RELAy1:OPeration?

CONFigure:RELAy2:CHannel {0|1|2}  
RELAy2:CHannel?

CONFigure:RELAy2:SETpoint <level>  
RELAy2:SETpoint?

CONFigure:RELAy2:OPeration {0|1}  
RELAy2:OPeration?

{RELAy1|RELAy2}:STATus?

**Alarm Configuration Commands and Queries**

*(see page 77 for more information)*

CONFigure:ALArm1:CHannel {0|1|2}  
ALArm1:CHannel?

CONFigure:ALArm1:SETPoint <level>  
ALArm1:SETPoint?

CONFigure:ALArm1:Operation {0|1}  
ALArm1:Operation?

CONFigure:ALArm2:CHannel {0|1|2}  
ALArm2:CHannel?

CONFigure:ALArm2:SETPoint <level>  
ALArm2:SETPoint?

CONFigure:ALArm2:Operation {0|1}  
ALArm2:Operation?

{ALArm1|ALArm2}:STATus?

ALArm:MUTE {0|1} *or* {NO|YES}  
ALArm:MUTE?

**Measurement Commands and Queries**

*(see page 79 for more information)*

MEASure:N2:LEVel?  
MEASure:N2:PERIod?

**Fill Control and Queries**

*(see page 79 for more information)*

CONFigure:FILL:CHannel {0|1|2}  
FILL:CHannel?

CONFigure:FILL:A <level>  
FILL:A?

CONFigure:FILL:B <level>

FILL:B?

CONFigure:INTerval:FILL <minutes>  
INTerval:FILL?

### Assignment Commands and Queries

*(see page 80 for more information)*

CONFigure:SOURCE:REC\_out {0|1|2}  
SOURCE:REC\_out?

CONFigure:SOURCE:CURrent\_LOOP {0|1|2}  
SOURCE:CURrent\_LOOP?

### N2 Channel Calibration Commands and Queries

*(see page 80 for more information)*

CONFigure:N2:LENGth <value>  
N2:LENGth?

MINCAL  
MINCAL?

MAXCAL  
MAXCAL?

NOSENSorCAL  
NOSENSorCAL?

APPROXMAXCAL <value>  
APPROXMAXCAL?

### Remote Units Commands and Queries

*(see page 81 for more information)*

CONFigure:N2:UNIT {0|1|2} or {PERCENT|INCH|CM}  
N2:UNIT?

---

## RS-232 CONFIGURATION

---

An RS-232 serial communication port is available as a 9-pin D-type connector on the rear panel of the instrument for serial communication function.

The Model 1700 uses the following *fixed* parameters related to the RS-232 interface:

- *Parity*: No Parity
- *Data Bits*: 8 Data Bits
- *Number of Start Bits*: 1 bit
- *Number of Stop Bits*: 1 bit
- *Flow Control*: None

The baud rate default is 115200, but is adjustable in the SYSTEM CONFIGURATION > RS232 SETUP menu (see page 41).

---

### SERIAL PORT CONNECTOR AND CABLING

A PC-compatible computer's serial port can be directly connected to the Model 1700 via a standard USB-to-serial cable. Refer to the computer's documentation to determine which ports are available on a computer and the required connector type.

The Model 1700 uses only three wires of the rear-panel DB9 connector: pin 2 (transmit), pin 3 (receive), and pin 7 (common). There are no software or hardware handshaking options. The Model 1700 is classified as a DCE (Data Communication Equipment) device since it transmits data on pin 3 and receives data on pin 2. The instrument to which the Model 1700 is attached must do the opposite, i.e., transmit on pin 2 and receive on pin 3 (the requirements for a DTE, or Data Terminal Equipment device).

---

### COMMAND/RETURN TERMINATION CHARACTERS

All commands are transmitted and received as ASCII values and are case insensitive. The Model 1700 can be configured to return <CR> (i.e. a *carriage return*) or <CR><LF> (i.e. a *carriage return* followed by a *linefeed*) at the end of a serial transmission (see page 41). The Model 1700 can accept <CR>, <LF>, <CR><LF>, or <LF><CR> as termination characters from an external computer.

---

## ETHERNET CONFIGURATION

---

The Model 1700 provides a 10/100Base-T Ethernet interface as a standard feature. It complies with the IEEE 802.3u 100Base-TX and 802.3 10Base-T standards.

The Model 1700 allows its IP address, subnet mask and gateway IP address to be assigned either statically or dynamically. To make these values static and assign them manually, set IP Address Assignment to Static (see page 40) and then set the values using the Edit parameter list. To enable the values to be dynamically assigned by a network DHCP server, set IP Address Assignment to DHCP.

The system name (also known as *host name* or *computer name*), can be set using the Network screen (see page 40) or via remote communications (either Ethernet or RS-232).

All network parameters (even those assigned by a DHCP server) can be viewed using the Network submenu (see page 40).

---

### ETHERNET CONNECTOR

The Model 1700 uses a standard RJ-45 jack for Ethernet communications. The Ethernet jack pinout is fully documented on page 96 in the *Appendix*.

---

### TERMINATION CHARACTERS

All commands and queries are transmitted and received as ASCII values and are case insensitive. The Model 1700 always transmits <CR><LF> (a *carriage return* followed by a *linefeed*) at the end of an Ethernet transmission. The Model 430 can accept <CR>, <LF>, <CR><LF>, or <LF><CR> as termination characters from an external computer.

---

### PORT ASSIGNMENT

The Model 1700 accepts remote connections to port 7180. Multiple connections to port 7180 are allowed.

---

## COMMAND REFERENCE

---

All commands sent to the Model 1700 are processed and the Model 1700 responds with a return value (if applicable) and termination. If the command is invalid, the Model 1700 will respond with an error code (see the *Error Codes* section). All return values including error codes are terminated with `<CR><LF>` (i.e. a *carriage return* followed by a *linefeed*). For those commands that do not return a value, the Model 1700 will return the `<CR><LF>` termination only.

The remote units settings are saved in non-volatile memory and are restored at power-up.

The Model 1700 instrument may be configured for reading liquid nitrogen, liquid helium, or both. Some commands will not be applicable if the instrument is not configured for certain level measurement.

---

### SYSTEM RELATED COMMANDS

- \*IDN?

Returns the identification string of the Model 1700. The identification string contains the manufacturer name, model number, serial number, and firmware revision code. Example output:

```
AMERICAN MAGNETICS INC.,MODEL 1700,1700-16-002,2.4.32015
```

- \*RST

Performs a Factory Restore if a restore file is available. All prior settings are lost!

- \*TST?

Returns a value incremented by "1" for each query to the requesting interface if unit is functioning. Return value does not indicate any operational status other than a functioning interface.

- N2?

Returns "0" if the instrument is not configured to read liquid nitrogen level, "1" if it is with the internal oscillator, and a "2" if it is with an external oscillator/transmitter.

- HE?

Returns "0" if the instrument is not configured to read liquid helium level.



- `SERial_NUMber?`

Returns the serial number of the instrument as a string, for example:  
1700-16-002.

- `DATE_MANUfacture?`

Returns the date of manufacture as a string in the form *yyyy-mm-dd*, for example: 2016-11-30.

- `HardWare_VERsion?`

Returns the hardware version of the instrument.

- `FirmWare_VERsion?`

Returns the firmware version of the instrument, for example:  
2.4.32015.

- `Scpi_VERsion?`

Returns the SCPI command version of the instrument.

- `SYStem:BEEPer:IMMediate <seconds>`

The receipt of this command with a valid argument causes an audible tone to be generated by the instrument. The duration time parameter is specified in seconds.

- `SYStem:BEEPer:STATe {0|1}`

Enables/disables the beeper. When `STATe "0"` (OFF) is selected, no instrument condition, except the `SYStem:BEEPer:IMMediate` command, shall cause an audible beep to be emitted. When the instrument restarts, this value is reset to a default of "1" (ON).

- `SYStem:BEEPer:STATe?`

Returns present beeper state. A value of "0" indicates the beeper is disabled. A value of "1" indicates it is enabled.

- `SYStem:KLOCK {0|1}`

This command locks the local controls of an instrument if set to a value of "1" (ON). This includes all front panel, keyboard, or other local interfaces. The state of the lock defaults to "0" (OFF) when the instrument is restarted.

- SYStem:KLOCK?

Returns the state of the instrument local controls lock. A return value of “0” indicates the lock is OFF. A return value of “1” indicates the local interfaces, such as the touch display, of the instrument are locked.

- SYStem:DATE <yyyy>,<mm>,<dd>

Sets the date per <yyyy> as a four-digit number; <mm> month with range of 1 to 12 inclusive; <dd> day of the month. Do not include the <> characters in the command.

- SYStem:DATE?

Returns the system date if the form <yyyy>,<mm>,<dd>.

- SYStem:TIME <hh>,<mm>,<ss>

This command is used to set the instrument's clock:

<hh> Range of hours is 0 to 23 inclusive.

<mm> Range of minutes is 0 to 59 inclusive.

<ss> Range of seconds is 0 to 60.

Do not include the <> characters in the command.

- SYStem:TIME?

Returns the system time if the form <hh>,<mm>,<ss>.

- SYStem:REBOOT

Reboots the instrument.

- SYStem:RESTORE

Reboots the instrument and sets all parameters back to factory defaults. All prior settings are lost!

- IPV4ADDR?

Returns the presently assigned TCP/IP version 4 address in the form:  
xxx.xxx.xxx.xxx

- MACADDR?

Returns the MAC address of the network interface.

- GATEWAY?

Returns the Gateway address of the network interface.

- HOSTname?

Returns the hostname of the network interface.

- SCREENCAP

Snaps a TGA-formatted (TARGA) image of the current local display screen of the instrument and places it in the “Log files” folder of the instrument’s network share.

- SYSTem:LOCal

Enables the front panel touchscreen. All front panels controls are enabled by default after a power-up or REBOOT command.

- SYSTem:REMOte

Disables front panel touchscreen.

---

## DISPLAY CONFIGURATION COMMANDS AND QUERIES

- DISPLAY:N2?

Returns a “0” if the instrument is not configured to display liquid nitrogen level on the home screen and a “1” if it is.

- CONFigure:NAME:SENSor:N2 "<string>"

Sets the displayed name of the nitrogen level sensor.

- NAME:SENSor:N2?

Returns the displayed name of the nitrogen level sensor.

- DISPLAY:HE?

Returns a “0” if the instrument is not configured to display liquid helium level on the home screen and a “1” if it is.

---

## RELAY CONFIGURATION COMMANDS AND QUERIES

- CONFigure:RELAy1:CHannel {0|1|2}

Assigns relay №1 as disabled (0), or to nitrogen (1) or helium (2).

- RELay1:CHannel?

Returns a "0" if relay №1 is disabled, a "1" if the relay is assigned to the nitrogen channel, and a "2" if the relay is assigned to the helium channel. An error return code is generated if attempting to assign the relay to a non-existent measurement channel in the present instrument configuration.

- CONFigure:RELay1:SETpoint <level>

Configures the relay №1 trip setpoint in the currently assigned channel's units.

- RELay1:SETpoint?

Returns the relay №1 setpoint in the current units.

- CONFigure:RELay1:OPeration {0|1}

Configures relay №1 such that it closes when the level is  $\leq$  the setpoint (0) or  $\geq$  the setpoint (1).

- RELay1:OPeration?

Returns a "0" if relay №1 closes when the level is  $\leq$  the setpoint and a "1" if the relay closes when the level is  $\geq$  the setpoint.

- CONFigure:RELay2:CHannel {0|1|2}

Assigns relay №2 to disabled (0), or to nitrogen (1) or helium (2). An error return code is generated if attempting to assign the relay to a non-existent measurement channel in the present instrument configuration.

- RELay2:CHannel?

Returns a "0" if relay №2 is disabled, a "1" if the relay is assigned to the nitrogen channel, and a "2" if the relay is assigned to the helium channel.

- CONFigure:RELay2:SETpoint <level>

Configures the relay №2 trip setpoint in the currently assigned channel's units.

- RELay2:SETpoint?

Returns the Relay №2 setpoint in the current units.

- `CONFigure:RELAy2:OPeration {0|1}`

Configures relay №2 such that it closes when the level is  $\leq$  the setpoint (0) or  $\geq$  the setpoint (1)

- `RELAy2:OPeration?`

Returns a “0” if relay №2 closes when the level is  $\leq$  the setpoint and a “1” if the relay closes when the level is  $\geq$  the setpoint.

- `{RELAy1|RELAy2}:STATus?`

Returns the present status of the specified relay. A return value of “0” indicates the relay is not closed. A value of “1” indicates the relay is closed.

---

## ALARM CONFIGURATION COMMANDS AND QUERIES

- `CONFigure:ALArml:CHannel {0|1|2}`

Assigns alarm №1 as disabled (0), or to nitrogen (1) or helium (2).

- `ALArml:CHannel?`

Returns a “0” if alarm №1 is disabled, a “1” if the alarm is assigned to the nitrogen channel, and a “2” if the alarm is assigned to the helium channel. An error return code is generated if attempting to assign the alarm to a non-existent measurement channel in the present instrument configuration.

- `CONFigure:ALArml:SETpoint <level>`

Configures the alarm №1 trip setpoint in the currently assigned channel's units.

- `ALArml:SETpoint?`

Returns the alarm №1 setpoint in the current units.

- `CONFigure:ALArml:OPeration {0|1}`

Configures alarm №1 such that it alarms when the level is  $\leq$  the setpoint (0) or  $\geq$  the setpoint (1).

- `ALArml:OPeration?`

Returns a “0” if alarm №1 alarms when the level is  $\leq$  the setpoint and a “1” if it alarms when the level is  $\geq$  the setpoint. Alarm №1 is considered as the “legacy” HI level alarm when the alarm condition is set to level  $\geq$  a setpoint.

- `CONFigure:ALArm2:CHannel {0|1|2}`

Assigns alarm №2 to disabled (0), or to nitrogen (1) or helium (2). An error return code is generated if attempting to assign the alarm to a non-existent measurement channel in the present instrument configuration.

- `ALArm2:CHannel?`

Returns a "0" if alarm №2 is disabled, a "1" if the alarm is assigned to the nitrogen channel, and a "2" if the alarm is assigned to the helium channel.

- `CONFigure:ALArm2:SEtpoint <level>`

Configures the alarm №2 trip setpoint in the currently assigned channel's units.

- `ALArm2:SEtpoint?`

Returns the alarm №2 setpoint in the current units.

- `CONFigure:ALArm2:OPeration {0|1}`

Configures alarm №2 such that it alarms when the level is  $\leq$  the setpoint (0) or  $\geq$  the setpoint (1)

- `ALArm2:OPeration?`

Returns a "0" if alarm №2 alarms when the level is  $\leq$  the setpoint and a "1" if it alarms when the level is  $\geq$  the setpoint. Alarm №2 is considered as the "legacy" LO level alarm when the alarm condition is set to level  $\leq$  a setpoint.

- `{ALArm1|ALArm2}:STATus?`

Returns the present status of the specified alarm. A return value of "0" indicates the alarm is inactive. A value of "1" indicates the alarm is active.

- `ALARm:MUTE {0|1} or {NO|YES}`

Mutes the audible alarm for any active alarm(s) when set to "1" or "YES". Unmutes the audible alarm if set to "0" or "NO".

- `ALARm:MUTE?`

Returns "0" if audible alarm is muted. Returns "1" if not muted.

## MEASUREMENT COMMANDS AND QUERIES

- `MEASure:N2:LEVel?`

Returns the liquid nitrogen level in the current units.

- `MEASure:N2:PERIod?`

Returns the liquid nitrogen level measurement period in microseconds.

## FILL CONTROL AND QUERIES

- `CONFigure:FILL:CHannel {0|1|2}`

Assigns the auto fill control relay to either no channel (disabled) (0), nitrogen (1), or helium (2). An error return code is generated if attempting to assign the relay to a non-existent measurement channel in the present instrument configuration.

- `FILL:CHannel?`

Returns a "0" if the auto fill relay is disabled, a "1" if the relay is assigned to the nitrogen channel, and a "2" if the relay is assigned to the helium channel.

- `CONFigure:FILL:A <level>`

Sets the A setpoint (control band upper limit) in the assigned channel's current units. The A setpoint (autofill stop) must be greater than the B setpoint (fill start) and must also be between 0% and 100%.

- `FILL:A?`

Returns the A setpoint (control band upper limit) in the assigned channel's current units.

- `CONFigure:FILL:B <level>`

Sets the B setpoint (control band lower limit) in the assigned channel's current units. The B setpoint (autofill start) must be less than the A setpoint (fill stop) and must also be between 0% and 100%.

- `FILL:B?`

Returns the B setpoint (control band lower limit) in the assigned channel's current units.

- `CONFigure:INTerval:FILL <minutes>`

Sets the fill timer in minutes. If the level does not reach or exceed the A setpoint within the fill time, the fill is terminated. Setting the value of FILL to "0" disables the fill timer function.

- `INTerval:FILL?`

Returns the fill timer setting in minutes.

---

## ASSIGNMENT COMMANDS AND QUERIES

- `CONFigure:SOURCE:REC_out {0|1|2}`

Configures the 0-10 V<sub>DC</sub> Recorder Output source to disabled (0), assigned to the nitrogen channel (1), or the helium channel (2). An error return code is generated if attempting to assign the output to a non-existent measurement channel in the present instrument configuration.

- `SOURCE:REC_out?`

Returns a “0” if the 0-10 V<sub>DC</sub> Recorder Output is disabled, a “1” if it is configured for the nitrogen channel, and a “2” if it is configured for the helium channel.

- `CONFigure:SOURCE:CURrent_LOOP {0|1|2}`

Configures the 4-20 mA Current Loop output source to disabled (0), assigned to the nitrogen channel (1), or the helium channel (2). An error return code is generated if attempting to assign the output to a non-existent measurement channel in the present instrument configuration.

- `SOURCE:CURrent_LOOP?`

Returns a “0” if the 4-20 mA Current Loop output is disabled, a “1” if it is configured for the nitrogen channel, and a “2” if it is configured for the helium channel.

---

## N2 CHANNEL CALIBRATION COMMANDS AND QUERIES

- `CONFigure:N2:LENGth <value>`

Configures the liquid nitrogen sensor active length in current units. Returns an error code if the current units are percent.

- `N2:LENGth?`

Returns the configured active length of the liquid nitrogen sensor in current units. Returns an error code if the current units are percent.

- `MINCAL`

Performs and saves a **MIN** calibration point calibration for the liquid nitrogen sensor.



- MINCAL?

Returns the last saved **MIN** calibration point in microseconds for the liquid nitrogen sensor.

- MAXCAL

Performs and saves a **MAX** calibration point calibration for the liquid nitrogen sensor.

- MAXCAL?

Returns the last saved **MAX** calibration point in microseconds for the liquid nitrogen sensor.

- NOSENSOrCAL

Calibrates and saves the **LOSS OF SENSOR** condition for the liquid nitrogen sensor.

- NOSENSOrCAL?

Returns the last saved **LOSS OF SENSOR** calibration point in microseconds for the liquid nitrogen sensor.

- APPROXMAXCAL <value>

Sets the **Approx Cal (Max) Multiplier** factor (see page 62) to a decimal value. The value is nominally 1.000 for no modification to the MAX calibration point.

- APPROXMAXCAL?

Returns the **Approx Cal (Max) Multiplier** factor (see page 62) as a decimal value.

---

## REMOTE UNITS COMMANDS AND QUERIES

- CONFigure:N2:UNIT {0|1|2} or {PERCENT|INCH|CM}

Sets the liquid nitrogen level units of measurement to percent (0 or PERCENT), inches (1 or INCH), or centimeters (2 or CM). The default is PERCENT units.

- N2:UNIT?

Returns the current liquid helium level units in use as either C, I, or %.

## ERROR CODES

The Model 1700 returns specific error codes for invalid commands and/or arguments. If an error condition is returned, the command is not processed and the configuration of the instrument is not modified. The table below provides a list of error codes, their meaning, and any associated limits.

Error Code	Meaning	Valid Range
-1	LO (or relay №2 <sup>a</sup> ) setpoint out of range	$0 \leq \text{LO (or relay №2}^a) \leq \text{LENGTH}$
-2	Fill B setpoint (fill start) out of range	$0 \leq B < A$
-3	Fill A setpoint (fill stop) out of range	$B < A \leq \text{LENGTH}$
-4	HI (or relay №1 <sup>a</sup> ) setpoint out of range	$0 \leq \text{HI (or relay №1}^a) \leq \text{LENGTH}$
-5	Attempted to set or query for LENGTH in PERCENT units mode	
-6	Invalid argument, value out of maximum calibration range	$1 \text{ cm} \leq \text{value} \leq 650 \text{ cm (LN2)}$
-7	INTERVAL setting out of range	$0 \leq \text{INTERVAL} \leq 999 \text{ min}^b \text{ (or)}$ $0 \leq \text{INTERVAL} \leq 99999 \text{ min}^c$
-8	Unrecognized command	
-9	Invalid argument, value was negative or non-numeric	
-10	Approximate calibration <i>factor</i> out of range	$0.1 \leq \text{factor} \leq 999.9$
-11	Command exceeds SCPI input buffer limit	256 characters, including spaces, etc.
-12	Command invalid for selected channel or interface in present configuration	

a. Applies to dual instrument configuration

b. For Version HW = 4712

c. For Version HW = 4900

## LEGACY COMMAND REFERENCE

The legacy command set is included for compatibility with existing Model 18x installations. New installations should consider using the modern command set that conforms more closely with the SCPI specification. The modern and legacy command sets can be mixed as desired.

All legacy command or queries sent to the Model 1700 are processed and the Model 1700 responds with a return value (if applicable) and termination. If the command is invalid, the Model 1700 will respond with an error code (see the *Error Codes* section). All return values including error codes are terminated with `<CR><LF>` (i.e. a *carriage return* followed by a *linefeed*). For those commands that do not return a value, the Model 1700 will return the `<CR><LF>` termination only.

### COMMANDS FOR CONTROLLING THE UNITS OF MEASUREMENT

Command:	CM	Function:	Sets the units of measurement to centimeters	Returns:	<code>&lt;CR&gt;&lt;LF&gt;</code>
Command:	INCH	Function:	Sets the units of measurement to inches	Returns:	<code>&lt;CR&gt;&lt;LF&gt;</code>
Command:	PERCENT	Function:	Sets the measurement to % of sensor length	Returns:	<code>&lt;CR&gt;&lt;LF&gt;</code>
Command:	UNIT	Function:	Returns the current units in use	Returns:	C, I, or % <code>&lt;CR&gt;&lt;LF&gt;</code>

The CM command sets the units of measurement to centimeters and the INCH command selects inches. The PERCENT command sets the units of measurement to the percentage of active sensor length that is immersed in liquid. **The units of measurement selected through the remote interface are controlled independently from the units shown on the front panel display.**

The remote units setting is *automatically* saved in permanent memory and is restored at power-up. The UNIT command returns a one character value (and termination) indicating the current units — C for centimeters, I for inches, or % for percentage.

## COMMANDS FOR CONFIGURING PERMANENT MEMORY

Command:	HI=<value>	Function:	Configures the HI setpoint limit	Returns:	<CR><LF>
Command:	LO=<value>	Function:	Configures the LO setpoint limit	Returns:	<CR><LF>
Command:	A=<value>	Function:	Configures the A setpoint (control band upper limit)	Returns:	<CR><LF>
Command:	B=<value>	Function:	Configures the B setpoint (control band lower limit)	Returns:	<CR><LF>
Command:	INTERVAL=<value>	Function:	Configures the fill timer in minutes	Returns:	<CR><LF>
Command:	LENGTH=<value>	Function:	Configures the active sensor length	Returns:	<CR><LF>
Command:	SAVE	Function:	None (for backward compatibility only).	Returns:	<CR><LF>

The HI and LO command configure the high and low setpoint limit values, respectively. For example, HI=90.0 would configure the high setpoint limit to 90.0 in whichever units of measurement last selected through the serial interface. The A and B commands configure the upper limit and lower limit of the control band, respectively. The HI, LO, A, and B commands are compatible with the percent units selection.

**NOTE** The HI setting is associated with hardware relay №1. The LO setting is associated with the hardware relay №2. This association is fixed and cannot be changed.

The LENGTH command configures the active sensor length setting in the current units. LENGTH=35.0 would configure the active sensor length to 35.0 units of centimeters or inches.

**NOTE** The LENGTH=<value> command will only function if CM or INCH are currently selected as the units of measurement. The LENGTH command does not configure the Model 1700 if the units of measurement are PERCENT.

The INTERVAL command sets the nitrogen fill timer in minutes as described in the *Operation* section on page 36. Setting the value of INTERVAL to 0 disables the fill timer function.

The HI, LO, A, B, INTERVAL, LENGTH, and current remote units settings are automatically saved to permanent memory. Saved settings are then recalled each time the power is turned off and then reapplied to the instrument on subsequent power on.

## QUERYING THE CONFIGURATION

Command:	HI	Function:	Returns the HI setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	LO	Function:	Returns the LO setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	A	Function:	Returns the A setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	B	Function:	Returns the B setpoint limit in the current units	Returns:	<value> <CR><LF>
Command:	INTERVAL	Function:	Returns the fill timer setting in minutes	Returns:	<value> <CR><LF>
Command:	LENGTH	Function:	Returns the active sensor length in the current units	Returns:	<value> <CR><LF>

The HI, LO, A, B, INTERVAL, and LENGTH queries return the current configuration of the instrument. Each return value is terminated with <CR><LF>.

## RETURNING A LEVEL MEASUREMENT

Command:	LEVEL	Function:	Returns the liquid level in the current units	Returns:	<value> <CR><LF>
----------	-------	-----------	---	----------	---------------------

The LEVEL query returns the liquid level in the current units selected through the communication interface.



# SERVICE AND REPAIR

## CLEANING

To prevent electrical shock, disconnect the instrument from AC mains power and disconnect all connected wiring before cleaning. Clean the outside of the instrument using a soft, lint-free, cloth slightly dampened with water.

Do not use detergent or solvents.

Do not attempt internal cleaning.

## USER REPLACEABLE PARTS

Replacement parts for the instrument are listed in the table below.

AMI Part Number	Description
HG0128	Instrument foot
SA 1045	Single Rack Mount Kit
SA 1046	Dual Rack Mount Kit
EF1700	Fuse, 3 A, 250 Vac, 5x20 mm, fast acting, UL/CSA recognized.
HG0005	Battery, 3V lithium, 20mm x 3.2 mm coin cell; CR2032.

---

## BATTERY REPLACEMENT

---

This section describes the procedure for replacing the battery on the instrument's main circuit board.

**WARNING** This procedure should only be performed by a technician who is familiar with electronic instrumentation and trained in electrical safety and ESD precautions. Always disconnect the power cord and any external wiring before removing the instrument cover.


Always disconnect all inputs, cords, and cables before disassembling the instrument.

---

### LOW BATTERY INDICATION


Battery status is displayed in the **SYSTEM DATE & TIME** dialog. To access the dialog:

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen


2. From the **MENU** screen, choose **SYSTEM**.



MENU Selection Screen




3. From the **SYSTEM CONFIGURATION** screen choose **SET TIME**.



SYSTEM CONFIGURATION Screen

4. The battery status is displayed in the **SYSTEM DATE & TIME** screen. If the battery is weak and should be replaced soon, the message will appear as shown at right.



Weak Battery Indication

5. If the battery is expired, the message will appear as shown at right.

**NOTE** The NOTIFY button state will cause a warning message to appear in the the HOME screen footer. NO NOTIFY supresses the HOME screen warning.



Expired Battery Indication

The following sections detail how to replace the battery.

## TOOLS REQUIRED

- Torx Plus (T.M.) size 10 driver (Wera 028034 or equivalent)
- Torx Plus (T.M.) size 15 driver (Wera 028035 or equivalent)
- Small, flat-blade screw driver (for prying)

## PROCEDURE

1. Unplug the instrument from the AC power source.
2. Using the T-15 driver, remove the four 8-32 machine screws on the sides of the instrument cover. Set these screws aside as they will be re-used.
3. Using the T-10 driver, remove the four 6-32 machine screws on the rear of the instrument cover. Set these screws aside as they will be re-used.
4. Lift the instrument cover off of the instrument chassis and set aside.

5. Using the small, flat-blade screwdriver, carefully pry the battery from the holder BH1.
6. Install the new battery into the battery holder BH1.
7. Replace the top cover and secure using the eight machine screws which were removed previously.

---

## FUSE REPLACEMENT

---

This section describes the procedure for replacing the two fuses on the instrument's main circuit board.

**WARNING** This procedure should only be performed by a technician who is familiar with electronic instrumentation and trained in electrical safety and ESD precautions. Always disconnect the power cord and any external wiring before removing the instrument cover.

Always disconnect all inputs, cords, and cables before disassembling the instrument.

---

### TOOLS REQUIRED

- Torx Plus (T.M.) size 10 driver (Wera 028034 or equivalent)
- Torx Plus (T.M.) size 15 driver (Wera 028035 or equivalent)
- Small, flat-blade screw driver (for prying)

---


### PROCEDURE

1. Unplug the instrument from the AC power source.
2. Using the T-15 driver, remove the four 8-32 machine screws on the sides of the instrument cover. Set these screws aside as they will be re-used.
3. Using the T-10 driver, remove the four 6-32 machine screws on the rear of the instrument cover. Set these screws aside as they will be re-used.
4. Lift the instrument cover off of the instrument chassis and set aside.
5. Using the small, flat-blade screwdriver, carefully pry the fuse(s) from the fuse holders F1 and/or F2.
6. Install the new fuse(s) into the fuse holder(s) F1 and/or F2.
7. Replace the top cover and secure using the eight machine screws which were removed previously.

## FIRMWARE UPGRADE VIA ETHERNET

The following procedure can be utilized by users to upgrade the Model 1700 by downloading firmware upgrade archives (.tar.gz) from the AMI support website.

1. Connect the instrument to a computer network via Ethernet cable and energize the instrument.
2. When the instrument has completed booting, note the IP Address: **MENU > NETWORK > Address**. There should be a message in green on the screen indicating "Network is up".




Network Configuration Screen


3. Determine the currently installed version of firmware. From the home page, press **MENU > SYSTEM > FW Version**.
4. Download the available new version of firmware and place on a PC on the same network as the instrument.

**NOTE** The following steps assume the computer is Windows-based.

5. Start Windows Explorer and enter the instrument's IP address noted from Step 2 above in the address field. Be sure to use the format such as \\192.168.1.101. Ensure the Explorer screen looks similar to the following.



- Open a second Windows Explorer and locate the previously downloaded firmware file:




- Drag the firmware file to the instrument's Application Update folder in the first Explorer screen.
- Install the new firmware by pressing **MENU > SYSTEM > UPDATE** on the instrument and choose the name of the firmware file copied to the instrument. It is likely the top choice of the files listed.
- Press **INSTALL** on the instrument when prompted.
- The instrument will install the new firmware and reboot when completed.

## UPGRADE VIA SCP

If the instrument's file shares are not visible on the local network, it is possible to copy the firmware upgrade directly to the instrument using the `root` login and password using SCP (Secure Copy) over Ethernet.

- Connect the instrument to the local computer network via Ethernet cable and energize the instrument.
- When the instrument has completed booting, note the IP Address: **MENU > NETWORK > Address**. There should be a message in green on the screen indicating "Network is up".



Network Configuration Screen

- Determine the currently installed version of firmware. From the home page, press **MENU > SYSTEM > FW Version**.
- Download the available new version of firmware and place on a Windows, Linux, or Mac computer on the same network as the instrument.

5. Open a command line application and set the current directory to the folder where the newly download firmware from the prior step resides.
6. Enter the following command where the example IP address of 192.168.1.101 is replaced by the IP of the instrument noted in step 2 above, and the example firmware file of Ver-2.22719.15.tar.gz is replaced by the actual file downloaded from AMI:

```
scp ./Ver-2.220719.15.tar.gz root@192.168.1.101:/update
```

Note that you will be prompted to enter the root account password to complete the file transfer. If you do not know the password, contact an Authorized AMI Technical Support Representative for assistance.


7. Install the new firmware by pressing **MENU > SYSTEM > UPDATE** on the instrument and choose the name of the firmware file copied to the instrument via SCP. It is likely the top choice of the files listed.
8. Press **INSTALL** on the instrument when prompted.

The instrument will install the new firmware and reboot when completed.



# APPENDIX

## SERIAL (RS-232) CONNECTOR




### Serial (RS-232) Pin Out

The RS-232 connector is a 9-pin D-sub female connector to connect standard DTE 9-pin D-sub male connector using a standard straight (not NULL) cable.

### Serial (RS-232) Pin Definitions

Pin	Mnemonic	Function
1	N/C	
2	TXD	Transmit Data
3	RXD	Receive Data
4	N/C	
5	GND	Signal Ground
6	N/C	
7	N/C	
8	N/C	
9	N/C	

## ETHERNET CONNECTOR



**Ethernet Connector Socket Pin out**

### **Ethernet RJ-45 Connector Pin Definitions**

Pin	Mnemonic	Function
1	TXD+	Transmit differential output +
2	TXD-	Transmit differential output -
3	RXD+	Transmit differential input +
4	not used	
5		
6	RXD	Transmit differential input -
7	not used	
8		




## AUX I/O CONNECTOR



**Aux I/O Connector**

The AUX I/O connector is a 9-pin D-sub female connector.

**Aux I/O Pin Definitions**

Pin	Function	Polarity
1	4-20 mA Current Loop	+
2		—
3	0-10 VDC Output	+
4		—
5	Relay № 1 Dry Contact	N/A
6		
7	Relay № 2 Dry Contact	N/A
8		
9	N/A	

---

## DIELECTRIC CONSTANTS FOR CRYOGENIC LIQUIDS

---

**Common Cryogenic Liquid Dielectric Constants<sup>a</sup>**

Liquid	Dielectric Constant
Argon (Ar)	1.5034
Carbon Dioxide (CO <sub>2</sub> )	2.5911
Hydrogen (H <sub>2</sub> )	1.2307
Methane (CH <sub>4</sub> )	1.6257
Nitrogen (N <sub>2</sub> )	1.4327
Propane (C <sub>3</sub> H <sub>8</sub> )	1.8029
Oxygen (O <sub>2</sub> )	1.4872

a. NIST RefProp version 9.1

---

## TROUBLESHOOTING

---

The following paragraphs serve as an aid to assist the user in troubleshooting a potential problem with the Model 1700 Instrument. If the user is not comfortable in troubleshooting the system, contact an AMI Technical Support.

If the cause of the problem cannot be located, contact an AMI Technical Support Representative at +1 (865) 482-1056 for assistance. The AMI technical support group may also be reached by internet e-mail at:

**support@americanmagnetics.com**

---

### INSTRUMENT DISPLAYS “LOSS OF SENSOR” CONDITION FOR LN<sub>2</sub> LEVEL MEASUREMENT

This message indicates one of two things: the signal from the external oscillator/transmitter has been lost (disconnection) or the capacitance measured by the instrument is significantly less than the minimum capacitance expected during level measurement<sup>1</sup>. Ensure the instrument has been calibrated (especially “Performing Loss of Sensor Calibration” on page 53) and is properly connected via coaxial cable(s) to the BNC connector on the back of the instrument.

---

### INSTRUMENT DISPLAYS “SENSOR SHORTED” CONDITION FOR LN<sub>2</sub> LEVEL MEASUREMENT

This message indicates that the sensor oscillator circuit has ceased oscillating. This is usually caused by a contaminated level sensor where a conductive substance (usually water) is shorting out the sensor tubes.<sup>2</sup> Dry or clean out the sensor and this condition is usually removed. Note that the sensor resistance measured across the BNC connector should be  $>10^7$  ohms.<sup>3</sup>

---

### SYSTEM TEST SCREENS

The instrument provides various system test functions that can be performed from the front panel touchscreen. The system test functions are not available to the web browser interface nor the SCPI remote interface over serial or Ethernet connections. To access the test functions:


- 
1. The threshold is the midpoint between the MINimum CALibration and the NO SENSOR CALibration.
  2. The threshold is approximately 200 kohms with a sensor capacitance of 2 nF.
  3. Any resistance between the sensor elements will cause the instrument to read higher than actual level.

1. From the home screen, choose **MENU**.



Menu Selection From Home Screen

2. From the **MENU** screen, choose **SYSTEM**.




MENU Selection Screen

3. From the **SYSTEM CONFIGURATION** screen choose **SYSTEM TEST**.



SYSTEM CONFIGURATION Screen


4. In the **SYSTEM TEST** screen, several selections are available to test the touchscreen colors. Pressing a pixel test button will cause the screen to show a scrolling test pattern. When touched again, the display will then return to normal.



SYSTEM TEST Pixel Tests


5. To calibrate or test the touch function of the screen, use the **TOUCH SCREEN CALIBRATE** or **TOUCH SCREEN TEST** selections in the **SYSTEM TEST** screen.

The calibrate function will display a screen with touch targets. Touch all the targets as they are displayed to complete the calibration. The instrument will then automatically store the new touch calibration.




SYSTEM TEST Touch Cal/Test

6. The **HARDWARE CONTROL** selection in the **SYSTEM TEST** screen allows the user to observe in real time and/or manually test various functions of the instrument hardware. Press the **HARDWARE CONTROL** button to move to the detailed test screen.



Hardware Test Selection

7. The **HARDWARE CONTROL** screen includes several controls to manually exercise the hardware. However, to do so the operator must press the Control button to toggle it from **AUTO** to **MANUAL** mode.




Hardware Control Screen

8. Once in **MANUAL** mode, the user can touch any of the buttons to toggle the state of the hardware. The **Analog Outputs** button is, however, not a toggle but will move the operator to the **ANALOG OUTPUTS TEST** screen.

**NOTE** In the **AUTO** mode, the Test screen displays the real time measurement or state of each labeled function.

**NOTE** When changing to **MANUAL** test mode, please be aware that the automatic functions of the instrument will be overridden. Functions such as autofill will be interrupted, or a fill initiated, by the manual actions. Be sure the system is in a *safe condition* before exercising manual control.


9. The **ANALOG OUTPUTS TEST** screen requires the operator press the Control button to toggle from **AUTO** to the **MANUAL** mode to exercise control over the output values. The output values range from 0 to 4095 where 0 is the minimum output value and 4095 is the maximum output value. Tap a value to enter a manual override.



Analog Outputs Tests Screen

**NOTE** The hardware control state will return to AUTO after a few seconds if the operator does not explicitly return to automatic control. However, the state of the hardware or autofill functions are *not guaranteed to be the same* as when the MANUAL control state was entered.

10. Finally, the **CURSOR ON/OFF** button in the **SYSTEM TEST** screen allows the operator to turn the cursor off and on. The cursor display can be useful for verifying the touchscreen response if it appears to function incorrectly.



CURSOR Display Selection


11. Use the **HOME** button to return to the main level display.

## SYSTEM LOGS

The Model 1700 provides a logging feature that saves detailed information regarding the liquid level and system events. The information can be very useful in troubleshooting suspected system errors.

The Model 17XX family of liquid level instruments keeps log files in a local Ethernet network accessible folder via Samba file shares. The log

folder is labeled “Log files” and can found on a network by opening a file explorer at the instrument IP address or hostname.



In the Log files folder, there will be a variety of logging files:

**He.log:** Contains a log of liquid helium measurement events, with each entry timestamped and a record of the hardware bits at that time.

**N2-1.log:** Contains a log of liquid nitrogen measurement events, with each entry timestamped and a record of the hardware bits at that time.

**operations.csv:** Contains a log of all operations/events, with each entry timestamped and appropriate data for the event in a comma separated file format.

## LOG FILE FORMATS

The text based logging format is in a comma separated value format. It consist of a unix timestamp, followed by the current level, and then the status bits encoded in a hexadecimal number.

The status bits contains information for both the channels, even if the channels are disabled.

The status bits are a summation of the hardware flags, each flag has a unique bit position, and so the status byte displayed is a hexadecimal representation of those bits.

Refer to the Status Bits table on page 108 for the bit values assigned for each hardware flag.

A log entry is made when either the level changes by 0.1% or if the status bits change from last entry. Levels and bits are checked every second.



A Unix timestamp is the time in seconds since Jan 1, 1970 in Coordinated Universal Time.

### DESCRIPTION OF THE HE.LOG FORMAT

Here is an example of a measurement, with the helium level changing from 26.3 to 45.1:

```

      timestamp, %level, status bits
#1  1626706702,  26.3,028
#2  1626706703,  45.1,038

```

**Entry #1:** timestamped at Mon Jul 19 2021 10:58:22 (Eastern Daylight Time) with a previously read level of 26.3. The status bits (028) indicated that the external/internal N2 sensor oscillator selection was internal, and that the helium sensor was active and dirty sensor mode was off. The level updated when the helium measurement cycle completed, as recorded in the next entry.

**NOTE** Note that even if an instrument is not configured for nitrogen measuring, there will be inactive status bits for that nitrogen channel in the status word.

**Entry #2:** timestamped at Mon Jul 19 2021 10:58:22 with a last read level of 45.1. The status bits (028) indicated that the external/internal N2 sensor oscillator selection was internal, and that the helium sensor was active and dirty sensor mode was off. The level had changed between entry #2 and #3, so when the unit stopped measuring (turning off the helium sensor) it updated the level, which changed to 45.1

Refer to the Status Bits table on page 108 for the bit values assigned for each hardware flag.

### DESCRIPTION OF THE N2-1.LOG FORMAT

Here is an example of measurements of changing level, with a loss of sensor event for the liquid nitrogen measurements:

```

      timestamp, %level, status bits
#1  1626708988,  88.2,000020
#2  1626708989,  77.8,000020
#3  1626708990,  54.1,000020
#4  1626708991,  41.3,000020
#5  1626708992,  34.4,000020
#6  1626708994,  16.1,000020
#7  1626708996,   0.0,000221

```

#8 1626708999, 16.1,000020

**Entry #1:** timestamped at Mon Jul 19 2021 11:36:28 (Eastern Daylight Time) with a read level of 88.2%. The status bits (000020) indicated that the external/internal N2 sensor oscillator selection was internal.

**Entries #2-#6:** shows the level changing over a 5 second span, settling at 16.1%

**Entry #7:** timestamped at Mon Jul 19 2021 11:36:36 showed a loss of sensor event. The level will be 0.0%, and the status bits show that the sensor was disconnected and the alarm bit was set.

**Entry #8:** timestamped at Mon Jul 19 2021 11:36:36 showed back to normal operation, with a reading of 16.1%.

Refer to the Log Status Bits table on page 108 for information on the status bits.

### DESCRIPTION OF THE OPERATIONS LOGGING FILE

The `operations.csv` file is a text file that records a log of all operations/events: remote, local, and SCPI.

This file contains operations log entries in a comma separated value format, as follows:

```
<unix timestamp>, <class>, <data field 1>, <data field 2>, <data field 3>
```

A data field can be empty, that is, no text.

An overview of the classes:

Class	Description of Event
AF	Master System Autofill
AL	Alarm(s)
BT	Battery Status
CE	System commands
N1	Liquid Nitrogen channel sensor
LE	GUI page load event
LF	GUI link (AML function link)
LL	GUI link (AML page link)
LS	GUI link (system variable)
RY	Hardware Relay Status
SC	SCPI command that modified the instrument
SF	State file operation
PU	Power Up
TG	GUI toggle

VM            Modification of a instrument variable

More detailed information on each class is available from your Authorized AMI Technical Representative.

Here is a simple example of operations logging file entries for a shutdown initiation:

```
      timestamp, class,data1,data2,data3
#1 1626712173,LL,MENU,menu.aml,
#2 1626712191,LL,SHUTDOWN,../shutdown.aml,
#3 1626712200,LL,SHUTDOWN ,../shutdown.aml,
#4 1626712200,SF,shmem_dump,checksum,6A7E48E4
#5 1626712200,CE,shutdown,,
```

**Line by line description:**

```
1626712173,LL,MENU,menu.aml,
```

The operator pressed the MENU button on the home page.

Menu.aml indicated the screen page loaded was menu.aml, that is, the MENU page.

The class LL indicates that this was a link label button event, that is, the label on the button was MENU, and that the linked page loaded was menu.aml.

```
1626712191,LL,SHUTDOWN,../shutdown.aml,
```

The operator pressed the SHUTDOWN menu button which then loaded the ../shutdown.aml page.

```
1626712200,LL,SHUTDOWN ,../shutdown.aml,
```

The operator touched the SHUTDOWN command button which then executed the shutdown procedure.

```
1626712200,SF,shmem_dump,checksum,6A7E48E4
```

The instrument saved the variables to the state file. Class SF indicates this is a state file operation.

```
1626712200,CE,shutdown,,
```

The command executed (CE) was the shutdown command. The data fields #2 and #3 in this line are empty as shown by the two commas after the shutdown command.

**LOG STATUS BITS TABLE****Status Bit Representation (in hexadecimal)**

STATUS_Alarm	0x00001
STATUS_relay1	0x00002
STATUS_relay2	0x00004
STATUS_Ext_1_osc	0x00020
STATUS_Valve_1_open	0x00040
STATUS_N2_1_LossOfSensor	0x00200
STATUS_N2_1_SensorShorted	0x00400
STATUS_PowerUp	0x00800
STATUS_He_1_Burnout_Timeout	0x04000
STATUS_He_1_LossOfSensor	0x08000

# GLOSSARY

## ABBREVIATIONS AND ACRONYMS

Term	Meaning
AC; ac	Alternating Current; strictly, electrical <i>current</i> that periodically reverses direction. Typically used also to describe an electrical power source in terms of the <i>voltage</i> . For example, 240 VAC.
ASCII	American Standard Code for Information Interchange; numerical representation of characters such as 'a' or '@' or an action (such as line-feed); 'plain' raw text with no formatting such as tabs, bold or underscoring
BNC connector	A miniature quick connect/disconnect RF connector used for coaxial cable, featuring two bayonet lugs on the female connector.
CR or <CR>	Text Carriage-Return character
Cryogen	A substance for obtaining low temperatures. In the case of use with the Model 1700 instrument, a cryogen is a liquefied gas such as liquid nitrogen or liquid helium.
D-Sub	Term referring to the family of connectors containing an odd number of pins in two parallel rows with a 1-pin difference in pins-per-row (DB9, DB15, and DB25 are most common)
DB9	Type of electrical connector containing 9 pins arranged in two parallel rows of 4 pins and 5 pins each
DCE	Data Circuit-terminating Equipment - a device that sits between the Data Terminal Equipment (DTE) and a data transmission circuit.
DHCP	Dynamic Host Configuration Protocol; a computer networking protocol which dynamically distributes the IP address to networked devices
dt	Rate of change
DTE	Data Terminal Equipment - an end instrument that converts user information into signals or reconverts received signals. A DTE device communicates with the Data Circuit-terminating Equipment (DCE).

Term	Meaning
ECL	Electrical Connection Lubricant - also known as Dielectric Tune-up Grease, a protective lubricant that prevents corrosion.
$E_o$	Power supply output voltage
i, I	Electrical current flow
$I_o$	Power supply output current
IEC	International Electrotechnical Commission
IEEE	Institute of Electrical and Electronics Engineers
I/O	Input/Output; The hardware and associated protocol that implement communication between information processing systems and/or devices. Inputs are the signals or data received by the system or device, and outputs are the signals or data sent from it.
IP	Internet Protocol; when used with “address”, refers to a numerical Internet address
kG	kilogauss: a magnetic field unit of measurement
LED	Light-Emitting Diode; a semiconductor device that emits light when energized - used for visual status indication
LF or <LF>	Text Linefeed character
LHe	Liquid Helium
LN2	Liquid Nitrogen
Max	Maximum
Min	Minimum
MSDS	Material Safety Data Sheet - provides workers and emergency personnel with procedures for handling or working with a specific substance in a safe manner and includes information such as physical data, toxicity, health effects, first aid, reactivity, storage, disposal, protective equipment, and spill-handling procedures.
RG-59/U	A specific type of coaxial cable, often used for low-power video and RF signal connections, with a characteristic impedance of 75 ohms.
$R_{lead}$	Electrical circuit lead or wiring resistance
RS-232	RS-232 is a long-established standard and protocol for relatively low speed serial data communication between computers and related devices; originally established for teletypewriter communication.
SCPI	Standard Commands for Programmable Instruments
V	Volts
VA	Volt-amperes ( $V \times I$ ); a unit of electrical reactive power
$V_{lead}$	Voltage ( $I \times R$ ) developed across circuit lead or wiring resistance due to current flow

Term	Meaning
$V_m$	Magnet voltage
$V_s$	Power supply voltage





**A**

abbreviations and acronyms 109

abnormal operation 42

Alarm Status screen 31

alarms 31

- alarm status screen 31
- configure source and level 33
- indication and acknowledgement 33
- multiple 32
- muting 31, 34
- remote configuration 77
- states 33

AMI contacts 100

analog outputs

- connection 18
- connector 98
- current loop diagram 19
- current loop output 18
- recorder output 18
- remote configuration 80
- simultaneous usage 18
- source selection 39
- specifications 8

approximate calibration 62–64

autofill

- AUTO-OFF 37
- AUTO-ON 37
- home screen indicator 22
- loss of sensor response 17
- manual overrides 37
- M-CLOSED 37
- M-OPEN 37
- output rating 9
- parameters 36
- relief valve caution 16
- remote configuration 79
- socket 9
- source selection 36
- standard 16
- start and stop levels 36
- states 37
- system diagram 15
- system setup 15
- timeout 9, 32, 36
- timeout reset 37

AUTO-OFF 37

AUTO-ON 37

**B**

battery backup 9

battery replacement 88

**C**

cables

- capacitance sensor coax 13
- long-distance 13
- serial 70

calibration

- approximate 48, 62–64
  - approx cal equation 63
- choosing a method 49
- closed dewar 57–61
  - completing closed dewar cal 59
  - increasing the range 60
  - preset MIN/MAX 57
- diagram 50
- dielectric constants 99
- liquid level sensors 47
- loss of sensor detection 53
- open dewar 55–56
- presetting MAX/MIN 57
- remote LN2 calibration 80
- stabilizing level during cal 59
- substitute liquid or partial length 62
- touch screen 102
- using analog outputs 59

capacitance sensor

- active length 47
- approximate calibration 62–64
- avoiding icing 12
- calibration diagram 50
- calibration methods 47
- cleaning 43
- closed dewar calibration 57–61
- coax cable specs 13
- configuration 30
- connecting 13
- dielectric variations 48
- ECL coating 12
- exposure to humidity 43
- installation 12
- loss of sensor 53, 100
- measurement method 48
- open dewar calibration 55–56
- sensor length 52
- sensor name 52
- shorted condition 100
- top vent hole 12
- vent holes 47

choosing units 38

closed dewar calibration 57, 57–61

- presetting MIN/MAX 57

command conventions 65

command summary 65–69

commands - see remote commands

communications

- commands - see remote commands
- error codes 82
- Ethernet 9
- Ethernet connector 9
- protocols 9
- remote, browser-based 4
- response format and termination 65
- serial 9
- serial connector 95

## configuration

- alarm source and level 33
- alarm state 33
- alarms and relays 31
- analog outputs 39
- autofill 36
- capacitance xiii
- custom xiii
- damping 38
- DST 46
- Ethernet 40
- external oscillator 1, 30
- helium xiii
- instrument options xii
- internal oscillator 1, 30
- level-based alarms 31
- line cord xii
- mounting xii
- muting alarms 34
- oscillator 51
- relay actuation 35
- relays source and level 34
- reset to defaults 43
- set date/time 45–46
- shutdown 44
- timezone 46

## connections

- capacitance sensor 13
- loss of sensor 53

## connector

- analog outputs 98
- Ethernet 71, 96
- serial port 70

## connectors

- RJ-45 - see Ethernet

## contact AMI 100

## cryogenic liquids viii

- first aid, first aid ix
- safety viii

## current loop

- connection diagram 19

## current loop specifications 8

## D

## dielectric constants 99

## dielectric effects 48

## dimensions 9

## display

- ALARM indicator/button 22
- autofill control 37
- AUTOFILL indicator/button 22
- BACK button 23
- field editing footer 22
- fill timeout indicators 32
- HELP button 22
- home screen 21
- home screen footer 22
- LN2 level 30
- MENU button 22, 23
- multiple alarm conditions 32
- remote configuration 75
- SAVE/CANCEL buttons 22
- size and resolution 3, 8
- units 38

## DST support 46

## E

## ECL use 12, 13

## environmental limits 10

## equipment

- recommended for safety x

## Ethernet

- configuration 40
- connector 9, 71
- connector pin out 96
- IP addressing 9
- ip parameters 71
- remote commands 72
- remote port 7180 71
- static ip parameters 40
- termination characters 71

## Ethernet connector pin out 96

## F

## fill timeout 36

- disable 36

## fill timeout for LN2 32

- reset 37

## footer

- editing 22
- home screen 22

## front panel layout 6

- power switch 21

## fuse replacement 90

**G**

glossary of terms 109  
DHCP 109

**H**

home screen 21  
  AUTOFILL button 37  
  AUTO-OFF button 37  
  AUTO-ON button 37  
  footer 22  
  M-CLOSED button 37  
  menu button 23  
  M-OPEN button 27  
  show LN2 level 30  
  units selection 38

**I**

ice formation 12

installation

  autofill system 15  
  capacitance sensors 12  
  mounting 11  
  power 18  
  unpacking 11

installing

  earth ground xi

instrument

  battery replacement 88  
  cleaning 87  
  fuse replacement 90  
  part number definition xii  
  replacement parts 87  
  troubleshooting 100

instrument description 1

ip address 71

**L**

level damping 39

level units 38

liquid level system terminology 109

logging 103

loss of sensor

  autofill 17  
  LN2 troubleshooting 100

**M**

M-CLOSED 37

measurement

  remote query 79

menu tree

  overview 24

menus 23

  descriptions 25  
  overall structure 24

method of measurement 48

moisture protection 12, 13, 43

M-OPEN 37

mounting

  coaxial cable 12  
  external oscillator 12  
  removing bottom feet 11  
  top vent hole 12

muting alarms 34

**N**

network

  address 40  
  configuration 40  
  connector 71  
  DHCP 71  
  DHCP definition 109  
  gateway 40  
  hostname 40  
  netmask 40  
  port 7180 71  
  static or DHCP mode 40

no sensor calibration 54

**O**

open dewar calibration 55–56

operation warnings x, 11, 18, 88, 90

**P**

part number definition xii

power on/off 21

power requirements 9, 18

power supply

  operating parameters 8

**R**

rear panel layout 7

recorder output specifications 8

relays 31

  actuation states 35  
  configure source and level 34  
  ratings 9  
  remote configuration 75

- type 9
- remote commands
  - alarms configuration 77
  - analog output assignment 80
  - display configuration 75
  - error codes 82
  - fill control 79
  - legacy commands 83
  - measurement 79
  - N2 calibration 80
  - relays configuration 75
  - system related 72
  - system-related commands 66
  - units 81
- remote control
  - browser-based 4
- remote error codes 82
- remote interface reference
  - command summary 68, 69
- remote interface reference - see remote commands
- reset to factory defaults 43
- RJ-45 connector - see Ethernet
- RS-232 configuration
  - parameters 70
- RS-232 connector 95
- RS-232 setup 41
- S**
- safety
  - cryogenic liquids viii
  - equipment x
  - legend x
  - relief valve 16
  - warnings xi
- screen descriptions
  - DAMPING 25
  - HISTORY 29
  - NETWORK 27
  - OUTPUTS
    - ALARMS 26
    - ANALOG OUTPUTS 27
    - AUTOFILL 27
    - RELAYS 26
  - SENSORS 25
    - CAL N2 25
    - SENSOR NAME 26
  - SHUTDOWN 29
  - SYSTEM 28
    - RS232 SETUP 28
    - SET TIME 28
  - SYSTST 29
  - UPDATE 28
- sensor
  - active length 47
  - installation 12
  - shorted LN2 troubleshooting 100
- sensor physical parameters
  - explanation 47
- serial communications
  - baud rate 9, 42
  - commands 72
  - connector 9
  - connector pin out 95
  - factory reset 74
  - interactive communication 65
  - set date 74
  - set time 74
  - system reboot 74
  - terminators 70
- serial port
  - baud rate 42
  - echo 42
  - function 42
  - setup 41
  - termination 42, 70
- serial port connector/cables 70
- service
  - AMI contact 100
- service - see instrument
- shutdown 44
- specifications
  - analog output 8
  - autofill 9
  - display 8
  - environmental limits 10
  - level measurements 8
  - physical 9
  - power 9
  - relays 9
  - standards conformance 10
- system
  - date and time 45–46
  - logging 103
  - remote configuration 72
  - shutdown 44
- system features 1
- system specifications 8
- system test
  - analog outputs 103
  - hardware controls 102
  - screen pixels 101
  - touch screen cal 102

**T**

TIMEZONE selection 46

troubleshooting common issues 100  
    system test screens 100

**U**

units 38

    remote configuration 81

user replaceable parts 87

**V**

vent holes 47

**W**

warnings

    equipment xi

weight 9

